

PARS-SEQ  
 SHAPE-SEQ  
 MDA  
 MALBAC  
 CHIRP-SEQ  
 RNA-SEQ  
 MAINE-SEQ  
 RRBS-SEQ  
 ICE  
 OS-SEQ  
 RIBO-SEQ  
 CLASH-SEQ  
 TIF-SEQ/PEAT  
 IN-SEQ  
 TAB-SEQ  
 CHIP-SEQ  
 5-C  
 SMMIP  
 OXBS-SEQ  
 RIP-SEQ  
 4-C  
 TC-SEQ  
 NET-SEQ  
 UMI  
 CAP-SEQ  
 FAIRE-SEQ  
 DUPLEX-SEQ  
 DNASE-SEQ  
 PAR-CLIP-SEQ  
 BS-SEQ  
 MEDIP-SEQ  
 GRO-SEQ  
 MERIP-SEQ  
 DIGITAL  
 CHIA-PET  
 ATAC-SEQ  
 CIP-TAP  
 ICLIP  
 MBDCAP-SEQ  
 HITS-CLIP  
 PARE-SEQ/GMUT  
 HI-C/3-C  
 TRAP-SEQ  
 RC-SEQ  
 FRAG-SEQ

# Sequencing Methods Review

A review of publications featuring Illumina® Technology

## TABLE OF CONTENTS

<b>Table of Contents</b>	<b>2</b>
<b>Introduction</b>	<b>4</b>
<b>RNA Transcription</b>	<b>5</b>
Chromatin Isolation by RNA Purification (ChIRP-Seq)	7
Global Run-on Sequencing (GRO-Seq)	9
Ribosome Profiling Sequencing (Ribo-Seq)/ARTseq™	12
RNA Immunoprecipitation Sequencing (RIP-Seq)	15
High-Throughput Sequencing of CLIP cDNA library (HITS-CLIP) or	17
Crosslinking and Immunoprecipitation Sequencing (CLIP-Seq)	17
Photoactivatable Ribonucleoside-Enhanced Crosslinking and Immunoprecipitation (PAR-CLIP)	19
Individual Nucleotide Resolution CLIP (iCLIP)	22
Native Elongating Transcript Sequencing (NET-Seq)	24
Targeted Purification of Polysomal mRNA (TRAP-Seq)	25
Crosslinking, Ligation, and Sequencing of Hybrids (CLASH-Seq)	26
Parallel Analysis of RNA Ends Sequencing (PARE-Seq) or	27
Genome-Wide Mapping of Uncapped Transcripts (GMUCT)	27
Transcript Isoform Sequencing (TIF-Seq) or	29
Paired-End Analysis of TSSs (PEAT)	29
<b>RNA Structure</b>	<b>30</b>
Selective 2'-Hydroxyl Acylation Analyzed by Primer Extension Sequencing (SHAPE-Seq)	31
Parallel Analysis of RNA Structure (PARS-Seq)	32
Fragmentation Sequencing (FRAG-Seq)	33
CXXC Affinity Purification Sequencing (CAP-Seq)	34
Alkaline Phosphatase, Calf Intestine-Tobacco Acid Pyrophosphatase Sequencing (CIP-TAP)	36
Inosine Chemical Erasing Sequencing (ICE)	38
m6A-Specific Methylated RNA Immunoprecipitation Sequencing (MeRIP-Seq)	39
<b>Low-Level RNA Detection</b>	<b>40</b>
Digital RNA Sequencing	42
Whole-Transcript Amplification for Single Cells (Quartz-Seq)	43
Designed Primer-Based RNA Sequencing (DP-Seq)	44
Switch Mechanism at the 5' End of RNA Templates (Smart-Seq)	45
Switch Mechanism at the 5' End of RNA Templates Version 2 (Smart-Seq2)	47
Unique Molecular Identifiers (UMI)	49
Cell Expression by Linear Amplification Sequencing (CEL-Seq)	51
Single-Cell Tagged Reverse Transcription Sequencing (STRT-Seq)	52
<b>Low-Level DNA Detection</b>	<b>53</b>
Single-Molecule Molecular Inversion Probes (smMIP)	55
Multiple Displacement Amplification (MDA)	56
Multiple Annealing and Looping-Based Amplification Cycles (MALBAC)	59
Oligonucleotide-Selective Sequencing (OS-Seq)	61
Duplex Sequencing (Duplex-Seq)	62
<b>DNA Methylation</b>	<b>63</b>
Bisulfite Sequencing (BS-Seq)	65
Post-Bisulfite Adapter Tagging (PBAT)	70
Tagmentation-Based Whole Genome Bisulfite Sequencing (T-WGBS)	72
Oxidative Bisulfite Sequencing (oxBS-Seq)	73

Tet-Assisted Bisulfite Sequencing (TAB-Seq)	74
Methylated DNA Immunoprecipitation Sequencing (MeDIP-Seq)	76
Methylation-Capture (MethylCap) Sequencing or	79
Methyl-Binding-Domain-Capture (MBDCap) Sequencing	79
Reduced-Representation Bisulfite Sequencing (RRBS-Seq)	81
<b>DNA-Protein Interactions</b>	<b>83</b>
DNase I Hypersensitive Sites Sequencing (DNase-Seq)	85
MNase-Assisted Isolation of Nucleosomes Sequencing (MAINE-Seq)	88
Chromatin Immunoprecipitation Sequencing (ChIP-Seq)	91
Formaldehyde-Assisted Isolation of Regulatory Elements (FAIRE-Seq)	94
Assay for Transposase-Accessible Chromatin Sequencing (ATAC-Seq)	96
Chromatin Interaction Analysis by Paired-End Tag Sequencing (ChIA-PET)	97
Chromatin Conformation Capture (Hi-C/3C-Seq)	99
Circular Chromatin Conformation Capture (4-C or 4C-Seq)	101
Chromatin Conformation Capture Carbon Copy (5-C)	104
<b>Sequence Rearrangements</b>	<b>105</b>
Retrotransposon Capture Sequencing (RC-Seq)	107
Transposon Sequencing (Tn-Seq) or Insertion Sequencing (INSeq)	109
Translocation-Capture Sequencing (TC-Seq)	111
<b>Bibliography</b>	<b>113</b>
<b>Appendix</b>	<b>131</b>
DNA/RNA Purification Kits	131
DNA Sequencing	133
RNA Sequencing	168
Arrays	190
PCR and Enzyme Solutions	194
Instruments	198

## INTRODUCTION

This collection of next-generation sequencing (NGS) sample preparation protocols was compiled from the scientific literature to demonstrate the wide range of scientific questions that can be addressed by Illumina's sequencing by synthesis technology. It is both a tribute to the creativity of the users and the versatility of the technology. We hope it will inspire researchers to use these methods or to develop new ones to address new scientific challenges.

These methods were developed by users, so readers should refer to the original publications for detailed descriptions and protocols.

Have we missed anything? Please contact us if you are aware of a protocol that should be listed.



## RNA TRANSCRIPTION

The regulation of RNA transcription and processing directly affects protein synthesis. Proteins, in turn, mediate cellular functions to establish the phenotype of the cell. Dysregulated RNAs are the cause for some diseases and cancers<sup>1,2</sup>. Sequencing RNA provides information about both the abundance and sequence of the RNA molecules. Careful analysis of the results, along with adaptation of the sample preparation protocols, can provide remarkable insight into all the various aspects of RNA processing and control of transcription. Examples of these measures include: post-translational modifications, RNA splicing, RNA bound to RNA binding proteins (RBP), RNA expressed at various stages, unique RNA isoforms, RNA degradation, and regulation of other RNA species<sup>3,4</sup>. Studies of RNA transcription and translation are leading to a better understanding of the implications of RNA production, processing, and regulation for cellular phenotype.



Scientists have discovered a link between long term memory and protein synthesis in brain<sup>5,6</sup>.

- 1 Kloosterman W. P. and Plasterk R. H. (2006) The diverse functions of microRNAs in animal development and disease. *Dev Cell* 11: 441-450
- 2 Castello A., Fischer B., Hentze M. W. and Preiss T. (2013) RNA-binding proteins in Mendelian disease. *Trends Genet* 29: 318-327
- 3 McGettigan P. A. (2013) Transcriptomics in the RNA-seq era. *Curr Opin Chem Biol* 17: 4-11
- 4 Feng H., Qin Z. and Zhang X. (2013) Opportunities and methods for studying alternative splicing in cancer with RNA-Seq. *Cancer Lett* 340: 179-191
- 5 Davis H. P. and Squire L. R. (1984) Protein synthesis and memory: a review. *Psychol Bull* 96: 518-559
- 6 Holt C. E. and Schuman E. M. (2013) The central dogma decentralized: new perspectives on RNA function and local translation in neurons. *Neuron* 80: 648-657

## Reviews

Castello A., Fischer B., Hentze M. W. and Preiss T. (2013) RNA-binding proteins in Mendelian disease. *Trends Genet* 29: 318-327

Feng H., Qin Z. and Zhang X. (2013) Opportunities and methods for studying alternative splicing in cancer with RNA-Seq. *Cancer Lett* 340: 179-191

Holt C. E. and Schuman E. M. (2013) The central dogma decentralized: new perspectives on RNA function and local translation in neurons. *Neuron* 80: 648-657

Law G. L., Korth M. J., Benecke A. G. and Katze M. G. (2013) Systems virology: host-directed approaches [to viral pathogenesis and drug targeting]. *Nat Rev Microbiol* 11: 455-466

Licatalosi D. D. and Darnell R. B. (2010) RNA processing and its regulation: global insights into biological networks. *Nat Rev Genet* 11: 75-87

## CHROMATIN ISOLATION BY RNA PURIFICATION (CHIRP-SEQ)

Chromatin isolation by RNA purification (ChIRP-Seq) is a protocol to detect the locations on the genome where non-coding RNAs (ncRNAs), such as long non-coding RNAs (lncRNAs), and their proteins are bound<sup>7</sup>. In this method, samples are first crosslinked and sonicated. Biotinylated tiling oligos are hybridized to the RNAs of interest, and the complexes are captured with streptavidin magnetic beads. After treatment with RNase H the DNA is extracted and sequenced. With deep sequencing the lncRNA/protein interaction site can be determined at single-base resolution.



### Pros

- Binding sites can be found anywhere on the genome
- New binding sites can be discovered
- Specific RNAs of interest can be selected

### Cons

- Nonspecific oligo interactions can lead to misinterpretation of binding sites
- Chromatin can be disrupted during the preparation stage
- The sequence of the RNA of interest must be known

## References

Li Z., Chao T. C., Chang K. Y., Lin N., Patil V. S., et al. (2014) The long noncoding RNA THRIL regulates TNFalpha expression through its interaction with hnRNPL. *Proc Natl Acad Sci U S A* 111: 1002-1007

The non-protein-coding parts of the mammalian genome encode thousands of large intergenic non-coding RNAs (lincRNAs). To identify lincRNAs associated with activation of the innate immune response, this study applied custom microarrays and Illumina RNA sequencing for THP1 macrophages. A panel of 159 lincRNAs was found to be differentially expressed following innate activation. Further analysis of the RNA-Seq data revealed that linc1992 was required for expression of many immune-response genes, including cytokines and regulators of TNF-alpha expression.

**Illumina Technology:** HiSeq 2000®

<sup>7</sup> Chu C., Qu K., Zhong F. L., Artandi S. E. and Chang H. Y. (2011) Genomic maps of long noncoding RNA occupancy reveal principles of RNA-chromatin interactions. *Mol Cell* 44: 667-678

Li W., Notani D., Ma Q., Tanasa B., Nunez E., et al. (2013) Functional roles of enhancer RNAs for oestrogen-dependent transcriptional activation. *Nature* 498: 516-520

Enhancers are regions of DNA with regulatory function. Through binding of transcription factors and cis-interactions with promoters, target gene expression may be increased. In addition, both lncRNAs and bidirectional ncRNAs may be transcribed on enhancers and are referred to as enhancer RNAs (eRNAs). This study examined eRNA expression in breast cancer cells using a combination of sequencing protocols on HiSeq 2000 (ChIRP-seq, GRO-seq, ChIP-Seq, 3C, 3D-DSL) to discover a global increase in eRNA transcription on enhancers adjacent to E2-upregulated coding genes. These data suggest that eRNAs may play an important role in transcriptional regulation.

**Illumina Technology:** HiSeq 2000

Chu C., Qu K., Zhong F. L., Artandi S. E. and Chang H. Y. (2011) Genomic maps of long noncoding RNA occupancy reveal principles of RNA-chromatin interactions. *Mol Cell* 44: 667-678

#### **Associated Kits**


ScriptSeq™ Complete Kit

TruSeq® RNA Sample Prep Kit

TruSeq® Small RNA Sample Prep Kit

## GLOBAL RUN-ON SEQUENCING (GRO-SEQ)

Global run-on sequencing (GRO-Seq) maps binding sites of transcriptionally active RNA polymerase II. In this method, active RNA polymerase II<sup>8</sup> is allowed to run on in the presence of Br-UTP. RNAs are hydrolyzed and purified using beads coated with Brd-UTP antibody. The eluted RNA undergoes cap removal and end repair prior to reverse transcription to cDNA. Deep sequencing of the cDNA provides sequences of RNAs that are actively transcribed by RNA polymerase II.



Pros	Cons
<ul style="list-style-type: none"> <li>• Maps position of transcriptionally-engaged RNA polymerases</li> <li>• Determines relative activity of transcription sites</li> <li>• Detects sense and antisense transcription</li> <li>• Detects transcription anywhere on the genome</li> <li>• No prior knowledge of transcription sites is needed</li> </ul>	<ul style="list-style-type: none"> <li>• The protocol is limited to cell cultures and other artificial systems due to the requirement for incubation in the presence of labeled nucleotides</li> <li>• Artifacts may be introduced during the preparation of the nuclei<sup>9</sup></li> <li>• New initiation events may occur during the run-on step</li> <li>• Physical impediments may block the polymerases</li> </ul>

## References

Heinz S., Romanoski C. E., Benner C., Allison K. A., Kaikkonen M. U., et al. (2013) Effect of natural genetic variation on enhancer selection and function. *Nature* 503: 487-492

Previous work in epigenetics has proposed a model where lineage-determining transcription factors (LDTF) collaboratively compete with nucleosomes to bind DNA in a cell type-specific manner. In order to determine the sequence variants that guide transcription factor binding, the authors of this paper tested this model in vivo by comparing the SNPs that disrupted transcription factor binding sites in two inbred mouse strains. The authors used GRO-seq in combination with ChIP-seq and RNA-Seq to determine expression and transcription factor binding. The SNPs of the two strains were then classified based on their ability to perturb transcription factor binding and the authors found substantial evidence to support the model.

**Illumina Technology:** TruSeq RNA Sample Prep Kit, HiSeq 2000

<sup>8</sup> Core L. J., Waterfall J. J. and Lis J. T. (2008) Nascent RNA sequencing reveals widespread pausing and divergent initiation at human promoters. *Science* 322: 1845-1848

<sup>9</sup> Adelman K. and Lis J. T. (2012) Promoter-proximal pausing of RNA polymerase II: emerging roles in metazoans. *Nat Rev Genet* 13: 720-731

Jin F., Li Y., Dixon J. R., Selvaraj S., Ye Z., et al. (2013) A high-resolution map of the three-dimensional chromatin interactome in human cells. *Nature* 503: 290-294

Cis-acting regulatory elements in the genome interact with their target gene promoter by transcription factors, bringing the two locations close together in the 3D conformation of the chromatin. In this study the chromosome conformation is examined by a genome-wide analysis method (Hi-C) using the Illumina HiSeq 2000 system. The authors determined over one million long-range chromatin interactions in human fibroblasts. In addition, they characterized the dynamics of promoter-enhancer contacts after TNF-alpha signaling and discovered pre-existing chromatin looping with the TNF-alpha-responsive enhancers, suggesting the three-dimensional chromatin conformation may be stable over time.

**Illumina Technology:** HiSeq 2000

Kaikkonen M. U., Spann N. J., Heinz S., Romanoski C. E., Allison K. A., et al. (2013) Remodeling of the enhancer landscape during macrophage activation is coupled to enhancer transcription. *Mol Cell* 51: 310-325

Enhancers have been shown to specifically bind lineage-determining transcription factors in a cell-type-specific manner. Toll-like receptor 4 (TLR4) signaling primarily regulates macrophage gene expression through a pre-existing enhancer landscape. In this study the authors used GRO-seq and ChIP-seq to discover that enhancer transcription precedes local mono- and dimethylation of histone H3 lysine 4 (H3K4).

**Illumina Technology:** Genome Analyzer<sup>ix</sup>

Kim Y. J., Greer C. B., Cecchini K. R., Harris L. N., Tuck D. P., et al. (2013) HDAC inhibitors induce transcriptional repression of high copy number genes in breast cancer through elongation blockade. *Oncogene* 32: 2828-2835

Histone deacetylase inhibitors (HDACI) are a promising class of cancer-repressing drugs. This study investigated the molecular mechanism of HDACI by using GRO-seq in combination with expression analysis. The authors show that HDACI preferentially represses transcription of highly expressed genes which, in cancers, are typically misregulated oncogenes supporting further development of HDACI as a general cancer inhibitor.

**Illumina Technology:** Genome Analyzer<sup>ix</sup>, Human Gene Expression—BeadArray; 35 bp reads

Li W., Notani D., Ma Q., Tanasa B., Nunez E., et al. (2013) Functional roles of enhancer RNAs for oestrogen-dependent transcriptional activation. *Nature* 498: 516-520

Enhancers are regions of DNA with regulatory function. Through binding of transcription factors and cis-interactions with promoters, target gene expression may be increased. In addition, both lncRNAs and bidirectional ncRNAs may be transcribed on enhancers and are referred to as enhancer RNAs (eRNAs). This study examined eRNA expression in breast cancer cells using a combination of sequencing protocols on HiSeq 2000 (ChIRP-seq, GRO-seq, ChIP-Seq, 3C, 3D-DSL) to discover a global increase in eRNA transcription on enhancers adjacent to E2-upregulated coding genes. These data suggest that eRNAs may play an important role in transcriptional regulation.

**Illumina Technology:** HiSeq 2000

Saunders A., Core L. J., Sutcliffe C., Lis J. T. and Ashe H. L. (2013) Extensive polymerase pausing during *Drosophila* axis patterning enables high-level and pliable transcription. *Genes Dev* 27: 1146-1158

*Drosophila* embryogenesis has been intensively studied for the expression patterns of genes corresponding to differentiation of embryonal tissue. In this study, gene regulation was examined using GRO-seq to map the details of RNA polymerase distribution over the genome during early embryogenesis. The authors found that certain groups of genes were more highly paused than others, and that bone morphogenetic protein (BMP) target gene expression requires the pause-inducing negative elongation factor complex (NELF).

**Illumina Technology:** Genome Analyzer<sub>lix</sub>

Ji X., Zhou Y., Pandit S., Huang J., Li H., et al. (2013) SR proteins collaborate with 7SK and promoter-associated nascent RNA to release paused polymerase. *Cell* 153: 855-868

Lam M. T., Cho H., Lesch H. P., Gosselin D., Heinz S., et al. (2013) Rev-Erbs repress macrophage gene expression by inhibiting enhancer-directed transcription. *Nature* 498: 511-515

Li P., Spann N. J., Kaikkonen M. U., Lu M., Oh da Y., et al. (2013) NCoR repression of LXRs restricts macrophage biosynthesis of insulin-sensitizing omega 3 fatty acids. *Cell* 155: 200-214

Chopra V. S., Hendrix D. A., Core L. J., Tsui C., Lis J. T., et al. (2011) The Polycomb Group Mutant *esc* Leads to Augmented Levels of Paused Pol II in the *Drosophila* Embryo. *Mol Cell* 42: 837-844

Hah N., Danko C. G., Core L., Waterfall J. J., Siepel A., et al. (2011) A rapid, extensive, and transient transcriptional response to estrogen signaling in breast cancer cells. *Cell* 145: 622-634

Larschan E., Bishop E. P., Kharchenko P. V., Core L. J., Lis J. T., et al. (2011) X chromosome dosage compensation via enhanced transcriptional elongation in *Drosophila*. *Nature* 471: 115-118

Wang D., Garcia-Bassets I., Benner C., Li W., Su X., et al. (2011) Reprogramming transcription by distinct classes of enhancers functionally defined by eRNA. *Nature* 474: 390-394

Core L. J., Waterfall J. J. and Lis J. T. (2008) Nascent RNA sequencing reveals widespread pausing and divergent initiation at human promoters. *Science* 322: 1845-1848

#### Associated Kits

ScriptSeq™ Complete Kit

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA® and Total RNA® Sample Preparation Kit

TruSeq Targeted RNA® Expression Kit

## RIBOSOME PROFILING SEQUENCING (RIBO-SEQ)/ARTSEQ™

Active mRNA Translation Sequencing (ARTseq), also called ribosome profiling (Ribo-Seq), isolates RNA that is being processed by the ribosome in order to monitor the translation process<sup>10</sup>. In this method ribosome-bound RNA first undergoes digestion. The RNA is then extracted and the rRNA is depleted. Extracted RNA is reverse-transcribed to cDNA. Deep sequencing of the cDNA provides the sequences of RNAs bound by ribosomes during translation. This method has been refined to improve the quality and quantitative nature of the results. Careful attention should be paid to: (1) generation of cell extracts in which ribosomes have been faithfully halted along the mRNA they are translating in vivo; (2) nuclease digestion of RNAs that are not protected by the ribosome followed by recovery of the ribosome-protected mRNA fragments; (3) quantitative conversion of the protected RNA fragments into a DNA library that can be analyzed by deep sequencing<sup>11</sup>. The addition of harringtonine (an alkaloid that inhibits protein biosynthesis) causes ribosomes to accumulate precisely at initiation codons and assists in their detection.



Pros	Cons
<ul style="list-style-type: none"> <li>• Reveals a snapshot with the precise location of ribosomes on the RNA</li> <li>• Ribosome profiling more closely reflects the rate of protein synthesis than mRNA levels</li> <li>• No prior knowledge of the RNA or ORFs is required</li> <li>• The whole genome is surveyed</li> <li>• Can be used to identify protein-coding regions</li> </ul>	<ul style="list-style-type: none"> <li>• Initiation from multiple sites within a single transcript makes it challenging to define all ORFs</li> <li>• Does not provide the kinetics of translational elongation</li> </ul>

### References

Becker A. H., Oh E., Weissman J. S., Kramer G. and Bukau B. (2013) Selective ribosome profiling as a tool for studying the interaction of chaperones and targeting factors with nascent polypeptide chains and ribosomes. *Nat Protoc* 8: 2212-2239

A plethora of factors is involved in the maturation of newly synthesized proteins, including chaperones, membrane targeting factors, and enzymes. This paper presents an assay for selective ribosome profiling (SeRP) to determine the interaction of factors with ribosome-nascent chain complexes (RNCs). The protocol is based on Illumina sequencing of ribosome-bound mRNA fragments combined with selection for RNCs associated with the factor of interest.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>10</sup> Ingolia N. T., Ghaemmaghani S., Newman J. R. and Weissman J. S. (2009) Genome-wide analysis in vivo of translation with nucleotide resolution using ribosome profiling. *Science* 324: 218-223

<sup>11</sup> Ingolia N. T., Lareau L. F. and Weissman J. S. (2011) Ribosome Profiling of Mouse Embryonic Stem Cells Reveals the Complexity and Dynamics of Mammalian Proteomes. *Cell* 147: 789-802



Lee M. T., Bonneau A. R., Takacs C. M., Bazzini A. A., DiVito K. R., et al. (2013) Nanog, Pou5f1 and SoxB1 activate zygotic gene expression during the maternal-to-zygotic transition. *Nature* 503: 360-364

In the developmental transition from egg to zygote, the fertilized egg must clear maternal mRNAs and initiate the zygote development program—the zygotic genome activation (ZGA). In this paper, the ZGA was studied in zebrafish using Illumina sequencing to determine the factors that activate the zygotic program. Using a combination of ribosome profiling and mRNA sequencing, the authors identified several hundred genes directly activated by maternal factors, constituting the first wave of zygotic transcription.

**Illumina Technology:** HiSeq 2000/2500

Stumpf C. R., Moreno M. V., Olshen A. B., Taylor B. S. and Ruggero D. (2013) The translational landscape of the Mammalian cell cycle. *Mol Cell* 52: 574-582

The regulation of gene expression accounts for the differences seen between different cell types and tissues that share the same genomic information. Regulation may vary over time, and the mechanism and extent is still poorly understood. This study applied Illumina HiSeq technology to sequence total mRNA and total ribosome-occupied mRNA throughout the cell cycle of synchronized HeLa cells to study the translational regulation by ribosome occupancy. The authors identified a large number of mRNAs that undergo significant changes in translation between phases of the cell cycle, and they found 112 mRNAs that were translationally regulated exclusively between specific phases of the cell cycle. The authors suggest translational regulation is a particularly well-suited mechanism for controlling dynamic processes, such as the cell cycle.

**Illumina Technology:** HiSeq 2000/2500

Wang T., Cui Y., Jin J., Guo J., Wang G., et al. (2013) Translating mRNAs strongly correlate to proteins in a multivariate manner and their translation ratios are phenotype specific. *Nucleic Acids Res* 41: 4743-4754

It is well known that the abundance of total mRNAs correlates poorly to protein levels. This study set out to analyze the relative abundances of mRNAs, ribosome-nascent chain complex (RNC)-mRNAs, and proteins on a genome-wide scale. A human lung cancer cell line and normal bronchial epithelial cells were analyzed with RNA-seq and the protein abundance measured. The authors created a multivariate linear model showing strong correlation of RNA and protein abundance by integrating the mRNA length as a key factor.

**Illumina Technology:** Genome Analyzer<sub>ix</sub> and HiSeq 2000

Liu B., Han Y. and Qian S. B. (2013) Cotranslational response to proteotoxic stress by elongation pausing of ribosomes. *Mol Cell* 49: 453-463

Liu X., Jiang H., Gu Z. and Roberts J. W. (2013) High-resolution view of bacteriophage lambda gene expression by ribosome profiling. *Proc Natl Acad Sci U S A* 110: 11928-11933

Cho J., Chang H., Kwon S. C., Kim B., Kim Y., et al. (2012) LIN28A is a suppressor of ER-associated translation in embryonic stem cells. *Cell* 151: 765-777

- Fritsch C., Herrmann A., Nothnagel M., Szafranski K., Huse K., et al. (2012) Genome-wide search for novel human uORFs and N-terminal protein extensions using ribosomal footprinting. *Genome Res* 22: 2208-2218
- Gerashchenko M. V., Lobanov A. V. and Gladyshev V. N. (2012) Genome-wide ribosome profiling reveals complex translational regulation in response to oxidative stress. *Proc Natl Acad Sci U S A* 109: 17394-17399
- Han Y., David A., Liu B., Magadan J. G., Binnik J. R., et al. (2012) Monitoring cotranslational protein folding in mammalian cells at codon resolution. *Proc Natl Acad Sci U S A* 109: 12467-12472
- Hsieh A. C., Liu Y., Edlind M. P., Ingolia N. T., Janes M. R., et al. (2012) The translational landscape of mTOR signalling steers cancer initiation and metastasis. *Nature* 485: 55-61
- Lee S., Liu B., Lee S., Huang S. X., Shen B., et al. (2012) Global mapping of translation initiation sites in mammalian cells at single-nucleotide resolution. *Proc Natl Acad Sci U S A* 109: E2424-2432
- Li G. W., Oh E. and Weissman J. S. (2012) The anti-Shine-Dalgarno sequence drives translational pausing and codon choice in bacteria. *Nature* 484: 538-541
- Stadler M., Artiles K., Pak J. and Fire A. (2012) Contributions of mRNA abundance, ribosome loading, and post- or peri-translational effects to temporal repression of *C. elegans* heterochronic miRNA targets. *Genome Res* 22: 2418-2426
- Darnell J. C., Van Driesche S. J., Zhang C., Hung K. Y., Mele A., et al. (2011) FMRP Stalls Ribosomal Translocation on mRNAs Linked to Synaptic Function and Autism. *Cell* 146: 247-261
- Ingolia N. T., Lareau L. F. and Weissman J. S. (2011) Ribosome Profiling of Mouse Embryonic Stem Cells Reveals the Complexity and Dynamics of Mammalian Proteomes. *Cell* 147: 789-802
- Oh E., Becker A. H., Sandikci A., Huber D., Chaba R., et al. (2011) Selective ribosome profiling reveals the cotranslational chaperone action of trigger factor in vivo. *Cell* 147: 1295-1308
- Han Y., David A., Liu B., Magadan J. G., Binnik J. R., et al. (2012) Monitoring cotranslational protein folding in mammalian cells at codon resolution. *Proc Natl Acad Sci U S A* 109: 12467-12472
- Ingolia N. T. (2010) Genome-wide translational profiling by ribosome footprinting. *Methods Enzymol* 470: 119-142


### Associated Kits

ARTseq™ Ribosome Profiling Kit

Ribo-Zero® Kit

## RNA IMMUNOPRECIPITATION SEQUENCING (RIP-SEQ)

RNA immunoprecipitation sequencing (RIP-Seq) maps the sites where proteins are bound to the RNA within RNA-protein complexes<sup>12</sup>. In this method, RNA-protein complexes are immunoprecipitated with antibodies targeted to the protein of interest. After RNase digestion, RNA covered by protein is extracted and reverse-transcribed to cDNA. The locations can then be mapped back to the genome. Deep sequencing of cDNA provides single-base resolution of bound RNA.



Pros	Cons
<ul style="list-style-type: none"> <li>• Maps specific protein-RNA complexes, such as polycomb-associated RNAs</li> <li>• Low background and higher resolution of binding site due to RNase digestion</li> <li>• No prior knowledge of the RNA is required</li> <li>• Genome-wide RNA screen</li> </ul>	<ul style="list-style-type: none"> <li>• Requires antibodies to the targeted proteins</li> <li>• Nonspecific antibodies will precipitate nonspecific complexes</li> <li>• Lack of crosslinking or stabilization of the complexes may lead to false negatives</li> <li>• RNase digestion must be carefully controlled</li> </ul>

### References

Kanematsu S., Tanimoto K., Suzuki Y. and Sugano S. (2014) Screening for possible miRNA-mRNA associations in a colon cancer cell line. *Gene* 533: 520-531

MicroRNAs (miRNAs) are small ncRNAs mediating the regulation of gene expression in various biological contexts, including carcinogenesis. This study examined the putative associations between miRNAs and mRNAs via Argonaute1 (Ago1) or Ago2 immunoprecipitation in a colon cancer cell line. The mRNA sequencing and RIP-seq was performed on an Illumina Genome Analyzer<sub>ix</sub> system. From this analysis the authors found specific associations of Ago1 with genes having constitutive cellular functions, whereas putative miRNA-mRNA associations detected with Ago2 IP appeared to be related to signal transduction genes.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Udan-Johns M., Bengoechea R., Bell S., Shao J., Diamond M. I., et al. (2014) Prion-like nuclear aggregation of TDP-43 during heat shock is regulated by HSP40/70 chaperones. *Hum Mol Genet* 23: 157-170

Aberrant aggregation of the protein TDP-43 is a key feature of the pathology of amyotrophic lateral sclerosis (ALS). Studying the mechanism of TDP-43 aggregation, this paper presents an analysis of gene expression and RNA-binding partners in human and mouse cell lines. The aggregation of TDP-43 was observed during heat shock and potential interaction partners were identified. The authors suggest TDP-43 shares properties with physiologic prions from yeast, requiring chaperone proteins for aggregation.

**Illumina Technology:** HiSeq 2000

<sup>12</sup> Zhao J., Ohsumi T. K., Kung J. T., Ogawa Y., Grau D. J., et al. (2010) Genome-wide identification of polycomb-associated RNAs by RIP-seq. *Mol Cell* 40: 939-953

Wang X., Lu Z., Gomez A., Hon G. C., Yue Y., et al. (2014) N<sup>6</sup>-methyladenosine-dependent regulation of messenger RNA stability. *Nature* 505: 117-120

N<sup>6</sup>-methyladenosine (m<sup>6</sup>A) is the most prevalent internal (non-cap) modification present in the messenger RNA of all higher eukaryotes. To understand the role of m<sup>6</sup>A modification in mammalian cells, the authors of this study applied Illumina sequencing to characterize the YTH domain family 2 (YTHDF2) reader protein regulation of mRNA degradation. The authors performed m<sup>6</sup>A-seq (MeRIP-Seq), RIP-seq, mRNA-Seq, photoactivatable-ribonucleoside-enhanced crosslinking and immunoprecipitation (PAR-CLIP), and ribosome profiling for HeLa cells on an Illumina HiSeq system with 100 bp single-end reads. They demonstrated that m<sup>6</sup>A is selectively recognized by YTHDF2, affecting the translation status and lifetime of mRNA.

**Illumina Technology:** HiSeq 2000; 100 bp single-end reads

Di Ruscio A., Ebralidze A. K., Benoukraf T., Amabile G., Goff L. A., et al. (2013) DNMT1-interacting RNAs block gene-specific DNA methylation. *Nature* 503: 371-376

DNA methylation is one of the many epigenetic factors that influence the regulation of gene expression. In this paper, the authors show that a novel RNA from the CEBPA gene locus is critical in regulating the local DNA methylation profile, and thus co-influences gene regulation. Using RIP-seq and RNA-Seq on Illumina platforms, the authors showed that this novel RNA binds DNA (cytosine-5)-methyltransferase 1 (DNMT1) and prevents methylation of the CEBPA gene locus.

**Illumina Technology:** Genome Analyzer<sub>15k</sub> and HiSeq 2000

Meyer K. D., Saletore Y., Zumbo P., Elemento O., Mason C. E., et al. (2012) Comprehensive analysis of mRNA methylation reveals enrichment in 3' UTRs and near stop codons. *Cell* 149: 1635-1646

Cernilogar F. M., Onorati M. C., Kothe G. O., Burroughs A. M., Parsi K. M., et al. (2011) Chromatin-associated RNA interference components contribute to transcriptional regulation in *Drosophila*. *Nature* 480: 391-395

Salton M., Elkon R., Borodina T., Davydov A., Yaspo M. L., et al. (2011) Matrin 3 binds and stabilizes mRNA. *PLoS One* 6: e23882

Zhao J., Ohsumi T. K., Kung J. T., Ogawa Y., Grau D. J., et al. (2010) Genome-wide identification of polycomb-associated RNAs by RIP-seq. *Mol Cell* 40: 939-953

#### Associated Kits

ARTseq™ Ribosome Profiling Kit

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit


TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

## HIGH-THROUGHPUT SEQUENCING OF CLIP CDNA LIBRARY (HITS-CLIP) OR CROSSLINKING AND IMMUNOPRECIPITATION SEQUENCING (CLIP-SEQ)

High-throughput sequencing of CLIP cDNA library (HITS-CLIP) or crosslinking and immunoprecipitation sequencing (CLIP-Seq) maps protein-RNA binding sites *in vivo*<sup>13</sup>. This approach is similar to RIP-Seq, but uses crosslinking to stabilize the protein-RNA complexes. In this method, RNA-protein complexes are UV crosslinked and immunoprecipitated. The protein-RNA complexes are treated with RNase followed by Proteinase K. RNA is extracted and reverse-transcribed to cDNA. Deep sequencing of cDNA provides single-base resolution mapping of protein binding to RNAs.



Pros	Cons
<ul style="list-style-type: none"> <li>• Crosslinking stabilizes the protein-target binding</li> <li>• UV crosslinking can be carried out <i>in vivo</i></li> <li>• Low background and higher resolution of binding site due to RNase digestion</li> <li>• No prior knowledge of the RNA is required</li> <li>• Genome-wide RNA screen</li> </ul>	<ul style="list-style-type: none"> <li>• Antibodies not specific to the target may precipitate nonspecific complexes</li> <li>• UV crosslinking is not very efficient and requires very close protein-RNA interactions</li> <li>• Artifacts may be introduced during the crosslinking process</li> </ul>

### References

Poulos M. G., Batra R., Li M., Yuan Y., Zhang C., et al. (2013) Progressive impairment of muscle regeneration in muscleblind-like 3 isoform knockout mice. *Hum Mol Genet* 22: 3547-3558

The human muscleblind-like (*MBNL*) genes encode alternative splicing factors essential for development of multiple tissues. In the neuromuscular disease myotonic dystrophy, C(C)UG repeats in RNA inhibit MBNL activity. This paper reports a study of the Mbnl3 protein isoform in a mouse model to determine the function of Mbnl3 in muscle regeneration and muscle function. The authors used an Illumina Genome Analyzer system for RNA-Seq and HITS-CLIP to determine Mbnl3-RNA interaction.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>13</sup> Chi SW, Zang JB, Mele A, Darnell RB; (2009) Argonaute HITS-CLIP decodes microRNA-mRNA interaction maps. *Nature* 460: 479-86

Xu D., Shen W., Guo R., Xue Y., Peng W., et al. (2013) Top3beta is an RNA topoisomerase that works with fragile X syndrome protein to promote synapse formation. *Nat Neurosci* 16: 1238-1247

Topoisomerases are crucial for solving DNA topological problems, but they have not previously been linked to RNA metabolism. In this study the human topoisomerase 3beta (Top3B), which is known to regulate the translation of mRNAs, was found to bind multiple mRNAs encoded by genes with neuronal functions linked to schizophrenia and autism.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Charizanis K., Lee K. Y., Batra R., Goodwin M., Zhang C., et al. (2012) Muscleblind-like 2-mediated alternative splicing in the developing brain and dysregulation in myotonic dystrophy. *Neuron* 75: 437-450

Chi S. W., Hannon G. J. and Darnell R. B. (2012) An alternative mode of microRNA target recognition. *Nat Struct Mol Biol* 19: 321-327

Riley K. J., Rabinowitz G. S., Yario T. A., Luna J. M., Darnell R. B., et al. (2012) EBV and human microRNAs co-target oncogenic and apoptotic viral and human genes during latency. *EMBO J* 31: 2207-2221

Vourekas A., Zheng Q., Alexiou P., Maragkakis M., Kirino Y., et al. (2012) Mili and Miwi target RNA repertoire reveals piRNA biogenesis and function of Miwi in spermiogenesis. *Nat Struct Mol Biol* 19: 773-781

Darnell J. C., Van Driesche S. J., Zhang C., Hung K. Y., Mele A., et al. (2011) FMRP Stalls Ribosomal Translocation on mRNAs Linked to Synaptic Function and Autism. *Cell* 146: 247-261

Polymenidou M., Lagier-Tourenne C., Hutt K. R., Huelga S. C., Moran J., et al. (2011) Long pre-mRNA depletion and RNA missplicing contribute to neuronal vulnerability from loss of TDP-43. *Nat Neurosci* 14: 459-468

Zhang C. and Darnell R. B. (2011) Mapping in vivo protein-RNA interactions at single-nucleotide resolution from HITS-CLIP data. *Nat Biotechnol* 29: 607-614

McKenna L. B., Schug J., Vourekas A., McKenna J. B., Bramswig N. C., et al. (2010) MicroRNAs control intestinal epithelial differentiation, architecture, and barrier function. *Gastroenterology* 139: 1654-1664, 1664 e1651

Yano M., Hayakawa-Yano Y., Mele A. and Darnell R. B. (2010) Nova2 regulates neuronal migration through an RNA switch in disabled-1 signaling. *Neuron* 66: 848-858

Zhang C., Frias M. A., Mele A., Ruggiu M., Eom T., et al. (2010) Integrative modeling defines the Nova splicing-regulatory network and its combinatorial controls. *Science* 329: 439-443

Chi S. W., Zang J. B., Mele A. and Darnell R. B. (2009) Argonaute HITS-CLIP decodes microRNA-mRNA interaction maps. *Nature* 460: 479-486

### Associated Kits

ARTseq™ Ribosome Profiling Kit

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit


TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

## PHOTOACTIVATABLE RIBONUCLEOSIDE-ENHANCED CROSSLINKING AND IMMUNOPRECIPITATION (PAR-CLIP)

Photoactivatable ribonucleoside-enhanced crosslinking and immunoprecipitation (PAR-CLIP) maps RNA-binding proteins (RBPs)<sup>14</sup>. This approach is similar to HITS-CLIP and CLIP-Seq, but uses much more efficient crosslinking to stabilize the protein-RNA complexes. The requirement to introduce a photoactivatable ribonucleoside limits this approach to cell culture and *in vitro* systems. In this method, 4-thiouridine (4-SU) and 6-thioguanosine (6-SG) are incorporated into transcripts of cultured cells. UV irradiation crosslinks 4-SU/6-SG-labeled transcripts to interacting RBPs. The targeted complexes are immunoprecipitated and digested with RNase T1, followed by Proteinase K, before RNA extraction. The RNA is reverse-transcribed to cDNA and sequenced. Deep sequencing of cDNA accurately maps RBPs interacting with labeled transcripts.



Pros	Cons
<ul style="list-style-type: none"> <li>Highly accurate mapping of RNA-protein interactions</li> <li>Labeling with 4-SU/6-SG improves crosslinking efficiency</li> </ul>	<ul style="list-style-type: none"> <li>Antibodies not specific to target may precipitate nonspecific complexes</li> <li>Limited to cell culture and <i>in vitro</i> systems</li> </ul>

### References

Kaneko S., Bonasio R., Saldana-Meyer R., Yoshida T., Son J., et al. (2014) Interactions between JARID2 and Noncoding RNAs Regulate PRC2 Recruitment to Chromatin. *Mol Cell* 53: 290-300

JARID2 is an accessory component of Polycomb repressive complex-2 (PRC2) required for the differentiation of embryonic stem cells (ESCs). In this study the molecular role of JARID2 in gene silencing was elucidated using RIP, ChIP, and PAR-CLIP combined with sequencing on an Illumina HiSeq 2000 system. The authors found that Meg3 and other lncRNAs from the Dlk1-Dio3 locus interact with PRC2 via JARID2. These findings suggest a more general mechanism by which lncRNAs contribute to PRC2 recruitment.

**Illumina Technology:** HiSeq 2000

14 Hafner M., Landgraf P., Ludwig J., Rice A., Ojo T., et al. (2008) Identification of microRNAs and other small regulatory RNAs using cDNA library sequencing. *Methods* 44: 3-12

Liu Y., Hu W., Murakawa Y., Yin J., Wang G., et al. (2013) Cold-induced RNA-binding proteins regulate circadian gene expression by controlling alternative polyadenylation. *Sci Rep* 3: 2054

In an effort to understand the concert of gene regulation by the circadian rhythm, the authors of this study used a mouse model with a fixed light/dark cycle, to determine genes regulated by variations in body temperature. The authors applied RNA-Seq and PAR-CLIP sequencing on an Illumina Genome Analyzer system to determine *Cirbp* and *Rbm3* as important regulators for the temperature entrained circadian gene expression. They discovered that these two proteins regulate the peripheral clocks by controlling the oscillation of alternative polyadenylation sites.

**Illumina Technology:** Genome Analyzer®; 76 bp single-end reads

Stoll G., Pietilainen O. P., Linder B., Suvisaari J., Brosi C., et al. (2013) Deletion of *TOP3beta*, a component of FMRP-containing mRNPs, contributes to neurodevelopmental disorders. *Nat Neurosci* 16: 1228-1237

Genetic studies, including studies of mRNA-binding proteins, have brought new light to the connection of mRNA metabolism to disease. In this study the authors found the deletion of the topoisomerase 3 $\beta$  (*TOP3 $\beta$* ) gene was associated with neurodevelopmental disorders in the Northern Finnish population. Combining genotyping with immunoprecipitation of mRNA-bound proteins (PAR-CLIP), the authors found that the recruitment of *TOP3 $\beta$*  to cytosolic messenger ribonucleoproteins (mRNPs) was coupled to the co-recruitment of FMRP, the disease gene involved in fragile X syndrome mental disorders.

**Illumina Technology:** Human Gene Expression—BeadArray, Human610-Quad (Infinium GT®), HumanHap300 (Duo/Duo+) (Infinium GT), HumanCNV370-Duo (Infinium GT)

Whisnant A. W., Bogerd H. P., Flores O., Ho P., Powers J. G., et al. (2013) In-depth analysis of the interaction of HIV-1 with cellular microRNA biogenesis and effector mechanisms. *MBio* 4: e000193

The question of how HIV-1 interfaces with cellular miRNA biogenesis and effector mechanisms has been highly controversial. In this paper, the authors used the Illumina HiSeq 2000 platform for deep sequencing of small RNAs in two different infected cell lines and two types of primary human cells. They unequivocally demonstrated that HIV-1 does not encode any viral miRNAs.

**Illumina Technology:** TruSeq RNA Sample Prep Kit, HiSeq 2000

Majoros W. H., Lekprasert P., Mukherjee N., Skalsky R. L., Corcoran D. L., et al. (2013) MicroRNA target site identification by integrating sequence and binding information. *Nat Methods* 10: 630-633

Mandal P. K., Ewing A. D., Hancks D. C. and Kazazian H. H., Jr. (2013) Enrichment of processed pseudogene transcripts in L1-ribonucleoprotein particles. *Hum Mol Genet* 22: 3730-3748

Hafner M., Lianoglou S., Tuschl T. and Betel D. (2012) Genome-wide identification of miRNA targets by PAR-CLIP. *Methods* 58: 94-105

Sievers C., Schlumpf T., Sawarkar R., Comoglio F. and Paro R. (2012) Mixture models and wavelet transforms reveal high confidence RNA-protein interaction sites in MOV10 PAR-CLIP data. *Nucleic Acids Res* 40: e160



Skalsky R. L., Corcoran D. L., Gottwein E., Frank C. L., Kang D., et al. (2012) The viral and cellular microRNA targetome in lymphoblastoid cell lines. *PLoS Pathog* 8: e1002484

Uniacke J., Holterman C. E., Lachance G., Franovic A., Jacob M. D., et al. (2012) An oxygen-regulated switch in the protein synthesis machinery. *Nature* 486: 126-129

Gottwein E., Corcoran D. L., Mukherjee N., Skalsky R. L., Hafner M., et al. (2011) Viral microRNA targetome of KSHV-infected primary effusion lymphoma cell lines. *Cell Host Microbe* 10: 515-526

Jungkamp A. C., Stoeckius M., Mecnas D., Grun D., Mastrobuoni G., et al. (2011) In vivo and transcriptome-wide identification of RNA binding protein target sites. *Mol Cell* 44: 828-840

Kishore S., Jaskiewicz L., Burger L., Hausser J., Khorshid M., et al. (2011) A quantitative analysis of CLIP methods for identifying binding sites of RNA-binding proteins. *Nat Methods* 8: 559-564

Lebedeva S., Jens M., Theil K., Schwanhauser B., Selbach M., et al. (2011) Transcriptome-wide analysis of regulatory interactions of the RNA-binding protein HuR. *Mol Cell* 43: 340-352

Mukherjee N., Corcoran D. L., Nusbaum J. D., Reid D. W., Georgiev S., et al. (2011) Integrative regulatory mapping indicates that the RNA-binding protein HuR couples pre-mRNA processing and mRNA stability. *Mol Cell* 43: 327-339

Hafner M., Landthaler M., Burger L., Khorshid M., Hausser J., et al. (2010) Transcriptome-wide identification of RNA-binding protein and microRNA target sites by PAR-CLIP. *Cell* 141: 129-141

Hafner M., Landthaler M., Burger L., Khorshid M., Hausser J., et al. (2010) PAR-CLIP--a method to identify transcriptome-wide the binding sites of RNA binding proteins. *J Vis Exp*

#### **Associated Kits**

ARTseq™ Ribosome Profiling Kit

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit


TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

## INDIVIDUAL NUCLEOTIDE RESOLUTION CLIP (ICLIP)

Individual nucleotide resolution CLIP (iCLIP) maps protein-RNA interactions similar to HITS-CLIP and PAR-CLIP<sup>15</sup>. This approach includes additional steps to digest the proteins after crosslinking and to map the crosslink sites with reverse transcriptase. In this method specific crosslinked RNA-protein complexes are immunoprecipitated. The complexes are then treated with proteinase K, as the protein crosslinked at the binding site remains undigested. Upon reverse transcription, cDNA truncates at the binding site and is circularized. These circularized fragments are then linearized and PCR-amplified. Deep sequencing of these amplified fragments provides nucleotide resolution of protein-binding site.



Pros	Cons
<ul style="list-style-type: none"> <li>• Nucleotide resolution of protein-binding site</li> <li>• Avoids the use of nucleases</li> <li>• Amplification allows the detection of rare events</li> </ul>	<ul style="list-style-type: none"> <li>• Antibodies not specific to target will precipitate nonspecific complexes</li> <li>• Non-linear PCR amplification can lead to biases affecting reproducibility</li> <li>• Artifacts may be introduced in the circularization step</li> </ul>

### References

Broughton J. P. and Pasquinelli A. E. (2013) Identifying Argonaute binding sites in *Caenorhabditis elegans* using iCLIP. *Methods* 63: 119-125

The identification of endogenous targets remains an important challenge in understanding miRNA function. New approaches include iCLIP-sequencing, using Illumina sequencing, for high-throughput detection of miRNA targets. In this study the iCLIP protocol was adapted for use in *Caenorhabditis elegans* to identify endogenous sites targeted by the worm Argonaute protein primarily responsible for miRNA function.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Zarnack K., König J., Tajnik M., Martincorena I., Eustermann S., et al. (2013) Direct competition between hnRNP C and U2AF65 protects the transcriptome from the exonization of Alu elements. *Cell* 152: 453-466

Alu elements are a certain type of repeat scattered all over the human genome. Interestingly, Alu elements may be found within gene regions and contain cryptic splice sites. This study investigated the mechanism by which the Alu splice sites are prevented from disrupting normal gene splicing and expression. By using CLIP with Illumina sequencing, the authors profiled mRNAs bound by protein and showed that heterogeneous nuclear riboprotein (hnRNP) C competes with the splicing factor at many genuine and cryptic splice sites. These results suggest hnRNP C acts as a genome-wide protection against transcription disruption by Alu elements.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>15</sup> König J., Zarnack K., Rot G., Curk T., Kayikci M., et al. (2010) iCLIP reveals the function of hnRNP particles in splicing at individual nucleotide resolution. *Nat Struct Mol Biol* 17: 909-915

Zund D., Gruber A. R., Zavolan M. and Muhlemann O. (2013) Translation-dependent displacement of UPF1 from coding sequences causes its enrichment in 3' UTRs. *Nat Struct Mol Biol* 20: 936-943

UPF1 is a factor involved in nonsense-mediated mRNA decay (NMD). The target binding sites and timing of the binding to target mRNAs has been investigated. In this report the binding sites of UPF1 were studied using transcriptome-wide mapping by CLIP-seq on an Illumina HiSeq 2000 system. The authors show how UPF1 binds RNA before translation and is displaced by translating ribosomes. This observation suggests that the triggering of NMD occurs after the binding of UPF1, presumably through aberrant translation termination.

**Illumina Technology:** HiSeq 2000

Rogelj B., Easton L. E., Bogu G. K., Stanton L. W., Rot G., et al. (2012) Widespread binding of FUS along nascent RNA regulates alternative splicing in the brain. *Sci Rep* 2: 603

Tollervey J. R., Curk T., Rogelj B., Briese M., Cereda M., et al. (2011) Characterizing the RNA targets and position-dependent splicing regulation by TDP-43. *Nat Neurosci* 14: 452-458

Konig J., Zarnack K., Rot G., Curk T., Kayikci M., et al. (2010) iCLIP reveals the function of hnRNP particles in splicing at individual nucleotide resolution. *Nat Struct Mol Biol* 17: 909-915

#### **Associated Kits**

ARTseq™ Ribosome Profiling Kit

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

## NATIVE ELONGATING TRANSCRIPT SEQUENCING (NET-SEQ)

Native elongating transcript sequencing (NET-Seq) maps transcription through the capture of 3' RNA<sup>16</sup>. In this method the RNA polymerase II elongation complex is immunoprecipitated, and RNA is extracted and reverse-transcribed to cDNA. Deep sequencing of the cDNA allows for 3'-end sequencing of nascent RNA, providing nucleotide resolution at transcription.



Pros	Cons
<ul style="list-style-type: none"> <li>• Mapping of nascent RNA-bound protein</li> <li>• Transcription is mapped at nucleotide resolution</li> </ul>	<ul style="list-style-type: none"> <li>• Antibodies not specific to target will precipitate nonspecific complexes</li> </ul>

### References

Larson M. H., Gilbert L. A., Wang X., Lim W. A., Weissman J. S., et al. (2013) CRISPR interference (CRISPRi) for sequence-specific control of gene expression. *Nat Protoc* 8: 2180-2196

This paper describes a protocol for selective gene repression based on clustered regularly interspaced palindromic repeats interference (CRISPRi). The protocol provides a simplified approach for rapid gene repression within 1-2 weeks. The method can also be adapted for high-throughput interrogation of genome-wide gene functions and genetic interactions, thus providing a complementary approach to standard RNA interference protocols.

**Illumina Technology:** HiSeq 2000


### Associated Kits

- ARTseq™ Ribosome Profiling Kit
- Ribo-Zero Kit
- TruSeq RNA Sample Prep Kit
- TruSeq Small RNA Sample Prep Kit
- TruSeq Stranded mRNA and Total RNA Sample Prep Kit
- TruSeq Targeted RNA Expression Kit

<sup>16</sup> Churchman L. S. and Weissman J. S. (2011) Nascent transcript sequencing visualizes transcription at nucleotide resolution. *Nature* 469: 368-373

## TARGETED PURIFICATION OF POLYSOMAL MRNA (TRAP-SEQ)

Targeted purification of polysomal mRNA (TRAP-Seq) maps translating mRNAs under various conditions<sup>17</sup>. In this method, tagged ribosomal proteins are expressed in cells. The tagged ribosomal proteins are then purified and the RNA isolated. RNAs are reverse-transcribed to cDNA. Deep sequencing of the cDNA provides single-base resolution of translating RNA.



Pros	Cons
<ul style="list-style-type: none"> <li>• Allows detection of translating RNAs</li> <li>• RNAs translated by specific targeted ribosomes can be assessed</li> <li>• No prior knowledge of the RNA is required</li> <li>• Genome-wide RNA screen</li> </ul>	<ul style="list-style-type: none"> <li>• Not as specific as more recently developed methods, such as Ribo-Seq</li> </ul>

### References

Mellen M., Ayata P., Dewell S., Kriaucionis S. and Heintz N. (2012) MeCP2 binds to 5hmC enriched within active genes and accessible chromatin in the nervous system. *Cell* 151: 1417-1430

Epigenetic markers, such as chromatin-binding factors and modifications to the DNA itself, are important for regulation of gene expression and differentiation. In this study, the DNA methylation 5-hydroxymethylcytosine (5hmC) was profiled in differentiated central nervous system cells *in vivo*. The authors found 5hmC enriched in active genes along with a strong depletion of the alternative methylation 5mC. The authors hypothesize that binding of 5hmC by methyl CpG binding protein 2 (MeCP2) plays a central role in the epigenetic regulation of neural chromatin and gene expression.

**illumina Technology:** TruSeq DNA Sample Prep Kit, HiSeq 2000

### Associated Kits

ARTseq™ Ribosome Profiling Kit

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

<sup>17</sup> Jiao Y. and Meyerowitz E. M. (2010) Cell-type specific analysis of translating RNAs in developing flowers reveals new levels of control. *Mol Syst Biol* 6: 419

## CROSSLINKING, LIGATION, AND SEQUENCING OF HYBRIDS (CLASH-SEQ)

Crosslinking, ligation, and sequencing of hybrids (CLASH-Seq) maps RNA-RNA interactions<sup>18</sup>. In this method RNA-protein complexes are UV crosslinked and affinity-purified. RNA-RNA hybrids are then ligated, isolated, and reverse-transcribed to cDNA. Deep sequencing of the cDNA provides high-resolution chimeric reads of RNA-RNA interactions.



Pros	Cons
<ul style="list-style-type: none"><li>• Maps RNA-RNA interactions</li><li>• Performed <i>in vivo</i></li></ul>	<ul style="list-style-type: none"><li>• Hybrid ligation may be difficult between short RNA fragments</li></ul>

### References

Kudla G., Granneman S., Hahn D., Beggs J. D. and Tollervey D. (2011) Cross-linking, ligation, and sequencing of hybrids reveals RNA-RNA interactions in yeast. *Proc Natl Acad Sci U S A* 108: 10010-10015

### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

<sup>18</sup> Kudla G., Granneman S., Hahn D., Beggs J. D. and Tollervey D. (2011) Cross-linking, ligation, and sequencing of hybrids reveals RNA-RNA interactions in yeast. *Proc Natl Acad Sci U S A* 108: 10010-10015

## PARALLEL ANALYSIS OF RNA ENDS SEQUENCING (PARE-SEQ) OR GENOME-WIDE MAPPING OF UNCAPPED TRANSCRIPTS (GMUCT)

Parallel analysis of RNA ends sequencing (PARE-Seq) or genome-wide mapping of uncapped transcripts (GMUCT) maps miRNA cleavage sites. Various RNA degradation processes impart characteristic sequence ends. By analyzing the cleavage sites, the degradation processes can be inferred<sup>19</sup>. In this method, degraded capped mRNA is adapter-ligated and reverse-transcribed. Fragments are then Mmel-digested, purified, 3'-adapter-ligated, and PCR-amplified. Deep sequencing of the cDNA provides information about uncapped transcripts that undergo degradation.



Pros	Cons
<ul style="list-style-type: none"> <li>• Maps degrading RNA</li> <li>• miRNA cleavage sites are identified</li> <li>• No prior knowledge of the target RNA sequence is required</li> </ul>	<ul style="list-style-type: none"> <li>• Non-linear PCR amplification can lead to biases, affecting reproducibility</li> <li>• Amplification errors caused by polymerases will be represented and sequenced incorrectly</li> </ul>

### References

Karlova R, van Haarst JC, Maliepaard C, van de Geest H, Bovy AG, Lammers M, Angenent GC, de Maagd RA; (2013) Identification of microRNA targets in tomato fruit development using high-throughput sequencing and degradome analysis. *J Exp Bot* 64: 1863-78

The biochemical and genetic processes of fruit development and ripening are of great interest for the food production industry. In this study, the involvement of miRNA in gene regulation was investigated for tomato plants to determine the fruit development processes regulated by miRNA. Using PARE-Seq, the authors identified a total of 119 target genes of miRNAs. Auxin response factors as well as two known ripening regulators were among the identified target genes, indicating an involvement of miRNAs in regulation of fruit ripening.

**Illumina Technology:** HiSeq 2000

Yang X, Wang L, Yuan D, Lindsey K, Zhang X; (2013) Small RNA and degradome sequencing reveal complex miRNA regulation during cotton somatic embryogenesis. *J Exp Bot* 64: 1521-36

The authors used PARE-seq to study miRNA expression during cotton somatic embryogenesis. They identified 25 novel miRNAs, as well as their target genes during development.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

<sup>19</sup> German M. A., Pillay M., Jeong D. H., Hetawal A., Luo S., et al. (2008) Global identification of microRNA-target RNA pairs by parallel analysis of RNA ends. *Nat Biotechnol* 26: 941-946

Shamimuzzaman M, Vodkin L; (2012) Identification of soybean seed developmental stage-specific and tissue-specific miRNA targets by degradome sequencing. *BMC Genomics* 13: 310

Bracken CP, Szubert JM, Mercer TR, Dinger ME, Thomson DW, Mattick JS, Michael MZ, Goodall GJ; (2011) Global analysis of the mammalian RNA degradome reveals widespread miRNA-dependent and miRNA-independent endonucleolytic cleavage. *Nucleic Acids Res* 39: 5658-68

Mercer TR, Neph S, Dinger ME, Crawford J, Smith MA, Shearwood AM, Haugen E, Bracken CP, Rackham O, Stamatoyannopoulos JA, Filipovska A, Mattick JS; (2011) The human mitochondrial transcriptome. *Cell* 146: 645-58

#### **Associated Kits**

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit



## TRANSCRIPT ISOFORM SEQUENCING (TIF-SEQ) OR PAIRED-END ANALYSIS OF TSS (PEAT)

Transcript isoform sequencing (TIF-Seq)<sup>20</sup> or paired-end analysis of transcription start sites (TSSs) (PEAT)<sup>21</sup> maps RNA isoforms. In this method, the 5' cap is removed with tobacco acid pyrophosphatase (TAP) treatment, then a "5'-oligocap" oligonucleotide is ligated and the RNA is reverse-transcribed. Biotinylated primers are incorporated and the circularized fragment is purified. Deep sequencing of the cDNA provides high-resolution information of the 5' and 3' ends of transcripts.



Pros	Cons
<ul style="list-style-type: none"> <li>• <i>Transcript isoforms are identified by 5' and 3' paired-end sequencing</i></li> </ul>	<ul style="list-style-type: none"> <li>• <i>Low-level transcripts may be missed or underrepresented</i></li> <li>• <i>Artifacts may be introduced during the circularization step</i></li> </ul>

### References

Pelechano V., Wei W. and Steinmetz L. M. (2013) Extensive transcriptional heterogeneity revealed by isoform profiling. Nature 497: 127-131

Identifying gene transcripts by sequencing allows high-throughput profiling of gene expression. However, methods that identify either 5' or 3' transcripts individually do not convey information about the occurrence of transcript isoforms. This paper presents TIF-Seq, a new assay for transcript isoform sequencing. By jointly determining both transcript ends for millions of RNA molecules, this method provides genome-wide detection and annotation of transcript isoforms. The authors demonstrate the TIF-Seq assay for yeast and note that over 26 major transcript isoforms per protein-coding gene were found to be expressed in yeast, suggesting a much higher genome expression repertoire than previously expected.

**Illumina Technology:** HiSeq 2000

Ni T., Corcoran D. L., Rach E. A., Song S., Spana E. P., et al. (2010) A paired-end sequencing strategy to map the complex landscape of transcription initiation. Nat Methods 7: 521-527

### Associated Kits

Ribo-Zero Kit

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

Enzyme Solutions:

Tobacco Acid Pyrophosphatase (TAP)

<sup>20</sup> Pelechano V., Wei W. and Steinmetz L. M. (2013) Extensive transcriptional heterogeneity revealed by isoform profiling. Nature 497: 127-131

<sup>21</sup> Ni T., Corcoran D. L., Rach E. A., Song S., Spana E. P., et al. (2010) A paired-end sequencing strategy to map the complex landscape of transcription initiation. Nat Methods 7: 521-527

## RNA STRUCTURE

RNA has the ability to form secondary structures that can either promote or inhibit RNA-protein or protein-protein interactions<sup>22,23</sup>. The most diverse secondary and tertiary structures are found in transfer RNAs (tRNAs) and are thought to play a major role in modulating protein translation. RNA structures were first studied in *Tetrahymena thermophila* using X-ray crystallography, but those studies are inherently cumbersome and limited<sup>24</sup>. Sequencing not only provides information on secondary structures, but it can also determine point mutation effects on RNA structures in a large number of samples. Recent studies have shown that sequencing is a powerful tool to identify RNA structures and determine their significance.



Paramecia species were one of the first model organisms used to study tRNA structure.

### Reviews

Lai D., Proctor J. R. and Meyer I. M. (2013) On the importance of cotranscriptional RNA structure formation. *RNA* 19: 1461-1473

Thapar R., Denmon A. P. and Nikonowicz E. P. (2014) Recognition modes of RNA tetraloops and tetraloop-like motifs by RNA-binding proteins. *Wiley Interdiscip Rev RNA* 5: 49-67

---



22 Osborne R. J. and Thornton C. A. (2006) RNA-dominant diseases. *Hum Mol Genet* 15 Spec No 2: R162-169

23 Thapar R., Denmon A. P. and Nikonowicz E. P. (2014) Recognition modes of RNA tetraloops and tetraloop-like motifs by RNA-binding proteins. *Wiley Interdiscip Rev RNA* 5: 49-67

24 Rich A. and RajBhandary U. L. (1976) Transfer RNA: molecular structure, sequence, and properties. *Annu Rev Biochem* 45: 805-860

## SELECTIVE 2'-HYDROXYL ACYLATION ANALYZED BY PRIMER EXTENSION SEQUENCING (SHAPE-SEQ)

Selective 2'-hydroxyl acylation analyzed by primer extension sequencing (SHAPE-Seq)<sup>25</sup> provides structural information about RNA. In this method, a unique barcode is first added to the 3' end of RNA, and the RNA is then allowed to fold under pre-established *in vitro* conditions. The barcoded and folded RNA is treated with a SHAPE reagent, 1M7, that blocks reverse transcription. The RNA is then reverse-transcribed to cDNA. Deep sequencing of the cDNA provides single-nucleotide sequence information for the positions occupied by 1M7. The structural information of the RNA can then be deduced.



Pros	Cons
<ul style="list-style-type: none"> <li>• Provides RNA structural information</li> <li>• Multiplexed analysis of barcoded RNAs provides information for multiple RNAs</li> <li>• Effect of point mutations on RNA structure can be assessed</li> <li>• Alternative to mass spectrometry, NMR, and crystallography</li> </ul>	<ul style="list-style-type: none"> <li>• Need positive and negative controls to account for transcriptase drop-off</li> <li>• Need pre-established conditions for RNA folding</li> <li>• The folding <i>in vitro</i> may not reflect actual folding <i>in vivo</i></li> </ul>

### References

Lucks J. B., Mortimer S. A., Trapnell C., Luo S., Aviran S., et al. (2011) Multiplexed RNA structure characterization with selective 2'-hydroxyl acylation analyzed by primer extension sequencing (SHAPE-Seq). *Proc Natl Acad Sci U S A* 108: 11063-11068


### Associated Kits

TruSeq Small RNA Sample Prep Kit

<sup>25</sup> Lucks J. B., Mortimer S. A., Trapnell C., Luo S., Aviran S., et al. (2011) Multiplexed RNA structure characterization with selective 2'-hydroxyl acylation analyzed by primer extension sequencing (SHAPE-Seq). *Proc Natl Acad Sci U S A* 108: 11063-11068

## PARALLEL ANALYSIS OF RNA STRUCTURE (PARS-SEQ)

Parallel analysis of RNA structure (PARS-Seq)<sup>26</sup> mapping gives information about the secondary and tertiary structure of RNA. In this method RNA is digested with RNases that are specific for double-stranded and single-stranded RNA, respectively. The resulting fragments are reverse-transcribed to cDNA. Deep sequencing of the cDNA provides high-resolution sequences of the RNA. The RNA structure can be deduced by comparing the digestion patterns of the various RNases.



Pros	Cons
<ul style="list-style-type: none"> <li>• Provides RNA structural information</li> <li>• Distinguishes between paired and unpaired bases</li> <li>• Alternative to mass spectrometry, NMR, and crystallography</li> </ul>	<ul style="list-style-type: none"> <li>• Enzyme digestion can be nonspecific</li> <li>• Digestion conditions must be carefully controlled</li> <li>• RNA can be overdigested</li> </ul>

### References

Wan Y, Qu K, Ouyang Z, Chang HY; (2013) Genome-wide mapping of RNA structure using nuclease digestion and high-throughput sequencing. Nat Protoc 8: 849-69

RNA structure is important for RNA function and regulation, and there is growing interest in determining the RNA structure of many transcripts. This is the first paper to describe the PARS protocol. In this method, enzymatic footprinting is coupled with high-throughput sequencing to retrieve information about secondary RNA structure for thousands of RNAs simultaneously.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

### Associated Kits


TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

26 Pelechano V., Wei W. and Steinmetz L. M. (2013) Extensive transcriptional heterogeneity revealed by isoform profiling. Nature 497: 127-131

## FRAGMENTATION SEQUENCING (FRAG-SEQ)

Fragmentation sequencing (FRAG-Seq)<sup>27</sup> is a method for probing RNA structure. In this method, RNA is digested using nuclease P1, followed by reverse transcription. Deep sequencing of the cDNA provides high-resolution single-stranded reads, which can be used to determine the structure of RNA by mapping P1 endonuclease digestion sites.



Pros	Cons
<ul style="list-style-type: none"> <li>• Simple and fast protocol compared to PARS-seq</li> <li>• High throughput</li> <li>• Alternative to mass spectrometry, NMR, and crystallography</li> </ul>	<ul style="list-style-type: none"> <li>• Need endogenous controls</li> <li>• Potential for contamination between samples and controls</li> </ul>

### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

<sup>27</sup> Underwood J. G., Uzirov A. V., Katzman S., Onodera C. S., Mainzer J. E., et al. (2010) FragSeq: transcriptome-wide RNA structure probing using high-throughput sequencing. Nat Methods 7: 995-1001

## CXXC AFFINITY PURIFICATION SEQUENCING (CAP-SEQ)

CXXC affinity purification sequencing (CAP-Seq)<sup>28</sup> maps the 5' end of RNAs anchored to RNA polymerase II. In this method, RNA transcripts are treated with a terminator, calf intestine alkaline phosphatase (CIP), and then tobacco acid pyrophosphatase (TAP), followed by linker ligation and reverse transcription to cDNA. Deep sequencing of the cDNA provides high-resolution sequences of RNA polymerase II transcripts.



### Pros

- Maps RNAs anchored to RNA polymerase II

### Cons

- Multiple steps and treatments can lead to loss of material

## References

Farcas A. M., Blackledge N. P., Sudbery I., Long H. K., McGouran J. F., et al. (2012) KDM2B links the Polycomb Repressive Complex 1 (PRC1) to recognition of CpG islands. *Elife* 1: e00205

DNA methylation occurs naturally throughout the genome, mostly at positions where cytosine is bonded to guanine to form a CpG dinucleotide. Many stretches of CpGs, also called CpG islands, contain a high proportion of unmethylated CpGs. In this study, the unmethylated CpG islands were studied for possible mechanisms favoring the unmethylated sites. Using ChIP-Seq experiments for various transcription factors, the authors showed that CpG islands are occupied by low levels of polycomb repressive complex 1 throughout the genome, potentially making the sites susceptible to polycomb-mediated silencing.

**Illumina Technology:** HiSeq 2000

28 Illingworth R. S., Gruenewald-Schneider U., Webb S., Kerr A. R., James K. D., et al. (2010) Orphan CpG islands identify numerous conserved promoters in the mammalian genome. *PLoS Genet* 6: e1001134

Gu W., Lee H. C., Chaves D., Youngman E. M., Pazour G. J., et al. (2012) CapSeq and CIP-TAP identify Pol II start sites and reveal capped small RNAs as *C. elegans* piRNA precursors. *Cell* 151: 1488-1500

Small RNA molecules account for many different functions in the cell. Piwi-interacting RNAs (piRNAs) represent one type of germline-expressed small RNAs linked to epigenetic programming. This study presents CAP-Seq, an assay developed to characterize the transcription of piRNAs in *C. elegans*. To their surprise, the authors found that likely piRNA precursors are capped small RNAs that initiate precisely 2 bp upstream of mature piRNAs. In addition, they identified a new class of piRNAs, further adding to the complexity of small RNA molecules.

**Illumina Technology:** Genome Analyzer<sub>IIx</sub>, HiSeq 2000

Clouaire T., Webb S., Skene P., Illingworth R., Kerr A., et al. (2012) Cfp1 integrates both CpG content and gene activity for accurate H3K4me3 deposition in embryonic stem cells. *Genes Dev* 26: 1714-1728

Gendrel A. V., Apedale A., Coker H., Termanis A., Zvetkova I., et al. (2012) Smcnd1-dependent and -independent pathways determine developmental dynamics of CpG island methylation on the inactive x chromosome. *Dev Cell* 23: 265-279

Matsushita H., Vesely M. D., Koboldt D. C., Rickert C. G., Uppaluri R., et al. (2012) Cancer exome analysis reveals a T-cell-dependent mechanism of cancer immunoediting. *Nature* 482: 400-404

Illingworth R. S., Gruenewald-Schneider U., Webb S., Kerr A. R., James K. D., et al. (2010) Orphan CpG islands identify numerous conserved promoters in the mammalian genome. *PLoS Genet* 6: e1001134

#### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA<sup>®</sup> Sample Preparation Kit

Enzyme Solutions:


Tobacco Acid Pyrophosphatase (TAP)

Calf Intestinal Phosphatase (CIP)

APex Heat-Labile Alkaline Phosphatase

## ALKALINE PHOSPHATASE, CALF INTESTINE-TOBACCO ACID PYROPHOSPHATASE SEQUENCING (CIP-TAP)

Alkaline phosphatase, calf intestine-tobacco acid pyrophosphatase sequencing (CIP-TAP) maps capped small RNAs<sup>29</sup>. In this method, RNA is treated with CIP followed by 3'-end linker ligation, then treated with TAP followed by 5'-end linker ligation. The fragments are then reverse-transcribed to cDNA, PCR-amplified, and sequenced. Deep sequencing provides single-nucleotide resolution reads of the capped small RNAs.



Pros	Cons
<ul style="list-style-type: none"> <li>Identifies capped small RNAs missed by CAP-Seq</li> <li>High throughput</li> </ul>	<ul style="list-style-type: none"> <li>Non-linear PCR amplification can lead to biases affecting reproducibility</li> <li>Amplification errors caused by polymerases</li> </ul>

### References

Yang L., Lin C., Jin C., Yang J. C., Tanasa B., et al. (2013) lncRNA-dependent mechanisms of androgen-receptor-regulated gene activation programs. *Nature* 500: 598-602

lncRNAs have recently been indicated to play a role in physiological aspects of cell-type determination and tissue homeostasis. In this paper, the authors applied three sequencing assays (GRO-Seq, ChIRP-Seq, and ChIP-Seq) using the Illumina HiSeq 2000 platform to study expression and epigenetic profiles of prostate cancer cells. The authors found two lncRNAs highly overexpressed and showed that they enhance androgen-receptor-mediated gene activation programs and proliferation of prostate cancer cells.

**Illumina Technology:** HiSeq 2000

29 Gu W., Lee H. C., Chaves D., Youngman E. M., Pazour G. J., et al. (2012) CapSeq and CIP-TAP identify Pol II start sites and reveal capped small RNAs as *C. elegans* piRNA precursors. *Cell* 151: 1488-1500



Gu W., Lee H. C., Chaves D., Youngman E. M., Pazour G. J., et al. (2012) CapSeq and CIP-TAP identify Pol II start sites and reveal capped small RNAs as *C. elegans* piRNA precursors. *Cell* 151: 1488-1500

Small RNA molecules account for many different functions in the cell. Piwi-interacting RNAs (piRNAs) represent one type of germline-expressed small RNAs linked to epigenetic programming. This study presents CAP-Seq, an assay developed to characterize the transcription of piRNAs in *C. elegans*. To their surprise, the authors found that likely piRNA precursors are capped small RNAs that initiate precisely 2 nt upstream of mature piRNAs. In addition, they identified a new class of piRNAs, further adding to the complexity of small RNA molecules.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

#### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

Enzyme Solutions:


Tobacco Acid Pyrophosphatase (TAP)

Calf Intestinal Phosphatase (CIP)

APex Heat-Labile Alkaline Phosphatase

## INOSINE CHEMICAL ERASING SEQUENCING (ICE)

Inosine chemical erasing (ICE)<sup>30</sup> identifies adenosine to inosine editing. In this method, RNA is treated with acrylonitrile, while control RNA is untreated. Control and treated RNAs are then reverse-transcribed and PCR-amplified. Inosines in RNA fragments treated with acrylonitrile cannot be reverse-transcribed. Deep sequencing of the cDNA of control and treated RNA provides high-resolution reads of inosines in RNA fragments.



Pros	Cons
<ul style="list-style-type: none"> <li>• Mapping of adenosine to inosine editing</li> <li>• Can be performed with limited material</li> </ul>	<ul style="list-style-type: none"> <li>• Non-linear PCR amplification can lead to biases, affecting reproducibility</li> <li>• Amplification errors caused by polymerases will be represented and sequenced incorrectly</li> </ul>

### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

<sup>30</sup> Sakurai M., Yano T., Kawabata H., Ueda H. and Suzuki T. (2010) Inosine cyanoethylation identifies A-to-I RNA editing sites in the human transcriptome. Nat Chem Biol 6: 733-740

## M<sup>6</sup>A-SPECIFIC METHYLATED RNA IMMUNOPRECIPITATION SEQUENCING (MERIP-SEQ)

m<sup>6</sup>A-specific methylated RNA immunoprecipitation with next generation sequencing (MeRIP-Seq)<sup>31</sup> maps m<sup>6</sup>A methylated RNA. In this method, m<sup>6</sup>A-specific antibodies are used to immunoprecipitate RNA. RNA is then reverse-transcribed to cDNA and sequenced. Deep sequencing provides high resolution reads of m<sup>6</sup>A-methylated RNA.



### Pros

- Maps m<sup>6</sup>A methylated RNA

### Cons

- Antibodies not specific to target will precipitate nonspecific RNA modifications

## References

Meyer K. D., Saletore Y., Zumbo P., Elemento O., Mason C. E., et al. (2012) Comprehensive analysis of mRNA methylation reveals enrichment in 3' UTRs and near stop codons. *Cell* 149: 1635-1646

In addition to DNA, RNA may also carry epigenetic modifications. Methylation of the N6 position of adenosine (m<sup>6</sup>A) has been implicated in the regulation of physiological processes. In this study, the authors apply MeRIP-Seq to determine mammalian genes containing m<sup>6</sup>A in their mRNA. The sites of m<sup>6</sup>A residues are enriched near stop codons and in 3'-untranslated regions (3'-UTRs), pointing to a non-random distribution and possibly functional relevance of methylated RNA transcripts.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000


## Associated Kits

- EpiGnome™ Methyl-Seq® Kit
- TruSeq RNA Sample Prep Kit
- TruSeq Small RNA Sample Prep Kit
- TruSeq Stranded mRNA and Total RNA Sample Preparation Kit
- TruSeq Targeted RNA Expression Kit

<sup>31</sup> Meyer K. D., Saletore Y., Zumbo P., Elemento O., Mason C. E., et al. (2012) Comprehensive analysis of mRNA methylation reveals enrichment in 3' UTRs and near stop codons. *Cell* 149: 1635-1646

## LOW-LEVEL RNA DETECTION

Low-level RNA detection refers to both detection of rare RNA molecules in a cell-free environment, such as circulating tumor RNA, or the expression patterns of single cells. Tissues consist of a multitude of different cell types, each with a distinctly different set of functions. Even within a single cell type, the transcriptomes are highly dynamic and reflect temporal, spatial, and cell cycle–dependent changes. Cell harvesting, handling, and technical issues with sensitivity and bias during amplification add an additional level of complexity. To resolve this multi-tiered complexity would require the analysis of many thousands of cells. The use of unique barcodes has greatly increased the number of samples that can be multiplexed and pooled, with little to no decrease in reads associated with each sample. Recent improvements in cell capture and sample preparation will provide more information, faster, and at lower cost<sup>32</sup>. This promises to fundamentally expand our understanding of cell function with significant implications for research and human health<sup>33</sup>.



Organs, such as the kidney depicted in this cross-section, consist of a myriad of phenotypically distinct cells. Single-cell transcriptomics can characterize the function of each of these cell types.

### Reviews

Blainey P. C. (2013) The future is now: single-cell genomics of bacteria and archaea. *FEMS Microbiol Rev* 37: 407-427

Lovett M. (2013) The applications of single-cell genomics. *Hum Mol Genet* 22: R22-26

Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

Spaethling J. M. and Eberwine J. H. (2013) Single-cell transcriptomics for drug target discovery. *Curr Opin Pharmacol* 13: 786-790

32 Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

33 Spaethling J. M. and Eberwine J. H. (2013) Single-cell transcriptomics for drug target discovery. *Curr Opin Pharmacol* 13: 786-790

## References

Shalek A. K., Satija R., Adiconis X., Gertner R. S., Gaublomme J. T., et al. (2013) Single-cell transcriptomics reveals bimodality in expression and splicing in immune cells. *Nature* 498: 236-240


Xue Z., Huang K., Cai C., Cai L., Jiang C. Y., et al. (2013) Genetic programs in human and mouse early embryos revealed by single-cell RNA sequencing. *Nature* 500: 593-597

Yan L., Yang M., Guo H., Yang L., Wu J., et al. (2013) Single-cell RNA-Seq profiling of human preimplantation embryos and embryonic stem cells. *Nat Struct Mol Biol* 20: 1131-1139

Goetz J. J. and Trimarchi J. M. (2012) Transcriptome sequencing of single cells with Smart-Seq. *Nat Biotechnol* 30: 763-765

## DIGITAL RNA SEQUENCING

Digital RNA sequencing is an approach to RNA-Seq that removes sequence-dependent PCR amplification biases by barcoding the RNA molecules before amplification<sup>34</sup>. RNA is reverse-transcribed to cDNA, then an excess of adapters, each with a unique barcode, is added to the preparation. This barcoded cDNA is then amplified and sequenced. Deep sequencing reads are compared, and barcodes are used to determine the actual ratio of RNA abundance.



Pros	Cons
<ul style="list-style-type: none"> <li>• Low amplification bias during PCR</li> <li>• Information about abundance of RNA</li> <li>• Detection of low-copy-number RNA</li> <li>• Single-copy resolution</li> </ul>	<ul style="list-style-type: none"> <li>• Some amplification bias still persists</li> <li>• Barcodes may miss targets during ligation</li> </ul>

### References

Shiroguchi K., Jia T. Z., Sims P. A. and Xie X. S. (2012) Digital RNA sequencing minimizes sequence-dependent bias and amplification noise with optimized single-molecule barcodes. *Proc Natl Acad Sci U S A* 109: 1347-1352

Experimental protocols that include PCR as an amplification step are subject to the sequence-dependent bias of the PCR. For RNA-Seq, this results in difficulties in quantifying expression levels, especially at very low copy numbers. In this study, digital RNA-Seq is introduced as an accurate method for quantitative measurements by appending unique barcode sequences to the pool of RNA fragments. The authors demonstrate how digital RNA-Seq allows transcriptome profiling of *Escherichia coli* with more accurate and reproducible quantification than conventional RNA-Seq. The efficacy of optimization was estimated by comparison to simulated data.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit


TruSeq Stranded mRNA and Total RNA Sample Preparation Kit

TruSeq Targeted RNA Expression Kit

<sup>34</sup> Shiroguchi K., Jia T. Z., Sims P. A. and Xie X. S. (2012) Digital RNA sequencing minimizes sequence-dependent bias and amplification noise with optimized single-molecule barcodes. *Proc Natl Acad Sci U S A* 109: 1347-1352

## WHOLE-TRANSCRIPT AMPLIFICATION FOR SINGLE CELLS (QUARTZ-SEQ)

The Quartz-Seq method optimizes whole-transcript amplification (WTA) of single cells<sup>35</sup>. In this method, a reverse-transcription (RT) primer with a T7 promoter and PCR target is first added to extracted mRNA. Reverse transcription synthesizes first-strand cDNA, after which the RT primer is digested by exonuclease I. A poly(A) tail is then added to the 3' ends of first-strand cDNA, along with a dT primer containing a PCR target. After second-strand generation, a blocking primer is added to ensure PCR enrichment in sufficient quantity for sequencing. Deep sequencing allows for accurate, high-resolution representation of the whole transcriptome of a single cell.



Pros	Cons
<ul style="list-style-type: none"> <li>• Single-tube reaction suitable for automation</li> <li>• Digestion of RT primers by exonuclease I eliminates amplification of byproducts</li> <li>• Short fragments and byproducts are suppressed during enrichment</li> </ul>	<ul style="list-style-type: none"> <li>• PCR biases can underrepresent GC-rich templates</li> <li>• Amplification errors caused by polymerases will be represented and sequenced incorrectly</li> <li>• Targets smaller than 500 bp are preferentially amplified by polymerases during PCR</li> </ul>

### References

Sasagawa Y., Nikaido I., Hayashi T., Danno H., Uno K. D., et al. (2013) Quartz-Seq: a highly reproducible and sensitive single-cell RNA sequencing method, reveals non-genetic gene-expression heterogeneity. *Genome Biol* 14: R31

Individual cells may exhibit variable gene expression even if they share the same genome. The analysis of single-cell variability in gene expression requires robust protocols with a minimum of bias. This paper presents a novel single-cell RNA-Seq method, Quartz-Seq, based on Illumina sequencing that has a simpler protocol and higher reproducibility and sensitivity than existing methods. The authors implemented improvements in three main areas: 1) they optimized the protocol for suppression of byproduct synthesis; 2) they identified a robust PCR enzyme to allow a single-tube reaction; and 3) they determined optimal conditions for RT and second-strand synthesis.

**Illumina Technology:** TruSeq RNA Sample Prep Kit, HiSeq 2000

### Associated Kits

TruSeq RNA Sample Prep Kit


TruSeq Small RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

<sup>35</sup> Sasagawa Y., Nikaido I., Hayashi T., Danno H., Uno K. D., et al. (2013) Quartz-Seq: a highly reproducible and sensitive single-cell RNA sequencing method, reveals non-genetic gene-expression heterogeneity. *Genome Biol* 14: R31

## DESIGNED PRIMER-BASED RNA SEQUENCING (DP-SEQ)

Designed Primer-based RNA sequencing (DP-Seq) is a method that amplifies mRNA from limited starting material, as low as 50 pg<sup>36</sup>. In this method, a specific set of heptamer primers are first designed. Enriched poly(A)-selected mRNA undergoes first-strand cDNA synthesis. Designed primers are then hybridized to first-strand cDNA, followed by second strand synthesis and PCR. Deep sequencing of amplified DNA allows for accurate detection of specific mRNA expression at the single-cell level.



Pros	Cons
<ul style="list-style-type: none"> <li>• As little as 50 pg of starting material can be used</li> <li>• Little transcript-length bias</li> </ul>	<ul style="list-style-type: none"> <li>• The sequences of the target areas must be known to design the heptamers</li> <li>• Exponential amplification during PCR can lead to primer-dimers and spurious PCR products<sup>37</sup></li> <li>• Some read-length bias</li> </ul>

### References

Bhargava V., Ko P., Willems E., Mercola M. and Subramaniam S. (2013) Quantitative transcriptomics using designed primer-based amplification. *Sci Rep* 3: 1740

Standard amplification of RNA transcripts before sequencing is prone to introduce bias. This paper presents a protocol for selecting a unique subset of primers to target the majority of expressed transcripts in mouse for amplification while preserving their relative abundance. This protocol was developed for Illumina sequencing platforms and the authors show how the protocol yielded high levels of amplification from as little as 50 pg of mRNA, while offering a dynamic range of over five orders of magnitude.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

### Associated Kits

TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit


<sup>36</sup> Sasagawa Y., Nikaido I., Hayashi T., Danno H., Uno K. D., et al. (2013) Quartz-Seq: a highly reproducible and sensitive single-cell RNA sequencing method, reveals non-genetic gene-expression heterogeneity. *Genome Biol* 14: R31

<sup>37</sup> Bhargava V., Head S. R., Ordoukhanian P., Mercola M. and Subramaniam S. (2014) Technical variations in low-input RNA-seq methodologies. *Sci Rep* 4: 3678



## SWITCH MECHANISM AT THE 5' END OF RNA TEMPLATES (SMART-SEQ)

Smart-Seq was developed as a single-cell sequencing protocol with improved read coverage across transcripts<sup>38</sup>. Complete coverage across the genome allows the detection of alternative transcript isoforms and single-nucleotide polymorphisms. In this protocol, cells are lysed and the RNA hybridized to an oligo(dT)-containing primer. The first strand is then created with the addition of a few untemplated C nucleotides. This poly(C) overhang is added exclusively to full-length transcripts. An oligonucleotide primer is then hybridized to the poly(C) overhang and used to synthesize the second strand. Full-length cDNAs are PCR-amplified to obtain nanogram amounts of DNA. The PCR products are purified for sequencing.



Pros	Cons
<ul style="list-style-type: none"> <li>• As little as 50 pg of starting material can be used</li> <li>• The sequence of the mRNA does not have to be known</li> <li>• Improved coverage across transcripts</li> <li>• High level of mappable reads</li> </ul>	<ul style="list-style-type: none"> <li>• Not strand-specific</li> <li>• No early multiplexing<sup>39</sup></li> <li>• Transcript length bias with inefficient transcription of reads over 4 Kb<sup>40</sup></li> <li>• Preferential amplification of high-abundance transcripts</li> <li>• The purification step may lead to loss of material</li> <li>• Could be subject to strand-invasion bias<sup>41</sup></li> </ul>

## References

Kadkhodaei B., Alvarsson A., Schintu N., Ramsköld D., Volakakis N., et al. (2013) Transcription factor Nurr1 maintains fiber integrity and nuclear-encoded mitochondrial gene expression in dopamine neurons. *Proc Natl Acad Sci U S A* 110: 2360-2365

Developmental transcription factors important in early neuron differentiation are often found expressed also in the adult brain. This study set out to investigate the development of ventral midbrain dopamine (DA) neurons by studying the transcriptional expression in a mouse model system. By using the Smart-Seq method, which allows sequencing from low amounts of total RNA, the authors could sequence RNA from laser-microdissected DA neurons. Their analysis showed transcriptional activation of the essential transcription factor Nurr1 and its key role in sustaining healthy DA cells.

**Illumina Technology:** HiSeq 2000, Genomic DNA Sample Prep Kit (FC-102-1001; Illumina)

38 Ramsköld D., Luo S., Wang Y. C., Li R., Deng Q., et al. (2012) Full-length mRNA-Seq from single-cell levels of RNA and individual circulating tumor cells. *Nat Biotechnol* 30: 777-782

39 Shapiro E., Blezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

40 Bhargava V., Head S. R., Ordoukhanian P., Mercola M. and Subramaniam S. (2014) Technical variations in low-input RNA-seq methodologies. *Sci Rep* 4: 3678

41 Tang D. T., Plessey C., Salmullah M., Suzuki A. M., Calligaris R., et al. (2013) Suppression of artifacts and barcode bias in high-throughput transcriptome analyses utilizing template switching. *Nucleic Acids Res* 41: e44

Marinov G. K., Williams B. A., McCue K., Schroth G. P., Gertz J., et al. (2014) From single-cell to cell-pool transcriptomes: Stochasticity in gene expression and RNA splicing. *Genome Res* 24: 496-510

Recent studies are increasingly discovering cell-to-cell variability in gene expression levels and transcriptional regulation. This study examined the lymphoblastoid cell line GM12878 using the Smart-Seq single-cell RNA-Seq protocol on the Illumina HiSeq 2000 platform to determine variation in transcription among individual cells. The authors determined, through careful quantification, that there are significant differences in expression among individual cells, over and above technical variation. In addition, they showed that the transcriptomes from small pools of 30-100 cells approach the information content and reproducibility of contemporary pooled RNA-Seq analysis from large amounts of input material.

**Illumina Technology:** Nextera DNA® Sample Prep Kit, HiSeq 2000

Shalek A. K., Satija R., Adiconis X., Gertner R. S., Gaublomme J. T., et al. (2013) Single-cell transcriptomics reveals bimodality in expression and splicing in immune cells. *Nature* 498: 236-240

Individual cells can exhibit substantial differences in gene expression, and only recently have genome profiling methods been developed to monitor the expression of single cells. This study applied the Smart-Seq single-cell RNA sequencing on the Illumina HiSeq 2000 platform to investigate heterogeneity in the response of mouse bone marrow-derived dendritic cells (BMDCs) to lipopolysaccharide. The authors found extensive bimodal variation in mRNA abundance and splicing patterns, which was subsequently validated using RNA fluorescence in situ hybridization for select transcripts.

**Illumina Technology:** HiSeq 2000

Yamaguchi S., Hong K., Liu R., Inoue A., Shen L., et al. (2013) Dynamics of 5-methylcytosine and 5-hydroxymethylcytosine during germ cell reprogramming. *Cell Res* 23: 329-339

Mouse primordial germ cells (PGCs) undergo genome-wide DNA methylation reprogramming to reset the epigenome for totipotency. In this study, the dynamics between 5-methylcytosine (5mC) and 5-hydroxymethylcytosine (5hmC) were characterized using immunostaining techniques and analyzed in combination with transcriptome profiles obtained with Illumina RNA sequencing. The study revealed that the dynamics of 5mC and 5hmC during PGC reprogramming support a model in which DNA demethylation in PGCs occurs through multiple steps, with both active and passive mechanisms. In addition, the transcriptome study suggests that PGC reprogramming may have an important role in the activation of a subset of meiotic and imprinted genes.

**Illumina Technology:** HiSeq 2000

Ramskold D., Luo S., Wang Y. C., Li R., Deng Q., et al. (2012) Full-length mRNA-Seq from single-cell levels of RNA and individual circulating tumor cells. *Nat Biotechnol* 30: 777-782

Yamaguchi S., Hong K., Liu R., Shen L., Inoue A., et al. (2012) Tet1 controls meiosis by regulating meiotic gene expression. *Nature* 492: 443-447

#### Associated Kits

Nextera DNA Sample Prep Kit



TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

## SWITCH MECHANISM AT THE 5' END OF RNA TEMPLATES VERSION 2 (SMART-SEQ2)

Smart-Seq2 includes several improvements over the original Smart-Seq protocol<sup>42,43</sup>. The new protocol includes a locked nucleic acid (LNA), an increased  $MgCl_2$  concentration, betaine, and elimination of the purification step to significantly improve the yield. In this protocol, single cells are lysed in a buffer that contains free dNTPs and oligo(dT)-tailed oligonucleotides with a universal 5'-anchor sequence. Reverse transcription is performed, which adds 2–5 untemplated nucleotides to the cDNA 3' end. A template-switching oligo (TSO) is added, carrying two riboguanosines and a modified guanosine to produce a LNA as the last base at the 3' end. After the first-strand reaction, the cDNA is amplified using a limited number of cycles. Tagmentation is then used to quickly and efficiently construct sequencing libraries from the amplified cDNA.



Betaine

Pros	Cons
<ul style="list-style-type: none"> <li>• The sequence of the mRNA does not have to be known</li> <li>• As little as 50 pg of starting material can be used</li> <li>• Improved coverage across transcripts</li> <li>• High level of mappable reads</li> </ul>	<ul style="list-style-type: none"> <li>• Not strand-specific</li> <li>• No early multiplexing</li> <li>• Applicable only to poly(A)+ RNA</li> </ul>

42 Picelli S., Bjorklund A. K., Faridani O. R., Sagasser S., Winberg G., et al. (2013) Smart-seq2 for sensitive full-length transcriptome profiling in single cells. Nat Methods 10: 1096-1098

43 Picelli S., Faridani O. R., Bjorklund Å. K., Winberg G., Sagasser S., et al. (2014) Full-length RNA-Seq from single cells using Smart-seq2. Nat. Protocols 9: 171-181

## References

Picelli S., Bjorklund A. K., Faridani O. R., Sagasser S., Winberg G., et al. (2013) Smart-seq2 for sensitive full-length transcriptome profiling in single cells. *Nat Methods* 10: 1096-1098

Single-cell gene expression analyses hold promise for characterizing cellular heterogeneity, but current methods compromise on the coverage, sensitivity, or throughput. This paper introduces Smart-Seq2 with improved reverse transcription, template switching, and preamplification to increase both yield and length of cDNA libraries generated from individual cells. The authors evaluated the efficacy of the Smart-Seq2 protocol using the Illumina HiSeq 2000 platform and concluded that Smart-Seq2 transcriptome libraries have improved detection, coverage, bias, and accuracy compared to Smart-Seq libraries. In addition, they are generated with off-the-shelf reagents at lower cost.

**Illumina Technology:** Nextera DNA Sample Prep Kit, HiSeq 2000


## Associated Kits

Nextera DNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

## UNIQUE MOLECULAR IDENTIFIERS (UMI)

Unique molecular identifiers (UMI) is a method that uses molecular tags to detect and quantify unique mRNA transcripts<sup>44</sup>. In this method, mRNA libraries are generated by fragmentation and then reverse-transcribed to cDNA. Oligo(dT) primers with specific sequencing linkers are added to cDNA. Another sequencing linker with a 10 bp random label and an index sequence is added to the 5' end of the template, which is amplified and sequenced. Sequencing allows for high-resolution reads, enabling accurate detection of true variants.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can sequence unique mRNA transcripts</li> <li>• Can be used to detect transcripts occurring at low frequencies</li> <li>• Transcripts can be quantified based on sequencing reads specific to each barcode</li> <li>• Can be applied to multiple platforms to karyotype chromosomes as well</li> </ul>	<ul style="list-style-type: none"> <li>• Targets smaller than 500 bp are preferentially amplified by polymerases during PCR</li> </ul>

## References

Islam S., Zeisel A., Joost S., La Manno G., Zajac P., et al. (2014) Quantitative single-cell RNA-seq with unique molecular identifiers. Nat Methods 11: 163-166

Gene expression varies among different tissues, in effect giving rise to different tissue types out of undifferentiated cells; however, expression also varies among different cells in the same tissue. Most assays for measuring gene expression depend on input material from multiple cells, but in this study a method for single-cell RNA sequencing is presented based on Illumina sequencing technology. This technology can be applied to characterize sources of transcriptional noise, or to study expression in early embryos and other sample types where the cell count is naturally limited. One attractive possibility is the application of single-cell sequencing to assess cell type diversity in complex tissues.

**Illumina Technology:** HiSeq 2000

<sup>44</sup> Kivioja T., Vaharautio A., Karlsson K., Bonke M., Enge M., et al. (2012) Counting absolute numbers of molecules using unique molecular identifiers. Nat Methods 9: 72-74

Murtaza M., Dawson S. J., Tsui D. W., Gale D., Forshew T., et al. (2013) Non-invasive analysis of acquired resistance to cancer therapy by sequencing of plasma DNA. *Nature* 497: 108-112

Recent studies have shown that genomic alterations in solid cancers can be characterized by sequencing of circulating cell-free tumor DNA released from cancer cells into plasma, representing a non-invasive liquid biopsy. This study describes how this approach was applied using Illumina HiSeq sequencing technology to track the genomic evolution of metastatic cancers in response to therapy. Six patients with breast, ovarian, and lung cancers were followed over 1–2 years. For two cases, synchronous biopsies were also analyzed, confirming genome-wide representation of the tumor genome in plasma and establishing the proof-of-principle of exome-wide analysis of circulating tumor DNA.

**Illumina Technology:** TruSeq Exome® Enrichment Kit, HiSeq 2000

Kivioja T., Vaharautio A., Karlsson K., Bonke M., Enge M., et al. (2012) Counting absolute numbers of molecules using unique molecular identifiers. *Nat Methods* 9: 72-74

This is the first paper to describe the UMI method and its utility as a tool for sequencing. The authors use UMIs, which make each molecule in a population distinct for genome-scale karyotyping and mRNA sequencing.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

#### Associated Kits


TruSeq RNA Sample Prep Kit

TruSeq Small RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

## CELL EXPRESSION BY LINEAR AMPLIFICATION SEQUENCING (CEL-SEQ)

Cell expression by linear amplification sequencing (CEL-Seq) is a method that utilizes barcoding and pooling of RNA to overcome challenges from low input<sup>45</sup>. In this method, each cell undergoes reverse transcription with a unique barcoded primer in its individual tube. After second-strand synthesis, cDNAs from all reaction tubes are pooled, and PCR-amplified. Paired-end deep sequencing of the PCR products allows for accurate detection of sequence derived from sequencing both strands.



Pros	Cons
<ul style="list-style-type: none"> <li>• <i>Barcoding and pooling allow for multiplexing and studying many different single cells at a time</i></li> <li>• <i>Cross-contamination is greatly reduced due to using one tube per cell</i></li> <li>• <i>Fewer steps than STRT-Seq</i></li> <li>• <i>Very little read-length bias<sup>46</sup></i></li> <li>• <i>Strand-specific</i></li> </ul>	<ul style="list-style-type: none"> <li>• <i>Strongly 3' biased<sup>47</sup></i></li> <li>• <i>Abundant transcripts are preferentially amplified</i></li> <li>• <i>Requires at least 400 pg of total RNA</i></li> </ul>

### References

Hashimshony T., Wagner F., Sher N. and Yanai I. (2012) CEL-Seq: single-cell RNA-Seq by multiplexed linear amplification. *Cell Rep* 2: 666-673

High-throughput sequencing has allowed for unprecedented detail in gene expression analyses, yet its efficient application to single cells is challenged by the small starting amounts of RNA. This paper presents the CEL-Seq protocol, which uses barcoding, pooling of samples, and linear amplification with one round of *in vitro* transcription. The assay is designed around a modified version of the Illumina directional RNA protocol and sequencing is done on the Illumina HiSeq 2000 system. The authors demonstrate their method by single-cell expression profiling of early *C. elegans* embryonic development.

**Illumina Technology:** HiSeq 2000

### Associated Kits

TruSeq RNA Sample Prep Kit


45 Hashimshony T., Wagner F., Sher N. and Yanai I. (2012) CEL-Seq: single-cell RNA-Seq by multiplexed linear amplification. *Cell Rep* 2: 666-673

46 Bhargava V., Head S. R., Ordoukhanian P., Mercola M. and Subramaniam S. (2014) Technical variations in low-input RNA-seq methodologies. *Sci Rep* 4: 3678

47 Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

## SINGLE-CELL TAGGED REVERSE TRANSCRIPTION SEQUENCING (STRT-SEQ)

Single-cell tagged reverse transcription sequencing (STRT-Seq) is a method similar to CEL-seq that involves unique barcoding and sample pooling to overcome the challenges of samples with limited material<sup>48</sup>. In this method, single cells are first picked in individual tubes, where first-strand cDNA synthesis occurs using an oligo(dT) primer with the addition of 3–6 cytosines. A helper oligo promotes template switching, which introduces the barcode on the cDNA. Barcoded cDNA is then amplified by single-primer PCR. Deep sequencing allows for accurate transcriptome sequencing of individual cells.



Pros	Cons
<ul style="list-style-type: none"> <li>• <i>Barcoding and pooling allows for multiplexing and studying many different single cells at a time</i></li> <li>• <i>Sample handling and the potential for cross-contamination are greatly reduced due to using one tube per cell</i></li> </ul>	<ul style="list-style-type: none"> <li>• <i>PCR biases can underrepresent GC-rich templates</i></li> <li>• <i>Non-linear PCR amplification can lead to biases affecting reproducibility</i></li> <li>• <i>Amplification errors caused by polymerases will be represented and sequenced incorrectly</i></li> <li>• <i>Loss of accuracy due to PCR bias</i></li> <li>• <i>Targets smaller than 500 bp are preferentially amplified by polymerases during PCR</i></li> </ul>

### References

Islam S., Kjallquist U., Moliner A., Zajac P., Fan J. B., et al. (2011) Characterization of the single-cell transcriptional landscape by highly multiplex RNA-seq. *Genome Res* 21: 1160-1167

Gene expression varies among different tissues, in effect giving rise to different tissue types out of undifferentiated cells; however, expression also varies among different cells in the same tissue. Most assays for measuring gene expression depend on input material from multiple cells, but in this study a method for single-cell RNA sequencing is presented based on Illumina sequencing technology. This technology can be applied to characterize sources of transcriptional noise, or to study expression in early embryos and other sample types where the cell count is naturally limited. One attractive possibility is the application of single-cell sequencing to assess cell type diversity in complex tissues.

**Illumina Technology:** HiSeq 2000

### Associated Kits

TruSeq RNA Sample Prep Kit


TruSeq Small RNA Sample Prep Kit

TruSeq Targeted RNA Expression Kit

48 Islam S., Kjallquist U., Moliner A., Zajac P., Fan J. B., et al. (2011) Characterization of the single-cell transcriptional landscape by highly multiplex RNA-seq. *Genome Res* 21: 1160-1167



Single-cell genomics can be used to identify and study circulating tumor cells, cell-free DNA, microbes, uncultured microbes, for preimplantation diagnosis, and to help us better understand tissue-specific cellular differentiation<sup>49, 50</sup>. DNA replication during cell division is not perfect; as a result, progressive generations of cells accumulate unique somatic mutations. Consequently, each cell in our body has a unique genomic signature, which allows the reconstruction of cell lineage trees with very high precision.<sup>51</sup> These cell lineage trees can predict the existence of small populations of stem cells. This information is important for fields as diverse as cancer development<sup>52, 53</sup> preimplantation, and genetic diagnosis.<sup>54, 55</sup>



Single-cell genomics can help characterize and identify circulating tumor cells as well as microbes.

### Reviews:

Blainey P. C. (2013) The future is now: single-cell genomics of bacteria and archaea. *FEMS Microbiol Rev* 37: 407-427

Lovett M. (2013) The applications of single-cell genomics. *Hum Mol Genet* 22: R22-26

Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

---

49 Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630

50 Blainey P. C. (2013) The future is now: single-cell genomics of bacteria and archaea. *FEMS Microbiol Rev* 37: 407-427

51 Frumkin D., Wasserstrom A., Kaplan S., Feige U. and Shapiro E. (2005) Genomic variability within an organism exposes its cell lineage tree. *PLoS Comput Biol* 1: e5

52 Navin N., Kendall J., Troge J., Andrews P., Rodgers L., et al. (2011) Tumour evolution inferred by single-cell sequencing. *Nature* 472: 90-94

53 Potter N. E., Ermini L., Papaemmanuil E., Cazzaniga G., Vijayaraghavan G., et al. (2013) Single-cell mutational profiling and clonal phylogeny in cancer. *Genome Res* 23: 2115-2125

54 Van der Aa N., Esteki M. Z., Vermeesch J. R. and Voet T. (2013) Preimplantation genetic diagnosis guided by single-cell genomics. *Genome Med* 5: 71

55 Hou Y., Fan W., Yan L., Li R., Lian Y., et al. (2013) Genome analyses of single human oocytes. *Cell* 155: 1492-1506

Baslan T., Kendall J., Rodgers L., Cox H., Riggs M., et al. (2012) Genome-wide copy number analysis of single cells. *Nat Protoc* 7: 1024-1041

Böttcher R., Amberg R., Ruzius F. P., Guryev V., Verhaegh W. F., et al. (2012) Using a priori knowledge to align sequencing reads to their exact genomic position. *Nucleic Acids Res* 40: e125

Kalisky T. and Quake S. R. (2011) Single-cell genomics. *Nat Methods* 8: 311-314

Navin N. and Hicks J. (2011) Future medical applications of single-cell sequencing in cancer. *Genome Med* 3: 31

Yilmaz S. and Singh A. K. (2011) Single cell genome sequencing. *Curr Opin Biotechnol* 23: 437-443


## References

Voet T., Kumar P., Van Loo P., Cooke S. L., Marshall J., et al. (2013) Single-cell paired-end genome sequencing reveals structural variation per cell cycle. *Nucleic Acids Res* 41: 6119-6138

Hou Y., Song L., Zhu P., Zhang B., Tao Y., et al. (2012) Single-cell exome sequencing and monoclonal evolution of a JAK2-negative myeloproliferative neoplasm. *Cell* 148: 873-885

## SINGLE-MOLECULE MOLECULAR INVERSION PROBES (SMIP)

The single-molecule molecular inversion probes (smMIP) method uses single-molecule tagging and molecular inversion probes to detect and quantify genetic variations occurring at very low frequencies<sup>56</sup>. In this method, probes are used to detect targets in genomic DNA. After the probed targets are copied, exonuclease digestion leaves the target with a tag, which undergoes PCR amplification and sequencing. Sequencing allows for high-resolution sequence reads of targets, while greater depth allows for better alignment for every unique molecular tag.



Pros	Cons
<ul style="list-style-type: none"> <li>• Detection of low-frequency targets</li> <li>• Can perform single-cell sequencing or sequencing for samples with very limited starting material</li> </ul>	<ul style="list-style-type: none"> <li>• PCR amplification errors</li> <li>• PCR biases can underrepresent GC-rich templates</li> <li>• Targets smaller than 500 bp are preferentially amplified by polymerases during PCR</li> </ul>

### References

Hiatt J. B., Pritchard C. C., Salipante S. J., O’Roak B. J. and Shendure J. (2013) Single molecule molecular inversion probes for targeted, high-accuracy detection of low-frequency variation. *Genome Res* 23: 843-854

This is the first paper to describe the smMIP assay, along with its practicality, ability for multiplexing, scaling, and compatibility with desktop sequencing for rapid data collection. The authors demonstrated the assay by resequencing 33 clinically informative cancer genes in 8 cell lines and 45 clinical cancer samples, retrieving accurate data.

**Illumina Technology:** MiSeq®, HiSeq 2000

### Associated Kits

TruSeq Nano DNA® Sample Prep Kit

TruSeq DNA PCR-Free® Sample Prep Kit

Nextera DNA Sample Prep Kit

Nextera XT DNA® Sample Prep Kit

Nextera Rapid Capture Exome/Custom® Enrichment Kit

<sup>56</sup> Hiatt J. B., Pritchard C. C., Salipante S. J., O’Roak B. J. and Shendure J. (2013) Single molecule molecular inversion probes for targeted, high-accuracy detection of low-frequency variation. *Genome Res* 23: 843-854

## MULTIPLE DISPLACEMENT AMPLIFICATION (MDA)

Multiple displacement amplification (MDA) is a method commonly used for sequencing microbial genomes due to its ability to amplify templates larger than 0.5 Mbp, but it can also be used to study genomes of other sizes<sup>57</sup>. In this method, 3'-blocked random hexamer primers are hybridized to the template, followed by synthesis with Phi 29 polymerase. Phi 29 performs strand-displacement DNA synthesis, allowing for efficient and rapid DNA amplification. Deep sequencing of the amplified DNA allows for accurate representation of reads, while sequencing depth provides better alignment and consensus for sequences.



Pros	Cons
<ul style="list-style-type: none"> <li>• Templates used for this method can be circular DNA (plasmids, bacterial DNA)</li> <li>• Can sequence large templates</li> <li>• Can perform single-cell sequencing or sequencing for samples with very limited starting material</li> </ul>	<ul style="list-style-type: none"> <li>• Strong amplification bias. Genome coverage as low as ~6%<sup>58</sup></li> <li>• PCR biases can underrepresent GC-rich templates</li> <li>• Contaminated reagents can impact results<sup>59</sup></li> </ul>

## References

Embree M., Nagarajan H., Movahedi N., Chitsaz H. and Zengler K. (2013) Single-cell genome and metatranscriptome sequencing reveal metabolic interactions of an alkane-degrading methanogenic community. *ISME J*

Microbial communities amass a wealth of biochemical processes, and metagenomics approaches are often unable to decipher the key functions of individual microorganisms. This study analyzed a microbial community by first determining the genome sequence of a dominant bacterial member of the genus *Smithella*, using a single-cell sequencing approach on the Illumina Genome Analyzer. After establishing a working draft genome of *Smithella*, the authors used low-input metatranscriptomics to determine which genes were active during alkane degradation. The authors then designed a genome-scale metabolic model to integrate the genomic and transcriptomic data.

**Illumina Technology:** Nextera DNA Sample Prep Kit, MiSeq, Genome Analyzer<sub>ix</sub>

<sup>57</sup> Dean F. B., Nelson J. R., Giesler T. L. and Lasken R. S. (2001) Rapid amplification of plasmid and phage DNA using Phi 29 DNA polymerase and multiply-primed rolling circle amplification. *Genome Res* 11: 1095-1099

<sup>58</sup> Navin N., Kendall J., Troge J., Andrews P., Rodgers L., et al. (2011) Tumour evolution inferred by single-cell sequencing. *Nature* 472: 90-94

<sup>59</sup> Woyke T., Sczyrba A., Lee J., Rinke C., Tighe D., et al. (2011) Decontamination of MDA reagents for single cell whole genome amplification. *PLoS ONE* 6: e26161

Hou Y., Fan W., Yan L., Li R., Lian Y., et al. (2013) Genome analyses of single human oocytes. *Cell* 155: 1492-1506

Chromosomal crossover occurs in the oocyte, producing unique combinations of the parent chromosomes in the fertilized egg. This paper presents a protocol for single-cell genome analysis of human oocytes. Using multiple annealing and looping-based amplification cycle (MALBAC)-based sequencing, the authors sequenced triads of the first and second polar bodies from oocyte pronuclei. These pronuclei were derived from the same female egg donors and the authors phased their genomes to determine crossover maps for the oocytes. This breakthrough assay makes important progress toward using whole-genome sequencing for meiosis research and embryo selection for *in vitro* fertilization.

**Illumina Technology:** HiSeq 2000

McLean J. S., Lombardo M. J., Ziegler M. G., Novotny M., Yee-Greenbaum J., et al. (2013) Genome of the pathogen *Porphyromonas gingivalis* recovered from a biofilm in a hospital sink using a high-throughput single-cell genomics platform. *Genome Res* 23: 867-877

Single-cell genomics is becoming an accepted method to capture novel genomes, primarily in marine and soil environments. This study shows, for the first time, that it also enables comparative genomic analysis of strain variation in a pathogen captured from complex biofilm samples in a healthcare facility. The authors present a nearly complete genome representing a novel strain of the periodontal pathogen *Porphyromonas gingivalis* using the single-cell assembly tool SPAdes.

**Illumina Technology:** Nextera DNA Sample Prep Kit, Genome Analyzer<sub>ix</sub>

Seth-Smith H. M., Harris S. R., Skilton R. J., Radebe F. M., Golparian D., et al. (2013) Whole-genome sequences of *Chlamydia trachomatis* directly from clinical samples without culture. *Genome Res* 23: 855-866

The use of whole-genome sequencing as a tool to study infectious bacteria is of growing clinical interest. Cultures of *Chlamydia trachomatis* have, until now, been a prerequisite to obtaining DNA for whole-genome sequencing. Unfortunately, culturing *C. trachomatis* is a technically demanding and time-consuming procedure. This paper presents IMS-MDA: a new approach combining immunomagnetic separation (IMS) and multiple-displacement amplification (MDA) for whole-genome sequencing of bacterial genomes directly from clinical samples.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

Dunowska M., Biggs P. J., Zheng T. and Perrott M. R. (2012) Identification of a novel nidovirus associated with a neurological disease of the Australian brushtail possum (*Trichosurus vulpecula*). *Vet Microbiol* 156: 418-424

Wobbly possum disease (WPD) is a fatal neurological disease of the Australian brushtail possum. In this study, the previously unconfirmed mechanism of disease transmission was identified as a novel virus. The identification utilized enrichment for viral DNA followed by sequencing on an Illumina Genome Analyzer.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Chitsaz H., Yee-Greenbaum J. L., Tesler G., Lombardo M. J., Dupont C. L., et al. (2011) Efficient de novo assembly of single-cell bacterial genomes from short-read data sets. *Nat Biotechnol* 29: 915-921

Woyke T., Tighe D., Mavromatis K., Clum A., Copeland A., et al. (2010) One bacterial cell, one complete genome. *PLoS ONE* 5: e10314

Valentim C. L., LoVerde P. T., Anderson T. J. and Criscione C. D. (2009) Efficient genotyping of *Schistosoma mansoni* miracidia following whole genome amplification. *Mol Biochem Parasitol* 166: 81-84

Jasmine F., Ahsan H., Andrusis I. L., John E. M., Chang-Claude J., et al. (2008) Whole-genome amplification enables accurate genotyping for microarray-based high-density single nucleotide polymorphism array. *Cancer Epidemiol Biomarkers Prev* 17: 3499-3508

### **Associated Kits**

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit


Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Rapid Capture Exome/Custom Enrichment Kit

## MULTIPLE ANNEALING AND LOOPING-BASED AMPLIFICATION CYCLES (MALBAC)

Multiple annealing and looping-based amplification cycles (MALBAC) is intended to address some of the shortcomings of MDA<sup>60</sup>. In this method, MALBAC primers randomly anneal to a DNA template. A polymerase with displacement activity at elevated temperatures amplifies the template, generating “semi-amplicons.” As the amplification and annealing process is repeated, the semi-amplicons are amplified into full amplicons that have a 3’ end complementary to the 5’ end. As a result, full-amplicon ends hybridize to form a looped structure, inhibiting further amplification of the looped amplicon, while only the semi-amplicons and genomic DNA undergo amplification. Deep sequencing of the full-amplicon sequences allows for accurate representation of reads, while sequencing depth provides improved alignment for consensus sequences.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can sequence large templates</li> <li>• Can perform single-cell sequencing or sequencing for samples with very limited starting material</li> <li>• Full-amplicon looping inhibits over-representation of templates, reducing PCR bias</li> <li>• Can amplify GC-rich regions</li> <li>• Uniform genome coverage</li> <li>• Lower allele drop-out rate compared to MDA</li> </ul>	<ul style="list-style-type: none"> <li>• Polymerase is relatively error prone compared to Phi 29</li> <li>• Temperature-sensitive protocol</li> <li>• Genome coverage up to ~90%,<sup>61</sup> but some regions of the genome are consistently underrepresented<sup>62</sup></li> </ul>

### References

Hou Y., Fan W., Yan L., Li R., Lian Y., et al. (2013) Genome analyses of single human oocytes. *Cell* 155: 1492-1506

Chromosomal crossover occurs in the oocyte, producing unique combinations of the parent chromosomes in the fertilized egg. This paper presents a protocol for single-cell genome analysis in human oocytes. Using multiple annealing and looping-based amplification cycle (MALBAC)-based sequencing, the authors sequenced triads of the first and second polar bodies from oocyte pronuclei. These pronuclei were derived from the same female egg donors and the authors phased their genomes to determine crossover maps for the oocytes. This breakthrough assay makes important progress toward using whole-genome sequencing for meiosis research and embryo selection for *in vitro* fertilization.

**Illumina Technology:** HiSeq 2000

60 Zong C., Lu S., Chapman A. R. and Xie X. S. (2012) Genome-wide detection of single-nucleotide and copy-number variations of a single human cell. *Science* 338: 1622-1626

61 Lovett M. (2013) The applications of single-cell genomics. *Hum Mol Genet* 22: R22-26

62 Lasken R. S. (2013) Single-cell sequencing in its prime. *Nat Biotechnol* 31: 211-212

Ni X., Zhuo M., Su Z., Duan J., Gao Y., et al. (2013) Reproducible copy number variation patterns among single circulating tumor cells of lung cancer patients. *Proc Natl Acad Sci U S A* 110: 21083-21088

There is a great deal of interest in identifying and studying circulating tumor cells (CTCs). Cells from primary tumors enter the bloodstream and can seed metastases. A major barrier to such analysis is low input amounts from single cells, leading to lower coverage. In this study the authors use MALBAC for whole-genome sequencing of single CTCs from patients with lung cancer. They identify copy-number variations that were consistent in patients with the same cancer subtype. Such information about cancers can help identify drug resistance and cancer subtypes, and offers potential for diagnostics, allowing for individualized treatment.

**Illumina Technology:** MiSeq, HiSeq 2000

Zong C., Lu S., Chapman A. R. and Xie X. S. (2012) Genome-wide detection of single-nucleotide and copy-number variations of a single human cell. *Science* 338: 1622-1626

This is the first paper that describes the MALBAC method, which the authors indicate has a higher detection efficiency than the traditional MDA method for single-cell studies. The authors show detection of copy-number variations and single-nucleotide variations of single cancer cells with no false positives.

**Illumina Technology:** HiSeq 2000

#### **Associated Kits**


TruSeq Nano DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit



## OLIGONUCLEOTIDE-SELECTIVE SEQUENCING (OS-SEQ)

Oligonucleotide-selective sequencing (OS-Seq)<sup>63</sup> was developed to improve targeted resequencing, by capturing and sequencing gene targets directly on the flow cell. In this method target sequences with adaptors are used to modify the flow cell primers. Targets in the template are captured onto the flow cell with the modified primers. Further extension, denaturation, and hybridization provide sequence reads for target genes. Deep sequencing provides accurate representation of reads.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can resequence multiple targets at a time</li> <li>• No gel excision or narrow size purification required</li> <li>• Very fast (single-day) protocol</li> <li>• Samples can be multiplexed</li> <li>• Reduced PCR bias due to removal of amplification steps</li> <li>• Avoids loss of material</li> </ul>	<ul style="list-style-type: none"> <li>• Primers may interact with similar target sequences, leading to sequence ambiguity</li> </ul>

### References

Myllykangas S., Buenrostro J. D., Natsoulis G., Bell J. M. and Ji H. P. (2011) Efficient targeted resequencing of human germline and cancer genomes by oligonucleotide-selective sequencing. *Nat Biotechnol* 29: 1024-1027

As a new method for targeted genome resequencing, the authors present OS-Seq. The method uses a modification of the immobilized lawn of oligonucleotide primers on the flow cell to function as both a capture and sequencing substrate. The method is demonstrated by targeted sequencing of tumor/normal tissue from colorectal cancer.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>


### Associated Kits

- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Rapid Capture Exome/Custom Enrichment Kit

63 Myllykangas S., Buenrostro J. D., Natsoulis G., Bell J. M. and Ji H. P. (2011) Efficient targeted resequencing of human germline and cancer genomes by oligonucleotide-selective sequencing. *Nat Biotechnol* 29: 1024-1027

## DUPLEX SEQUENCING (DUPLEX-SEQ)

Duplex sequencing is a tag-based error correction method to improve sequencing accuracy<sup>64</sup>. In this method, adapters (with primer sequences and random 12 bp indices) are ligated onto the template and amplified using PCR. Deep sequencing provides consensus sequence information from every unique molecular tag. Based on molecular tags and sequencing primers, duplex sequences are aligned, determining the true sequence on each DNA strand.



Pros	Cons
<ul style="list-style-type: none"> <li>• Very low error rate due to duplex tagging system</li> <li>• PCR amplification errors can be detected and removed from analysis</li> <li>• No additional library preparation steps after addition of adaptors</li> </ul>	<ul style="list-style-type: none"> <li>• PCR amplification errors</li> <li>• Non-linear PCR amplification can lead to biases affecting reproducibility</li> <li>• PCR biases can underrepresent GC-rich templates</li> </ul>

### References

Kennedy S. R., Salk J. J., Schmitt M. W. and Loeb L. A. (2013) Ultra-sensitive sequencing reveals an age-related increase in somatic mitochondrial mutations that are inconsistent with oxidative damage. *PLoS Genet* 9: e1003794

Studies of mitochondrial DNA (mtDNA) mutations have been limited due to technical limitations of the protocols applied. In this paper, the authors present a highly sensitive Duplex-Seq method, based on the HiSeq platform, which can detect a single mutation among >107 wild-type molecules. The authors applied the method to study the accumulation of mutations in mtDNA over the course of 80 years of life. Their results show that the mutation spectra of brain tissue of old compared to young individuals are dominated by transition mutations and not G to T mutations, which are the characteristic mutations caused by oxidative damage.

**Illumina Technology:** HiSeq 2000/2500; 101 bp paired-end reads

Schmitt M. W., Kennedy S. R., Salk J. J., Fox E. J., Hiatt J. B., et al. (2012) Detection of ultra-rare mutations by next-generation sequencing. *Proc Natl Acad Sci U S A* 109: 14508-14513

The authors propose a tag-based error correction method to improve sequencing accuracy, especially in heterogeneous samples. The method allows double-stranded DNA sequence read collection, proving mutation status on both strands. The method is demonstrated by sequencing M13mp2 DNA. This method is proposed to be useful for assessing mutations due to DNA damage, as well as the determining the mutational status of genes on both DNA strands.

**Illumina Technology:** HiSeq 2000

### Associated Kits

- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Rapid Capture Exome/Custom Enrichment Kit

64 Schmitt M. W., Kennedy S. R., Salk J. J., Fox E. J., Hiatt J. B., et al. (2012) Detection of ultra-rare mutations by next-generation sequencing. *Proc Natl Acad Sci U S A* 109: 14508-14513

## DNA METHYLATION

DNA methylation and hydroxymethylation are involved in development, X-chromosome inactivation, cell differentiation, tissue-specific gene expression, plant epigenetic variation, imprinting, cancers, and diseases<sup>65,66,67,68</sup>. Methylation usually occurs at the 5' position of cytosines and plays a crucial role in gene regulation and chromatin remodeling.



**The active agouti gene in mice codes for yellow coat color. When pregnant mice with the active agouti gene are fed a diet rich in methyl donors, the offspring are born with the agouti gene turned off<sup>69</sup>. This effect has been used as an epigenetic biosensor for nutritional and environmental alterations on the fetal epigenome<sup>70</sup>.**

Most cytosine methylation occurs on cytosines located near guanines, called CpG sites. These CpG sites are often located upstream of promoters, or within the gene body. CpG islands are defined as regions that are greater than 500 bp in length with greater than 55% GC and an expected/observed CpG ratio of  $> 0.65$ .

While cytosine methylation (5mC) is known as a silencing mark that represses genes, cytosine hydroxymethylation (5hmC) is shown to be an activating mark that promotes gene expression and is a proposed intermediate in the DNA demethylation pathway<sup>1,4,6</sup>. Similar to 5mC, 5hmC is involved during development, cancers, cell differentiation, and diseases<sup>71</sup>.

5mC and/or 5hmC can be a diagnostic tool to help identify the effects of nutrition, carcinogens<sup>72</sup>, and environmental factors in relation to diseases. The impact of these modifications on gene regulation depends on their locations within the genome. It is therefore important to determine the exact position of the modified bases.

65 Smith Z. D. and Meissner A. (2013) DNA methylation: roles in mammalian development. *Nat Rev Genet* 14: 204-220

66 Jullien P. E. and Berger F. (2010) DNA methylation reprogramming during plant sexual reproduction? *Trends Genet* 26: 394-399

67 Schmitz R. J., He Y., Valdes-Lopez O., Khan S. M., Joshi T., et al. (2013) Epigenome-wide inheritance of cytosine methylation variants in a recombinant inbred population. *Genome Res* 23: 1663-1674

68 Koh K. P. and Rao A. (2013) DNA methylation and methylcytosine oxidation in cell fate decisions. *Curr Opin Cell Biol* 25: 152-161

69 Dolinoy D. C., Weidman J. R., Waterland R. A. and Jirtle R. L. (2006) Maternal genistein alters coat color and protects Avy mouse offspring from obesity by modifying the fetal epigenome. *Environ Health Perspect* 114: 567-572

70 Dolinoy D. C. (2008) The agouti mouse model: an epigenetic biosensor for nutritional and environmental alterations on the fetal epigenome. *Nutr Rev* 66 Suppl 1: S7-11, Dolinoy D. C. and Faulk C. (2012) Introduction: The use of animals models to advance epigenetic science. *ILAR J* 53: 227-231

71 Pfeifer G. P., Kadam S. and Jin S. G. (2013) 5-hydroxymethylcytosine and its potential roles in development and cancer. *Epigenetics Chromatin* 6: 10

72 Thomson J. P., Lempainen H., Hackett J. A., Nestor C. E., Muller A., et al. (2012) Non-genotoxic carcinogen exposure induces defined changes in the 5-hydroxymethylome. *Genome Biol* 13: R93

Base	Sequence	BS Sequence	oxBS Sequence	TAB Sequence	RRBS Sequence
C	C	T	T	T	T
5mC	C	C	C	T	C
5hmC	C	C	T	C	C

Sequencing reads created by various methods

## Reviews

- Koh K. P. and Rao A. (2013) DNA methylation and methylcytosine oxidation in cell fate decisions. *Curr Opin Cell Biol* 25: 152-161
- Lister R., Mukamel E. A., Nery J. R., Urich M., Puddifoot C. A., et al. (2013) Global epigenomic reconfiguration during mammalian brain development. *Science* 341: 1237905
- Pfeifer G. P., Kadam S. and Jin S. G. (2013) 5-hydroxymethylcytosine and its potential roles in development and cancer. *Epigenetics Chromatin* 6: 10
- Piccolo F. M. and Fisher A. G. (2014) Getting rid of DNA methylation. *Trends Cell Biol* 24: 136-143
- Rivera C. M. and Ren B. (2013) Mapping human epigenomes. *Cell* 155: 39-55
- Schweiger M. R., Barmeyer C. and Timmermann B. (2013) Genomics and epigenomics: new promises of personalized medicine for cancer patients. *Brief Funct Genomics* 12: 411-421
- Smith Z. D. and Meissner A. (2013) DNA methylation: roles in mammalian development. *Nat Rev Genet* 14: 204-220
- Telese F., Gamliel A., Skowronska-Krawczyk D., Garcia-Bassets I. and Rosenfeld M. G. (2013) "Seq-ing" insights into the epigenetics of neuronal gene regulation. *Neuron* 77: 606-623
- Veluchamy A., Lin X., Maumus F., Rivarola M., Bhavsar J., et al. (2013) Insights into the role of DNA methylation in diatoms by genome-wide profiling in *Phaeodactylum tricornutum*. *Nat Commun* 4: 2091
- Vidaki A., Daniel B. and Court D. S. (2013) Forensic DNA methylation profiling--Potential opportunities and challenges. *Forensic Sci Int Genet* 7: 499-507


## References

- Meaburn E. and Schulz R. (2012) Next generation sequencing in epigenetics: insights and challenges. *Semin Cell Dev Biol* 23: 192-199
- Thomson J. P., Lempiainen H., Hackett J. A., Nestor C. E., Muller A., et al. (2012) Non-genotoxic carcinogen exposure induces defined changes in the 5-hydroxymethylome. *Genome Biol* 13: R93
- Jin S. G., Kadam S. and Pfeifer G. P. (2010) Examination of the specificity of DNA methylation profiling techniques towards 5-methylcytosine and 5-hydroxymethylcytosine. *Nucleic Acids Res* 38: e125
- Dolinoy D. C., Weidman J. R., Waterland R. A. and Jirtle R. L. (2006) Maternal genistein alters coat color and protects Avy mouse offspring from obesity by modifying the fetal epigenome. *Environ Health Perspect* 114: 567-572

## BISULFITE SEQUENCING (BS-SEQ)

Bisulfite sequencing (BS-Seq) or whole-genome bisulfite sequencing (WGBS) is a well-established protocol to detect methylated cytosines in genomic DNA<sup>73</sup>. In this method, genomic DNA is treated with sodium bisulfite and then sequenced, providing single-base resolution of methylated cytosines in the genome. Upon bisulfite treatment, unmethylated cytosines are deaminated to uracils which, upon sequencing, are converted to thymidines. Simultaneously, methylated cytosines resist deamination and are read as cytosines. The location of the methylated cytosines can then be determined by comparing treated and untreated sequences. Bisulfite treatment of DNA converts unmethylated cytosines to thymidines, leading to reduced sequence complexity. Very accurate deep sequencing serves to mitigate this loss of complexity

The EpiGnome™ Kit uses a unique library construction method that incorporates bisulfite conversion as the first step. The EpiGnome method retains sample diversity while providing uniform coverage.



### Pros

### Cons

#### BS-Seq or WGBS

- CpG and non-CpG methylation throughout the genome is covered at single-base resolution
- 5mC in dense, less dense, and repeat regions are covered

- Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments
- NPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion
- Bisulfite conversion does not distinguish between 5mC and 5hmC

#### EpiGnome

- Pre-library bisulfite conversion
- Low input gDNA (50 ng)
- Uniform CpG, CHG, and CHH coverage
- No fragmentation and no methylated adapters
- Retention of sample diversity

- Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments
- SNPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion
- Bisulfite conversion does not distinguish between 5mC and 5hmC
- Higher duplicate percentage

73 Feil R., Charlton J., Bird A. P., Walter J. and Reik W. (1994) Methylation analysis on individual chromosomes: improved protocol for bisulphite genomic sequencing. Nucleic Acids Res 22: 695-696

## References

Gustems M., Woellmer A., Rothbauer U., Eck S. H., Wieland T., et al. (2013) c-Jun/c-Fos heterodimers regulate cellular genes via a newly identified class of methylated DNA sequence motifs. *Nucleic Acids Res*

Transcription factors bind with specificity to their preferred DNA sequence motif. However, a virus-encoded transcription factor Zta was the first example of a sequence-specific transcription factor binding selectively and preferentially to methylated CpG residues. In this study the authors present their finding of a novel AP-1 binding site, termed meAP-1, which contains a CpG nucleotide. Using ChIP-Seq with Illumina sequencing, they show how the methylation state of this nucleotide affects binding by c-Jun/c-Fos in vitro and in vivo.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

Habibi E., Brinkman A. B., Arand J., Kroeze L. I., Kerstens H. H., et al. (2013) Whole-genome bisulfite sequencing of two distinct interconvertible DNA methylomes of mouse embryonic stem cells. *Cell Stem Cell* 13: 360-369

Mouse embryonic stem cells (ESCs) provide an excellent model system for studying mammalian cell differentiation on the molecular level. This study uses two kinase inhibitors (2i) to derive mouse ESCs in the pluripotent ground state to study the deposition and loss of DNA methylation during differentiation. The epigenetic state and expression of the cells were monitored using ChIP-Seq and RNA-Seq on the Illumina HiSeq platform.

**Illumina Technology:** HiSeq 2000, MiSeq

Hussain S., Sajini A. A., Blanco S., Dietmann S., Lombard P., et al. (2013) NSun2-mediated cytosine-5 methylation of vault noncoding RNA determines its processing into regulatory small RNAs. *Cell Rep* 4: 255-261

This paper presents miCLIP: a new technique for identifying RNA methylation sites in transcriptomes. The authors use the miCLIP method with Illumina sequencing to determine site-specific methylation in tRNAs and additional messenger and noncoding RNAs. As a case study, the authors studied the methyltransferase NSun2 and showed that loss of cytosine-5 methylation in vault RNAs causes aberrant processing that may interrupt processing of small RNA fragments, such as microRNAs.

**Illumina Technology:** TruSeq RNA Kit, Genome Analyzer<sub>ix</sub>

Kozlenkov A., Roussos P., Timashpolsky A., Barbu M., Rudchenko S., et al. (2014) Differences in DNA methylation between human neuronal and glial cells are concentrated in enhancers and non-CpG sites. *Nucleic Acids Res* 42: 109-127

Epigenetic regulation by DNA methylation varies among different cell types. In this study, the authors compared the methylation status of neuronal and non-neuronal nuclei using Illumina Human Methylation450k arrays. They classified the differentially methylated (DM) sites into those undermethylated in the neuronal cell type, and those that were undermethylated in non-neuronal cells. Using this approach, they identified sets of cell type-specific patterns and characterized these by their genomic locations.

**Illumina Technology:** HumanMethylation450 BeadChip, HumanOmni1-Quad (Infinium GT), HiSeq 2000

Lun F. M., Chiu R. W., Sun K., Leung T. Y., Jiang P., et al. (2013) Noninvasive prenatal methylomic analysis by genomewide bisulfite sequencing of maternal plasma DNA. *Clin Chem* 59: 1583-1594

The presence of fetal DNA in maternal plasma opens up possibilities for non-invasive prenatal DNA testing of the fetus through blood samples from the mother. Using SNP differences between mother and fetus to identify fetal molecules, this study inspected the genome-wide methylome of the unborn child by bisulfite sequencing. The authors determined the methylation density over each 1 Mbp region of the genome for samples taken in each trimester and after delivery to show how the fetal methylome is established gradually throughout pregnancy.

**Illumina Technology:** HiSeq 2000, HumanMethylation450 BeadChip

Regulski M., Lu Z., Kendall J., Donoghue M. T., Reinders J., et al. (2013) The maize methylome influences mRNA splice sites and reveals widespread paramutation-like switches guided by small RNA. *Genome Res* 23: 1651-1662

The maize genome encompasses a widely unexplored landscape for epigenetic mechanisms of paramutation and imprinting. In this study whole-exome bisulfite sequencing was applied to map the cytosine methylation profile of two maize inbred lines. The analysis revealed that frequent methylation switches, guided by siRNA, may persist for up to eight generations, suggesting that epigenetic inheritance resembling paramutation is much more common than previously supposed.

**Illumina Technology:** HiSeq 2000, Genome Analyzer<sub>ix</sub>

Schmitz R. J., He Y., Valdes-Lopez O., Khan S. M., Joshi T., et al. (2013) Epigenome-wide inheritance of cytosine methylation variants in a recombinant inbred population. *Genome Res* 23: 1663-1674

In an effort to elucidate the mammalian DNA methylome, this study applied whole-genome bisulfite sequencing using the Illumina HiSeq platform and gene expression analysis to define functional classes of hypomethylated regions (HMRs). Comparing HMR profiles in embryonic stem and primary blood cells, the authors showed that the HMRs in intergenic space (iHMRs) mark an exclusive subset of active DNase hypersensitive sites. The authors went on to compare primate-specific and human population variation at iHMRs, and they derived models of the cellular timelines for DHS and iHMR establishment.

**Illumina Technology:** HiSeq 2000

Schlesinger F., Smith A. D., Gingeras T. R., Hannon G. J. and Hodges E. (2013) De novo DNA demethylation and noncoding transcription define active intergenic regulatory elements. *Genome Res* 23: 1601-1614

In an effort to elucidate the mammalian DNA methylome, this study applied whole-genome bisulfite sequencing using the Illumina HiSeq platform and gene expression analysis to define functional classes of hypomethylated regions (HMRs). Comparing HMR profiles in embryonic stem and primary blood cells, the authors showed that the HMRs in intergenic space (iHMRs) mark an exclusive subset of active DNase hypersensitive sites. The authors went on to compare primate-specific and human population variation at iHMRs, and they derived models of the cellular timelines for DHS and iHMR establishment.

**Illumina Technology:** HiSeq 2000

Xie W., Schultz M. D., Lister R., Hou Z., Rajagopal N., et al. (2013) Epigenomic analysis of multilineage differentiation of human embryonic stem cells. *Cell* 153: 1134-1148

The authors studied the differentiation of hESCs into four cell types: trophoblast-like cells, mesendoderm, neural progenitor cells, and mesenchymal stem cells. DNA methylation (WGBS) and histone modifications were examined for each cell type. The study provides insight into the dynamic changes that accompany lineage-specific cell differentiation in hESCs.

**Illumina Technology:** HiSeq 2000

Yamaguchi S., Shen L., Liu Y., Sandler D. and Zhang Y. (2013) Role of Tet1 in erasure of genomic imprinting. *Nature* 504: 460-464

Genomic imprinting is the cellular mechanism for switching off one of two alleles by DNA methylation. This allele-specific gene expression system is very important for mammalian development and function. In this study, the Tet1 protein was studied for its function in primordial germ cells, the phase of development where the imprinting methylation mark of the parent is erased. Using ChIP-Seq and bisulfite sequencing on the Illumina HiSeq platform, the authors showed that Tet1 knockout males exhibited aberrant hypermethylation in the paternal allele of differential methylated regions.

**Illumina Technology:** HiSeq 2500®

Blaschke K., Ebata K. T., Karimi M. M., Zepeda-Martinez J. A., Goyal P., et al. (2013) Vitamin C induces Tet-dependent DNA demethylation and a blastocyst-like state in ES cells. *Nature* 500: 222-226

Potok M. E., Nix D. A., Parnell T. J. and Cairns B. R. (2013) Reprogramming the maternal zebrafish genome after fertilization to match the paternal methylation pattern. *Cell* 153: 759-772

Rodrigues J. A., Ruan R., Nishimura T., Sharma M. K., Sharma R., et al. (2013) Imprinted expression of genes and small RNA is associated with localized hypomethylation of the maternal genome in rice endosperm. *Proc Natl Acad Sci U S A* 110: 7934-7939

Shirane K., Toh H., Kobayashi H., Miura F., Chiba H., et al. (2013) Mouse oocyte methylomes at base resolution reveal genome-wide accumulation of non-CpG methylation and role of DNA methyltransferases. *PLoS Genet* 9: e1003439

Warden C. D., Lee H., Tompkins J. D., Li X., Wang C., et al. (2013) COHCAP: an integrative genomic pipeline for single-nucleotide resolution DNA methylation analysis. *Nucleic Acids Res* 41: e117

Adey A. and Shendure J. (2012) Ultra-low-input, tagmentation-based whole-genome bisulfite sequencing. *Genome Res* 22: 1139-1143

Diep D., Plongthongkum N., Gore A., Fung H. L., Shoemaker R., et al. (2012) Library-free methylation sequencing with bisulfite padlock probes. *Nat Methods* 9: 270-272

Seisenberger S., Andrews S., Krueger F., Arand J., Walter J., et al. (2012) The dynamics of genome-wide DNA methylation reprogramming in mouse primordial germ cells. *Mol Cell* 48: 849-862

Feng S., Cokus S. J., Zhang X., Chen P. Y., Bostick M., et al. (2010) Conservation and divergence of methylation patterning in plants and animals. *Proc Natl Acad Sci U S A* 107: 8689-8694

Li N., Ye M., Li Y., Yan Z., Butcher L. M., et al. (2010) Whole genome DNA methylation analysis based on high throughput sequencing technology. *Methods* 52: 203-212

Lyko F., Foret S., Kucharski R., Wolf S., Falckenhayn C., et al. (2010) The honey bee epigenomes: differential methylation of brain DNA in queens and workers. *PLoS Biol* 8: e1000506



- Ball M. P., Li J. B., Gao Y., Lee J. H., LeProust E. M., et al. (2009) Targeted and genome-scale strategies reveal gene-body methylation signatures in human cells. *Nat Biotechnol* 27: 361-368
- Gehring M., Bubb K. L. and Henikoff S. (2009) Extensive demethylation of repetitive elements during seed development underlies gene imprinting. *Science* 324: 1447-1451
- Hodges E., Smith A. D., Kendall J., Xuan Z., Ravi K., et al. (2009) High definition profiling of mammalian DNA methylation by array capture and single molecule bisulfite sequencing. *Genome Res* 19: 1593-1605
- Hsieh T. F., Ibarra C. A., Silva P., Zemach A., Eshed-Williams L., et al. (2009) Genome-wide demethylation of Arabidopsis endosperm. *Science* 324: 1451-1454
- Jacob Y., Feng S., Leblanc C. A., Bernatavichute Y. V., Stroud H., et al. (2009) ATXR5 and ATXR6 are H3K27 monomethyltransferases required for chromatin structure and gene silencing. *Nat Struct Mol Biol* 16: 763-768
- Cokus S. J., Feng S., Zhang X., Chen Z., Merriman B., et al. (2008) Shotgun bisulphite sequencing of the Arabidopsis genome reveals DNA methylation patterning. *Nature* 452: 215-219
- He Y., Vogelstein B., Velculescu V. E., Papadopoulos N. and Kinzler K. W. (2008) The antisense transcriptomes of human cells. *Science* 322: 1855-1857
- Meissner A., Mikkelsen T. S., Gu H., Wernig M., Hanna J., et al. (2008) Genome-scale DNA methylation maps of pluripotent and differentiated cells. *Nature* 454: 766-770

#### **Associated Kits**

EpiGnome™ Methyl-Seq® Kit

Infinium HumanMethylation450® Arrays

## POST-BISULFITE ADAPTER TAGGING (PBAT)

To avoid the bisulfite-induced loss of intact sequencing templates, in post-bisulfite adaptor tagging (PBAT)<sup>74</sup> bisulfite treatment is followed by adaptor tagging and two rounds of random primer extension. This procedure generates a substantial number of unamplified reads from as little as subnanogram quantities of DNA.



Pros	Cons
<ul style="list-style-type: none"><li>• Requires only 100 ng of DNA for amplification-free WGBS of mammalian genomes</li></ul>	<ul style="list-style-type: none"><li>• Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments</li><li>• SNPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion</li><li>• Bisulfite conversion does not distinguish between 5mC and 5hmC</li></ul>

### References

Kobayashi H., Sakurai T., Miura F., Imai M., Mochiduki K., et al. (2013) High-resolution DNA methylome analysis of primordial germ cells identifies gender-specific reprogramming in mice. *Genome Res* 23: 616-627

Dynamic epigenetic reprogramming occurs during mammalian germ cell development. One of these processes is DNA methylation and demethylation, which is commonly studied using bisulfite sequencing. This study used an Illumina HiSeq 2000 system for WGBS to characterize the DNA methylation profiles of male and female mouse primordial germ cells (PGCs) at different stages of embryonic development. The authors found sex- and chromosome-specific differences in genome-wide CpG and CGI methylation during early- to late-stage PGC development. They also obtained high-resolution details of DNA methylation changes, for instance, that LINE/LTR retrotransposons were resistant to DNA methylation at high CpG densities.

**Illumina Technology:** HiSeq 2000

74 Miura F., Enomoto Y., Dairiki R. and Ito T. (2012) Amplification-free whole-genome bisulfite sequencing by post-bisulfite adaptor tagging. *Nucleic Acids Res* 40: e136

Shirane K., Toh H., Kobayashi H., Miura F., Chiba H., et al. (2013) Mouse oocyte methylomes at base resolution reveal genome-wide accumulation of non-CpG methylation and role of DNA methyltransferases. *PLoS Genet* 9: e1003439

DNA methylation is an epigenetic modification that plays a crucial role in normal mammalian development, retrotransposon silencing, and cellular reprogramming. Using amplification-free WGBS, the authors constructed the base-resolution methylome maps of germinal vesicle oocytes (GVOs), non-growing oocytes (NGOs), and mutant GVOs lacking the DNA methyltransferases Dnmt1, Dnmt3a, Dnmt3b, or Dnmt3L. They found that nearly two-thirds of all methylcytosines occur in a non-CG context in GVOs. The distribution of non-CG methylation closely resembled that of CG methylation throughout the genome and showed clear enrichment in gene bodies.

**Illumina Technology:** HiSeq 2000

Kobayashi H., Sakurai T., Imai M., Takahashi N., Fukuda A., et al. (2012) Contribution of Intragenic DNA Methylation in Mouse Gametic DNA Methylomes to Establish Oocyte-Specific Heritable Marks. *PLoS Genet* 8: e1002440

Miura F., Enomoto Y., Dairiki R. and Ito T. (2012) Amplification-free whole-genome bisulfite sequencing by post-bisulfite adaptor tagging. *Nucleic Acids Res* 40: e136


#### **Associated Kits**

EpiGenome™ Methyl-Seq Kit

Infinium HumanMethylation450 Arrays

## TAGMENTATION-BASED WHOLE GENOME BISULFITE SEQUENCING (T-WGBS)

Tagmentation-based whole-genome bisulfite sequencing (T-WGBS) is a protocol that utilizes the Epicentre® Tn5 transposome and bisulfite conversion to study 5mC<sup>75</sup>. In this method, DNA is incubated with Tn5 transposome containing methylated primers, which fragments the DNA and ligates adapters. Tagged DNA first undergoes oligo displacement, followed by methylated oligo replacement and gap repair, assuring methylated adapter addition to tagmented DNA. DNA is then treated with sodium bisulfite, PCR-amplified, and sequenced. Deep sequencing provides single-base resolution of 5mC in the genome.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can sequence samples with very limited starting material (~20 ng)</li> <li>• Fast protocol with few steps</li> <li>• Elimination of multiple steps prevents loss of DNA</li> </ul>	<ul style="list-style-type: none"> <li>• Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments</li> <li>• SNPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion</li> <li>• Bisulfite conversion does not distinguish between 5mC and 5hmC</li> </ul>

### References

Wang Q., Gu L., Adey A., Radlwimmer B., Wang W., et al. (2013) Tagmentation-based whole-genome bisulfite sequencing. Nat Protoc 8: 2022-2032

Scaling up bisulfite sequencing to genome-wide analysis has been hindered by the requirements for large amounts of DNA and high sequencing costs. This paper presents a protocol for T-WGBS with sequencing on the Illumina HiSeq 2000 system. The authors demonstrate the robustness of the protocol in comparison with conventional WGBS. T-WGBS requires not more than 20 ng of input DNA; hence, the protocol allows the comprehensive methylome analysis of limited amounts of DNA isolated from precious biological specimens.

**Illumina Technology:** Nextera DNA Sample Prep Kit, HiSeq 2000; 101 bp paired-end reads

### Associated Kits

EpiGnome™ Methyl-Seq Kit

Infinium HumanMethylation450 Arrays

Nextera DNA Sample Prep Kit


Nextera XT DNA Sample Prep Kit

Nextera Rapid Capture Exome/Custom Enrichment Kit

<sup>75</sup> Wang Q., Gu L., Adey A., Radlwimmer B., Wang W., et al. (2013) Tagmentation-based whole-genome bisulfite sequencing. Nat Protoc 8: 2022-2032

## OXIDATIVE BISULFITE SEQUENCING (OXBS-SEQ)

Oxidative bisulfite sequencing (oxBS-Seq) differentiates between 5mC and 5hmC<sup>76</sup>. With oxBS, 5hmC is oxidized to 5formylcytosine (5fC) with an oxidizing agent, while 5mC remains unchanged. Sodium bisulfite treatment of oxidized 5hmC results in its deamination to uracil which, upon sequencing, is read as a thymidine. Deep sequencing of oxBS-treated DNA and sequence comparison of treated vs. untreated can identify 5mC locations at base resolution.



Pros	Cons
<ul style="list-style-type: none"> <li>• CpG and non-CpG methylation throughout the genome is covered at single-base resolution</li> <li>• 5mC dense and less dense in repeat regions are covered</li> <li>• Method clearly differentiates between 5mC and 5hmC, precisely identifying 5mC</li> </ul>	<ul style="list-style-type: none"> <li>• Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments</li> <li>• SNPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion</li> </ul>

### References

Booth M. J., Ost T. W., Beraldi D., Bell N. M., Branco M. R., et al. (2013) Oxidative bisulfite sequencing of 5-methylcytosine and 5-hydroxymethylcytosine. Nat Protoc 8: 1841-1851

This is the first paper to report a method combining chemical treatment of DNA with the well-established bisulfite protocol, highlighting Illumina's TruSeq kit and calling for the use of MiSeq or HiSeq platforms. The OxBS-Seq protocol helps distinguish between 5mC and 5hmC, while standard bisulfite sequencing is incapable of distinguishing between 5mC and 5hmC. Genomic DNA is first treated with an oxidizing agent that reacts with 5hmC, promoting its deamination to uracil, while the 5mC modification remains unchanged and is read as cytosine. Using Illumina technology, this method allows base resolution of the exact location of 5hmC and 5mC modifications.

**Illumina Technology:** TruSeq DNA Sample Prep Kit, MiSeq, HiSeq 2000

### Associated Kits

EpiGnome™ Methyl-Seq Kit

TruSeq DNA Sample Prep Kit


TruSeq DNA PCR-Free Sample Prep Kit

TruSeq Nano DNA Sample Prep Kit

76 Booth M. J., Branco M. R., Ficiz G., Oxley D., Krueger F., et al. (2012) Quantitative sequencing of 5-methylcytosine and 5-hydroxymethylcytosine at single-base resolution. Science 336: 934-937

## TET-ASSISTED BISULFITE SEQUENCING (TAB-SEQ)

TAB-Seq is a novel method that uses bisulfite conversion and Tet proteins to study 5hmC<sup>77</sup>. In this protocol, 5hmC is first protected with a glucose moiety that allows selective interaction and subsequent oxidation of 5mC with the Tet proteins. The oxidized genomic DNA is then treated with bisulfite, where 5hmC remains unchanged and is read as a cytosine, while 5mC and unmethylated cytosines are deaminated to uracil and read as thymidines upon sequencing. Deep sequencing of TAB-treated DNA compared with untreated DNA provides accurate representation of 5hmC localization in the genome.



Pros	Cons
<ul style="list-style-type: none"> <li>• CpG and non-CpG hydroxymethylation throughout the genome is covered at single-base resolution</li> <li>• Dense, less dense, and 5hmC in repeat regions are covered</li> <li>• Method clearly differentiates between 5hmC and 5mC, specifically identifying 5hmC</li> </ul>	<ul style="list-style-type: none"> <li>• Bisulfite converts unmethylated cytosines to thymidines, reducing sequence complexity, which can make it difficult to create alignments</li> <li>• SNPs where a cytosine is converted to thymidine will be missed upon bisulfite conversion</li> <li>• Requires deep sequencing to provide sufficient depth to cover the entire genome and accurately map the low amounts 5hmC<sup>78</sup></li> </ul>

## References

Kim M., Park Y. K., Kang T. W., Lee S. H., Rhee Y. H., et al. (2013) Dynamic changes in DNA methylation and hydroxymethylation when hES cells undergo differentiation toward a neuronal lineage. *Hum Mol Genet* 23: 657-667

Epigenetic markers on chromatin include the methylation of DNA. Several forms of DNA methylation exist and their function and interaction is a field of intensive study. This paper describes how an in vitro model system of gradual differentiation of hESCs underwent dramatic genome-wide changes in 5mC and 5hmC methylation patterns during lineage commitment. The authors used Illumina BeadArray for expression profiling and Genome Analyzer hMeDIP-sequencing to study the correlation between gene expression and DNA methylation.

**Illumina Technology:** Human-6 Whole-Genome Expression BeadChip, Genome Analyzer<sub>ix</sub>, HiScanSQ<sup>®</sup> Scanner, Infinium HumanMethylation 450 BeadChip

<sup>77</sup> Yu M., Hon G. C., Szulwach K. E., Song C. X., Zhang L., et al. (2012) Base-resolution analysis of 5-hydroxymethylcytosine in the Mammalian genome. *Cell* 149: 1368-1380

<sup>78</sup> Thomson J. P., Hunter J. M., Nestor C. E., Dunican D. S., Terranova R., et al. (2013) Comparative analysis of affinity-based 5-hydroxymethylation enrichment techniques. *Nucleic Acids Res* 41: e206

Lister R., Mukamel E. A., Nery J. R., Urich M., Puddifoot C. A., et al. (2013) Global epigenomic reconfiguration during mammalian brain development. *Science* 341: 1237905

DNA methylation is implicated in mammalian brain development and plasticity underlying learning and memory. This paper reports the genome-wide composition, patterning, cell specificity, and dynamics of DNA methylation at single-base resolution in human and mouse frontal cortex throughout their lifespan. The extensive methylome profiling was performed with CHIP-Seq on an Illumina HiSeq sequencer at single-base resolution.

**Illumina Technology:** TruSeq RNA Sample Prep Kit, TruSeq DNA Sample Prep Kit, HiSeq 2000

Wang T., Wu H., Li Y., Szulwach K. E., Lin L., et al. (2013) Subtelomeric hotspots of aberrant 5-hydroxymethylcytosine-mediated epigenetic modifications during reprogramming to pluripotency. *Nat Cell Biol* 15: 700-711

The transcriptional reprogramming that allows mammalian somatic cells to be reprogrammed into pluripotent stem cells (iPSCs) includes a complete reconfiguration of the epigenetic marks in the genome. This study examined the levels of 5hmC in hESCs during reprogramming to iPSCs. The authors found reprogramming hotspots in subtelomeric regions, most of which featured incomplete hydroxymethylation at CG sites.

**Illumina Technology:** HiSeq 2000, HiScanSQ, MiSeq

Jiang L., Zhang J., Wang J. J., Wang L., Zhang L., et al. (2013) Sperm, but not oocyte, DNA methylome is inherited by zebrafish early embryos. *Cell* 153: 773-784

Song C. X., Szulwach K. E., Dai Q., Fu Y., Mao S. Q., et al. (2013) Genome-wide profiling of 5-formylcytosine reveals its roles in epigenetic priming. *Cell* 153: 678-691

Yu M., Hon G. C., Szulwach K. E., Song C. X., Jin P., et al. (2012) Tet-assisted bisulfite sequencing of 5-hydroxymethylcytosine. *Nat Protoc* 7: 2159-2170

Yu M., Hon G. C., Szulwach K. E., Song C. X., Zhang L., et al. (2012) Base-resolution analysis of 5-hydroxymethylcytosine in the mammalian genome. *Cell* 149: 1368-1380


#### Associated Kits

EpiGnome™ Methyl-Seq Kit

Infinium HumanMethylation450 Arrays

## METHYLATED DNA IMMUNOPRECIPITATION SEQUENCING (MEDIP-SEQ)

Methylated DNA immunoprecipitation sequencing (MeDIP-Seq) is commonly used to study 5mC or 5hmC modification<sup>79</sup>. Specific antibodies can be used to study cytosine modifications. If using 5mC-specific antibodies, methylated DNA is isolated from genomic DNA via immunoprecipitation. Anti-5mC antibodies are incubated with fragmented genomic DNA and precipitated, followed by DNA purification and sequencing. Deep sequencing provides greater genome coverage, representing the majority of immunoprecipitated methylated DNA.



Pros	Cons
<ul style="list-style-type: none"> <li>Covers CpG and non-CpG 5mC throughout the genome</li> <li>5mC in dense, less dense, and repeat regions are covered</li> <li>Antibody-based selection is independent of sequence and does not enrich for 5hmC due to antibody specificity</li> </ul>	<ul style="list-style-type: none"> <li>Base-pair resolution is lower (~150 bp) as opposed to single base resolution</li> <li>Antibody specificity and selectivity must be tested to avoid nonspecific interaction</li> <li>Antibody-based selection is biased towards hypermethylated regions</li> </ul>

### References

Puszyk W., Down T., Grimwade D., Chomienne C., Oakey R. J., et al. (2013) The epigenetic regulator PLZF represses L1 retrotransposition in germ and progenitor cells. *EMBO J* 32: 1941-1952

Each transcription factor in the human cell may regulate a large number of target genes through specific chromatin interactions. Promyelocytic leukemia zinc finger protein (PLZF) acts as an epigenetic regulator of stem cell maintenance in germ cells and hematopoietic stem cells. In this study, L1 retrotransposons were identified as the primary targets of PLZF. Using ChIP-Seq and MeDIP-Seq on Illumina Genome Analyzer, the authors identified how PLZF-mediated DNA methylation induces silencing of L1 and inhibits L1 retrotransposition.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Shen H., Qiu C., Li J., Tian Q. and Deng H. W. (2013) Characterization of the DNA methylome and its interindividual variation in human peripheral blood monocytes. *Epigenomics* 5: 255-269

Peripheral blood monocytes (PBMs) play multiple and critical roles in the immune response, and abnormalities in PBMs have been linked to a variety of human disorders. In this study, the epigenome-wide DNA methylation profiles of purified PBMs were identified using MeDIP-Seq on an Illumina Genome Analyzer. Interestingly, the authors observed substantial interindividual variation in DNA methylation across the individual PBM methylomes.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>79</sup> Weber M., Davies J. J., Wittig D., Oakeley E. J., Haase M., et al. (2005) Chromosome-wide and promoter-specific analyses identify sites of differential DNA methylation in normal and transformed human cells. *Nat Genet* 37: 853-862



Tan L., Xiong L., Xu W., Wu F., Huang N., et al. (2013) Genome-wide comparison of DNA hydroxymethylation in mouse embryonic stem cells and neural progenitor cells by a new comparative hMeDIP-seq method. *Nucleic Acids Res* 41: e84

The genome-wide distribution patterns of the “sixth base” 5hmC in many tissues and cells have recently been revealed by hydroxymethylated DNA immunoprecipitation (hMeDIP) followed by high throughput sequencing or tiling arrays. This paper presents a new comparative hMeDIP-seq method which allows for direct genome-wide comparison of DNA hydroxymethylation across multiplesamples. The authors demonstrate the new method by profiling DNA hydroxymethylation and gene expression during neural differentiation.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Saied M. H., Marzec J., Khalid S., Smith P., Down T. A., et al. (2012) Genome wide analysis of acute myeloid leukemia reveal leukemia specific methylome and subtype specific hypomethylation of repeats. *PLoS One* 7: e33213

Epigenetic modifications in the form of DNA methylation are part of the regulatory machinery of the cell. By studying the patterns of DNA methylation in disease tissue, we may characterize disease mechanisms. In this study, bone marrow samples from 12 patients with acute myeloid leukemia (AML) were analyzed with MeDIP-Seq and compared to normal bone marrow. The investigators found considerable cytogenetic subtype specificity in the methylomes affecting different genomic features.

**Illumina Technology:** HumanMethylation27 arrays, Genome Analyzer<sub>ix</sub>

Taiwo O., Wilson G. A., Morris T., Seisenberger S., Reik W., et al. (2012) Methylome analysis using MeDIP-seq with low DNA concentrations. *Nat Protoc* 7: 617-636

DNA methylation can be assayed at high throughput using MeDIP-Seq, but the application has been limited to samples where the amount of DNA was sufficient for the assay (5–20 µg). This study presents a new optimized protocol for MeDIP-Seq, requiring as little as 50 ng of starting DNA.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Bian C. and Yu X. (2013) PGC7 suppresses TET3 for protecting DNA methylation. *Nucleic Acids Res*

Colquitt B. M., Allen W. E., Barnea G. and Lomvardas S. (2013) Alteration of genic 5-hydroxymethylcytosine patterning in olfactory neurons correlates with changes in gene expression and cell identity. *Proc Natl Acad Sci U S A* 110: 14682-14687

Neri F., Krepelova A., Incarnato D., Maldotti M., Parlato C., et al. (2013) Dnmt3L Antagonizes DNA Methylation at Bivalent Promoters and Favors DNA Methylation at Gene Bodies in ESCs. *Cell* 155: 121-134

Stevens M., Cheng J. B., Li D., Xie M., Hong C., et al. (2013) Estimating absolute methylation levels at single-CpG resolution from methylation enrichment and restriction enzyme sequencing methods. *Genome Res* 23: 1541-1553

Zhang B., Zhou Y., Lin N., Lowdon R. F., Hong C., et al. (2013) Functional DNA methylation differences between tissues, cell types, and across individuals discovered using the M&M algorithm. *Genome Res* 23: 1522-1540

Zilbauer M., Rayner T. F., Clark C., Coffey A. J., Joyce C. J., et al. (2013) Genome-wide methylation analyses of primary human leukocyte subsets identifies functionally important cell-type-specific hypomethylated regions. *Blood* 122: e52-60


- Sati S., Tanwar V. S., Kumar K. A., Patowary A., Jain V., et al. (2012) High resolution methylome map of rat indicates role of intragenic DNA methylation in identification of coding region. *PLoS One* 7: e31621
- Gao Q., Steine E. J., Barrasa M. I., Hockemeyer D., Pawlak M., et al. (2011) Deletion of the de novo DNA methyltransferase Dnmt3a promotes lung tumor progression. *Proc Natl Acad Sci U S A* 108: 18061-18066
- Bock C., Tomazou E. M., Brinkman A. B., Muller F., Simmer F., et al. (2010) Quantitative comparison of genome-wide DNA methylation mapping technologies. *Nat Biotechnol* 28: 1106-1114
- Chavez L., Jozefczuk J., Grimm C., Dietrich J., Timmermann B., et al. (2010) Computational analysis of genome-wide DNA methylation during the differentiation of human embryonic stem cells along the endodermal lineage. *Genome Res* 20: 1441-1450
- Harris R. A., Wang T., Coarfa C., Nagarajan R. P., Hong C., et al. (2010) Comparison of sequencing-based methods to profile DNA methylation and identification of monoallelic epigenetic modifications. *Nat Biotechnol* 28: 1097-1105
- Li N., Ye M., Li Y., Yan Z., Butcher L. M., et al. (2010) Whole genome DNA methylation analysis based on high throughput sequencing technology. *Methods* 52: 203-212
- Maunakea A. K., Nagarajan R. P., Bilenky M., Ballinger T. J., D'Souza C., et al. (2010) Conserved role of intragenic DNA methylation in regulating alternative promoters. *Nature* 466: 253-257
- Ruike Y., Imanaka Y., Sato F., Shimizu K. and Tsujimoto G. (2010) Genome-wide analysis of aberrant methylation in human breast cancer cells using methyl-DNA immunoprecipitation combined with high-throughput sequencing. *BMC Genomics* 11: 137
- Hammoud S. S., Nix D. A., Zhang H., Purwar J., Carrell D. T., et al. (2009) Distinctive chromatin in human sperm packages genes for embryo development. *Nature* 460: 473-478
- Pomraning K. R., Smith K. M. and Freitag M. (2009) Genome-wide high throughput analysis of DNA methylation in eukaryotes. *Methods* 47: 142-150
- Down T. A., Rakyan V. K., Turner D. J., Flicek P., Li H., et al. (2008) A Bayesian deconvolution strategy for immunoprecipitation-based DNA methylome analysis. *Nat Biotechnol* 26: 779-785

### **Associated Kits**

- Infinium HumanMethylation450 Arrays
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Rapid Capture Exome/Custom Enrichment Kit

## METHYLATION-CAPTURE (METHYLCAP) SEQUENCING OR METHYL-BINDING-DOMAIN-CAPTURE (MBDCAP) SEQUENCING

MethylCap<sup>80, 81</sup> or MBDCap<sup>82, 83</sup> uses proteins to capture methylated DNA in the genome. Genomic DNA is first sonicated and incubated with tagged MBD proteins that can bind methylated cytosines. The protein-DNA complex is then precipitated with antibody-conjugated beads that are specific to the protein tag. Deep sequencing provides greater genome coverage, representing the majority of MBD-bound methylated DNA.



Pros	Cons
<ul style="list-style-type: none"> <li>• Genome-wide coverage of 5mC in dense CpG areas and repeat regions</li> <li>• MBD proteins do not interact with 5hmC</li> </ul>	<ul style="list-style-type: none"> <li>• Genome-wide CpGs and non-CpG methylation is not covered Areas with less dense 5mC are also missed</li> <li>• Base-pair resolution is lower (~150 bp) as opposed to single base resolution</li> <li>• Protein-based selection is biased towards hypermethylated regions</li> </ul>

### References

Kim M., Park Y. K., Kang T. W., Lee S. H., Rhee Y. H., et al. (2013) Dynamic changes in DNA methylation and hydroxymethylation when hES cells undergo differentiation toward a neuronal lineage. *Hum Mol Genet* 23: 657-667

Epigenetic markers on chromatin include the methylation of DNA. Several forms of DNA methylation exist and their function and interaction is a field of intensive study. This paper describes how an in vitro model system of gradual differentiation of hESCs underwent dramatic genome-wide changes in 5mC and 5hmC methylation patterns during lineage commitment. The authors used Illumina BeadArray for expression profiling and Genome Analyzer hMedIP-sequencing to study the correlation between gene expression and DNA methylation.

**Illumina Technology:** Human-6 Whole-Genome Expression BeadChip, Genome Analyzer<sub>ix</sub>, HiScanSQ Scanner, Infinium HumanMethylation 450 BeadChip

80 Bock C., Tomazou E. M., Brinkman A. B., Muller F., Simmer F., et al. (2010) Quantitative comparison of genome-wide DNA methylation mapping technologies. *Nat Biotechnol* 28: 1106-1114

81 Brinkman A. B., Simmer F., Ma K., Kaan A., Zhu J., et al. (2010) Whole-genome DNA methylation profiling using MethylCap-seq. *Methods* 52: 232-236

82 Rauch T. A., Zhong X., Wu X., Wang M., Kernstine K. H., et al. (2008) High-resolution mapping of DNA hypermethylation and hypomethylation in lung cancer. *Proc Natl Acad Sci U S A* 105: 252-257

83 Rauch T. A. and Pfeifer G. P. (2009) The MIRA method for DNA methylation analysis. *Methods Mol Biol* 507: 65-75

Huang T. T., Gonzales C. B., Gu F., Hsu Y. T., Jadhav R. R., et al. (2013) Epigenetic deregulation of the anaplastic lymphoma kinase gene modulates mesenchymal characteristics of oral squamous cell carcinomas. *Carcinogenesis* 34: 1717-1727

Promoter methylation is associated with silencing tumor suppressor genes in oral squamous cell carcinomas (OSCCs). The authors used MBDCap-Seq to study methylation in OSCC cell lines, sequencing on the Illumina HiSeq platform, and identifying differentially methylated regions. The authors note the ALK gene was susceptible to epigenetic silencing during oral tumorigenesis.

**Illumina Technology:** HiSeq 2000

Zhao Y., Guo S., Sun J., Huang Z., Zhu T., et al. (2012) Methylcap-seq reveals novel DNA methylation markers for the diagnosis and recurrence prediction of bladder cancer in a Chinese population. *PLoS ONE* 7: e35175

Bladder cancer (BC) has a high mortality rate and is the sixth most common cancer in the world. For successfully treated BCs, the relapse rate is 60-70% within the first 5 years, necessitating the development of efficient diagnostics and biomarkers for monitoring disease progression. The presence of cells in the urine allow for noninvasive genetic screening directly from urine. In this study, the authors identify and validate nine DNA methylation markers through genome-wide profiling of DNA methylation from clinical urine samples.

**Illumina Technology:** Genome Analyzer<sub>IIx</sub>

Brinkman A. B., Gu H., Bartels S. J., Zhang Y., Matarese F., et al. (2012) Sequential ChIP-bisulfite sequencing enables direct genome-scale investigation of chromatin and DNA methylation cross-talk. *Genome Res* 22: 1128-1138

Rodriguez B. A., Frankhouser D., Murphy M., Trimarchi M., Tam H. H., et al. (2012) Methods for high-throughput MethylCap-Seq data analysis. *BMC Genomics* 13 Suppl 6: S14

Yu W., Jin C., Lou X., Han X., Li L., et al. (2011) Global analysis of DNA methylation by Methyl-Capture sequencing reveals epigenetic control of cisplatin resistance in ovarian cancer cell. *PLoS One* 6: e29450

Bock C., Tomazou E. M., Brinkman A. B., Muller F., Simmer F., et al. (2010) Quantitative comparison of genome-wide DNA methylation mapping technologies. *Nat Biotechnol* 28: 1106-1114

Brinkman A. B., Simmer F., Ma K., Kaan A., Zhu J., et al. (2010) Whole-genome DNA methylation profiling using MethylCap-seq. *Methods* 52: 232-236

Serre D., Lee B. H. and Ting A. H. (2010) MBD-isolated Genome Sequencing provides a high-throughput and comprehensive survey of DNA methylation in the human genome. *Nucleic Acids Res* 38: 391-399

#### Associated Kits

Infinium HumanMethylation450 Arrays


Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Rapid Capture Exome/Custom Enrichment

## REDUCED-REPRESENTATION BISULFITE SEQUENCING (RRBS-SEQ)

Reduced-representation bisulfite sequencing (RRBS-Seq) is a protocol that uses one or multiple restriction enzymes on the genomic DNA to produce sequence-specific fragmentation<sup>84</sup>. The fragmented genomic DNA is then treated with bisulfite and sequenced. This is the method of choice to study specific regions of interest. It is particularly effective where methylation is high, such as in promoters and repeat regions.



Pros	Cons
<ul style="list-style-type: none"> <li>• Genome-wide coverage of CpGs in islands at single-base resolution</li> <li>• Areas dense in CpG methylation are covered</li> </ul>	<ul style="list-style-type: none"> <li>• Restriction enzymes cut at specific sites, providing biased sequence selection</li> <li>• Method measures 10-15% of all CpGs in genome</li> <li>• Cannot distinguish between 5mC and 5hmC</li> <li>• Does not cover non-CpG areas, genome-wide CpGs, and CpGs in areas without the enzyme restriction site</li> </ul>

### References

Kozlenkov A., Roussos P., Timashpolsky A., Barbu M., Rudchenko S., et al. (2014) Differences in DNA methylation between human neuronal and glial cells are concentrated in enhancers and non-CpG sites. *Nucleic Acids Res* 42: 109-127

Epigenetic regulation by DNA methylation varies among different cell types. In this study, the authors compared the methylation status of neuronal and non-neuronal nuclei using Illumina Human Methylation450k arrays. They classified the differentially methylated (DM) sites into those undermethylated in the neuronal cell type, and those that were undermethylated in non-neuronal cells. Using this approach, they identified sets of cell type-specific patterns and characterized these by their genomic locations.

**Illumina Technology:** HumanMethylation450 BeadChip, HumanOmni1-Quad (Infinium GT), HiSeq 2000

Schillebeeckx M., Schrade A., Lobs A. K., Pihlajoki M., Wilson D. B., et al. (2013) Laser capture microdissection-reduced representation bisulfite sequencing (LCM-RRBS) maps changes in DNA methylation associated with gonadectomy-induced adrenocortical neoplasia in the mouse. *Nucleic Acids Res* 41: e116

DNA methylation profiling by sequencing is challenging due to inaccurate cell enrichment methods and low DNA yields. This proof-of-concept study presents a new method for genome-wide DNA methylation profiling using down to 1 ng of input DNA. The method—laser-capture microdissection reduced-representation bisulfite sequencing (LCM-RRBS)—combines Illumina HiSeq sequencing with customized methylated adapter sequences and bisulfite-PCR. The protocol allows for base-pair resolution of methylated sites.

**Illumina Technology:** HiSeq 2000, MiSeq

<sup>84</sup> Meissner A., Gnirke A., Bell G. W., Ramsahoye B., Lander E. S., et al. (2005) Reduced representation bisulfite sequencing for comparative high-resolution DNA methylation analysis. *Nucleic Acids Res* 33: 5868-5877

Stevens M., Cheng J. B., Li D., Xie M., Hong C., et al. (2013) Estimating absolute methylation levels at single-CpG resolution from methylation enrichment and restriction enzyme sequencing methods. *Genome Res* 23: 1541-1553

Current methods for sequencing-based DNA methylation profiling are continuously improving, but each common method, on its own, is insufficient in providing a genome-wide single-CpG resolution of DNA methylation at a low cost. In this paper the authors present a novel algorithm, methylCRF, which enables integration of data from MeDIP-Seq and MRE-Seq to provide single-CpG classification of methylation state. The method provides similar or higher accuracy than any array or sequencing method on its own. The authors demonstrate the algorithm on whole-genome bisulfite sequencing on Illumina HiSeq 2000 systems and Methylation450 arrays.

**Illumina Technology:** HumanMethylation450 BeadChip, HumanOmni1-Quad (Infinium GT), HiSeq 2000

Will B., Vogler T. O., Bartholdy B., Garrett-Bakelman F., Mayer J., et al. (2013) *Satb1* regulates the self-renewal of hematopoietic stem cells by promoting quiescence and repressing differentiation commitment. *Nat Immunol* 14: 437-445

This study evaluated genome-wide DNA cytosine methylation by enhanced reduced-representation bisulfite sequencing (ERRBS). DNA was digested with *MspI*, then end-repaired and ligated to paired-end Illumina sequencing adapters. This was followed by size selection and bisulfite treatment, clean-up, and PCR prior to sequencing.

**Illumina Technology:** HiSeq 2000

Xi Y., Bock C., Muller F., Sun D., Meissner A., et al. (2012) RRBSMAP: a fast, accurate and user-friendly alignment tool for reduced representation bisulfite sequencing. *Bioinformatics* 28: 430-432

Bock C., Kiskinis E., Verstappen G., Gu H., Boulting G., et al. (2011) Reference Maps of human ES and iPS cell variation enable high-throughput characterization of pluripotent cell lines. *Cell* 144: 439-452

Gu H., Smith Z. D., Bock C., Boyle P., Gnirke A., et al. (2011) Preparation of reduced representation bisulfite sequencing libraries for genome-scale DNA methylation profiling. *Nat Protoc* 6: 468-481

Gertz J., Varley K. E., Reddy T. E., Bowling K. M., Pauli F., et al. (2011) Analysis of DNA methylation in a three-generation family reveals widespread genetic influence on epigenetic regulation. *PLoS Genet* 7: e1002228

Smallwood S. A., Tomizawa S., Krueger F., Ruf N., Carli N., et al. (2011) Dynamic CpG island methylation landscape in oocytes and preimplantation embryos. *Nat Genet* 43: 811-814

Ziller M. J., Muller F., Liao J., Zhang Y., Gu H., et al. (2011) Genomic distribution and inter-sample variation of non-CpG methylation across human cell types. *PLoS Genet* 7: e1002389

Bock C., Tomazou E. M., Brinkman A. B., Muller F., Simmer F., et al. (2010) Quantitative comparison of genome-wide DNA methylation mapping technologies. *Nat Biotechnol* 28: 1106-1114

Smith Z. D., Gu H., Bock C., Gnirke A. and Meissner A. (2009) High-throughput bisulfite sequencing in mammalian genomes. *Methods* 48: 226-232

#### Associated Kits

EpiGenome™ Methyl-Seq Kit

Infinium HumanMethylation450 Arrays

Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Rapid Capture Exome/Custom Enrichment Kit

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit

Chromatin remodeling is a dynamic process driven by factors that change DNA-protein interactions. These epigenetic factors can involve protein modifications, such as histone methylation, acetylation, phosphorylation, and ubiquitination<sup>85</sup>. Histone modifications determine gene activation by recruiting regulatory factors and maintaining an open or closed chromatin state. Epigenetic factors play roles in tissue development<sup>86</sup>, embryogenesis, cell fate, immune response, and diseases such as cancer<sup>87</sup>. Bacterial pathogens can elicit transcriptional repression of immune genes by chromatin remodeling<sup>88</sup>. The study of protein-DNA interactions has also demonstrated that chromatin remodeling can respond to external factors such as excessive alcohol-seeking behaviors<sup>89</sup>, cigarette smoking<sup>90</sup>, and clinical drugs.



**Cigarette smoking disrupts DNA-protein interactions leading to the development of cancers or pulmonary diseases.**

### Reviews

Capell B. C. and Berger S. L. (2013) Genome-wide epigenetics. *J Invest Dermatol* 133: e9

Jakopovic M., Thomas A., Balasubramaniam S., Schrupp D., Giaccone G., et al. (2013) Targeting the Epigenome in Lung Cancer: Expanding Approaches to Epigenetic Therapy. *Front Oncol* 3: 261

Kahrstrom C. T. (2013) Epigenetics: Legionella makes its mark on histones. *Nat Rev Genet* 14: 370 and local translation in neurons. *Neuron* 80: 648-657

Lee T. I. and Young R. A. (2013) Transcriptional regulation and its misregulation in disease. *Cell* 152: 1237-1251

Rivera C. M. and Ren B. (2013) Mapping human epigenomes. *Cell* 155: 39-55

---

85 Rivera C. M. and Ren B. (2013) Mapping human epigenomes. *Cell* 155: 39-55

86 Pinello L., Xu J., Orkin S. H. and Yuan G. C. (2014) Analysis of chromatin-state plasticity identifies cell-type-specific regulators of H3K27me3 patterns. *Proc Natl Acad Sci U S A* 111: E344-353

87 Yeh H. H., Young D., Gelovani J. G., Robinson A., Davidson Y., et al. (2013) Histone deacetylase class II and acetylated core histone immunohistochemistry in human brains with Huntington's disease. *Brain Res* 1504: 16-24

88 Kahrstrom C. T. (2013) Epigenetics: Legionella makes its mark on histones. *Nat Rev Genet* 14: 370

89 Warnault V., Darco E., Levine A., Barak S. and Ron D. (2013) Chromatin remodeling--a novel strategy to control excessive alcohol drinking. *Transl Psychiatry* 3: e231

90 Marwick J. A., Kirkham P. A., Stevenson C. S., Danahay H., Giddings J., et al. (2004) Cigarette smoke alters chromatin remodeling and induces proinflammatory genes in rat lungs. *Am J Respir Cell Mol Biol* 31: 633-642

Rocha P. P., Chaumeil J. and Skok J. A. (2013) Molecular biology. Finding the right partner in a 3D genome. *Science* 342: 1333-1334

Ronan J. L., Wu W. and Crabtree G. R. (2013) From neural development to cognition: unexpected roles for chromatin. *Nat Rev Genet* 14: 347-359

Telese F., Gamliel A., Skowronska-Krawczyk D., Garcia-Bassets I. and Rosenfeld M. G. (2013) "Seq-ing" insights into the epigenetics of neuronal gene regulation. *Neuron* 77: 606-623

Yadon A. N. and Tsukiyama T. (2013) DNA looping-dependent targeting of a chromatin remodeling factor. *Cell Cycle* 12: 1809-1810

de Wit E. and de Laat W. (2012) A decade of 3C technologies: insights into nuclear organization. *Genes Dev* 26: 11-24


Sajan S. A. and Hawkins R. D. (2012) Methods for identifying higher-order chromatin structure. *Annu Rev Genomics Hum Genet* 13: 59-82

Zentner G. E. and Henikoff S. (2012) Surveying the epigenomic landscape, one base at a time. *Genome Biol* 13: 250



## DNASE I HYPERSENSITIVE SITES SEQUENCING (DNASE-SEQ)

DNase I hypersensitive sites sequencing (DNase-Seq) is based on a well-established DNase I footprinting protocol<sup>91</sup> that was optimized for sequencing<sup>92</sup>. In this method, DNA-protein complexes are treated with DNase I, and the DNA is then extracted and sequenced. Sequences bound by regulatory proteins are protected from DNase I digestion. Deep sequencing provides accurate representation of the location of regulatory proteins in genome. In a variation on this approach, the DNA-protein complexes are stabilized by formaldehyde crosslinking before DNase I digestion. The crosslinking is reversed before DNA purification. In an alternative modification, called GeF-Seq, both the crosslinking and the DNase I digestion are carried out *in vivo*, within permeabilized cells<sup>93</sup>.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can detect "open" chromatin<sup>94</sup></li> <li>• No prior knowledge of the sequence or binding protein is required</li> <li>• Compared to FAIRE-Seq, has greater sensitivity at promoters<sup>95</sup></li> </ul>	<ul style="list-style-type: none"> <li>• DNase I is sequence-specific and hypersensitive sites might not account for the entire genome</li> <li>• Integration of DNase I with ChIP data is necessary to identify and differentiate similar protein-binding sites</li> </ul>

### References

Chumsakul O., Nakamura K., Kurata T., Sakamoto T., Hobman J. L., et al. (2013) High-resolution mapping of *in vivo* genomic transcription factor binding sites using *in situ* DNase I footprinting and ChIP-seq. *DNA Res* 20: 325-338

This study describes an improvement and combination of DNase-Seq with ChIP-Seq, called genome footprinting by high throughput sequencing (GeF-Seq). The authors claim GeF-seq provides better alignment due to shorter reads, resulting in higher resolution of DNA-binding factor recognition sites.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

91 Galas D. J. and Schmitz A. (1978) DNase footprinting: a simple method for the detection of protein-DNA binding specificity. *Nucleic Acids Res* 5: 3157-3170

92 Anderson S. (1981) Shotgun DNA sequencing using cloned DNase I-generated fragments. *Nucleic Acids Res* 9: 3015-3027

93 Chumsakul O., Nakamura K., Kurata T., Sakamoto T., Hobman J. L., et al. (2013) High-resolution mapping of *in vivo* genomic transcription factor binding sites using *in situ* DNase I footprinting and ChIP-seq. *DNA Res* 20: 325-338

94 Zentner G. E. and Henikoff S. (2012) Surveying the epigenomic landscape, one base at a time. *Genome Biol* 13: 250

95 Kumar V., Muratani M., Rayan N. A., Kraus P., Lufkin T., et al. (2013) Uniform, optimal signal processing of mapped deep-sequencing data. *Nat Biotechnol* 31: 615-622

Deng T., Zhu Z. I., Zhang S., Leng F., Cherukuri S., et al. (2013) HMGN1 modulates nucleosome occupancy and DNase I hypersensitivity at the CpG island promoters of embryonic stem cells. *Mol Cell Biol* 33: 3377-3389

The authors use mouse ESCs and NPCs to study the interplay between histone H1 variants and high-mobility group (HMG) proteins in chromatin remodeling. They use ChIP-Seq and DNase-Seq to elucidate the role of HMGN1 (a HMG protein) in affecting chromatin structure at transcription start sites of promoters.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Iwata M., Sandstrom R. S., Delrow J. J., Stamatoyannopoulos J. A. and Torok-Storb B. (2013) Functionally and Phenotypically Distinct Subpopulations of Marrow Stromal Cells Are Fibroblast in Origin and Induce Different Fates in Peripheral Blood Monocytes. *Stem Cells Dev*

Individual cell growth and differentiation is under constant influence by the surrounding tissue and nearby cell types. This study examined marrow stromal cells (MSCs) and their gene expression profiles in comparison to monocyte-derived macrophages that often exist in close proximity to MSCs. Using Illumina sequencing for DNase 1 hypersensitivity mapping, the authors showed a lineage association between two types of MSCs (CD146+, CD146-) and marrow fibroblasts. Subpopulations of CD146+ MSCs were found to increase the expression of genes relevant to hematopoietic regulation upon contact with monocytes, indicating an interaction of fibroblast-macrophage expression.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Ballare C., Castellano G., Gaveglia L., Althammer S., Gonzalez-Vallinas J., et al. (2013) Nucleosome-driven transcription factor binding and gene regulation. *Mol Cell* 49: 67-79

Bertucci P. Y., Nacht A. S., Allo M., Rocha-Viegas L., Ballare C., et al. (2013) Progesterone receptor induces bcl-x expression through intragenic binding sites favoring RNA polymerase II elongation. *Nucleic Acids Res* 41: 6072-6086

Lazarovici A., Zhou T., Shafer A., Dantas Machado A. C., Riley T. R., et al. (2013) Probing DNA shape and methylation state on a genomic scale with DNase I. *Proc Natl Acad Sci U S A* 110: 6376-6381

Lo K. A., Labadorf A., Kennedy N. J., Han M. S., Yap Y. S., et al. (2013) Analysis of in vitro insulin-resistance models and their physiological relevance to in vivo diet-induced adipose insulin resistance. *Cell Rep* 5: 259-270

Degner J. F., Pai A. A., Pique-Regi R., Veyrieras J. B., Gaffney D. J., et al. (2012) DNase I sensitivity QTLs are a major determinant of human expression variation. *Nature* 482: 390-394

Dunowska M., Biggs P. J., Zheng T. and Perrott M. R. (2012) Identification of a novel nidovirus associated with a neurological disease of the Australian brushtail possum (*Trichosurus vulpecula*). *Vet Microbiol* 156: 418-424

He H. H., Meyer C. A., Chen M. W., Jordan V. C., Brown M., et al. (2012) Differential DNase I hypersensitivity reveals factor-dependent chromatin dynamics. *Genome Res* 22: 1015-1025

Lassen K. S., Schultz H., Heegaard N. H. and He M. (2012) A novel DNaseq program for enhanced analysis of Illumina GATK data: a case study on antibody complementarity-determining regions. *N Biotechnol* 29: 271-278

Liu M., Li C. L., Stamatoyannopoulos G., Dorschner M. O., Humbert R., et al. (2012) Gammaretroviral Vector Integration Occurs Overwhelmingly Within and Near DNase Hypersensitive Sites. *Hum Gene Ther* 23: 231-237

Maurano M. T., Humbert R., Rynes E., Thurman R. E., Haugen E., et al. (2012) Systematic localization of common disease-associated variation in regulatory DNA. *Science* 337: 1190-1195


- Neph S., Vierstra J., Stergachis A. B., Reynolds A. P., Haugen E., et al. (2012) An expansive human regulatory lexicon encoded in transcription factor footprints. *Nature* 489: 83-90
- Rosseeel T., Scheuch M., Hoper D., De Regge N., Cajj A. B., et al. (2012) DNase SISPA-next generation sequencing confirms Schmallenberg virus in Belgian field samples and identifies genetic variation in Europe. *PLoS One* 7: e41967
- Wang Y. M., Zhou P., Wang L. Y., Li Z. H., Zhang Y. N., et al. (2012) Correlation between DNase I hypersensitive site distribution and gene expression in HeLa S3 cells. *PLoS One* 7: e42414
- Zhang W., Wu Y., Schnable J. C., Zeng Z., Freeling M., et al. (2012) High-resolution mapping of open chromatin in the rice genome. *Genome Res* 22: 151-162
- Zhang W., Zhang T., Wu Y. and Jiang J. (2012) Genome-wide identification of regulatory DNA elements and protein-binding footprints using signatures of open chromatin in Arabidopsis. *Plant Cell* 24: 2719-2731
- Boyle A. P., Song L., Lee B. K., London D., Keefe D., et al. (2011) High-resolution genome-wide in vivo footprinting of diverse transcription factors in human cells. *Genome Res* 21: 456-464
- Stadler M. B., Murr R., Burger L., Ivanek R., Lienert F., et al. (2011) DNA-binding factors shape the mouse methylome at distal regulatory regions. *Nature* 480: 490-495
- McDaniell R., Lee B. K., Song L., Liu Z., Boyle A. P., et al. (2010) Heritable individual-specific and allele-specific chromatin signatures in humans. *Science* 328: 235-239
- Turnbaugh P. J., Quince C., Faith J. J., McHardy A. C., Yatsunenkov T., et al. (2010) Organismal, genetic, and transcriptional variation in the deeply sequenced gut microbiomes of identical twins. *Proc Natl Acad Sci U S A* 107: 7503-7508
- Audit B., Zaghoul L., Vaillant C., Chevereau G., d'Aubenton-Carafa Y., et al. (2009) Open chromatin encoded in DNA sequence is the signature of 'master' replication origins in human cells. *Nucleic Acids Res* 37: 6064-6075
- Turnbaugh P. J., Ridaura V. K., Faith J. J., Rey F. E., Knight R., et al. (2009) The effect of diet on the human gut microbiome: a metagenomic analysis in humanized gnotobiotic mice. *Sci Transl Med* 1: 6ra14

### Associated Kits

- TruSeq ChIP-Seq kit
- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Preparation Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit

## MNASE-ASSISTED ISOLATION OF NUCLEOSOMES SEQUENCING (MAINE-SEQ)

Micrococcal nuclease (MNase)-assisted isolation of nucleosomes sequencing (MAINE-Seq)<sup>96, 97</sup>, is a variation on the well-established use of MNase digestion to map nucleosome positions (MNase-Seq)<sup>98</sup>. It is estimated that almost half the genome contains regularly spaced arrays of nucleosomes, which are enriched in active chromatin domains<sup>99</sup>. In MAINE-Seq, genomic DNA is treated with MNase. The DNA from the DNA-protein complexes is then extracted and sequenced. Sequences bound by regulatory proteins are protected from MNase digestion. Deep sequencing provides accurate representation of the location of regulatory proteins in the genome<sup>100</sup>. To identify the regulatory proteins, MNase-Seq can be followed by ChIP (NChIP)<sup>101</sup>.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can map nucleosomes and other DNA-binding proteins<sup>102</sup></li> <li>• Identifies location of various regulatory proteins in the genome</li> <li>• Covers broad range of regulatory sites</li> </ul>	<ul style="list-style-type: none"> <li>• MNase sites might not account for the entire genome</li> <li>• Does not provide much information about the kind of regulatory elements</li> <li>• Integration of MNase with ChIP data is necessary to identify and differentiate similar protein-binding sites</li> </ul>

## References

Ballare C., Castellano G., Gaveglia L., Althammer S., Gonzalez-Vallinas J., et al. (2013) Nucleosome-driven transcription factor binding and gene regulation. *Mol Cell* 49: 67-79

This study combines DNase, ChIP, and MAINE sequencing to understand the effects of chromatin remodeling at hormone-responsive regions and thereby the access of hormone receptors to hormone-responsive elements. The authors report nucleosomal involvement in progesterone receptor binding and hormonal gene regulation.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

96 Cusick M. E., Herman T. M., DePamphilis M. L. and Wassarman P. M. (1981) Structure of chromatin at deoxyribonucleic acid replication forks: prenucleosomal deoxyribonucleic acid is rapidly excised from replicating simian virus 40 chromosomes by micrococcal nuclease. *Biochemistry* 20: 6648-6658

97 Ponts N., Harris E. Y., Prudhomme J., Wick I., Eckhardt-Ludka C., et al. (2010) Nucleosome landscape and control of transcription in the human malaria parasite. *Genome Res* 20: 228-238

98 Schlesinger F., Smith A. D., Gingeras T. R., Hannon G. J. and Hodges E. (2013) De novo DNA demethylation and noncoding transcription define active intergenic regulatory elements. *Genome Res* 23: 1601-1614

99 Gaffney D. J., McVicker G., Pai A. A., Fondufe-Mittendorf Y. N., Lewellen N., et al. (2012) Controls of nucleosome positioning in the human genome. *PLoS Genet* 8: e1003036

100 Schones D. E., Cui K., Cuddapah S., Roh T. Y., Barski A., et al. (2008) Dynamic regulation of nucleosome positioning in the human genome. *Cell* 132: 887-898

101 Boyd-Kirkup J. D., Green C. D., Wu G., Wang D. and Han J. D. (2013) Epigenomics and the regulation of aging. *Epigenomics* 5: 205-227

102 Zentner G. E. and Henikoff S. (2012) Surveying the epigenomic landscape, one base at a time. *Genome Biol* 13: 250

Deng T., Zhu Z. I., Zhang S., Leng F., Cherukuri S., et al. (2013) HMGN1 modulates nucleosome occupancy and DNase I hypersensitivity at the CpG island promoters of embryonic stem cells. *Mol Cell Biol* 33: 3377-3389

Chromatin structure and the interaction of DNA with epigenetic factors and chromatin-remodeling complexes play key roles in regulating gene expression and embryonic differentiation. In this study, the authors applied ChIP-Seq, DNase I-Seq, and MNase-Seq on an Illumina Genome Analyzer to analyze the organization of nucleosomes in relation to DNase I hypersensitivity and transcription in mouse ESCs. They found that loss of HMG protein HMGN1 affects two important aspects of chromatin organization: altering the nucleosome positioning at the TSS and reducing the number of DNase I hypersensitivity sites.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Chai X., Nagarajan S., Kim K., Lee K. and Choi J. K. (2013) Regulation of the boundaries of accessible chromatin. *PLoS Genet* 9: e1003778

Grøntved L., John S., Baek S., Liu Y., Buckley J. R., et al. (2013) C/EBP maintains chromatin accessibility in liver and facilitates glucocorticoid receptor recruitment to steroid response elements. *EMBO J* 32: 1568-1583

Maruyama H., Harwood J. C., Moore K. M., Paszkiewicz K., Durley S. C., et al. (2013) An alternative beads-on-a-string chromatin architecture in *Thermococcus kodakarensis*. *EMBO Rep* 14: 711-717

Nagarajavel V., Iben J. R., Howard B. H., Marais R. J. and Clark D. J. (2013) Global 'bootprinting' reveals the elastic architecture of the yeast TFIIB-TFIIC transcription complex in vivo. *Nucleic Acids Res* 41: 8135-8143

Nishida H., Katayama T., Suzuki Y., Kondo S. and Horiuchi H. (2013) Base composition and nucleosome density in exonic and intronic regions in genes of the filamentous ascomycetes *Aspergillus nidulans* and *Aspergillus oryzae*. *Gene* 525: 10-May

Tolstorukov M. Y., Sansam C. G., Lu P., Koellhoffer E. C., Helming K. C., et al. (2013) Swi/Snf chromatin remodeling/tumor suppressor complex establishes nucleosome occupancy at target promoters. *Proc Natl Acad Sci U S A* 110: 10165-10170

Allan J., Fraser R. M., Owen-Hughes T. and Keszenman-Pereyra D. (2012) Micrococcal nuclease does not substantially bias nucleosome mapping. *J Mol Biol* 417: 152-164

Gaffney D. J., McVicker G., Pai A. A., Fondufe-Mittendorf Y. N., Lewellen N., et al. (2012) Controls of nucleosome positioning in the human genome. *PLoS Genet* 8: e1003036

Guertin M. J., Martins A. L., Siepel A. and Lis J. T. (2012) Accurate prediction of inducible transcription factor binding intensities in vivo. *PLoS Genet* 8: e1002610

Zhang X., Robertson G., Woo S., Hoffman B. G. and Gottardo R. (2012) Probabilistic Inference for Nucleosome Positioning with MNase-Based or Sonicated Short-Read Data. *PLoS ONE* 7: e32095

Kent N. A., Adams S., Moorhouse A. and Paszkiewicz K. (2011) Chromatin particle spectrum analysis: a method for comparative chromatin structure analysis using paired-end mode next-generation DNA sequencing. *Nucleic Acids Res* 39: e26

Xi Y., Yao J., Chen R., Li W. and He X. (2011) Nucleosome fragility reveals novel functional states of chromatin and poises genes for activation. *Genome Res* 21: 718-724

Ponts N., Harris E. Y., Prudhomme J., Wick I., Eckhardt-Ludka C., et al. (2010) Nucleosome landscape and control of transcription in the human malaria parasite. *Genome Res* 20: 228-238

**Associated Kits**

TruSeq ChIP-Seq® Kit

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Preparation Kit


TruSeq DNA PCR-Free Sample Prep Kit

Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

## CHROMATIN IMMUNOPRECIPITATION SEQUENCING (CHIP-SEQ)

Chromatin immunoprecipitation sequencing (ChIP-Seq) is a well-established method to map specific protein-binding sites<sup>103</sup>. In this method, DNA-protein complexes are crosslinked *in vivo*. Samples are then fragmented and treated with an exonuclease to trim unbound oligonucleotides. Protein-specific antibodies are used to immunoprecipitate the DNA-protein complex. The DNA is extracted and sequenced, giving high-resolution sequences of the protein-binding sites.



Pros	Cons
<ul style="list-style-type: none"> <li>• Base-pair resolution of protein-binding site</li> <li>• Specific regulatory factors or proteins can be mapped</li> <li>• The use of exonuclease eliminates contamination by unbound DNA<sup>104</sup></li> </ul>	<ul style="list-style-type: none"> <li>• Nonspecific antibodies can dilute the pool of DNA-protein complexes of interest</li> <li>• The target protein must be known and able to raise an antibody</li> </ul>

### References

Berkseth M., Ikegami K., Arur S., Lieb J. D. and Zarkower D. (2013) TRA-1 ChIP-seq reveals regulators of sexual differentiation and multilevel feedback in nematode sex determination. *Proc Natl Acad Sci U S A* 110: 16033-16038

In an effort to identify targets of the nematode global sexual regulator Transformer 1 (TRA-1), this study applied Illumina sequencing for genome-wide ChIP-Seq analysis of TRA-1 binding sites. The authors identified DNA-binding sites driving male-specific expression patterns and TRA-1 binding sites adjacent to a number of regulatory genes, some of which drive male-specific expression. Overall, the results suggest that TRA-1 mediates sex-specific expression.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

Bowman S. K., Simon M. D., Deaton A. M., Tolstorukov M., Borowsky M. L., et al. (2013) Multiplexed Illumina sequencing libraries from picogram quantities of DNA. *BMC Genomics* 14: 466

This study reports a simple and fast library construction method from sub-nanogram quantities of DNA. This protocol yields conventional libraries with barcodes suitable for multiplexed sample analysis on the Illumina platform. The authors demonstrate the method by constructing a ChIP-Seq library from 100 pg of ChIP DNA that shows equivalent coverage of target regions to a library produced from a larger-scale experiment.

**Illumina Technology:** HiSeq 2000

<sup>103</sup> Solomon M. J., Larsen P. L. and Varshavsky A. (1988) Mapping protein-DNA interactions in vivo with formaldehyde: evidence that histone H4 is retained on a highly transcribed gene. *Cell* 53: 937-947

<sup>104</sup> Zentner G. E. and Henikoff S. (2012) Surveying the epigenomic landscape, one base at a time. *Genome Biol* 13: 250

Kumar V., Muratani M., Rayan N. A., Kraus P., Lufkin T., et al. (2013) Uniform, optimal signal processing of mapped deep-sequencing data. *Nat Biotechnol* 31: 615-622

ChIP-Seq experiments are used to determine the occupation of chromatin by DNA-binding proteins. Data analysis requires detection of binding signals above the background noise, and a common secondary analysis is the prediction of an effect, e.g., expression, from the level of the ChIP-Seq signal. This paper presents algorithms adapted from signal processing theory to solve the two general problems of signal detection and signal estimation from ChIP-Seq data. Using existing data and a new ChIP-Seq data set from an Illumina Genome Analyzer, the two tools DFilter and EFilter are shown to outperform the most commonly used methods in the field, including MACS and Quest.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Lesch B. J., Dokshin G. A., Young R. A., McCarrey J. R. and Page D. C. (2013) A set of genes critical to development is epigenetically poised in mouse germ cells from fetal stages through completion of meiosis. *Proc Natl Acad Sci U S A* 110: 16061-16066

At conception the zygote is totipotent: incorporating the potential to differentiate into any specialized cell in the body. This study used gene expression profiling and epigenetic regulatory marks (H3K4me3 and H3K37me3) to examine how germ cells change as they progress from differentiated cell to totipotent zygote. The authors used ChIP-Seq and RNA-Seq on the Illumina HiSeq platform for both male and female germ cells at three time points surrounding sex differentiation, meiosis, and post-meiosis. Their results indicate central regulatory genes are maintained in an epigenetically poised state, permitting establishment of totipotency following fertilization.

**Illumina Technology:** HiSeq2000

Schauer T., Schwalie P. C., Handley A., Margulies C. E., Flicek P., et al. (2013) CAST-ChIP maps cell-type-specific chromatin states in the *Drosophila* central nervous system. *Cell Rep* 5: 271-282

Accurate assays for epigenetic markers have been limited by the amount of input material required. This study presents a new assay (CAST-ChIP), based on Illumina sequencing, that allows for characterization of chromatin-associated proteins from specific cell types in complex tissues. The study validates the assay by profiling PolII and H2A.Z across both glia and neurons in *Drosophila* brain tissue.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Koldamova R., Schug J., Lefterova M., Cronican A. A., Fitz N. F., et al. (2014) Genome-wide approaches reveal EGR1-controlled regulatory networks associated with neurodegeneration. *Neurobiol Dis* 63: 107-114

McMullen P. D., Bhattacharya S., Woods C. G., Sun B., Yarborough K., et al. (2014) A map of the PPARalpha transcription regulatory network for primary human hepatocytes. *Chem Biol Interact* 209: 14-24

Waszak S. M., Kilpinen H., Gschwind A. R., Orioli A., Raghav S. K., et al. (2014) Identification and removal of low-complexity sites in allele-specific analysis of ChIP-seq data. *Bioinformatics* 30: 165-171

Wolchinsky Z., Shvitiel S., Kouwenhoven E. N., Putin D., Sprecher E., et al. (2014) Angiomodulin is required for cardiogenesis of embryonic stem cells and is maintained by a feedback loop network of p63 and Activin-A. *Stem Cell Res* 12: 49-59

Biancolella M., B K. F., Tring S., Plummer S. J., Mendoza-Fandino G. A., et al. (2013) Identification and characterization of functional risk variants for colorectal cancer mapping to chromosome 11q23.1. *Hum Mol Genet*

Chang C. Y., Pasolli H. A., Giannopoulou E. G., Guasch G., Gronostajski R. M., et al. (2013) NFIB is a governor of epithelial-melanocyte stem cell behaviour in a shared niche. *Nature* 495: 98-102



- Chauhan C., Zraly C. B. and Dingwall A. K. (2013) The Drosophila COMPASS-like Cmi-Trr coactivator complex regulates dpp/BMP signaling in pattern formation. *Dev Biol* 380: 185-198
- Jain A., Bacolla A., Del Mundo I. M., Zhao J., Wang G., et al. (2013) DHX9 helicase is involved in preventing genomic instability induced by alternatively structured DNA in human cells. *Nucleic Acids Res* 41: 10345-10357
- Lai C. F., Flach K. D., Alexi X., Fox S. P., Ottaviani S., et al. (2013) Co-regulated gene expression by oestrogen receptor alpha and liver receptor homolog-1 is a feature of the oestrogen response in breast cancer cells. *Nucleic Acids Res* 41: 10228-10240
- Lo K. A., Labadorf A., Kennedy N. J., Han M. S., Yap Y. S., et al. (2013) Analysis of in vitro insulin-resistance models and their physiological relevance to in vivo diet-induced adipose insulin resistance. *Cell Rep* 5: 259-270
- Maehara K., Odawara J., Harada A., Yoshimi T., Nagao K., et al. (2013) A co-localization model of paired ChIP-seq data using a large ENCODE data set enables comparison of multiple samples. *Nucleic Acids Res* 41: 54-62
- Mazzoni E. O., Mahony S., Closser M., Morrison C. A., Nedelec S., et al. (2013) Synergistic binding of transcription factors to cell-specific enhancers programs motor neuron identity. *Nat Neurosci* 16: 1219-1227
- Mendoza-Parra M. A., Van Gool W., Mohamed Saleem M. A., Ceschin D. G. and Gronemeyer H. (2013) A quality control system for profiles obtained by ChIP sequencing. *Nucleic Acids Res* 41: e196
- Neyret-Kahn H., Benhamed M., Ye T., Le Gras S., Cossec J. C., et al. (2013) Sumoylation at chromatin governs coordinated repression of a transcriptional program essential for cell growth and proliferation. *Genome Res* 23: 1563-1579
- Olovnikov I., Ryazansky S., Shpiz S., Lavrov S., Abramov Y., et al. (2013) De novo piRNA cluster formation in the Drosophila germ line triggered by transgenes containing a transcribed transposon fragment. *Nucleic Acids Res* 41: 5757-5768
- Prickett A. R., Barkas N., McCole R. B., Hughes S., Amante S. M., et al. (2013) Genome-wide and parental allele-specific analysis of CTCF and cohesin DNA binding in mouse brain reveals a tissue-specific binding pattern and an association with imprinted differentially methylated regions. *Genome Res* 23: 1624-1635
- Schmolka N., Serre K., Grosso A. R., Rei M., Pennington D. J., et al. (2013) Epigenetic and transcriptional signatures of stable versus plastic differentiation of proinflammatory gammadelta T cell subsets. *Nat Immunol* 14: 1093-1100
- Sakabe N. J., Aneas I., Shen T., Shokri L., Park S. Y., et al. (2012) Dual transcriptional activator and repressor roles of TBX20 regulate adult cardiac structure and function. *Hum Mol Genet* 21: 2194-2204
- Tanaka Y., Joshi A., Wilson N. K., Kinston S., Nishikawa S., et al. (2012) The transcriptional programme controlled by Runx1 during early embryonic blood development. *Dev Biol* 366: 404-419
- Tang B., Becanovic K., Desplats P. A., Spencer B., Hill A. M., et al. (2012) Forkhead box protein p1 is a transcriptional repressor of immune signaling in the CNS: implications for transcriptional dysregulation in Huntington disease. *Hum Mol Genet* 21: 3097-3111
- Zhan Q., Fang Y., He Y., Liu H. X., Fang J., et al. (2012) Function annotation of hepatic retinoid x receptor alpha based on genome-wide DNA binding and transcriptome profiling. *PLoS One* 7: e50013

### Associated Kits

- TruSeq ChIP-Seq Kit
- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Preparation Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit

## FORMALDEHYDE-ASSISTED ISOLATION OF REGULATORY ELEMENTS (FAIRE-SEQ)

Formaldehyde-assisted isolation of regulatory elements (FAIRE-Seq)<sup>105,106</sup> is based on differences in crosslinking efficiencies between DNA and nucleosomes or sequence-specific DNA-binding proteins. In this method, DNA-protein complexes are briefly crosslinked *in vivo* using formaldehyde. The sample is then lysed and sonicated. After phenol/chloroform extraction, the DNA in the aqueous phase is purified and sequenced. Sequencing provides information for regions of DNA that are not occupied by histones.



Pros	Cons
<ul style="list-style-type: none"> <li>• Simple and highly reproducible protocol</li> <li>• Does not require antibodies</li> <li>• Does not require enzymes, such as DNase or MNase, avoiding the optimization and extra steps necessary for enzymatic processing</li> <li>• Does not require a single-cell suspension or nuclear isolation, so it is easily adapted for use on tissue samples<sup>107</sup></li> </ul>	<ul style="list-style-type: none"> <li>• Cannot identify regulatory proteins bound to DNA</li> <li>• DNase-Seq may be better at identifying nucleosome-depleted promoters of highly expressed genes<sup>108</sup></li> </ul>

### References

Hilton I. B., Simon J. M., Lieb J. D., Davis I. J., Damania B., et al. (2013) The open chromatin landscape of Kaposi's sarcoma-associated herpesvirus. *J Virol* 87: 11831-11842

Kaposi's sarcoma-associated herpesvirus (KSHV) is a gammaherpesvirus that, upon infection, remains in a latent state. Histone modifications occupy inactive regions of the latent viral genome. The authors use FAIRE-Seq on the Illumina HiSeq 2000 system to study open chromatin regions in the KSHV genome, allowing them to identify regions of open chromatin in the latent virus. By integrating data on histone modifications, they were able to generate a genome-wide KSHV landscape, which indicated localization of active histone modifications near nucleosome-depleted sites.

**Illumina Technology:** TruSeq Sample Prep Kit, HiSeq 2000

105 Giresi P. G. and Lieb J. D. (2009) Isolation of active regulatory elements from eukaryotic chromatin using FAIRE (Formaldehyde Assisted Isolation of Regulatory Elements). *Methods* 48: 233-239

106 Hogan G. J., Lee C. K. and Lieb J. D. (2006) Cell cycle-specified fluctuation of nucleosome occupancy at gene promoters. *PLoS Genet* 2: e158

107 Simon J. M., Giresi P. G., Davis I. J. and Lieb J. D. (2012) Using formaldehyde-assisted isolation of regulatory elements (FAIRE) to isolate active regulatory DNA. *Nat Protoc* 7: 256-267

108 Song L., Zhang Z., Grasfeder L. L., Boyle A. P., Giresi P. G., et al. (2011) Open chromatin defined by DNaseI and FAIRE identifies regulatory elements that shape cell-type identity. *Genome Res* 21: 1757-1767

Meredith D. M., Borromeo M. D., Deering T. G., Casey B. H., Savage T. K., et al. (2013) Program specificity for Ptf1a in pancreas versus neural tube development correlates with distinct collaborating cofactors and chromatin accessibility. *Mol Cell Biol* 33: 3166-3179

Transcription factors (TFs) are the drivers of cell development and differentiation. The combined regulatory effects of different TFs allow any factor to play key roles in the different pathways of cell differentiation. This study examined how pancreas-specific transcription factor 1a (Ptf1a) is a critical driver for development of both the pancreas and nervous system. Using Illumina sequencing to perform ChIP-Seq for Ptf1a, FAIRE-Seq to detect open chromatin, and RNA-Seq for expression profiling, the authors characterized Fox and Sox factors as potential lineage-specific modifiers of Ptf1a binding.

**Illumina Technology:** HiSeq 2000, Genome Analyzer<sub>IIx</sub>

Paul D. S., Albers C. A., Rendon A., Voss K., Stephens J., et al. (2013) Maps of open chromatin highlight cell type-restricted patterns of regulatory sequence variation at hematological trait loci. *Genome Res* 23: 1130-1141

Genome-wide association studies (GWAS) have discovered many non-protein-coding loci associated with complex traits. However, due to the low resolution of GWAS, the exact location of the causative variant is often not known. In this study, the authors combined GWAS results with FAIRE-Seq to link complex hematopoietic traits to specific functional loci. They found that the majority of candidate functional variants coincided with binding sites of five transcription factors key to regulating megakaryopoiesis, and further found that 76.9% of the candidate regulatory variants affected protein binding at these sites. In conclusion, the combination of GWAS data with high-resolution epigenetic profiling by sequencing is a powerful assay for mapping complex genetic variants.

**Illumina Technology:** HiSeq 2000, Genome Analyzer<sub>IIx</sub>, Human Gene Expression—BeadArray

Chai X., Nagarajan S., Kim K., Lee K. and Choi J. K. (2013) Regulation of the boundaries of accessible chromatin. *PLoS Genet* 9: e1003778

Calabrese J. M., Sun W., Song L., Mugford J. W., Williams L., et al. (2012) Site-specific silencing of regulatory elements as a mechanism of X inactivation. *Cell* 151: 951-963

Simon J. M., Giresi P. G., Davis I. J. and Lieb J. D. (2012) Using formaldehyde-assisted isolation of regulatory elements (FAIRE) to isolate active regulatory DNA. *Nat Protoc* 7: 256-267

Paul D. S., Nisbet J. P., Yang T. P., Meacham S., Rendon A., et al. (2011) Maps of open chromatin guide the functional follow-up of genome-wide association signals: application to hematological traits. *PLoS Genet* 7: e1002139

Ponts N., Harris E. Y., Prudhomme J., Wick I., Eckhardt-Ludka C., et al. (2010) Nucleosome landscape and control of transcription in the human malaria parasite. *Genome Res* 20: 228-238

Auerbach R. K., Euskirchen G., Rozowsky J., Lamarre-Vincent N., Moqtaderi Z., et al. (2009) Mapping accessible chromatin regions using Sono-Seq. *Proc Natl Acad Sci U S A* 106: 14926-14931

#### Associated Kits

TruSeq ChIP-Seq Kit

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Prep Kit


TruSeq DNA PCR-Free Sample Prep Kit

Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

## ASSAY FOR TRANSPOSASE-ACCESSIBLE CHROMATIN SEQUENCING (ATAC-SEQ)

Assay for transposase-accessible chromatin using sequencing (ATAC-Seq) is a protocol that utilizes the Epicentre Tn5 transposome<sup>109</sup>. In this method, DNA is incubated with Tn5 transposome, which performs adaptor ligation and fragmentation of open chromatin regions. Deep sequencing of the purified regions provides base-pair resolution of nucleosome-free regions in the genome.




Pros	Cons
<ul style="list-style-type: none"> <li>• Two-step protocol with no adaptor ligation steps, gel purification, or crosslink reversal</li> <li>• Very high signal to noise ratio compared to FAIRE-Seq</li> </ul>	<ul style="list-style-type: none"> <li>• During mechanical sample processing, bound chromatin regions might open and be tagged by the transposome</li> </ul>

### References

Buenrostro J. D., Giresi P. G., Zaba L. C., Chang H. Y. and Greenleaf W. J. (2013) Transposition of native chromatin for fast and sensitive epigenomic profiling of open chromatin, DNA-binding proteins and nucleosome position. *Nat Methods* 10: 1213-1218

This is the first paper to describe ATAC-seq as a protocol to study regions of open chromatin. The authors identify the location of DNA-binding proteins in a B-cell line. They demonstrate that the protocol can analyze an individual's T-cell epigenome on a timescale compatible with clinical decision-making.

**Illumina Technology:** MiSeq, HiSeq 2000



ATAC-Seq enables real-time personal epigenomics.


### Associated Kits

- EpiGnome™ Methyl-Seq Kit
- TruSeq ChIP-Seq Kit
- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit

<sup>109</sup> Buenrostro J. D., Giresi P. G., Zaba L. C., Chang H. Y. and Greenleaf W. J. (2013) Transposition of native chromatin for fast and sensitive epigenomic profiling of open chromatin, DNA-binding proteins and nucleosome position. *Nat Methods* 10: 1213-1218

## CHROMATIN INTERACTION ANALYSIS BY PAIRED-END TAG SEQUENCING (CHIA-PET)

Chromatin interaction analysis by paired-end tag sequencing (ChIA-PET) is a variation of Hi-C that features an immunoprecipitation step to map long-range DNA interactions<sup>110, 111</sup>. In this method, DNA-protein complexes are crosslinked and fragmented. Specific antibodies are used to immunoprecipitate proteins of interest. Specific linkers are ligated to the DNA fragments, which ligate when in proximity. Linkers are then precipitated and digested with an enzyme and the DNA is sequenced. Deep sequencing provides base-pair resolution of ligated fragments. Hi-C and ChIA-PET currently provide the best balance of resolution and reasonable coverage in the human genome to map long-range interactions<sup>112</sup>



Pros	Cons
<ul style="list-style-type: none"> <li>• Suitable for detecting a large number of both long-range and short range chromatin interactions globally<sup>113</sup></li> <li>• Studies the interactions made by specific proteins or protein complexes</li> <li>• Provides information about DNA interactions aided by regulatory elements</li> <li>• Removes background generated during traditional ChIP assays</li> <li>• The immunoprecipitation step reduces data complexity<sup>113</sup></li> </ul>	<ul style="list-style-type: none"> <li>• Nonspecific antibodies can pull down unwanted protein complexes and contaminate the pool</li> <li>• Linkers can self-ligate, generating ambiguity about true DNA interactions</li> <li>• Limited sensitivity; may detect as little as 10% of interactions<sup>113</sup></li> </ul>

### References

DeMare L. E., Leng J., Cotney J., Reilly S. K., Yin J., et al. (2013) The genomic landscape of cohesin-associated chromatin interactions. *Genome Res* 23: 1224-1234

Knockdown of cohesin in ESCs results in aberrant gene expression and loss of pluripotency. Cohesin works to stabilize DNA by forming loops between distant-acting enhancers and their target promoters. The authors studied cohesin interaction in the developing limb using ChIA-PET, RNA-Seq, and ChIP-Seq analysis performed on a HiSeq 2000 system. They report tissue-specific enhancer-promoter interactions involving cohesin and the insulator protein CTCF. They also identified interactions that are maintained for tissue-specific activation or repression during development.

**Illumina Technology:** TruSeq Sample Prep Kit, HiSeq 2000

110 Li G., Fullwood M. J., Xu H., Mulawadi F. H., Velkov S., et al. (2010) ChIA-PET tool for comprehensive chromatin interaction analysis with paired-end tag sequencing. *Genome Biol* 11: R22

111 Fullwood M. J., Liu M. H., Pan Y. F., Liu J., Xu H., et al. (2009) An oestrogen-receptor-alpha-bound human chromatin interactome. *Nature* 462: 58-64

112 Dekker J., Marti-Renom M. A. and Mirny L. A. (2013) Exploring the three-dimensional organization of genomes: interpreting chromatin interaction data. *Nat Rev Genet* 14: 390-403

113 Sajjan S. A. and Hawkins R. D. (2012) Methods for identifying higher-order chromatin structure. *Annu Rev Genomics Hum Genet* 13: 59-82

Stadhouders R., Kolovos P., Brouwer R., Zuin J., van den Heuvel A., et al. (2013) Multiplexed chromosome conformation capture sequencing for rapid genome-scale high-resolution detection of long-range chromatin interactions. *Nat Protoc* 8: 509-524

This paper presents an assay for multiplexed chromosome conformation capture sequencing (3C-Seq) using an Illumina HiSeq 2000 system. This high-throughput assay outperforms PCR-based methods for ease of multiplexing, and outperforms 5C and Hi-C methods in terms of cost and ease of analysis. The preparation of multiplexed 3C-Seq libraries can be performed by any investigator with basic skills in molecular biology techniques, and the data analysis requires only basic expertise in bioinformatics.

**Illumina Technology:** HiSeq 2000

Li G., Ruan X., Auerbach R. K., Sandhu K. S., Zheng M., et al. (2012) Extensive promoter-centered chromatin interactions provide a topological basis for transcription regulation. *Cell* 148: 84-98

Zhang J., Poh H. M., Peh S. Q., Sia Y. Y., Li G., et al. (2012) ChIA-PET analysis of transcriptional chromatin interactions. *Methods* 58: 289-299

Tan S. K., Lin Z. H., Chang C. W., Varang V., Chng K. R., et al. (2011) AP-2gamma regulates oestrogen receptor-mediated long-range chromatin interaction and gene transcription. *EMBO J* 30: 2569-2581

Fullwood M. J., Han Y., Wei C. L., Ruan X. and Ruan Y. (2010) Chromatin interaction analysis using paired-end tag sequencing. *Curr Protoc Mol Biol* Chapter 21: Unit 21 15 21-25

Li G., Fullwood M. J., Xu H., Mulawadi F. H., Velkov S., et al. (2010) ChIA-PET tool for comprehensive chromatin interaction analysis with paired-end tag sequencing. *Genome Biol* 11: R22

#### **Associated Kits**

TruSeq ChIP-Seq Kit

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit


Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Mate® Pair Kit

## CHROMATIN CONFORMATION CAPTURE (HI-C/3C-SEQ)

Chromatin conformation capture sequencing (Hi-C)<sup>114</sup> or 3C-Seq<sup>115</sup> is used to analyze chromatin interactions. In this method, DNA-protein complexes are crosslinked using formaldehyde. The sample is fragmented and DNA ligated and digested. The resulting DNA fragments are PCR-amplified and sequenced. Deep sequencing provides base-pair resolution of ligated fragments.



Pros	Cons
<ul style="list-style-type: none"> <li>• Allows detection of long-range DNA interactions</li> <li>• High-throughput method</li> </ul>	<ul style="list-style-type: none"> <li>• Detection may result from random chromosomal collisions</li> <li>• 3C PCR is difficult and requires careful controls and experimental design</li> <li>• Needs further confirmation of interaction</li> <li>• Due to multiple steps, the method requires large amounts of starting material</li> </ul>

### References

Burton J. N., Adey A., Patwardhan R. P., Qiu R., Kitzman J. O., et al. (2013) Chromosome-scale scaffolding of de novo genome assemblies based on chromatin interactions. *Nat Biotechnol* 31: 1119-1125

The authors integrate shotgun fragment and short insert mate-pair sequences with Hi-C data to generate assemblies for human, mouse, and Drosophila genomes. The paper reports a bioinformatics tool used to compute the assemblies: ligating adjacent chromatin enables scaffolding in situ (LACHESIS).

**Illumina Technology:** HiSeq 2000

Jin F., Li Y., Dixon J. R., Selvaraj S., Ye Z., et al. (2013) A high-resolution map of the three-dimensional chromatin interactome in human cells. *Nature* 503: 290-294

Cis-acting regulatory elements in the genome interact with their target gene promoter by transcription factors bringing the two locations close in the three-dimensional conformation of the chromatin. In this study, the chromosome conformation is studied by a genome-wide analysis method (Hi-C) using the Illumina HiSeq 2000 system. The authors determined over one million long-range chromatin interactions in human fibroblasts. In addition, they characterized the dynamics of promoter-enhancer contacts after TNF-alpha signaling and discovered pre-existing chromatin looping with TNF-alpha-responsive enhancers, suggesting the three-dimensional chromatin conformation may be stable over time.

**Illumina Technology:** HiSeq 2000

114 Lieberman-Aiden E., van Berkum N. L., Williams L., Imakaev M., Ragoczy T., et al. (2009) Comprehensive mapping of long-range interactions reveals folding principles of the human genome. *Science* 326: 289-293

115 Duan Z., Andronescu M., Schutz K., Lee C., Shendure J., et al. (2012) A genome-wide 3C-method for characterizing the three-dimensional architectures of genomes. *Methods* 58: 277-288

- Nagano T., Lubling Y., Stevens T. J., Schoenfelder S., Yaffe E., et al. (2013) Single-cell Hi-C reveals cell-to-cell variability in chromosome structure. *Nature* 502: 59-64
- Belton J. M., McCord R. P., Gibcus J. H., Naumova N., Zhan Y., et al. (2012) Hi-C: a comprehensive technique to capture the conformation of genomes. *Methods* 58: 268-276
- Demichelis F., Setlur S. R., Banerjee S., Chakravarty D., Chen J. Y., et al. (2012) Identification of functionally active, low frequency copy number variants at 15q21.3 and 12q21.31 associated with prostate cancer risk. *Proc Natl Acad Sci U S A* 109: 6686-6691
- Dixon J. R., Selvaraj S., Yue F., Kim A., Li Y., et al. (2012) Topological domains in mammalian genomes identified by analysis of chromatin interactions. *Nature* 485: 376-380
- Imakaev M, Fudenberg G, McCord RP, Naumova N, Goloborodko A, Lajoie BR, Dekker J, Mirny LA; (2012) Iterative correction of Hi-C data reveals hallmarks of chromosome organization. *Nat Methods* 9: 999-1003
- Imakaev M., Fudenberg G., McCord R. P., Naumova N., Goloborodko A., et al. (2012) Iterative correction of Hi-C data reveals hallmarks of chromosome organization. *Nat Methods* 9: 999-1003
- Irimia M., Tena J. J., Alexis M. S., Fernandez-Minan A., Maeso I., et al. (2012) Extensive conservation of ancient microsynteny across metazoans due to cis-regulatory constraints. *Genome Res* 22: 2356-2367
- Lan X., Witt H., Katsumura K., Ye Z., Wang Q., et al. (2012) Integration of Hi-C and ChIP-seq data reveals distinct types of chromatin linkages. *Nucleic Acids Res* 40: 7690-7704
- Verma-Gaur J., Torkamani A., Schaffer L., Head S. R., Schork N. J., et al. (2012) Noncoding transcription within the Igh distal V(H) region at PAIR elements affects the 3D structure of the Igh locus in pro-B cells. *Proc Natl Acad Sci U S A* 109: 17004-17009
- Zhang Y., McCord R. P., Ho Y. J., Lajoie B. R., Hildebrand D. G., et al. (2012) Spatial organization of the mouse genome and its role in recurrent chromosomal translocations. *Cell* 148: 908-921
- Yaffe E. and Tanay A. (2011) Probabilistic modeling of Hi-C contact maps eliminates systematic biases to characterize global chromosomal architecture. *Nat Genet* 43: 1059-1065
- Yaffe E, Tanay A; (2011) Probabilistic modeling of Hi-C contact maps eliminates systematic biases to characterize global chromosomal architecture. *Nat Genet* 43: 1059-65
- Lieberman-Aiden E., van Berkum N. L., Williams L., Imakaev M., Ragozy T., et al. (2009) Comprehensive mapping of long-range interactions reveals folding principles of the human genome. *Science* 326: 289-293


#### **Associated Kits**

- TruSeq ChIP-Seq Kit
- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Mate Pair Kit



## CIRCULAR CHROMATIN CONFORMATION CAPTURE (4-C OR 4C-SEQ)

Circular chromatin conformation capture (4-C)<sup>116</sup>, also called 4C-Seq, is a method similar to 3-C and is sometimes called circular 3C. It allows the unbiased detection of all genomic regions that interact with a particular region of interest<sup>117</sup>. In this method, DNA-protein complexes are crosslinked using formaldehyde. The sample is fragmented, and the DNA is ligated and digested. The resulting DNA fragments self-circularize, followed by reverse PCR and sequencing. Deep sequencing provides base-pair resolution of ligated fragments.



Pros	Cons
<ul style="list-style-type: none"> <li>• 4C is the preferred strategy to assess the DNA contact profile of individual genomic sites.</li> <li>• Highly reproducible data</li> </ul>	<ul style="list-style-type: none"> <li>• Will miss local interactions (&lt; 50 kb) from the region of interest</li> <li>• The large circles do not PCR efficiently</li> </ul>

### References

de Wit E., Bouwman B. A., Zhu Y., Klous P., Splinter E., et al. (2013) The pluripotent genome in three dimensions is shaped around pluripotency factors. *Nature* 501: 227-231

Transcriptional regulation is influenced by the availability of specific transcription factors, but the evidence is increasing for the substantial importance of chromatin conformation within the nucleus. In this study, Illumina sequencing is used to analyze chromatin conformation by a genome-wide assay (4-C) demonstrating, along with ChIP-Seq data, that inactive chromatin is disorganized in PSC nuclei. In contrast to inactive chromatin, promoters are seen to engage in contacts between topological domains in a tissue-dependent manner, while enhancers have a more tissue-restricted interaction. The authors hypothesize that the chromatin interactions enhance the robustness of the pluripotent state.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

116 Zhao Z., Tavoosidana G., Sjolinder M., Gondor A., Mariano P., et al. (2006) Circular chromosome conformation capture (4C) uncovers extensive networks of epigenetically regulated intra- and interchromosomal interactions. *Nat Genet* 38: 1341-1347

117 Sajan S. A. and Hawkins R. D. (2012) Methods for identifying higher-order chromatin structure. *Annu Rev Genomics Hum Genet* 13: 59-82

Holwerda S. J., van de Werken H. J., Ribeiro de Almeida C., Bergen I. M., de Bruijn M. J., et al. (2013) Allelic exclusion of the immunoglobulin heavy chain locus is independent of its nuclear localization in mature B cells. *Nucleic Acids Res* 41: 6905-6916

Chromatin conformation is one of many mechanisms for regulating gene expression. In developing B cells, the immunoglobulin heavy chain (IgH) locus undergoes a scheduled genomic rearrangement of the V, D, and J gene segments. In this study, an allele-specific chromosome conformation capture sequencing technique (4C-Seq) was applied to unambiguously follow the individual IgH alleles in mature B lymphocytes. The authors found that IgH adopts a lymphoid-specific nuclear location, and in mature B cells the distal VH regions of both IgH alleles position themselves away from active chromatin.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq 2000

Wei Z., Gao F., Kim S., Yang H., Lyu J., et al. (2013) Klf4 organizes long-range chromosomal interactions with the oct4 locus in reprogramming and pluripotency. *Cell Stem Cell* 13: 36-47

PSCs are capable of differentiation into diverse cell types. The maintenance of pluripotency and the induction of differentiation are both highly regulated processes. This study examined the epigenetic mechanisms underlying reprogramming of PSCs. Using circular chromosome conformation capture with Illumina HiSeq sequencing technology (4C-Seq), the authors profiled the PSC-specific long-range chromosomal interactions during reprogramming to induced PSCs. The high-resolution genome-wide interaction map and a well-designed experimental setup allowed the authors to show evidence for a functional role of Kruppel-like factor 4 (Klf4) in facilitating long-range interactions.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>, HiSeq2000

Chaumeil J., Micsinai M., Ntziachristos P., Roth D. B., Aifantis I., et al. (2013) The RAG2 C-terminus and ATM protect genome integrity by controlling antigen receptor gene cleavage. *Nat Commun* 4: 2231

Delpretti S., Montavon T., Leleu M., Joye E., Tzika A., et al. (2013) Multiple Enhancers Regulate Hoxd Genes and the Hotdog LncRNA during Cecum Budding. *Cell Rep* 5: 137-150

Denholtz M., Bonora G., Chronis C., Splinter E., de Laat W., et al. (2013) Long-Range Chromatin Contacts in Embryonic Stem Cells Reveal a Role for Pluripotency Factors and Polycomb Proteins in Genome Organization. *Cell Stem Cell* 13: 602-616

Jankovic M., Feldhahn N., Oliveira T. Y., Silva I. T., Kieffer-Kwon K. R., et al. (2013) 53BP1 alters the landscape of DNA rearrangements and suppresses AID-induced B cell lymphoma. *Mol Cell* 49: 623-631

Medvedovic J., Ebert A., Tagoh H., Tamir I. M., Schwickert T. A., et al. (2013) Flexible long-range loops in the VH gene region of the Igh locus facilitate the generation of a diverse antibody repertoire. *Immunity* 39: 229-244

Yamane A., Robbiani D. F., Resch W., Bothmer A., Nakahashi H., et al. (2013) RPA Accumulation during Class Switch Recombination Represents 5'-3' DNA-End Resection during the S-G2/M Phase of the Cell Cycle. *Cell Rep* 3: 138-147

Hakim O., Resch W., Yamane A., Klein I., Kieffer-Kwon K. R., et al. (2012) DNA damage defines sites of recurrent chromosomal translocations in B lymphocytes. *Nature* 484: 69-74

Kovalchuk A. L., Ansarah-Sobrinho C., Hakim O., Resch W., Tolarova H., et al. (2012) Mouse model of endemic Burkitt translocations reveals the long-range boundaries of Ig-mediated oncogene deregulation. *Proc Natl Acad Sci U S A* 109: 10972-10977

Rocha P. P., Micsinai M., Kim J. R., Hewitt S. L., Souza P. P., et al. (2012) Close proximity to Igh is a contributing factor to AID-mediated translocations. *Mol Cell* 47: 873-885

Sexton T., Kurukuti S., Mitchell J. A., Umlauf D., Nagano T., et al. (2012) Sensitive detection of chromatin coassociations using enhanced chromosome conformation capture on chip. *Nat Protoc* 7: 1335-1350

Splinter E., de Wit E., van de Werken H. J., Klous P. and de Laat W. (2012) Determining long-range chromatin interactions for selected genomic sites using 4C-seq technology: from fixation to computation. *Methods* 58: 221-230

van de Werken H. J., Landan G., Holwerda S. J., Hoichman M., Klous P., et al. (2012) Robust 4C-seq data analysis to screen for regulatory DNA interactions. *Nat Methods* 9: 969-972

Hakim O., Sung M. H., Voss T. C., Splinter E., John S., et al. (2011) Diverse gene reprogramming events occur in the same spatial clusters of distal regulatory elements. *Genome Res* 21: 697-706

Oliveira T., Resch W., Jankovic M., Casellas R., Nussenzweig M. C., et al. (2011) Translocation capture sequencing: A method for high throughput mapping of chromosomal rearrangements. *J Immunol Methods* 375: 176-181

Robyr D., Friedli M., Gehrig C., Arcangeli M., Marin M., et al. (2011) Chromosome conformation capture uncovers potential genome-wide interactions between human conserved non-coding sequences. *PLoS ONE* 6: e17634

Rodley C. D., Pai D. A., Mills T. A., Engelke D. R. and O'Sullivan J. M. (2011) tRNA gene identity affects nuclear positioning. *PLoS One* 6: e29267

#### **Associated Kits**

TruSeq ChIP-Seq Kit

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit


Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Mate Pair Kit

## CHROMATIN CONFORMATION CAPTURE CARBON COPY (5-C)

Chromatin conformation capture carbon copy (5-C)<sup>118</sup> allows concurrent determination of interactions between multiple sequences and is a high-throughput version of 3-C<sup>119</sup>. In this method, DNA-protein complexes are crosslinked using formaldehyde. The sample is fragmented and the DNA ligated and digested. The resulting DNA fragments are amplified using ligation-mediated PCR and sequenced. Deep sequencing provides base-pair resolution of ligated fragments.



Pros	Cons
<ul style="list-style-type: none"> <li>• Different from 4-C, 5C provides a matrix of interaction frequencies for many pairs of sites</li> <li>• Can be used to reconstruct the (average) 3D conformation of larger genomic regions<sup>120</sup></li> </ul>	<ul style="list-style-type: none"> <li>• Detection may not necessarily mean an interaction, resulting from random chromosomal collisions</li> <li>• Needs further confirmation of interaction</li> <li>• Cannot scale to genome-wide studies that would require large amount of primers</li> </ul>

### References

Nora E. P., Lajoie B. R., Schulz E. G., Giorgetti L., Okamoto I., et al. (2012) Spatial partitioning of the regulatory landscape of the X-inactivation centre. *Nature* 485: 381-385

The authors use 5-C to analyze regulation of Xist, a non-protein coding transcript that is controlled by X-inactivation center (Xic) to initiate X chromosome inactivation in mouse. They identify a regulatory region of Xist antisense unit that produces a long overriding RNA.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

Sanyal A., Lajoie B. R., Jain G. and Dekker J. (2012) The long-range interaction landscape of gene promoters. *Nature* 489: 109-113

### Associated Kits

- TruSeq ChIP-Seq Kit
- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Mate Pair Kit

118 Dostie J. and Dekker J. (2007) Mapping networks of physical interactions between genomic elements using 5C technology. *Nat Protoc* 2: 988-1002

119 Sajan S. A. and Hawkins R. D. (2012) Methods for identifying higher-order chromatin structure. *Annu Rev Genomics Hum Genet* 13: 59-82

120 de Wit E. and de Laat W. (2012) A decade of 3C technologies: insights into nuclear organization. *Genes Dev* 26: 11-24

## SEQUENCE REARRANGEMENTS

A growing body of evidence suggests that somatic genomic rearrangements, such as retrotransposition and copy number variants (CNVs), are relatively common in healthy individuals<sup>121,122,123</sup>. Cancer genomes are also known to contain numerous complex rearrangements<sup>124</sup>. While many of these rearrangements can be detected during routine next-generation sequencing, specific techniques are available to study rearrangements such as transposable elements.

Transposable genetic elements (TEs) comprise a vast array of DNA sequences with the ability to move to new sites in genomes either directly by a cut-and-paste mechanism (transposons) or indirectly through an RNA intermediate (retrotransposons)<sup>125</sup>. TEs make up about 66-69% of the human genome<sup>126</sup> and play roles in ageing, cancers, brain, development, embryogenesis, and phenotypic variation in populations and evolution. TEs played a major role in dynamic arrangement of the sex determining region over evolution, giving us distinct X and Y chromosomes<sup>127</sup>.

Along with sequence rearrangements by TEs, chromosome and centromere rearrangements can also lead to multiple diseases and disorders<sup>128</sup>. Prenatal diagnostics to study rearrangements predict genetic abnormalities in the fetus. The role of specific TEs and the primary mechanism of chromosome and centromere rearrangements have yet to be elucidated; studying them will help understand their roles.



Transposable elements involved in the evolution of sex chromosomes.

121 O'Huallachain M., Weissman S. M. and Snyder M. P. (2013) The variable somatic genome. *Cell Cycle* 12: 5-6

122 Macosko E. Z. and McCarroll S. A. (2013) Genetics. Our fallen genomes. *Science* 342: 564-565

123 McConnell M. J., Lindberg M. R., Brennand K. J., Piper J. C., Voet T., et al. (2013) Mosaic copy number variation in human neurons. *Science* 342: 632-637

124 Malhotra A., Lindberg M., Faust G. G., Leibowitz M. L., Clark R. A., et al. (2013) Breakpoint profiling of 64 cancer genomes reveals numerous complex rearrangements spawned by homology-independent mechanisms. *Genome Res* 23: 762-776

125 Fedoroff N. V. (2012) Presidential address. Transposable elements, epigenetics, and genome evolution. *Science* 338: 758-767

126 Grandi F. C. and An W. (2013) Non-LTR retrotransposons and microsatellites: Partners in genomic variation. *Mob Genet Elements* 3: e25674

127 Gschwend A. R., Weingartner L. A., Moore R. C. and Ming R. (2012) The sex-specific region of sex chromosomes in animals and plants. *Chromosome Res* 20: 57-69

128 Chiarle R. (2013) Translocations in normal B cells and cancers: insights from new technical approaches. *Adv Immunol* 117: 39-71

## Reviews

Bunting S. F. and Nussenzweig A. (2013) End-joining, translocations and cancer. *Nat Rev Cancer* 13: 443-454

Chiarle R. (2013) Translocations in normal B cells and cancers: insights from new technical approaches. *Adv Immunol* 117: 39-71

Gifford W. D., Pfaff S. L. and Macfarlan T. S. (2013) Transposable elements as genetic regulatory substrates in early development. *Trends Cell Biol* 23: 218-226

Grandi F. C. and An W. (2013) Non-LTR retrotransposons and microsatellites: Partners in genomic variation. *Mob Genet Elements* 3: e25674

van Opijnen T. and Camilli A. (2013) Transposon insertion sequencing: a new tool for systems-level analysis of microorganisms. *Nat Rev Microbiol* 11: 435-442

Burns K. H. and Boeke J. D. (2012) Human transposon tectonics. *Cell* 149: 740-752

Fedoroff N. V. (2012) Presidential address. Transposable elements, epigenetics, and genome evolution. *Science* 338: 758-767


Gschwend A. R., Weingartner L. A., Moore R. C. and Ming R. (2012) The sex-specific region of sex chromosomes in animals and plants. *Chromosome Res* 20: 57-69

Hancks D. C. and Kazazian H. H., Jr. (2012) Active human retrotransposons: variation and disease. *Curr Opin Genet Dev* 22: 191-203

Febrer M., McLay K., Caccamo M., Twomey K. B. and Ryan R. P. (2011) Advances in bacterial transcriptome and transposon insertion-site profiling using second-generation sequencing. *Trends Biotechnol* 29: 586-594

## RETROTRANSPOSON CAPTURE SEQUENCING (RC-SEQ)

Retrotransposon capture sequencing (RC-Seq) is a high-throughput protocol to map and study retrotransposon insertions<sup>129</sup>. In this method, after genomic DNA is fractionated, retrotransposon binding sites on DNA hybridize to transposon binding sites on a microarray. Deep sequencing provides accurate information that can be aligned to a reference sequence to discover novel retrotransposition events.



Pros	Cons
<ul style="list-style-type: none"> <li>• Ability to clearly identify and detect novel retrotransposition events</li> <li>• Can specifically study transposon binding sites of interest</li> <li>• High-throughput protocol</li> </ul>	<ul style="list-style-type: none"> <li>• Different types of MEI require separate PCR experiments with different primers<sup>130</sup></li> <li>• Hybridization errors can lead to sequencing unwanted DNA fragments</li> <li>• PCR biases can underrepresent GC-rich templates</li> <li>• Similar transposition binding sites can lead to sequence ambiguity and detection for a transposition event</li> </ul>

### References

Shukla R., Upton K. R., Munoz-Lopez M., Gerhardt D. J., Fisher M. E., et al. (2013) Endogenous retrotransposition activates oncogenic pathways in hepatocellular carcinoma. *Cell* 153: 101-111

LINE-1 (L1) retrotransposons are mobile genetic elements comprising ~17% of the human genome. To investigate the significance of novel L1 insertions in cancer, this study used RC-Seq on an Illumina HiSeq 2000 system for 19 hepatocellular carcinoma (HCC) and colorectal cancers (MCC). From these data, the authors identified novel L1 insertion events: each individual genome contained on average 244 non-reference L1 insertions. Forty-five non-reference insertions were annotated as tumor-specific and three of these insertions coincided with strong inhibition of the tumor suppressor MCC. These data provide substantial evidence for L1-mediated retrotransposition playing a role in HCC development.

**Illumina Technology:** HiSeq 2000

129 Baillie J. K., Barnett M. W., Upton K. R., Gerhardt D. J., Richmond T. A., et al. (2011) Somatic retrotransposition alters the genetic landscape of the human brain. *Nature* 479: 534-537

130 Xing J., Witherspoon D. J. and Jorde L. B. (2013) Mobile element biology: new possibilities with high-throughput sequencing. *Trends Genet* 29: 280-289

Solyom S., Ewing A. D., Rahrmann E. P., Doucet T., Nelson H. H., et al. (2012) Extensive somatic L1 retrotransposition in colorectal tumors. *Genome Res* 22: 2328-2338

Baillie J. K., Barnett M. W., Upton K. R., Gerhardt D. J., Richmond T. A., et al. (2011) Somatic retrotransposition alters the genetic landscape of the human brain. *Nature* 479: 534-537

### **Associated Kits**

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit

Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Mate Pair Kit

Nextera Rapid Capture Exome/Custom Enrichment Kit



## TRANSPOSON SEQUENCING (TN-SEQ) OR INSERTION SEQUENCING (INSEQ)

Transposon sequencing (Tn-Seq) or insertion sequencing (INSeq) accurately determines quantitative genetic interactions<sup>131</sup>. In this method, a transposon with flanking MmeI digestion sites is transposed into bacteria which, after culturing, can help detect the frequency of mutations within the transposon. After MmeI digestion and subsequent adapter ligation, PCR amplification and sequencing can provide information about the transposon insertion sites.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can study mutational frequency of transposons</li> <li>• Method can be used to deduce fitness of genes within microorganisms</li> <li>• Protocol is robust, reproducible, and sensitive</li> </ul>	<ul style="list-style-type: none"> <li>• Limited to bacterial studies</li> <li>• Errors during PCR amplification can lead to inaccurate sequence reads</li> </ul>

### References

Dong T. G., Ho B. T., Yoder-Himes D. R. and Mekalanos J. J. (2013) Identification of T6SS-dependent effector and immunity proteins by Tn-seq in *Vibrio cholerae*. *Proc Natl Acad Sci U S A* 110: 2623-2628.

T6SS is an important protein for bacterial competition; however, T6SS-dependent effector and immunity proteins have not yet been determined. In this study, the authors use Tn-Seq to identify these proteins in *Vibrio cholerae*.

**Illumina Technology:** HiSeq 2000

Troy E. B., Lin T., Gao L., Lazinski D. W., Camilli A., et al. (2013) Understanding barriers to *Borrelia burgdorferi* dissemination during infection using massively parallel sequencing. *Infect Immun* 81: 2347-2357

Infection by *Borrelia burgdorferi* can cause chronic infections of skin, heart, joints, and the central nervous system of infected mammalian hosts. In this study, the authors characterized the population dynamics of mixed populations of *B. burgdorferi* during infection in a mouse model. Using Tn-Seq based on Illumina technology, they mapped the compositions of *B. burgdorferi* at both the injection site and in distal tissues. The authors found that the infection site was a population bottleneck that significantly altered the composition of the population; however, no such bottleneck was observed in colonization of distal tissues.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>131</sup> van Opijnen T., Bodi K. L. and Camilli A. (2009) Tn-seq: high-throughput parallel sequencing for fitness and genetic interaction studies in microorganisms. *Nat Methods* 6: 767-772


- Murray S. M., Panis G., Fumeaux C., Viollier P. H. and Howard M. (2013) Computational and genetic reduction of a cell cycle to its simplest, primordial components. *PLoS Biol* 11: e1001749
- Sarmiento F., Mrázek J. and Whitman W. B. (2013) Genome-scale analysis of gene function in the hydrogenotrophic methanogenic archaeon *Methanococcus maripaludis*. *Proc Natl Acad Sci U S A* 110: 4726-4731
- Skurnik D., Roux D., Cattoir V., Danilchanka O., Lu X., et al. (2013) Enhanced in vivo fitness of carbapenem-resistant *oprD* mutants of *Pseudomonas aeruginosa* revealed through high-throughput sequencing. *Proc Natl Acad Sci U S A* 110: 20747-20752
- Mann B., van Opijnen T., Wang J., Obert C., Wang Y. D., et al. (2012) Control of virulence by small RNAs in *Streptococcus pneumoniae*. *PLoS Pathog* 8: e1002788
- Qi X., Daily K., Nguyen K., Wang H., Mayhew D., et al. (2012) Retrotransposon profiling of RNA polymerase III initiation sites. *Genome Res* 22: 681-692
- van Opijnen T. and Camilli A. (2012) A fine scale phenotype-genotype virulence map of a bacterial pathogen. *Genome Res* 22: 2541-2551
- Zomer A., Burghout P., Bootsma H. J., Hermans P. W. and van Hijum S. A. (2012) ESSENTIALS: Software for Rapid Analysis of High Throughput Transposon Insertion Sequencing Data. *PLoS ONE* 7: e43012
- Goodman A. L., Wu M. and Gordon J. I. (2011) Identifying microbial fitness determinants by insertion sequencing using genome-wide transposon mutant libraries. *Nat Protoc* 6: 1969-1980

#### **Associated Kits**

- TruSeq Nano DNA Sample Prep Kit
- TruSeq DNA Sample Prep Kit
- TruSeq DNA PCR-Free Sample Prep Kit
- Nextera DNA Sample Prep Kit
- Nextera XT DNA Sample Prep Kit
- Nextera Mate Pair Kit
- Nextera Rapid Capture Exome/Custom Enrichment Kit

## TRANSLOCATION-CAPTURE SEQUENCING (TC-SEQ)

Translocation-capture sequencing (TC-Seq) is a method developed to study chromosomal rearrangements and translocations<sup>132</sup>. In this method, cells are infected with retrovirus expressing I-SceI sites in cells with and without activation-induced cytidine deaminase (AICDA or AID) protein. Genomic DNA from cells is sonicated, linker-ligated, purified, and amplified via semi-nested LM-PCR. The linker is then cleaved and the DNA is sequenced. Any AID-dependent chromosomal rearrangement will be amplified by LM-PCR, while AID-independent translocations will be discarded.



Pros	Cons
<ul style="list-style-type: none"> <li>• Can study chromosomal translocations within a given model or environment</li> <li>• Random sonication generates unique linker ligation points, and deep sequencing allows reading through rearrangement breakpoints</li> </ul>	<ul style="list-style-type: none"> <li>• PCR amplification errors</li> <li>• Non-linear PCR amplification can lead to biases affecting reproducibility</li> <li>• PCR biases can underrepresent GC-rich templates</li> </ul>

### References

Jankovic M., Feldhahn N., Oliveira T. Y., Silva I. T., Kieffer-Kwon K. R., et al. (2013) 53BP1 alters the landscape of DNA rearrangements and suppresses AID-induced B cell lymphoma. *Mol Cell* 49: 623-631

Programmed DNA rearrangement in lymphocytes is initiated by AID protein. The overexpression of AID is associated with cancer, but overexpression of AID alone is insufficient to produce malignancy. This study examines the roles of AID and tumor suppressor p53-binding protein 1 (53BP1) in combination. The results show that the combination of 53BP1 deficiency and AID deregulation increases the rate of rearrangements and results in B cell lymphoma in a mouse model. The rate of rearrangements and CNVs are studied using the Illumina Genome Analyzer.

**Illumina Technology:** Genome Analyzer<sub>ix</sub>

<sup>132</sup> Klein I. A., Resch W., Jankovic M., Oliveira T., Yamane A., et al. (2011) Translocation-capture sequencing reveals the extent and nature of chromosomal rearrangements in B lymphocytes. *Cell* 147: 95-106

Kovalchuk A. L., Ansarah-Sobrinho C., Hakim O., Resch W., Tolarova H., et al. (2012) Mouse model of endemic Burkitt translocations reveals the long-range boundaries of Ig-mediated oncogene deregulation. *Proc Natl Acad Sci U S A* 109: 10972-10977

Oliveira T. Y., Resch W., Jankovic M., Casellas R., Nussenzweig M. C., et al. (2012) Translocation capture sequencing: a method for high throughput mapping of chromosomal rearrangements. *J Immunol Methods* 375: 176-181

Qi X., Daily K., Nguyen K., Wang H., Mayhew D., et al. (2012) Retrotransposon profiling of RNA polymerase III initiation sites. *Genome Res* 22: 681-692

Rocha P. P., Micsinai M., Kim J. R., Hewitt S. L., Souza P. P., et al. (2012) Close proximity to Igh is a contributing factor to AID-mediated translocations. *Mol Cell* 47: 873-885

Klein I. A., Resch W., Jankovic M., Oliveira T., Yamane A., et al. (2011) Translocation-capture sequencing reveals the extent and nature of chromosomal rearrangements in B lymphocytes. *Cell* 147: 95-106

### **Associated Kits**

TruSeq Nano DNA Sample Prep Kit

TruSeq DNA Sample Prep Kit

TruSeq DNA PCR-Free Sample Prep Kit

Nextera DNA Sample Prep Kit

Nextera XT DNA Sample Prep Kit

Nextera Mate Pair Kit

Nextera Rapid Capture Exome/Custom Enrichment Kit

## BIBLIOGRAPHY

- Adelman K. and Lis J. T. (2012) Promoter-proximal pausing of RNA polymerase II: emerging roles in metazoans. *Nat Rev Genet* 13: 720-731
- Adey A. and Shendure J. (2012) Ultra-low input, tagmentation-based whole-genome bisulfite sequencing. *Genome Res* 22: 1139-1143
- Allan J., Fraser R. M., Owen-Hughes T. and Keszenman-Pereyra D. (2012) Micrococcal nuclease does not substantially bias nucleosome mapping. *J Mol Biol* 417: 152-164
- Anderson S. (1981) Shotgun DNA sequencing using cloned DNase I-generated fragments. *Nucleic Acids Res* 9: 3015-3027
- Audit B., Zaghloul L., Vaillant C., Chevereau G., d'Aubenton-Carafa Y., et al. (2009) Open chromatin encoded in DNA sequence is the signature of 'master' replication origins in human cells. *Nucleic Acids Res* 37: 6064-6075
- Auerbach R. K., Euskirchen G., Rozowsky J., Lamarre-Vincent N., Moqtaderi Z., et al. (2009) Mapping accessible chromatin regions using Sono-Seq. *Proc Natl Acad Sci U S A* 106: 14926-14931
- Baillie J. K., Barnett M. W., Upton K. R., Gerhardt D. J., Richmond T. A., et al. (2011) Somatic retrotransposition alters the genetic landscape of the human brain. *Nature* 479: 534-537
- Ball M. P., Li J. B., Gao Y., Lee J. H., LeProust E. M., et al. (2009) Targeted and genome-scale strategies reveal gene-body methylation signatures in human cells. *Nat Biotechnol* 27: 361-368
- Ballare C., Castellano G., Gaveglia L., Althammer S., Gonzalez-Vallinas J., et al. (2013) Nucleosome-driven transcription factor binding and gene regulation. *Mol Cell* 49: 67-79
- Baslan T., Kendall J., Rodgers L., Cox H., Riggs M., et al. (2012) Genome-wide copy number analysis of single cells. *Nat Protoc* 7: 1024-1041
- Becker A. H., Oh E., Weissman J. S., Kramer G. and Bukau B. (2013) Selective ribosome profiling as a tool for studying the interaction of chaperones and targeting factors with nascent polypeptide chains and ribosomes. *Nat Protoc* 8: 2212-2239
- Belton J. M., McCord R. P., Gibcus J. H., Naumova N., Zhan Y., et al. (2012) Hi-C: a comprehensive technique to capture the conformation of genomes. *Methods* 58: 268-276
- Berkseth M., Ikegami K., Arur S., Lieb J. D. and Zarkower D. (2013) TRA-1 ChIP-seq reveals regulators of sexual differentiation and multilevel feedback in nematode sex determination. *Proc Natl Acad Sci U S A* 110: 16033-16038
- Bertucci P. Y., Nacht A. S., Allo M., Rocha-Viegas L., Ballare C., et al. (2013) Progesterone receptor induces bcl-x expression through intragenic binding sites favoring RNA polymerase II elongation. *Nucleic Acids Res* 41: 6072-6086
- Bhargava V., Ko P., Willems E., Mercola M. and Subramaniam S. (2013) Quantitative transcriptomics using designed primer-based amplification. *Sci Rep* 3: 1740
- Bhargava V., Head S. R., Ordoukhanian P., Mercola M. and Subramaniam S. (2014) Technical variations in low-input RNA-seq methodologies. *Sci Rep* 4: 3678
- Bian C. and Yu X. (2013) PGC7 suppresses TET3 for protecting DNA methylation. *Nucleic Acids Res* Biancolella M., B K. F., Tring S., Plummer S. J., Mendoza-Fandino G. A., et al. (2013) Identification and characterization of functional risk variants for colorectal cancer mapping to chromosome 11q23.1. *Hum Mol Genet*
- Blainey P. C. (2013) The future is now: single-cell genomics of bacteria and archaea. *FEMS Microbiol Rev* 37: 407-427
- Blaschke K., Ebata K. T., Karimi M. M., Zepeda-Martinez J. A., Goyal P., et al. (2013) Vitamin C induces Tet-dependent DNA demethylation and a blastocyst-like state in ES cells. *Nature* 500: 222-226
- Bock C., Tomazou E. M., Brinkman A. B., Muller F., Simmer F., et al. (2010) Quantitative comparison of genome-wide DNA methylation mapping technologies. *Nat Biotechnol* 28: 1106-1114
- Bock C., Kiskinis E., Verstappen G., Gu H., Boulting G., et al. (2011) Reference Maps of human ES and iPS cell variation enable high-throughput characterization of pluripotent cell lines. *Cell* 144: 439-452
- Booth M. J., Branco M. R., Ficiz G., Oxley D., Krueger F., et al. (2012) Quantitative sequencing of 5-methylcytosine and 5-hydroxymethylcytosine at single-base resolution. *Science* 336: 934-937
- Booth M. J., Ost T. W., Beraldi D., Bell N. M., Branco M. R., et al. (2013) Oxidative bisulfite sequencing of 5-methylcytosine and 5-hydroxymethylcytosine. *Nat Protoc* 8: 1841-1851

- Böttcher R., Amberg R., Ruzius F. P., Guryev V., Verhaegh W. F., et al. (2012) Using a priori knowledge to align sequencing reads to their exact genomic position. *Nucleic Acids Res* 40: e125
- Bowman S. K., Simon M. D., Deaton A. M., Tolstorukov M., Borowsky M. L., et al. (2013) Multiplexed Illumina sequencing libraries from picogram quantities of DNA. *BMC Genomics* 14: 466
- Boyd-Kirkup J. D., Green C. D., Wu G., Wang D. and Han J. D. (2013) Epigenomics and the regulation of aging. *Epigenomics* 5: 205-227
- Boyle A. P., Song L., Lee B. K., London D., Keefe D., et al. (2011) High-resolution genome-wide in vivo footprinting of diverse transcription factors in human cells. *Genome Res* 21: 456-464
- Brinkman A. B., Simmer F., Ma K., Kaan A., Zhu J., et al. (2010) Whole-genome DNA methylation profiling using MethylCap-seq. *Methods* 52: 232-236
- Brinkman A. B., Gu H., Bartels S. J., Zhang Y., Matarese F., et al. (2012) Sequential ChIP-bisulfite sequencing enables direct genome-scale investigation of chromatin and DNA methylation cross-talk. *Genome Res* 22: 1128-1138
- Broughton J. P. and Pasquinelli A. E. (2013) Identifying Argonaute binding sites in *Caenorhabditis elegans* using iCLIP. *Methods* 63: 119-125
- Buenrostro J. D., Giresi P. G., Zaba L. C., Chang H. Y. and Greenleaf W. J. (2013) Transposition of native chromatin for fast and sensitive epigenomic profiling of open chromatin, DNA-binding proteins and nucleosome position. *Nat Methods* 10: 1213-1218
- Bunting S. F. and Nussenzweig A. (2013) End-joining, translocations and cancer. *Nat Rev Cancer* 13: 443-454
- Burns K. H. and Boeke J. D. (2012) Human transposon tectonics. *Cell* 149: 740-752
- Burton J. N., Adey A., Patwardhan R. P., Qiu R., Kitzman J. O., et al. (2013) Chromosome-scale scaffolding of de novo genome assemblies based on chromatin interactions. *Nat Biotechnol* 31: 1119-1125
- Calabrese J. M., Sun W., Song L., Mugford J. W., Williams L., et al. (2012) Site-specific silencing of regulatory elements as a mechanism of X inactivation. *Cell* 151: 951-963
- Capell B. C. and Berger S. L. (2013) Genome-wide epigenetics. *J Invest Dermatol* 133: e9
- Castello A., Fischer B., Hentze M. W. and Preiss T. (2013) RNA-binding proteins in Mendelian disease. *Trends Genet* 29: 318-327
- Cernilogar F. M., Onorati M. C., Kothe G. O., Burroughs A. M., Parsi K. M., et al. (2011) Chromatin-associated RNA interference components contribute to transcriptional regulation in *Drosophila*. *Nature* 480: 391-395
- Chai X., Nagarajan S., Kim K., Lee K. and Choi J. K. (2013) Regulation of the boundaries of accessible chromatin. *PLoS Genet* 9: e1003778
- Chang C. Y., Pasolli H. A., Giannopoulou E. G., Guasch G., Gronostajski R. M., et al. (2013) NFIB is a governor of epithelial-melanocyte stem cell behaviour in a shared niche. *Nature* 495: 98-102
- Charizanis K., Lee K. Y., Batra R., Goodwin M., Zhang C., et al. (2012) Muscleblind-like 2-mediated alternative splicing in the developing brain and dysregulation in myotonic dystrophy. *Neuron* 75: 437-450
- Chauhan C., Zraly C. B. and Dingwall A. K. (2013) The *Drosophila* COMPASS-like Cmi-Trr coactivator complex regulates dpp/BMP signaling in pattern formation. *Dev Biol* 380: 185-198
- Chaumeil J., Micsinai M., Ntziachristos P., Roth D. B., Aifantis I., et al. (2013) The RAG2 C-terminus and ATM protect genome integrity by controlling antigen receptor gene cleavage. *Nat Commun* 4: 2231
- Chavez L., Jozefczuk J., Grimm C., Dietrich J., Timmermann B., et al. (2010) Computational analysis of genome-wide DNA methylation during the differentiation of human embryonic stem cells along the endodermal lineage. *Genome Res* 20: 1441-1450
- Chen Y., Negre N., Li Q., Mieczkowska J. O., Slattery M., et al. (2012) Systematic evaluation of factors influencing ChIP-seq fidelity. *Nat Methods* 9: 609-614
- Chi S. W., Zang J. B., Mele A. and Darnell R. B. (2009) Argonaute HITS-CLIP decodes microRNA-mRNA interaction maps. *Nature* 460: 479-486
- Chi S. W., Hannon G. J. and Darnell R. B. (2012) An alternative mode of microRNA target recognition. *Nat Struct Mol Biol* 19: 321-327
- Chiarle R. (2013) Translocations in normal B cells and cancers: insights from new technical approaches. *Adv Immunol* 117: 39-71

- Chitsaz H., Yee-Greenbaum J. L., Tesler G., Lombardo M. J., Dupont C. L., et al. (2011) Efficient de novo assembly of single-cell bacterial genomes from short-read data sets. *Nat Biotechnol* 29: 915-921
- Cho J., Chang H., Kwon S. C., Kim B., Kim Y., et al. (2012) LIN28A is a suppressor of ER-associated translation in embryonic stem cells. *Cell* 151: 765-777
- Chopra V. S., Hendrix D. A., Core L. J., Tsui C., Lis J. T., et al. (2011) The Polycomb Group Mutant esc Leads to Augmented Levels of Paused Pol II in the Drosophila Embryo. *Mol Cell* 42: 837-844
- Chu C., Qu K., Zhong F. L., Artandi S. E. and Chang H. Y. (2011) Genomic maps of long noncoding RNA occupancy reveal principles of RNA-chromatin interactions. *Mol Cell* 44: 667-678
- Chumsakul O., Nakamura K., Kurata T., Sakamoto T., Hobman J. L., et al. (2013) High-resolution mapping of in vivo genomic transcription factor binding sites using in situ DNase I footprinting and ChIP-seq. *DNA Res* 20: 325-338
- Churchman L. S. and Weissman J. S. (2011) Nascent transcript sequencing visualizes transcription at nucleotide resolution. *Nature* 469: 368-373
- Clouaire T., Webb S., Skene P., Illingworth R., Kerr A., et al. (2012) Cfp1 integrates both CpG content and gene activity for accurate H3K4me3 deposition in embryonic stem cells. *Genes Dev* 26: 1714-1728
- Cokus S. J., Feng S., Zhang X., Chen Z., Merriman B., et al. (2008) Shotgun bisulphite sequencing of the Arabidopsis genome reveals DNA methylation patterning. *Nature* 452: 215-219
- Colquitt B. M., Allen W. E., Barnea G. and Lomvardas S. (2013) Alteration of genic 5-hydroxymethylcytosine patterning in olfactory neurons correlates with changes in gene expression and cell identity. *Proc Natl Acad Sci U S A* 110: 14682-14687
- Core L. J., Waterfall J. J. and Lis J. T. (2008) Nascent RNA sequencing reveals widespread pausing and divergent initiation at human promoters. *Science* 322: 1845-1848
- Core L. J., Waterfall J. J., Gilchrist D. A., Fargo D. C., Kwak H., et al. (2012) Defining the status of RNA polymerase at promoters. *Cell Rep* 2: 1025-1035
- Cusick M. E., Herman T. M., DePamphilis M. L. and Wassarman P. M. (1981) Structure of chromatin at deoxyribonucleic acid replication forks: prenucleosomal deoxyribonucleic acid is rapidly excised from replicating simian virus 40 chromosomes by micrococcal nuclease. *Biochemistry* 20: 6648-6658
- Darnell J. C., Van Driesche S. J., Zhang C., Hung K. Y., Mele A., et al. (2011) FMRP Stalls Ribosomal Translocation on mRNAs Linked to Synaptic Function and Autism. *Cell* 146: 247-261
- Davis H. P. and Squire L. R. (1984) Protein synthesis and memory: a review. *Psychol Bull* 96: 518-559
- de Wit E. and de Laat W. (2012) A decade of 3C technologies: insights into nuclear organization. *Genes Dev* 26: 11-24
- de Wit E., Bouwman B. A., Zhu Y., Klous P., Splinter E., et al. (2013) The pluripotent genome in three dimensions is shaped around pluripotency factors. *Nature* 501: 227-231
- Dean F. B., Nelson J. R., Giesler T. L. and Lasken R. S. (2001) Rapid amplification of plasmid and phage DNA using Phi 29 DNA polymerase and multiply-primed rolling circle amplification. *Genome Res* 11: 1095-1099
- Degner J. F., Pai A. A., Pique-Regi R., Veyrieras J. B., Gaffney D. J., et al. (2012) DNase I sensitivity QTLs are a major determinant of human expression variation. *Nature* 482: 390-394
- Dekker J., Marti-Renom M. A. and Mirny L. A. (2013) Exploring the three-dimensional organization of genomes: interpreting chromatin interaction data. *Nat Rev Genet* 14: 390-403
- Delpretti S., Montavon T., Leleu M., Joye E., Tzika A., et al. (2013) Multiple Enhancers Regulate Hoxd Genes and the Hotdog LncRNA during Cecum Budding. *Cell Rep* 5: 137-150
- DeMare L. E., Leng J., Cotney J., Reilly S. K., Yin J., et al. (2013) The genomic landscape of cohesin-associated chromatin interactions. *Genome Res* 23: 1224-1234
- Demichelis F., Setlur S. R., Banerjee S., Chakravarty D., Chen J. Y., et al. (2012) Identification of functionally active, low frequency copy number variants at 15q21.3 and 12q21.31 associated with prostate cancer risk. *Proc Natl Acad Sci U S A* 109: 6686-6691
- Deng T., Zhu Z. I., Zhang S., Leng F., Cherukuri S., et al. (2013) HMG1 modulates nucleosome occupancy and DNase I hypersensitivity at the CpG island promoters of embryonic stem cells. *Mol Cell Biol* 33: 3377-3389
- Denholtz M., Bonora G., Chronis C., Splinter E., de Laat W., et al. (2013) Long-Range Chromatin Contacts in Embryonic Stem Cells Reveal a Role for Pluripotency Factors and Polycomb Proteins in Genome Organization. *Cell Stem Cell* 13: 602-616
- Di Ruscio A., Ebralidze A. K., Benoukraf T., Amabile G., Goff L. A., et al. (2013) DNMT1-interacting RNAs block gene-specific DNA methylation. *Nature* 503: 371-376

- Diep D., Plongthongkum N., Gore A., Fung H. L., Shoemaker R., et al. (2012) Library-free methylation sequencing with bisulfite padlock probes. *Nat Methods* 9: 270-272
- Dixon J. R., Selvaraj S., Yue F., Kim A., Li Y., et al. (2012) Topological domains in mammalian genomes identified by analysis of chromatin interactions. *Nature* 485: 376-380
- Dolinoy D. C., Weidman J. R., Waterland R. A. and Jirtle R. L. (2006) Maternal genistein alters coat color and protects Avy mouse offspring from obesity by modifying the fetal epigenome. *Environ Health Perspect* 114: 567-572
- Dolinoy D. C. (2008) The agouti mouse model: an epigenetic biosensor for nutritional and environmental alterations on the fetal epigenome. *Nutr Rev* 66 Suppl 1: S7-11
- Dolinoy D. C. and Faulk C. (2012) Introduction: The use of animal models to advance epigenetic science. *ILAR J* 53: 227-231
- Dong T. G., Ho B. T., Yoder-Himes D. R. and Mekalanos J. J. (2013) Identification of T6SS-dependent effector and immunity proteins by Tn-seq in *Vibrio cholerae*. *Proc Natl Acad Sci U S A* 110: 2623-2628
- Dostie J. and Dekker J. (2007) Mapping networks of physical interactions between genomic elements using 5C technology. *Nat Protoc* 2: 988-1002
- Down T. A., Rakyen V. K., Turner D. J., Flicek P., Li H., et al. (2008) A Bayesian deconvolution strategy for immunoprecipitation-based DNA methylome analysis. *Nat Biotechnol* 26: 779-785
- Duan Z., Andronescu M., Schutz K., Lee C., Shendure J., et al. (2012) A genome-wide 3C-method for characterizing the three-dimensional architectures of genomes. *Methods* 58: 277-288
- Dunowska M., Biggs P. J., Zheng T. and Perrott M. R. (2012) Identification of a novel nidovirus associated with a neurological disease of the Australian brushtail possum (*Trichosurus vulpecula*). *Vet Microbiol* 156: 418-424
- Embree M., Nagarajan H., Movahedi N., Chitsaz H. and Zengler K. (2013) Single-cell genome and metatranscriptome sequencing reveal metabolic interactions of an alkane-degrading methanogenic community. *ISME J*
- Farcas A. M., Blackledge N. P., Sudbery I., Long H. K., McGouran J. F., et al. (2012) KDM2B links the Polycomb Repressive Complex 1 (PRC1) to recognition of CpG islands. *Elife* 1: e00205
- Febrer M., McLay K., Caccamo M., Twomey K. B. and Ryan R. P. (2011) Advances in bacterial transcriptome and transposon insertion-site profiling using second-generation sequencing. *Trends Biotechnol* 29: 586-594
- Fedoroff N. V. (2012) Presidential address. Transposable elements, epigenetics, and genome evolution. *Science* 338: 758-767
- Feil R., Charlton J., Bird A. P., Walter J. and Reik W. (1994) Methylation analysis on individual chromosomes: improved protocol for bisulphite genomic sequencing. *Nucleic Acids Res* 22: 695-696
- Feng H., Qin Z. and Zhang X. (2013) Opportunities and methods for studying alternative splicing in cancer with RNA-Seq. *Cancer Lett* 340: 179-191
- Feng S., Cokus S. J., Zhang X., Chen P. Y., Bostick M., et al. (2010) Conservation and divergence of methylation patterning in plants and animals. *Proc Natl Acad Sci U S A* 107: 8689-8694
- Fritsch C., Herrmann A., Nothnagel M., Szafranski K., Huse K., et al. (2012) Genome-wide search for novel human uORFs and N-terminal protein extensions using ribosomal footprinting. *Genome Res* 22: 2208-2218
- Frumkin D., Wasserstrom A., Kaplan S., Feige U. and Shapiro E. (2005) Genomic variability within an organism exposes its cell lineage tree. *PLoS Comput Biol* 1: e50
- Fullwood M. J., Liu M. H., Pan Y. F., Liu J., Xu H., et al. (2009) An oestrogen-receptor-alpha-bound human chromatin interactome. *Nature* 462: 58-64
- Fullwood M. J., Han Y., Wei C. L., Ruan X. and Ruan Y. (2010) Chromatin interaction analysis using paired-end tag sequencing. *Curr Protoc Mol Biol Chapter 21: Unit 21 15 21-25*
- Gaffney D. J., McVicker G., Pai A. A., Fondufe-Mittendorf Y. N., Lewellen N., et al. (2012) Controls of nucleosome positioning in the human genome. *PLoS Genet* 8: e1003036
- Galas D. J. and Schmitz A. (1978) DNase footprinting: a simple method for the detection of protein-DNA binding specificity. *Nucleic Acids Res* 5: 3157-3170
- Gao Q., Steine E. J., Barrasa M. I., Hockemeyer D., Pawlak M., et al. (2011) Deletion of the de novo DNA methyltransferase Dnmt3a promotes lung tumor progression. *Proc Natl Acad Sci U S A* 108: 18061-18066
- Gehring M., Bubb K. L. and Henikoff S. (2009) Extensive demethylation of repetitive elements during seed development underlies gene imprinting. *Science* 324: 1447-1451



- Gendrel A. V., Apedaile A., Coker H., Termanis A., Zvetkova I., et al. (2012) Smchd1-dependent and -independent pathways determine developmental dynamics of CpG island methylation on the inactive x chromosome. *Dev Cell* 23: 265-279
- Gerashchenko M. V., Lobanov A. V. and Gladyshev V. N. (2012) Genome-wide ribosome profiling reveals complex translational regulation in response to oxidative stress. *Proc Natl Acad Sci U S A* 109: 17394-17399
- German M. A., Pillay M., Jeong D. H., Hetawal A., Luo S., et al. (2008) Global identification of microRNA-target RNA pairs by parallel analysis of RNA ends. *Nat Biotechnol* 26: 941-946
- Gertz J., Varley K. E., Reddy T. E., Bowling K. M., Pauli F., et al. (2011) Analysis of DNA methylation in a three-generation family reveals widespread genetic influence on epigenetic regulation. *PLoS Genet* 7: e1002228
- Gifford W. D., Pfaff S. L. and Macfarlan T. S. (2013) Transposable elements as genetic regulatory substrates in early development. *Trends Cell Biol* 23: 218-226
- Giresi P. G. and Lieb J. D. (2009) Isolation of active regulatory elements from eukaryotic chromatin using FAIRE (Formaldehyde Assisted Isolation of Regulatory Elements). *Methods* 48: 233-239
- Goetz J. J. and Trimarchi J. M. (2012) Transcriptome sequencing of single cells with Smart-Seq. *Nat Biotechnol* 30: 763-765
- Gontan C., Achame E. M., Demmers J., Barakat T. S., Rentmeester E., et al. (2012) RNF12 initiates X-chromosome inactivation by targeting REX1 for degradation. *Nature* 485: 386-390
- Goodman A. L., Wu M. and Gordon J. I. (2011) Identifying microbial fitness determinants by insertion sequencing using genome-wide transposon mutant libraries. *Nat Protoc* 6: 1969-1980
- Gottwein E., Corcoran D. L., Mukherjee N., Skalsky R. L., Hafner M., et al. (2011) Viral microRNA targetome of KSHV-infected primary effusion lymphoma cell lines. *Cell Host Microbe* 10: 515-526
- Grandi F. C. and An W. (2013) Non-LTR retrotransposons and microsatellites: Partners in genomic variation. *Mob Genet Elements* 3: e25674
- Grøntved L., John S., Baek S., Liu Y., Buckley J. R., et al. (2013) C/EBP maintains chromatin accessibility in liver and facilitates glucocorticoid receptor recruitment to steroid response elements. *EMBO J* 32: 1568-1583
- Gschwend A. R., Weingartner L. A., Moore R. C. and Ming R. (2012) The sex-specific region of sex chromosomes in animals and plants. *Chromosome Res* 20: 57-69
- Gu H., Smith Z. D., Bock C., Boyle P., Gnirke A., et al. (2011) Preparation of reduced representation bisulfite sequencing libraries for genome-scale DNA methylation profiling. *Nat Protoc* 6: 468-481
- Gu W., Lee H. C., Chaves D., Youngman E. M., Pazour G. J., et al. (2012) CapSeq and CIP-TAP identify Pol II start sites and reveal capped small RNAs as *C. elegans* piRNA precursors. *Cell* 151: 1488-1500
- Guertin M. J., Martins A. L., Siepel A. and Lis J. T. (2012) Accurate prediction of inducible transcription factor binding intensities in vivo. *PLoS Genet* 8: e1002610
- Gustems M., Woellmer A., Rothbauer U., Eck S. H., Wieland T., et al. (2013) c-Jun/c-Fos heterodimers regulate cellular genes via a newly identified class of methylated DNA sequence motifs. *Nucleic Acids Res*
- Habibi E., Brinkman A. B., Arand J., Kroeze L. I., Kerstens H. H., et al. (2013) Whole-genome bisulfite sequencing of two distinct interconvertible DNA methylomes of mouse embryonic stem cells. *Cell Stem Cell* 13: 360-369
- Hafner M., Landgraf P., Ludwig J., Rice A., Ojo T., et al. (2008) Identification of microRNAs and other small regulatory RNAs using cDNA library sequencing. *Methods* 44: 3-12
- Hafner M., Landthaler M., Burger L., Khorshid M., Hausser J., et al. (2010) Transcriptome-wide identification of RNA-binding protein and microRNA target sites by PAR-CLIP. *Cell* 141: 129-141
- Hafner M., Landthaler M., Burger L., Khorshid M., Hausser J., et al. (2010) PAR-CLIP--a method to identify transcriptome-wide the binding sites of RNA binding proteins. *J Vis Exp*
- Hafner M., Lianoglou S., Tuschl T. and Betel D. (2012) Genome-wide identification of miRNA targets by PAR-CLIP. *Methods* 58: 94-105
- Hah N., Danko C. G., Core L., Waterfall J. J., Siepel A., et al. (2011) A rapid, extensive, and transient transcriptional response to estrogen signaling in breast cancer cells. *Cell* 145: 622-634
- Hakim O., Sung M. H., Voss T. C., Splinter E., John S., et al. (2011) Diverse gene reprogramming events occur in the same spatial clusters of distal regulatory elements. *Genome Res* 21: 697-706

- Hakim O., Resch W., Yamane A., Klein I., Kieffer-Kwon K. R., et al. (2012) DNA damage defines sites of recurrent chromosomal translocations in B lymphocytes. *Nature* 484: 69-74
- Hammoud S. S., Nix D. A., Zhang H., Purwar J., Carrell D. T., et al. (2009) Distinctive chromatin in human sperm packages genes for embryo development. *Nature* 460: 473-478
- Han Y., David A., Liu B., Magadan J. G., Bennink J. R., et al. (2012) Monitoring cotranslational protein folding in mammalian cells at codon resolution. *Proc Natl Acad Sci U S A* 109: 12467-12472
- Hancks D. C. and Kazazian H. H., Jr. (2012) Active human retrotransposons: variation and disease. *Curr Opin Genet Dev* 22: 191-203
- Harris R. A., Wang T., Coarfa C., Nagarajan R. P., Hong C., et al. (2010) Comparison of sequencing-based methods to profile DNA methylation and identification of monoallelic epigenetic modifications. *Nat Biotechnol* 28: 1097-1105
- Hashimshony T., Wagner F., Sher N. and Yanai I. (2012) CEL-Seq: single-cell RNA-Seq by multiplexed linear amplification. *Cell Rep* 2: 666-673
- He H. H., Meyer C. A., Chen M. W., Jordan V. C., Brown M., et al. (2012) Differential DNase I hypersensitivity reveals factor-dependent chromatin dynamics. *Genome Res* 22: 1015-1025
- He Y., Vogelstein B., Velculescu V. E., Papadopoulos N. and Kinzler K. W. (2008) The antisense transcriptomes of human cells. *Science* 322: 1855-1857
- Heinz S., Romanoski C. E., Benner C., Allison K. A., Kaikkonen M. U., et al. (2013) Effect of natural genetic variation on enhancer selection and function. *Nature* 503: 487-492
- Heisig J., Weber D., Englberger E., Winkler A., Kneitz S., et al. (2012) Target gene analysis by microarrays and chromatin immunoprecipitation identifies HEY proteins as highly redundant bHLH repressors. *PLoS Genet* 8: e1002728
- Hiatt J. B., Pritchard C. C., Salipante S. J., O'Roak B. J. and Shendure J. (2013) Single molecule molecular inversion probes for targeted, high-accuracy detection of low-frequency variation. *Genome Res* 23: 843-854
- Hilton I. B., Simon J. M., Lieb J. D., Davis I. J., Damania B., et al. (2013) The open chromatin landscape of Kaposi's sarcoma-associated herpesvirus. *J Virol* 87: 11831-11842
- Hodges E., Smith A. D., Kendall J., Xuan Z., Ravi K., et al. (2009) High definition profiling of mammalian DNA methylation by array capture and single molecule bisulfite sequencing. *Genome Res* 19: 1593-1605
- Hogan G. J., Lee C. K. and Lieb J. D. (2006) Cell cycle-specified fluctuation of nucleosome occupancy at gene promoters. *PLoS Genet* 2: e158
- Holt C. E. and Schuman E. M. (2013) The central dogma decentralized: new perspectives on RNA function and local translation in neurons. *Neuron* 80: 648-657
- Holwerda S. J., van de Werken H. J., Ribeiro de Almeida C., Bergen I. M., de Bruijn M. J., et al. (2013) Allelic exclusion of the immunoglobulin heavy chain locus is independent of its nuclear localization in mature B cells. *Nucleic Acids Res* 41: 6905-6916
- Hou Y., Song L., Zhu P., Zhang B., Tao Y., et al. (2012) Single-cell exome sequencing and monoclonal evolution of a JAK2-negative myeloproliferative neoplasm. *Cell* 148: 873-885
- Hou Y., Fan W., Yan L., Li R., Lian Y., et al. (2013) Genome analyses of single human oocytes. *Cell* 155: 1492-1506
- Hsieh A. C., Liu Y., Edlind M. P., Ingolia N. T., Janes M. R., et al. (2012) The translational landscape of mTOR signalling steers cancer initiation and metastasis. *Nature* 485: 55-61
- Hsieh T. F., Ibarra C. A., Silva P., Zemach A., Eshed-Williams L., et al. (2009) Genome-wide demethylation of Arabidopsis endosperm. *Science* 324: 1451-1454
- Huang T. T., Gonzales C. B., Gu F., Hsu Y. T., Jadhav R. R., et al. (2013) Epigenetic deregulation of the anaplastic lymphoma kinase gene modulates mesenchymal characteristics of oral squamous cell carcinomas. *Carcinogenesis* 34: 1717-1727
- Hussain S., Sajini A. A., Blanco S., Dietmann S., Lombard P., et al. (2013) NSun2-mediated cytosine-5 methylation of vault noncoding RNA determines its processing into regulatory small RNAs. *Cell Rep* 4: 255-261
- Illingworth R. S., Gruenewald-Schneider U., Webb S., Kerr A. R., James K. D., et al. (2010) Orphan CpG islands identify numerous conserved promoters in the mammalian genome. *PLoS Genet* 6: e1001134
- Imakaev M., Fudenberg G., McCord R. P., Naumova N., Goloborodko A., et al. (2012) Iterative correction of Hi-C data reveals hallmarks of chromosome organization. *Nat Methods* 9: 999-1003
- Ingolia N. T., Ghaemmaghami S., Newman J. R. and Weissman J. S. (2009) Genome-wide analysis in vivo of translation with nucleotide resolution using ribosome profiling. *Science* 324: 218-223

- Ingolia N. T. (2010) Genome-wide translational profiling by ribosome footprinting. *Methods Enzymol* 470: 119-142
- Ingolia N. T., Lareau L. F. and Weissman J. S. (2011) Ribosome Profiling of Mouse Embryonic Stem Cells Reveals the Complexity and Dynamics of Mammalian Proteomes. *Cell* 147: 789-802
- Irimia M., Tena J. J., Alexis M. S., Fernandez-Minan A., Maeso I., et al. (2012) Extensive conservation of ancient microsynteny across metazoans due to cis-regulatory constraints. *Genome Res* 22: 2356-2367
- Islam S., Kjallquist U., Moliner A., Zajac P., Fan J. B., et al. (2011) Characterization of the single-cell transcriptional landscape by highly multiplex RNA-seq. *Genome Res* 21: 1160-1167
- Islam S., Zeisel A., Joost S., La Manno G., Zajac P., et al. (2014) Quantitative single-cell RNA-seq with unique molecular identifiers. *Nat Methods* 11: 163-166
- Iwata M., Sandstrom R. S., Delrow J. J., Stamatoyannopoulos J. A. and Torok-Storb B. (2013) Functionally and Phenotypically Distinct Subpopulations of Marrow Stromal Cells Are Fibroblast in Origin and Induce Different Fates in Peripheral Blood Monocytes. *Stem Cells Dev*
- Jacob Y., Feng S., Leblanc C. A., Bernatavichute Y. V., Stroud H., et al. (2009) ATXR5 and ATXR6 are H3K27 monomethyltransferases required for chromatin structure and gene silencing. *Nat Struct Mol Biol* 16: 763-768
- Jain A., Bacolla A., Del Mundo I. M., Zhao J., Wang G., et al. (2013) DHX9 helicase is involved in preventing genomic instability induced by alternatively structured DNA in human cells. *Nucleic Acids Res* 41: 10345-10357
- Jakopovic M., Thomas A., Balasubramaniam S., Schrupp D., Giaccone G., et al. (2013) Targeting the Epigenome in Lung Cancer: Expanding Approaches to Epigenetic Therapy. *Front Oncol* 3: 261
- Jankovic M., Feldhahn N., Oliveira T. Y., Silva I. T., Kieffer-Kwon K. R., et al. (2013) 53BP1 alters the landscape of DNA rearrangements and suppresses AID-induced B cell lymphoma. *Mol Cell* 49: 623-631
- Jasmine F., Ahsan H., Andrusis I. L., John E. M., Chang-Claude J., et al. (2008) Whole-genome amplification enables accurate genotyping for microarray-based high-density single nucleotide polymorphism array. *Cancer Epidemiol Biomarkers Prev* 17: 3499-3508
- Ji X., Zhou Y., Pandit S., Huang J., Li H., et al. (2013) SR proteins collaborate with 7SK and promoter-associated nascent RNA to release paused polymerase. *Cell* 153: 855-868
- Jiang L., Zhang J., Wang J. J., Wang L., Zhang L., et al. (2013) Sperm, but not oocyte, DNA methylome is inherited by zebrafish early embryos. *Cell* 153: 773-784
- Jiao Y. and Meyerowitz E. M. (2010) Cell-type specific analysis of translating RNAs in developing flowers reveals new levels of control. *Mol Syst Biol* 6: 419
- Jin F., Li Y., Dixon J. R., Selvaraj S., Ye Z., et al. (2013) A high-resolution map of the three-dimensional chromatin interactome in human cells. *Nature* 503: 290-294
- Jin S. G., Kadam S. and Pfeifer G. P. (2010) Examination of the specificity of DNA methylation profiling techniques towards 5-methylcytosine and 5-hydroxymethylcytosine. *Nucleic Acids Res* 38: e125
- Jullien P. E. and Berger F. (2010) DNA methylation reprogramming during plant sexual reproduction? *Trends Genet* 26: 394-399
- Jumlongras D., Lachke S. A., O'Connell D. J., Aboukhalil A., Li X., et al. (2012) An evolutionarily conserved enhancer regulates Bmp4 expression in developing incisor and limb bud. *PLoS ONE* 7: e38568
- Jungkamp A. C., Stoeckius M., Mecnas D., Grun D., Mastrobuoni G., et al. (2011) In vivo and transcriptome-wide identification of RNA binding protein target sites. *Mol Cell* 44: 828-840
- Kadkhodaei B., Alvarsson A., Schintu N., Ramsköld D., Volakakis N., et al. (2013) Transcription factor Nurr1 maintains fiber integrity and nuclear-encoded mitochondrial gene expression in dopamine neurons. *Proc Natl Acad Sci U S A* 110: 2360-2365
- Kahrstrom C. T. (2013) Epigenetics: Legionella makes its mark on histones. *Nat Rev Genet* 14: 370
- Kaikkonen M. U., Spann N. J., Heinz S., Romanoski C. E., Allison K. A., et al. (2013) Remodeling of the enhancer landscape during macrophage activation is coupled to enhancer transcription. *Mol Cell* 51: 310-325
- Kalisky T. and Quake S. R. (2011) Single-cell genomics. *Nat Methods* 8: 311-314
- Kaneko S., Bonasio R., Saldana-Meyer R., Yoshida T., Son J., et al. (2014) Interactions between JARID2 and Noncoding RNAs Regulate PRC2 Recruitment to Chromatin. *Mol Cell* 53: 290-300
- Kanematsu S., Tanimoto K., Suzuki Y. and Sugano S. (2014) Screening for possible miRNA-mRNA associations in a colon cancer cell line. *Gene* 533: 520-531
- Kennedy S. R., Salk J. J., Schmitt M. W. and Loeb L. A. (2013) Ultra-sensitive sequencing reveals an age-related increase in somatic mitochondrial mutations that are inconsistent with oxidative damage. *PLoS Genet* 9: e1003794

- Kent N. A., Adams S., Moorhouse A. and Paszkiewicz K. (2011) Chromatin particle spectrum analysis: a method for comparative chromatin structure analysis using paired-end mode next-generation DNA sequencing. *Nucleic Acids Res* 39: e26
- Kim M., Park Y. K., Kang T. W., Lee S. H., Rhee Y. H., et al. (2013) Dynamic changes in DNA methylation and hydroxymethylation when hES cells undergo differentiation toward a neuronal lineage. *Hum Mol Genet* 23: 657-667
- Kim Y. J., Greer C. B., Cecchini K. R., Harris L. N., Tuck D. P., et al. (2013) HDAC inhibitors induce transcriptional repression of high copy number genes in breast cancer through elongation blockade. *Oncogene* 32: 2828-2835
- Kishore S., Jaskiewicz L., Burger L., Hausser J., Khorshid M., et al. (2011) A quantitative analysis of CLIP methods for identifying binding sites of RNA-binding proteins. *Nat Methods* 8: 559-564
- Kivioja T., Vaharautio A., Karlsson K., Bonke M., Enge M., et al. (2012) Counting absolute numbers of molecules using unique molecular identifiers. *Nat Methods* 9: 72-74
- Klein I. A., Resch W., Jankovic M., Oliveira T., Yamane A., et al. (2011) Translocation-capture sequencing reveals the extent and nature of chromosomal rearrangements in B lymphocytes. *Cell* 147: 95-106
- Kloosterman W. P. and Plasterk R. H. (2006) The diverse functions of microRNAs in animal development and disease. *Dev Cell* 11: 441-450
- Kobayashi H., Sakurai T., Imai M., Takahashi N., Fukuda A., et al. (2012) Contribution of Intragenic DNA Methylation in Mouse Gametic DNA Methylomes to Establish Oocyte-Specific Heritable Marks. *PLoS Genet* 8: e1002440
- Kobayashi H., Sakurai T., Miura F., Imai M., Mochiduki K., et al. (2013) High-resolution DNA methylome analysis of primordial germ cells identifies gender-specific reprogramming in mice. *Genome Res* 23: 616-627
- Koh K. P. and Rao A. (2013) DNA methylation and methylcytosine oxidation in cell fate decisions. *Curr Opin Cell Biol* 25: 152-16
- Koldamova R., Schug J., Lefterova M., Cronican A. A., Fitz N. F., et al. (2014) Genome-wide approaches reveal EGR1-controlled regulatory networks associated with neurodegeneration. *Neurobiol Dis* 63: 107-114
- Konig J., Zarnack K., Rot G., Curk T., Kayikci M., et al. (2010) iCLIP reveals the function of hnRNP particles in splicing at individual nucleotide resolution. *Nat Struct Mol Biol* 17: 909-915
- Kovalchuk A. L., Ansarah-Sobrinho C., Hakim O., Resch W., Tolarova H., et al. (2012) Mouse model of endemic Burkitt translocations reveals the long-range boundaries of Ig-mediated oncogene deregulation. *Proc Natl Acad Sci U S A* 109: 10972-10977
- Kozlenkov A., Roussos P., Timashpolsky A., Barbu M., Rudchenko S., et al. (2014) Differences in DNA methylation between human neuronal and glial cells are concentrated in enhancers and non-CpG sites. *Nucleic Acids Res* 42: 109-127
- Kudla G., Granneman S., Hahn D., Beggs J. D. and Tollervey D. (2011) Cross-linking, ligation, and sequencing of hybrids reveals RNA-RNA interactions in yeast. *Proc Natl Acad Sci U S A* 108: 10010-10015
- Kumar V., Muratani M., Rayan N. A., Kraus P., Lufkin T., et al. (2013) Uniform, optimal signal processing of mapped deep-sequencing data. *Nat Biotechnol* 31: 615-622
- Kuo T., Lew M. J., Mayba O., Harris C. A., Speed T. P., et al. (2012) Genome-wide analysis of glucocorticoid receptor-binding sites in myotubes identifies gene networks modulating insulin signaling. *Proc Natl Acad Sci U S A* 109: 11160-11165
- Lai C. F., Flach K. D., Alexi X., Fox S. P., Ottaviani S., et al. (2013) Co-regulated gene expression by oestrogen receptor alpha and liver receptor homolog-1 is a feature of the oestrogen response in breast cancer cells. *Nucleic Acids Res* 41: 10228-10240
- Lai D., Proctor J. R. and Meyer I. M. (2013) On the importance of cotranscriptional RNA structure formation. *RNA* 19: 1461-1473
- Lam M. T., Cho H., Lesch H. P., Gosselin D., Heinz S., et al. (2013) Rev-Erbs repress macrophage gene expression by inhibiting enhancer-directed transcription. *Nature* 498: 511-515
- Lan X., Witt H., Katsumura K., Ye Z., Wang Q., et al. (2012) Integration of Hi-C and ChIP-seq data reveals distinct types of chromatin linkages. *Nucleic Acids Res* 40: 7690-7704
- Larschan E., Bishop E. P., Kharchenko P. V., Core L. J., Lis J. T., et al. (2011) X chromosome dosage compensation via enhanced transcriptional elongation in *Drosophila*. *Nature* 471: 115-118
- Larson M. H., Gilbert L. A., Wang X., Lim W. A., Weissman J. S., et al. (2013) CRISPR interference (CRISPRi) for sequence-specific control of gene expression. *Nat Protoc* 8: 2180-2196
- Lasken R. S. (2013) Single-cell sequencing in its prime. *Nat Biotechnol* 31: 211-212
- Lassen K. S., Schultz H., Heegaard N. H. and He M. (2012) A novel DNaseq program for enhanced analysis of Illumina GAT data: a case study on antibody complementarity-determining regions. *N Biotechnol* 29: 271-278

- Law G. L., Korth M. J., Benecke A. G. and Katze M. G. (2013) Systems virology: host-directed approaches to viral pathogenesis and drug targeting. *Nat Rev Microbiol* 11: 455-466
- Lazarovici A., Zhou T., Shafer A., Dantas Machado A. C., Riley T. R., et al. (2013) Probing DNA shape and methylation state on a genomic scale with DNase I. *Proc Natl Acad Sci U S A* 110: 6376-6381
- Lebedeva S., Jens M., Theil K., Schwanhaussner B., Selbach M., et al. (2011) Transcriptome-wide analysis of regulatory interactions of the RNA-binding protein HuR. *Mol Cell* 43: 340-352
- Lee J., Seok S., Yu P., Kim K., Smith Z., et al. (2012) Genomic analysis of hepatic farnesoid X receptor binding sites reveals altered binding in obesity and direct gene repression by farnesoid X receptor in mice. *Hepatology* 56: 108-117
- Lee M. T., Bonneau A. R., Takacs C. M., Bazzini A. A., DiVito K. R., et al. (2013) Nanog, Pou5f1 and SoxB1 activate zygotic gene expression during the maternal-to-zygotic transition. *Nature* 503: 360-364
- Lee S., Liu B., Lee S., Huang S. X., Shen B., et al. (2012) Global mapping of translation initiation sites in mammalian cells at single-nucleotide resolution. *Proc Natl Acad Sci U S A* 109: E2424-2432
- Lee T. I. and Young R. A. (2013) Transcriptional regulation and its misregulation in disease. *Cell* 152: 1237-1251
- Lesch B. J., Dokshin G. A., Young R. A., McCarrey J. R. and Page D. C. (2013) A set of genes critical to development is epigenetically poised in mouse germ cells from fetal stages through completion of meiosis. *Proc Natl Acad Sci U S A* 110: 16061-16066
- Li G., Fullwood M. J., Xu H., Mulawadi F. H., Velkov S., et al. (2010) ChIA-PET tool for comprehensive chromatin interaction analysis with paired-end tag sequencing. *Genome Biol* 11: R22
- Li G., Ruan X., Auerbach R. K., Sandhu K. S., Zheng M., et al. (2012) Extensive promoter-centered chromatin interactions provide a topological basis for transcription regulation. *Cell* 148: 84-98
- Li G. W., Oh E. and Weissman J. S. (2012) The anti-Shine-Dalgarno sequence drives translational pausing and codon choice in bacteria. *Nature* 484: 538-541
- Li N., Ye M., Li Y., Yan Z., Butcher L. M., et al. (2010) Whole genome DNA methylation analysis based on high throughput sequencing technology. *Methods* 52: 203-212
- Li P., Spann N. J., Kaikkonen M. U., Lu M., Oh da Y., et al. (2013) NCoR repression of LXRs restricts macrophage biosynthesis of insulin-sensitizing omega 3 fatty acids. *Cell* 155: 200-214
- Li W., Notani D., Ma Q., Tanasa B., Nunez E., et al. (2013) Functional roles of enhancer RNAs for oestrogen-dependent transcriptional activation. *Nature* 498: 516-520
- Li Z., Chao T. C., Chang K. Y., Lin N., Patil V. S., et al. (2014) The long noncoding RNA THRIL regulates TNFalpha expression through its interaction with hnRNPL. *Proc Natl Acad Sci U S A* 111: 1002-1007
- R. B. (2010) RNA processing and its regulation: global insights into biological networks. *Nat Rev Genet* 11: 75-87
- Lieberman-Aiden E., van Berkum N. L., Williams L., Imakaev M., Ragoczy T., et al. (2009) Comprehensive mapping of long-range interactions reveals folding principles of the human genome. *Science* 326: 289-293
- Lister R., Mukamel E. A., Nery J. R., Urich M., Puddifoot C. A., et al. (2013) Global epigenomic reconfiguration during mammalian brain development. *Science* 341: 1237905
- Liu B., Han Y. and Qian S. B. (2013) Cotranslational response to proteotoxic stress by elongation pausing of ribosomes. *Mol Cell* 49: 453-463
- Liu M., Li C. L., Stamatoyannopoulos G., Dorschner M. O., Humbert R., et al. (2012) Gammaretroviral Vector Integration Occurs Overwhelmingly Within and Near DNase Hypersensitive Sites. *Hum Gene Ther* 23: 231-237
- Liu X., Jiang H., Gu Z. and Roberts J. W. (2013) High-resolution view of bacteriophage lambda gene expression by ribosome profiling. *Proc Natl Acad Sci U S A* 110: 11928-11933
- Liu Y., Hu W., Murakawa Y., Yin J., Wang G., et al. (2013) Cold-induced RNA-binding proteins regulate circadian gene expression by controlling alternative polyadenylation. *Sci Rep* 3: 2054
- Lo K. A., Labadorf A., Kennedy N. J., Han M. S., Yap Y. S., et al. (2013) Analysis of in vitro insulin-resistance models and their physiological relevance to in vivo diet-induced adipose insulin resistance. *Cell Rep* 5: 259-270
- Lovett M. (2013) The applications of single-cell genomics. *Hum Mol Genet* 22: R22-26
- Lucks J. B., Mortimer S. A., Trapnell C., Luo S., Aviran S., et al. (2011) Multiplexed RNA structure characterization with selective 2'-hydroxyl acylation analyzed by primer extension sequencing (SHAPE-Seq). *Proc Natl Acad Sci U S A* 108: 11063-11068
- Lun F. M., Chiu R. W., Sun K., Leung T. Y., Jiang P., et al. (2013) Noninvasive prenatal methylomic analysis by genomewide bisulfite sequencing of maternal plasma DNA. *Clin Chem* 59: 1583-1594



- Lyko F., Foret S., Kucharski R., Wolf S., Falckenhayn C., et al. (2010) The honey bee epigenomes: differential methylation of brain DNA in queens and workers. *PLoS Biol* 8: e1000506
- Macosko E. Z. and McCarroll S. A. (2013) Genetics. Our fallen genomes. *Science* 342: 564-565
- Maehara K., Odawara J., Harada A., Yoshimi T., Nagao K., et al. (2013) A co-localization model of paired ChIP-seq data using a large ENCODE data set enables comparison of multiple samples. *Nucleic Acids Res* 41: 54-62
- Majoros W. H., Lekprasert P., Mukherjee N., Skalsky R. L., Corcoran D. L., et al. (2013) MicroRNA target site identification by integrating sequence and binding information. *Nat Methods* 10: 630-633
- Malhotra A., Lindberg M., Faust G. G., Leibowitz M. L., Clark R. A., et al. (2013) Breakpoint profiling of 64 cancer genomes reveals numerous complex rearrangements spawned by homology-independent mechanisms. *Genome Res* 23: 762-776
- Mandal P. K., Ewing A. D., Hancks D. C. and Kazazian H. H., Jr. (2013) Enrichment of processed pseudogene transcripts in L1-ribonucleoprotein particles. *Hum Mol Genet* 22: 3730-3748
- Mann B., van Opijnen T., Wang J., Obert C., Wang Y. D., et al. (2012) Control of virulence by small RNAs in *Streptococcus pneumoniae*. *PLoS Pathog* 8: e1002788
- Marinov G. K., Williams B. A., McCue K., Schroth G. P., Gertz J., et al. (2014) From single-cell to cell-pool transcriptomes: Stochasticity in gene expression and RNA splicing. *Genome Res* 24: 496-510
- Maruyama H., Harwood J. C., Moore K. M., Paszkiewicz K., Durley S. C., et al. (2013) An alternative beads-on-a-string chromatin architecture in *Thermococcus kodakarensis*. *EMBO Rep* 14: 711-717
- Marwick J. A., Kirkham P. A., Stevenson C. S., Danahay H., Giddings J., et al. (2004) Cigarette smoke alters chromatin remodeling and induces proinflammatory genes in rat lungs. *Am J Respir Cell Mol Biol* 31: 633-642
- Matsushita H., Vesely M. D., Koboldt D. C., Rickert C. G., Uppaluri R., et al. (2012) Cancer exome analysis reveals a T-cell-dependent mechanism of cancer immunoediting. *Nature* 482: 400-404
- Maunakea A. K., Nagarajan R. P., Bilenky M., Ballinger T. J., D'Souza C., et al. (2010) Conserved role of intragenic DNA methylation in regulating alternative promoters. *Nature* 466: 253-257
- Maurano M. T., Humbert R., Rynes E., Thurman R. E., Haugen E., et al. (2012) Systematic localization of common disease-associated variation in regulatory DNA. *Science* 337: 1190-1195
- Mazzoni E. O., Mahony S., Closser M., Morrison C. A., Nedelec S., et al. (2013) Synergistic binding of transcription factors to cell-specific enhancers programs motor neuron identity. *Nat Neurosci* 16: 1219-1227
- McConnell M. J., Lindberg M. R., Brennand K. J., Piper J. C., Voet T., et al. (2013) Mosaic copy number variation in human neurons. *Science* 342: 632-637
- McDaniell R., Lee B. K., Song L., Liu Z., Boyle A. P., et al. (2010) Heritable individual-specific and allele-specific chromatin signatures in humans. *Science* 328: 235-239
- McGettigan P. A. (2013) Transcriptomics in the RNA-seq era. *Curr Opin Chem Biol* 17: 4-11
- McKenna L. B., Schug J., Vourekas A., McKenna J. B., Bramswig N. C., et al. (2010) MicroRNAs control intestinal epithelial differentiation, architecture, and barrier function. *Gastroenterology* 139: 1654-1664, 1664 e1651
- McLean J. S., Lombardo M. J., Ziegler M. G., Novotny M., Yee-Greenbaum J., et al. (2013) Genome of the pathogen *Porphyromonas gingivalis* recovered from a biofilm in a hospital sink using a high-throughput single-cell genomics platform. *Genome Res* 23: 867-877
- McMullen P. D., Bhattacharya S., Woods C. G., Sun B., Yarborough K., et al. (2014) A map of the PPARalpha transcription regulatory network for primary human hepatocytes. *Chem Biol Interact* 209: 14-24
- Meaburn E. and Schulz R. (2012) Next generation sequencing in epigenetics: insights and challenges. *Semin Cell Dev Biol* 23: 192-199
- Medvedovic J., Ebert A., Tagoh H., Tamir I. M., Schwickert T. A., et al. (2013) Flexible long-range loops in the VH gene region of the Igh locus facilitate the generation of a diverse antibody repertoire. *Immunity* 39: 229-244
- Meissner A., Gnirke A., Bell G. W., Ramsahoye B., Lander E. S., et al. (2005) Reduced representation bisulfite sequencing for comparative high-resolution DNA methylation analysis. *Nucleic Acids Res* 33: 5868-5877
- Meissner A., Mikkelsen T. S., Gu H., Wernig M., Hanna J., et al. (2008) Genome-scale DNA methylation maps of pluripotent and differentiated cells. *Nature* 454: 766-770
- Mellen M., Ayata P., Dewell S., Kriaucionis S. and Heintz N. (2012) MeCP2 binds to 5hmC enriched within active genes and accessible chromatin in the nervous system. *Cell* 151: 1417-1430

- Mendoza-Parra M. A., Van Gool W., Mohamed Saleem M. A., Ceschin D. G. and Gronemeyer H. (2013) A quality control system for profiles obtained by ChIP sequencing. *Nucleic Acids Res* 41: e196
- Meredith D. M., Borromeo M. D., Deering T. G., Casey B. H., Savage T. K., et al. (2013) Program specificity for Ptf1a in pancreas versus neural tube development correlates with distinct collaborating cofactors and chromatin accessibility. *Mol Cell Biol* 33: 3166-3179
- Meyer K. D., Saletore Y., Zumbo P., Elemento O., Mason C. E., et al. (2012) Comprehensive analysis of mRNA methylation reveals enrichment in 3' UTRs and near stop codons. *Cell* 149: 1635-1646
- Mitchell J. A., Clay I., Umlauf D., Chen C. Y., Moir C. A., et al. (2012) Nuclear RNA sequencing of the mouse erythroid cell transcriptome. *PLoS One* 7: e49274
- Miura F., Enomoto Y., Dairiki R. and Ito T. (2012) Amplification-free whole-genome bisulfite sequencing by post-bisulfite adaptor tagging. *Nucleic Acids Res* 40: e136
- Mukherjee N., Corcoran D. L., Nusbaum J. D., Reid D. W., Georgiev S., et al. (2011) Integrative regulatory mapping indicates that the RNA-binding protein HuR couples pre-mRNA processing and mRNA stability. *Mol Cell* 43: 327-339
- Murray S. M., Panis G., Fumeaux C., Viollier P. H. and Howard M. (2013) Computational and genetic reduction of a cell cycle to its simplest, primordial components. *PLoS Biol* 11: e1001749
- Murtaza M., Dawson S. J., Tsui D. W., Gale D., Forshew T., et al. (2013) Non-invasive analysis of acquired resistance to cancer therapy by sequencing of plasma DNA. *Nature* 497: 108-112
- Mylykangas S., Buenrostro J. D., Natsoulis G., Bell J. M. and Ji H. P. (2011) Efficient targeted resequencing of human germline and cancer genomes by oligonucleotide-selective sequencing. *Nat Biotechnol* 29: 1024-1027
- Nagano T., Lubling Y., Stevens T. J., Schoenfelder S., Yaffe E., et al. (2013) Single-cell Hi-C reveals cell-to-cell variability in chromosome structure. *Nature* 502: 59-64
- Nagarajavel V., Iben J. R., Howard B. H., Maraia R. J. and Clark D. J. (2013) Global 'bootprinting' reveals the elastic architecture of the yeast TFIIB-TFIIC transcription complex in vivo. *Nucleic Acids Res* 41: 8135-8143
- Navin N. and Hicks J. (2011) Future medical applications of single-cell sequencing in cancer. *Genome Med* 3: 31
- Navin N., Kendall J., Troge J., Andrews P., Rodgers L., et al. (2011) Tumour evolution inferred by single-cell sequencing. *Nature* 472: 90-94
- Neph S., Vierstra J., Stergachis A. B., Reynolds A. P., Haugen E., et al. (2012) An expansive human regulatory lexicon encoded in transcription factor footprints. *Nature* 489: 83-90
- Neri F., Krepelova A., Incarnato D., Maldotti M., Parlato C., et al. (2013) Dnmt3L Antagonizes DNA Methylation at Bivalent Promoters and Favors DNA Methylation at Gene Bodies in ESCs. *Cell* 155: 121-134
- Neyret-Kahn H., Benhamed M., Ye T., Le Gras S., Cossec J. C., et al. (2013) Sumoylation at chromatin governs coordinated repression of a transcriptional program essential for cell growth and proliferation. *Genome Res* 23: 1563-1579
- Ni T., Corcoran D. L., Rach E. A., Song S., Spana E. P., et al. (2010) A paired-end sequencing strategy to map the complex landscape of transcription initiation. *Nat Methods* 7: 521-527
- Ni X., Zhuo M., Su Z., Duan J., Gao Y., et al. (2013) Reproducible copy number variation patterns among single circulating tumor cells of lung cancer patients. *Proc Natl Acad Sci U S A* 110: 21083-21088
- Nishida H., Katayama T., Suzuki Y., Kondo S. and Horiuchi H. (2013) Base composition and nucleosome density in exonic and intronic regions in genes of the filamentous ascomycetes *Aspergillus nidulans* and *Aspergillus oryzae*. *Gene* 525: 10-May
- Nora E. P., Lajoie B. R., Schulz E. G., Giorgetti L., Okamoto I., et al. (2012) Spatial partitioning of the regulatory landscape of the X-inactivation centre. *Nature* 485: 381-385
- O'Huallachain M., Weissman S. M. and Snyder M. P. (2013) The variable somatic genome. *Cell Cycle* 12: 5-6
- Oh E., Becker A. H., Sandikci A., Huber D., Chaba R., et al. (2011) Selective ribosome profiling reveals the cotranslational chaperone action of trigger factor in vivo. *Cell* 147: 1295-1308
- Oliveira T., Resch W., Jankovic M., Casellas R., Nussenzweig M. C., et al. (2011) Translocation capture sequencing: A method for high throughput mapping of chromosomal rearrangements. *J Immunol Methods* 375: 176-181
- Oliveira T. Y., Resch W., Jankovic M., Casellas R., Nussenzweig M. C., et al. (2012) Translocation capture sequencing: a method for high throughput mapping of chromosomal rearrangements. *J Immunol Methods* 375: 176-181
- Olovnikov I., Ryazansky S., Shpiz S., Lavrov S., Abramov Y., et al. (2013) De novo piRNA cluster formation in the *Drosophila* germ line triggered by transgenes containing a transcribed transposon fragment. *Nucleic Acids Res* 41: 5757-5768

- Osborne R. J. and Thornton C. A. (2006) RNA-dominant diseases. *Hum Mol Genet* 15 Spec No 2: R162-169
- Paul D. S., Nisbet J. P., Yang T. P., Meacham S., Rendon A., et al. (2011) Maps of open chromatin guide the functional follow-up of genome-wide association signals: application to hematological traits. *PLoS Genet* 7: e1002139
- Paul D. S., Albers C. A., Rendon A., Voss K., Stephens J., et al. (2013) Maps of open chromatin highlight cell type-restricted patterns of regulatory sequence variation at hematological trait loci. *Genome Res* 23: 1130-1141
- Pelechano V., Wei W. and Steinmetz L. M. (2013) Extensive transcriptional heterogeneity revealed by isoform profiling. *Nature* 497: 127-131
- Pfeifer G. P., Kadam S. and Jin S. G. (2013) 5-hydroxymethylcytosine and its potential roles in development and cancer. *Epigenetics Chromatin* 6: 10
- Piccolo F. M. and Fisher A. G. (2014) Getting rid of DNA methylation. *Trends Cell Biol* 24: 136-143
- Picelli S., Bjorklund A. K., Faridani O. R., Sagasser S., Winberg G., et al. (2013) Smart-seq2 for sensitive full-length transcriptome profiling in single cells. *Nat Methods* 10: 1096-1098
- Pinello L., Xu J., Orkin S. H. and Yuan G. C. (2014) Analysis of chromatin-state plasticity identifies cell-type-specific regulators of H3K27me3 patterns. *Proc Natl Acad Sci U S A* 111: E344-353
- Polymenidou M., Lagier-Tourenne C., Hutt K. R., Huelga S. C., Moran J., et al. (2011) Long pre-mRNA depletion and RNA missplicing contribute to neuronal vulnerability from loss of TDP-43. *Nat Neurosci* 14: 459-468
- Pomraning K. R., Smith K. M. and Freitag M. (2009) Genome-wide high throughput analysis of DNA methylation in eukaryotes. *Methods* 47: 142-150
- Ponts N., Harris E. Y., Prudhomme J., Wick I., Eckhardt-Ludka C., et al. (2010) Nucleosome landscape and control of transcription in the human malaria parasite. *Genome Res* 20: 228-238
- Potok M. E., Nix D. A., Parnell T. J. and Cairns B. R. (2013) Reprogramming the maternal zebrafish genome after fertilization to match the paternal methylation pattern. *Cell* 153: 759-772
- Potter N. E., Ermini L., Papaemmanuil E., Cazzaniga G., Vijayaraghavan G., et al. (2013) Single-cell mutational profiling and clonal phylogeny in cancer. *Genome Res* 23: 2115-2125
- Poulos M. G., Batra R., Li M., Yuan Y., Zhang C., et al. (2013) Progressive impairment of muscle regeneration in muscleblind-like 3 isoform knockout mice. *Hum Mol Genet* 22: 3547-3558
- Prickett A. R., Barkas N., McCole R. B., Hughes S., Amante S. M., et al. (2013) Genome-wide and parental allele-specific analysis of CTCF and cohesin DNA binding in mouse brain reveals a tissue-specific binding pattern and an association with imprinted differentially methylated regions. *Genome Res* 23: 1624-1635
- Puszyk W., Down T., Grimwade D., Chomienne C., Oakey R. J., et al. (2013) The epigenetic regulator PLZF represses L1 retrotransposition in germ and progenitor cells. *EMBO J* 32: 1941-1952
- Qi L. S., Larson M. H., Gilbert L. A., Doudna J. A., Weissman J. S., et al. (2013) Repurposing CRISPR as an RNA-guided platform for sequence-specific control of gene expression. *Cell* 152: 1173-1183
- Qi X., Daily K., Nguyen K., Wang H., Mayhew D., et al. (2012) Retrotransposon profiling of RNA polymerase III initiation sites. *Genome Res* 22: 681-692
- Ramasubramanian S., Kanhere A., Osborn K., Flower K., Jenner R. G., et al. (2012) Genome-wide analyses of Zta binding to the Epstein-Barr virus genome reveals interactions in both early and late lytic cycles and an epigenetic switch leading to an altered binding profile. *J Virol* 86: 12494-12502
- Ramskold D., Luo S., Wang Y. C., Li R., Deng Q., et al. (2012) Full-length mRNA-Seq from single-cell levels of RNA and individual circulating tumor cells. *Nat Biotechnol* 30: 777-782
- Rauch T. A., Zhong X., Wu X., Wang M., Kernstine K. H., et al. (2008) High-resolution mapping of DNA hypermethylation and hypomethylation in lung cancer. *Proc Natl Acad Sci U S A* 105: 252-257
- Rauch T. A. and Pfeifer G. P. (2009) The MIRA method for DNA methylation analysis. *Methods Mol Biol* 507: 65-75
- Regulski M., Lu Z., Kendall J., Donoghue M. T., Reinders J., et al. (2013) The maize methylome influences mRNA splice sites and reveals widespread paramutation-like switches guided by small RNA. *Genome Res* 23: 1651-1662
- Rich A. and RajBhandary U. L. (1976) Transfer RNA: molecular structure, sequence, and properties. *Annu Rev Biochem* 45: 805-860
- Riley K. J., Rabinowitz G. S., Yario T. A., Luna J. M., Darnell R. B., et al. (2012) EBV and human microRNAs co-target oncogenic and apoptotic viral and human genes during latency. *EMBO J* 31: 2207-2221
- Rivera C. M. and Ren B. (2013) Mapping human epigenomes. *Cell* 155: 39-55



- Robyr D., Friedli M., Gehrig C., Arcangeli M., Marin M., et al. (2011) Chromosome conformation capture uncovers potential genome-wide interactions between human conserved non-coding sequences. *PLoS ONE* 6: e17634
- Rocha P. P., Micsinai M., Kim J. R., Hewitt S. L., Souza P. P., et al. (2012) Close proximity to Igh is a contributing factor to AID-mediated translocations. *Mol Cell* 47: 873-885
- Rocha P. P., Chaumeil J. and Skok J. A. (2013) Molecular biology. Finding the right partner in a 3D genome. *Science* 342: 1333-1334
- Rodley C. D., Pai D. A., Mills T. A., Engelke D. R. and O'Sullivan J. M. (2011) tRNA gene identity affects nuclear positioning. *PLoS One* 6: e29267
- Rodrigues J. A., Ruan R., Nishimura T., Sharma M. K., Sharma R., et al. (2013) Imprinted expression of genes and small RNA is associated with localized hypomethylation of the maternal genome in rice endosperm. *Proc Natl Acad Sci U S A* 110: 7934-7939
- Rodriguez B. A., Frankhouser D., Murphy M., Trimarchi M., Tam H. H., et al. (2012) Methods for high-throughput MethylCap-Seq data analysis. *BMC Genomics* 13 Suppl 6: S14
- Rogelj B., Easton L. E., Bogu G. K., Stanton L. W., Rot G., et al. (2012) Widespread binding of FUS along nascent RNA regulates alternative splicing in the brain. *Sci Rep* 2: 603
- Ronan J. L., Wu W. and Crabtree G. R. (2013) From neural development to cognition: unexpected roles for chromatin. *Nat Rev Genet* 14: 347-359
- Rosseeel T., Scheuch M., Hoper D., De Regge N., Caij A. B., et al. (2012) DNase SISPA-next generation sequencing confirms Schmallenberg virus in Belgian field samples and identifies genetic variation in Europe. *PLoS One* 7: e41967
- Ruike Y., Imanaka Y., Sato F., Shimizu K. and Tsujimoto G. (2010) Genome-wide analysis of aberrant methylation in human breast cancer cells using methyl-DNA immunoprecipitation combined with high-throughput sequencing. *BMC Genomics* 11: 137
- Saied M. H., Marzec J., Khalid S., Smith P., Down T. A., et al. (2012) Genome wide analysis of acute myeloid leukemia reveal leukemia specific methylome and subtype specific hypomethylation of repeats. *PLoS One* 7: e33213
- Sajan S. A. and Hawkins R. D. (2012) Methods for identifying higher-order chromatin structure. *Annu Rev Genomics Hum Genet* 13: 59-82
- Sakabe N. J., Aneas I., Shen T., Shokri L., Park S. Y., et al. (2012) Dual transcriptional activator and repressor roles of TBX20 regulate adult cardiac structure and function. *Hum Mol Genet* 21: 2194-2204
- Sakurai M., Yano T., Kawabata H., Ueda H. and Suzuki T. (2010) Inosine cyanoethylation identifies A-to-I RNA editing sites in the human transcriptome. *Nat Chem Biol* 6: 733-740
- Salton M., Elkon R., Borodina T., Davydov A., Yaspo M. L., et al. (2011) Matrin 3 binds and stabilizes mRNA. *PLoS One* 6: e23882
- Sanyal A., Lajoie B. R., Jain G. and Dekker J. (2012) The long-range interaction landscape of gene promoters. *Nature* 489: 109-113
- Sarmiento F., Mrázek J. and Whitman W. B. (2013) Genome-scale analysis of gene function in the hydrogenotrophic methanogenic archaeon *Methanococcus maripaludis*. *Proc Natl Acad Sci U S A* 110: 4726-4731
- Sasagawa Y., Nikaido I., Hayashi T., Danno H., Uno K. D., et al. (2013) Quartz-Seq: a highly reproducible and sensitive single-cell RNA sequencing method, reveals non-genetic gene-expression heterogeneity. *Genome Biol* 14: R31
- Sati S., Tanwar V. S., Kumar K. A., Patowary A., Jain V., et al. (2012) High resolution methylome map of rat indicates role of intragenic DNA methylation in identification of coding region. *PLoS One* 7: e31621
- Saunders A., Core L. J., Sutcliffe C., Lis J. T. and Ashe H. L. (2013) Extensive polymerase pausing during *Drosophila* axis patterning enables high-level and pliable transcription. *Genes Dev* 27: 1146-1158
- Schauer T., Schwalie P. C., Handley A., Margulies C. E., Flicek P., et al. (2013) CAST-ChIP maps cell-type-specific chromatin states in the *Drosophila* central nervous system. *Cell Rep* 5: 271-282
- Schillebeeckx M., Schrade A., Lobs A. K., Pihlajoki M., Wilson D. B., et al. (2013) Laser capture microdissection-reduced representation bisulfite sequencing (LCM-RRBS) maps changes in DNA methylation associated with gonadectomy-induced adrenocortical neoplasia in the mouse. *Nucleic Acids Res* 41: e116
- Schlesinger F., Smith A. D., Gingeras T. R., Hannon G. J. and Hodges E. (2013) De novo DNA demethylation and noncoding transcription define active intergenic regulatory elements. *Genome Res* 23: 1601-1614
- Schmitt M. W., Kennedy S. R., Salk J. J., Fox E. J., Hiatt J. B., et al. (2012) Detection of ultra-rare mutations by next-generation sequencing. *Proc Natl Acad Sci U S A* 109: 14508-14513
- Schmitz R. J., He Y., Valdes-Lopez O., Khan S. M., Joshi T., et al. (2013) Epigenome-wide inheritance of cytosine methylation variants in a recombinant inbred population. *Genome Res* 23: 1663-1674
- Schmolka N., Serre K., Grosso A. R., Rei M., Pennington D. J., et al. (2013) Epigenetic and transcriptional signatures of stable versus plastic differentiation of proinflammatory gamma delta T cell subsets. *Nat Immunol* 14: 1093-1100

- Schones D. E., Cui K., Cuddapah S., Roh T. Y., Barski A., et al. (2008) Dynamic regulation of nucleosome positioning in the human genome. *Cell* 132: 887-898
- Schweiger M. R., Barmeyer C. and Timmermann B. (2013) Genomics and epigenomics: new promises of personalized medicine for cancer patients. *Brief Funct Genomics* 12: 411-421
- Seisenberger S., Andrews S., Krueger F., Arand J., Walter J., et al. (2012) The dynamics of genome-wide DNA methylation reprogramming in mouse primordial germ cells. *Mol Cell* 48: 849-862
- Serre D., Lee B. H. and Ting A. H. (2010) MBD-isolated Genome Sequencing provides a high-throughput and comprehensive survey of DNA methylation in the human genome. *Nucleic Acids Res* 38: 391-399
- Seth-Smith H. M., Harris S. R., Skilton R. J., Radebe F. M., Golparian D., et al. (2013) Whole-genome sequences of *Chlamydia trachomatis* directly from clinical samples without culture. *Genome Res* 23: 855-866
- Sexton T., Kurukuti S., Mitchell J. A., Umlauf D., Nagano T., et al. (2012) Sensitive detection of chromatin coassociations using enhanced chromosome conformation capture on chip. *Nat Protoc* 7: 1335-1350
- Shalek A. K., Satija R., Adiconis X., Gertner R. S., Gaublomme J. T., et al. (2013) Single-cell transcriptomics reveals bimodality in expression and splicing in immune cells. *Nature* 498: 236-240
- Shang W. H., Hori T., Martins N. M., Toyoda A., Misu S., et al. (2013) Chromosome engineering allows the efficient isolation of vertebrate neocentromeres. *Dev Cell* 24: 635-648
- Shapiro E., Biezuner T. and Linnarsson S. (2013) Single-cell sequencing-based technologies will revolutionize whole-organism science. *Nat Rev Genet* 14: 618-630
- Shen H., Qiu C., Li J., Tian Q. and Deng H. W. (2013) Characterization of the DNA methylome and its interindividual variation in human peripheral blood monocytes. *Epigenomics* 5: 255-269
- Shirane K., Toh H., Kobayashi H., Miura F., Chiba H., et al. (2013) Mouse oocyte methylomes at base resolution reveal genome-wide accumulation of non-CpG methylation and role of DNA methyltransferases. *PLoS Genet* 9: e1003439
- Shiroguchi K., Jia T. Z., Sims P. A. and Xie X. S. (2012) Digital RNA sequencing minimizes sequence-dependent bias and amplification noise with optimized single-molecule barcodes. *Proc Natl Acad Sci U S A* 109: 1347-1352
- Shukla R., Upton K. R., Munoz-Lopez M., Gerhardt D. J., Fisher M. E., et al. (2013) Endogenous retrotransposition activates oncogenic pathways in hepatocellular carcinoma. *Cell* 153: 101-111
- Sievers C., Schlumpf T., Sawarkar R., Comoglio F. and Paro R. (2012) Mixture models and wavelet transforms reveal high confidence RNA-protein interaction sites in MOV10 PAR-CLIP data. *Nucleic Acids Res* 40: e160
- Simon J. M., Giresi P. G., Davis I. J. and Lieb J. D. (2012) Using formaldehyde-assisted isolation of regulatory elements (FAIRE) to isolate active regulatory DNA. *Nat Protoc* 7: 256-267
- Skalsky R. L., Corcoran D. L., Gottwein E., Frank C. L., Kang D., et al. (2012) The viral and cellular microRNA targetome in lymphoblastoid cell lines. *PLoS Pathog* 8: e1002484
- Skurnik D., Roux D., Cattoir V., Danilchanka O., Lu X., et al. (2013) Enhanced in vivo fitness of carbapenem-resistant oprD mutants of *Pseudomonas aeruginosa* revealed through high-throughput sequencing. *Proc Natl Acad Sci U S A* 110: 20747-20752
- Smallwood S. A., Tomizawa S., Krueger F., Ruf N., Carli N., et al. (2011) Dynamic CpG island methylation landscape in oocytes and preimplantation embryos. *Nat Genet* 43: 811-814
- Smith Z. D., Gu H., Bock C., Gnirke A. and Meissner A. (2009) High-throughput bisulfite sequencing in mammalian genomes. *Methods* 48: 226-232
- Smith Z. D. and Meissner A. (2013) DNA methylation: roles in mammalian development. *Nat Rev Genet* 14: 204-220
- Soleimani V. D., Palidwor G. A., Ramachandran P., Perkins T. J. and Rudnicki M. A. (2013) Chromatin tandem affinity purification sequencing. *Nat Protoc* 8: 1525-1534
- Solomon M. J., Larsen P. L. and Varshavsky A. (1988) Mapping protein-DNA interactions in vivo with formaldehyde: evidence that histone H4 is retained on a highly transcribed gene. *Cell* 53: 937-947
- Solyom S., Ewing A. D., Rahrmann E. P., Doucet T., Nelson H. H., et al. (2012) Extensive somatic L1 retrotransposition in colorectal tumors. *Genome Res* 22: 2328-2338
- Song C. X., Szulwach K. E., Dai Q., Fu Y., Mao S. Q., et al. (2013) Genome-wide profiling of 5-formylcytosine reveals its roles in epigenetic priming. *Cell* 153: 678-691

- Song L., Zhang Z., Grassefer L. L., Boyle A. P., Giresi P. G., et al. (2011) Open chromatin defined by DNaseI and FAIRE identifies regulatory elements that shape cell-type identity. *Genome Res* 21: 1757-1767
- Spaethling J. M. and Eberwine J. H. (2013) Single-cell transcriptomics for drug target discovery. *Curr Opin Pharmacol* 13: 786-790
- Splinter E., de Wit E., van de Werken H. J., Klous P. and de Laat W. (2012) Determining long-range chromatin interactions for selected genomic sites using 4C-seq technology: from fixation to computation. *Methods* 58: 221-230
- Stadhouders R., Kolovos P., Brouwer R., Zuin J., van den Heuvel A., et al. (2013) Multiplexed chromosome conformation capture sequencing for rapid genome-scale high-resolution detection of long-range chromatin interactions. *Nat Protoc* 8: 509-524
- Stadler M., Artilles K., Pak J. and Fire A. (2012) Contributions of mRNA abundance, ribosome loading, and post- or peri-translational effects to temporal repression of *C. elegans* heterochronic miRNA targets. *Genome Res* 22: 2418-2426
- Stadler M. B., Murr R., Burger L., Ivanek R., Lienert F., et al. (2011) DNA-binding factors shape the mouse methylome at distal regulatory regions. *Nature* 480: 490-495
- Stevens M., Cheng J. B., Li D., Xie M., Hong C., et al. (2013) Estimating absolute methylation levels at single-CpG resolution from methylation enrichment and restriction enzyme sequencing methods. *Genome Res* 23: 1541-1553
- Stevens M., Cheng J. B., Li D., Xie M., Hong C., et al. (2013) Estimating absolute methylation levels at single-CpG resolution from methylation enrichment and restriction enzyme sequencing methods. *Genome Res* 23: 1541-1553
- Stoll G., Pietilainen O. P., Linder B., Suvisaari J., Brosi C., et al. (2013) Deletion of TOP3beta, a component of FMRP-containing mRNPs, contributes to neurodevelopmental disorders. *Nat Neurosci* 16: 1228-1237
- Stumpf C. R., Moreno M. V., Olshen A. B., Taylor B. S. and Ruggero D. (2013) The translational landscape of the Mammalian cell cycle. *Mol Cell* 52: 574-582
- Taiwo O., Wilson G. A., Morris T., Seisenberger S., Reik W., et al. (2012) Methylome analysis using MeDIP-seq with low DNA concentrations. *Nat Protoc* 7: 617-636
- Tan L., Xiong L., Xu W., Wu F., Huang N., et al. (2013) Genome-wide comparison of DNA hydroxymethylation in mouse embryonic stem cells and neural progenitor cells by a new comparative hMeDIP-seq method. *Nucleic Acids Res* 41: e84
- Tan S. K., Lin Z. H., Chang C. W., Varang V., Chng K. R., et al. (2011) AP-2gamma regulates oestrogen receptor-mediated long-range chromatin interaction and gene transcription. *EMBO J* 30: 2569-2581
- Tanaka Y., Joshi A., Wilson N. K., Kinston S., Nishikawa S., et al. (2012) The transcriptional programme controlled by Runx1 during early embryonic blood development. *Dev Biol* 366: 404-419
- Tang B., Becanovic K., Desplats P. A., Spencer B., Hill A. M., et al. (2012) Forkhead box protein p1 is a transcriptional repressor of immune signaling in the CNS: implications for transcriptional dysregulation in Huntington disease. *Hum Mol Genet* 21: 3097-3111
- Tang D. T., Plessy C., Salimullah M., Suzuki A. M., Calligaris R., et al. (2013) Suppression of artifacts and barcode bias in high-throughput transcriptome analyses utilizing template switching. *Nucleic Acids Res* 41: e44
- Telese F., Gamiel A., Skowronska-Krawczyk D., Garcia-Bassets I. and Rosenfeld M. G. (2013) "Seq-ing" insights into the epigenetics of neuronal gene regulation. *Neuron* 77: 606-623
- Thapar R., Denmon A. P. and Nikonowicz E. P. (2014) Recognition modes of RNA tetraloops and tetraloop-like motifs by RNA-binding proteins. *Wiley Interdiscip Rev RNA* 5: 49-67
- Thomson J. P., Lempiainen H., Hackett J. A., Nestor C. E., Muller A., et al. (2012) Non-genotoxic carcinogen exposure induces defined changes in the 5-hydroxymethylome. *Genome Biol* 13: R93
- Thomson J. P., Hunter J. M., Nestor C. E., Dunican D. S., Terranova R., et al. (2013) Comparative analysis of affinity-based 5-hydroxymethylation enrichment techniques. *Nucleic Acids Res* 41: e206-
- Tollervy J. R., Curk T., Rogelj B., Briese M., Cereda M., et al. (2011) Characterizing the RNA targets and position-dependent splicing regulation by TDP-43. *Nat Neurosci* 14: 452-458
- Tolstorukov M. Y., Sansam C. G., Lu P., Koellhoffer E. C., Helming K. C., et al. (2013) Swi/Snf chromatin remodeling/tumor suppressor complex establishes nucleosome occupancy at target promoters. *Proc Natl Acad Sci U S A* 110: 10165-10170
- Troy E. B., Lin T., Gao L., Lazinski D. W., Camilli A., et al. (2013) Understanding barriers to *Borrelia burgdorferi* dissemination during infection using massively parallel sequencing. *Infect Immun* 81: 2347-2357
- Turnbaugh P. J., Ridaura V. K., Faith J. J., Rey F. E., Knight R., et al. (2009) The effect of diet on the human gut microbiome: a metagenomic analysis in humanized gnotobiotic mice. *Sci Transl Med* 1: 6ra14

- Turnbaugh P. J., Quince C., Faith J. J., McHardy A. C., Yatsunencko T., et al. (2010) Organismal, genetic, and transcriptional variation in the deeply sequenced gut microbiomes of identical twins. *Proc Natl Acad Sci U S A* 107: 7503-7508
- Turner J. R., Ray R., Lee B., Everett L., Xiang J., et al. (2013) Evidence from mouse and man for a role of neuregulin 3 in nicotine dependence. *Mol Psychiatry*
- Twiggs S. R., Vorgia E., McGowan S. J., Peraki I., Fenwick A. L., et al. (2013) Reduced dosage of ERF causes complex craniosynostosis in humans and mice and links ERK1/2 signaling to regulation of osteogenesis. *Nat Genet* 45: 308-313
- Udan-Johns M., Bengoechea R., Bell S., Shao J., Diamond M. I., et al. (2014) Prion-like nuclear aggregation of TDP-43 during heat shock is regulated by HSP40/70 chaperones. *Hum Mol Genet* 23: 157-170
- Underwood J. G., Uzilov A. V., Katzman S., Onodera C. S., Mainzer J. E., et al. (2010) FragSeq: transcriptome-wide RNA structure probing using high-throughput sequencing. *Nat Methods* 7: 995-1001
- Uniacke J., Holterman C. E., Lachance G., Franovic A., Jacob M. D., et al. (2012) An oxygen-regulated switch in the protein synthesis machinery. *Nature* 486: 126-129
- Valentim C. L., LoVerde P. T., Anderson T. J. and Criscione C. D. (2009) Efficient genotyping of *Schistosoma mansoni* miracidia following whole genome amplification. *Mol Biochem Parasitol* 166: 81-84
- van de Werken H. J., Landan G., Holwerda S. J., Hoichman M., Klous P., et al. (2012) Robust 4C-seq data analysis to screen for regulatory DNA interactions. *Nat Methods* 9: 969-972
- Van der Aa N., Esteki M. Z., Vermeesch J. R. and Voet T. (2013) Preimplantation genetic diagnosis guided by single-cell genomics. *Genome Med* 5: 71
- van Opijnen T., Bodi K. L. and Camilli A. (2009) Tn-seq: high-throughput parallel sequencing for fitness and genetic interaction studies in microorganisms. *Nat Methods* 6: 767-772
- van Opijnen T. and Camilli A. (2012) A fine scale phenotype-genotype virulence map of a bacterial pathogen. *Genome Res* 22: 2541-2551
- van Opijnen T. and Camilli A. (2013) Transposon insertion sequencing: a new tool for systems-level analysis of microorganisms. *Nat Rev Microbiol* 11: 435-442
- Veluchamy A., Lin X., Maumus F., Rivarola M., Bhavsar J., et al. (2013) Insights into the role of DNA methylation in diatoms by genome-wide profiling in *Phaeodactylum tricornutum*. *Nat Commun* 4: 2091
- Verma-Gaur J., Torkamani A., Schaffer L., Head S. R., Schork N. J., et al. (2012) Noncoding transcription within the Igh distal V(H) region at PAIR elements affects the 3D structure of the Igh locus in pro-B cells. *Proc Natl Acad Sci U S A* 109: 17004-17009
- Vidaki A., Daniel B. and Court D. S. (2013) Forensic DNA methylation profiling--Potential opportunities and challenges. *Forensic Sci Int Genet* 7: 499-507
- Voet T., Kumar P., Van Loo P., Cooke S. L., Marshall J., et al. (2013) Single-cell paired-end genome sequencing reveals structural variation per cell cycle. *Nucleic Acids Res* 41: 6119-6138
- Vourekas A., Zheng Q., Alexiou P., Maragkakis M., Kirino Y., et al. (2012) Mili and Miwi target RNA repertoire reveals piRNA biogenesis and function of Miwi in spermiogenesis. *Nat Struct Mol Biol* 19: 773-781
- Wan Y., Qu K., Ouyang Z. and Chang H. Y. (2013) Genome-wide mapping of RNA structure using nuclease digestion and high-throughput sequencing. *Nat Protoc* 8: 849-869
- Wang D., Garcia-Bassets I., Benner C., Li W., Su X., et al. (2011) Reprogramming transcription by distinct classes of enhancers functionally defined by eRNA. *Nature* 474: 390-394
- Wang Q., Gu L., Adey A., Radlwimmer B., Wang W., et al. (2013) Tagmentation-based whole-genome bisulfite sequencing. *Nat Protoc* 8: 2022-2032
- Wang T., Cui Y., Jin J., Guo J., Wang G., et al. (2013) Translating mRNAs strongly correlate to proteins in a multivariate manner and their translation ratios are phenotype specific. *Nucleic Acids Res* 41: 4743-4754
- Wang T., Wu H., Li Y., Szulwach K. E., Lin L., et al. (2013) Subtelomeric hotspots of aberrant 5-hydroxymethylcytosine-mediated epigenetic modifications during reprogramming to pluripotency. *Nat Cell Biol* 15: 700-711
- Wang X., Lu Z., Gomez A., Hon G. C., Yue Y., et al. (2014) N6-methyladenosine-dependent regulation of messenger RNA stability. *Nature* 505: 117-120
- Wang Y. M., Zhou P., Wang L. Y., Li Z. H., Zhang Y. N., et al. (2012) Correlation between DNase I hypersensitive site distribution and gene expression in HeLa S3 cells. *PLoS One* 7: e42414
- Warden C. D., Lee H., Tompkins J. D., Li X., Wang C., et al. (2013) COHCAP: an integrative genomic pipeline for single-nucleotide resolution DNA methylation analysis. *Nucleic Acids Res* 41: e117
- Warnault V., Darcq E., Levine A., Barak S. and Ron D. (2013) Chromatin remodeling--a novel strategy to control excessive alcohol drinking. *Transl Psychiatry* 3: e231

- Waszak S. M., Kilpinen H., Gschwind A. R., Orioli A., Raghav S. K., et al. (2014) Identification and removal of low-complexity sites in allele-specific analysis of ChIP-seq data. *Bioinformatics* 30: 165-171
- Weber M., Davies J. J., Wittig D., Oakeley E. J., Haase M., et al. (2005) Chromosome-wide and promoter-specific analyses identify sites of differential DNA methylation in normal and transformed human cells. *Nat Genet* 37: 853-862
- Wei Y., Chen R., Dimicoli S., Bueso-Ramos C., Neuberger D., et al. (2013) Global H3K4me3 genome mapping reveals alterations of innate immunity signaling and overexpression of JMJD3 in human myelodysplastic syndrome CD34+ cells. *Leukemia* 27: 2177-2186
- Wei Z., Gao F., Kim S., Yang H., Lyu J., et al. (2013) Klf4 organizes long-range chromosomal interactions with the oct4 locus in reprogramming and pluripotency. *Cell Stem Cell* 13: 36-47
- Whisnant A. W., Bogerd H. P., Flores O., Ho P., Powers J. G., et al. (2013) In-depth analysis of the interaction of HIV-1 with cellular microRNA biogenesis and effector mechanisms. *MBio* 4: e000193
- Will B., Vogler T. O., Bartholdy B., Garrett-Bakelman F., Mayer J., et al. (2013) Satb1 regulates the self-renewal of hematopoietic stem cells by promoting quiescence and repressing differentiation commitment. *Nat Immunol* 14: 437-445
- Wolchinsky Z., Shvritiel S., Kouwenhoven E. N., Putin D., Sprecher E., et al. (2014) Angiomodulin is required for cardiogenesis of embryonic stem cells and is maintained by a feedback loop network of p63 and Activin-A. *Stem Cell Res* 12: 49-59
- Woyke T., Tighe D., Mavromatis K., Clum A., Copeland A., et al. (2010) One bacterial cell, one complete genome. *PLoS ONE* 5: e10314
- Woyke T., Sczyrba A., Lee J., Rinke C., Tighe D., et al. (2011) Decontamination of MDA reagents for single cell whole genome amplification. *PLoS ONE* 6: e26161
- Xi Y., Yao J., Chen R., Li W. and He X. (2011) Nucleosome fragility reveals novel functional states of chromatin and poises genes for activation. *Genome Res* 21: 718-724
- Xi Y., Bock C., Muller F., Sun D., Meissner A., et al. (2012) RRBSMAP: a fast, accurate and user-friendly alignment tool for reduced representation bisulfite sequencing. *Bioinformatics* 28: 430-432
- Xie W., Schultz M. D., Lister R., Hou Z., Rajagopal N., et al. (2013) Epigenomic analysis of multilineage differentiation of human embryonic stem cells. *Cell* 153: 1134-1148
- Xing J., Witherspoon D. J. and Jorde L. B. (2013) Mobile element biology: new possibilities with high-throughput sequencing. *Trends Genet* 29: 280-289
- Xu D., Shen W., Guo R., Xue Y., Peng W., et al. (2013) Top3beta is an RNA topoisomerase that works with fragile X syndrome protein to promote synapse formation. *Nat Neurosci* 16: 1238-1247
- Xue Z., Huang K., Cai C., Cai L., Jiang C. Y., et al. (2013) Genetic programs in human and mouse early embryos revealed by single-cell RNA sequencing. *Nature* 500: 593-597
- Yadav S. P., Hao H., Yang H. J., Kautzmann M. A., Brooks M., et al. (2013) The transcription-splicing protein NonO/p54nrb and three NonO-interacting proteins bind to distal enhancer region and augment rhodopsin expression. *Hum Mol Genet*
- Yadon A. N. and Tsukiyama T. (2013) DNA looping-dependent targeting of a chromatin remodeling factor. *Cell Cycle* 12: 1809-1810
- Yaffe E. and Tanay A. (2011) Probabilistic modeling of Hi-C contact maps eliminates systematic biases to characterize global chromosomal architecture. *Nat Genet* 43: 1059-1065
- Yamaguchi S., Hong K., Liu R., Shen L., Inoue A., et al. (2012) Tet1 controls meiosis by regulating meiotic gene expression. *Nature* 492: 443-447
- Yamaguchi S., Hong K., Liu R., Inoue A., Shen L., et al. (2013) Dynamics of 5-methylcytosine and 5-hydroxymethylcytosine during germ cell reprogramming. *Cell Res* 23: 329-339
- Yamaguchi S., Shen L., Liu Y., Sandler D. and Zhang Y. (2013) Role of Tet1 in erasure of genomic imprinting. *Nature* 504: 460-464
- Yamane A., Robbiani D. F., Resch W., Bothmer A., Nakahashi H., et al. (2013) RPA Accumulation during Class Switch Recombination Represents 5'-3' DNA-End Resection during the S-G2/M Phase of the Cell Cycle. *Cell Rep* 3: 138-147
- Yan L., Yang M., Guo H., Yang L., Wu J., et al. (2013) Single-cell RNA-Seq profiling of human preimplantation embryos and embryonic stem cells. *Nat Struct Mol Biol* 20: 1131-1139
- Yano M., Hayakawa-Yano Y., Mele A. and Darnell R. B. (2010) Nova2 regulates neuronal migration through an RNA switch in disabled-1 signaling. *Neuron* 66: 848-858
- Yeh H. H., Young D., Gelovani J. G., Robinson A., Davidson Y., et al. (2013) Histone deacetylase class II and acetylated core histone immunohistochemistry in human brains with Huntington's disease. *Brain Res* 1504: 16-24
- Yilmaz S. and Singh A. K. (2011) Single cell genome sequencing. *Curr Opin Biotechnol* 23: 437-443



- Yokobayashi S., Liang C. Y., Kohler H., Nestorov P., Liu Z., et al. (2013) PRC1 coordinates timing of sexual differentiation of female primordial germ cells. *Nature* 495: 236-240
- Yu M., Hon G. C., Szulwach K. E., Song C. X., Jin P., et al. (2012) Tet-assisted bisulfite sequencing of 5-hydroxymethylcytosine. *Nat Protoc* 7: 2159-2170
- Yu M., Hon G. C., Szulwach K. E., Song C. X., Zhang L., et al. (2012) Base-resolution analysis of 5-hydroxymethylcytosine in the mammalian genome. *Cell* 149: 1368-1380
- Yu M., Hon G. C., Szulwach K. E., Song C. X., Zhang L., et al. (2012) Base-resolution analysis of 5-hydroxymethylcytosine in the mammalian genome. *Cell* 149: 1368-1380
- Yu W., Jin C., Lou X., Han X., Li L., et al. (2011) Global analysis of DNA methylation by Methyl-Capture sequencing reveals epigenetic control of cisplatin resistance in ovarian cancer cell. *PLoS One* 6: e29450
- Zarnack K., Konig J., Tajnik M., Martincorena I., Eustermann S., et al. (2013) Direct competition between hnRNP C and U2AF65 protects the transcriptome from the exonization of Alu elements. *Cell* 152: 453-466
- Zentner G. E. and Henikoff S. (2012) Surveying the epigenomic landscape, one base at a time. *Genome Biol* 13: 250
- Zhan Q., Fang Y., He Y., Liu H. X., Fang J., et al. (2012) Function annotation of hepatic retinoid x receptor alpha based on genome-wide DNA binding and transcriptome profiling. *PLoS One* 7: e50013
- Zhang B., Zhou Y., Lin N., Lowdon R. F., Hong C., et al. (2013) Functional DNA methylation differences between tissues, cell types, and across individuals discovered using the M&M algorithm. *Genome Res* 23: 1522-1540
- Zhang C., Frias M. A., Mele A., Ruggiu M., Eom T., et al. (2010) Integrative modeling defines the Nova splicing-regulatory network and its combinatorial controls. *Science* 329: 439-443
- Zhang C. and Darnell R. B. (2011) Mapping in vivo protein-RNA interactions at single-nucleotide resolution from HITS-CLIP data. *Nat Biotechnol* 29: 607-614
- Zhang J., Poh H. M., Peh S. Q., Sia Y. Y., Li G., et al. (2012) ChIA-PET analysis of transcriptional chromatin interactions. *Methods* 58: 289-299
- Zhang W., Wu Y., Schnable J. C., Zeng Z., Freeling M., et al. (2012) High-resolution mapping of open chromatin in the rice genome. *Genome Res* 22: 151-162
- Zhang W., Zhang T., Wu Y. and Jiang J. (2012) Genome-wide identification of regulatory DNA elements and protein-binding footprints using signatures of open chromatin in Arabidopsis. *Plant Cell* 24: 2719-2731
- Zhang X., Robertson G., Woo S., Hoffman B. G. and Gottardo R. (2012) Probabilistic Inference for Nucleosome Positioning with MNase-Based or Sonicated Short-Read Data. *PLoS ONE* 7: e32095
- Zhang Y., McCord R. P., Ho Y. J., Lajoie B. R., Hildebrand D. G., et al. (2012) Spatial organization of the mouse genome and its role in recurrent chromosomal translocations. *Cell* 148: 908-921
- Zhao J., Ohsumi T. K., Kung J. T., Ogawa Y., Grau D. J., et al. (2010) Genome-wide identification of polycomb-associated RNAs by RIP-seq. *Mol Cell* 40: 939-953
- Zhao Y., Guo S., Sun J., Huang Z., Zhu T., et al. (2012) Methylcap-seq reveals novel DNA methylation markers for the diagnosis and recurrence prediction of bladder cancer in a Chinese population. *PLoS ONE* 7: e35175
- Zhao Z., Tavoosidana G., Sjolinder M., Gondor A., Mariano P., et al. (2006) Circular chromosome conformation capture (4C) uncovers extensive networks of epigenetically regulated intra- and interchromosomal interactions. *Nat Genet* 38: 1341-1347
- Zilbauer M., Rayner T. F., Clark C., Coffey A. J., Joyce C. J., et al. (2013) Genome-wide methylation analyses of primary human leukocyte subsets identifies functionally important cell-type-specific hypomethylated regions. *Blood* 122: e52-60
- Ziller M. J., Muller F., Liao J., Zhang Y., Gu H., et al. (2011) Genomic distribution and inter-sample variation of non-CpG methylation across human cell types. *PLoS Genet* 7: e1002389
- Zomer A., Burghout P., Bootsma H. J., Hermans P. W. and van Hijum S. A. (2012) ESSENTIALS: Software for Rapid Analysis of High Throughput Transposon Insertion Sequencing Data. *PLoS ONE* 7: e43012
- Zong C., Lu S., Chapman A. R. and Xie X. S. (2012) Genome-wide detection of single-nucleotide and copy-number variations of a single human cell. *Science* 338: 1622-1626
- Zund D., Gruber A. R., Zavolan M. and Muhlemann O. (2013) Translation-dependent displacement of UPF1 from coding sequences causes its enrichment in 3' UTRs. *Nat Struct Mol Biol* 20: 936-943

DNA/RNA PURIFICATION KITS

MasterPure™ Complete DNA and RNA Purification Kit



# Pure Nucleic Acids

## MasterPure™ Complete DNA and RNA Purification kit

- ▶ Extract and purify total nucleic acids (TNA), DNA or RNA
- ▶ Pure for sequencing, qPCR and other molecular biology applications
  
- ▶ Scalable reactions
- ▶ High purity and yield
- ▶ Non-Toxic



The MasterPure™ Complete Kit purifies high yields of intact total nucleic acid, DNA, or RNA. MasterPure is suitable for every type of biological material.

MasterPure is optimized for use with:

- ▶ illumina® sequencing
- ▶ qPCR
- ▶ PCR
- ▶ Molecular biology applications

Table 1. Purify any sample.

Sample	Sample Size	TNA µg	DNA µg	RNA µg
HeLa/HL60 cells	1 X 10 <sup>6</sup> cells	10-30	3-12	7-15
Liver	5 mg	33-42	5-10	13-25
Brain	5 mg	9-13	6-9	4-11
Heart	5 mg	6-10	4-7	4-5
Blood	200 µl	3-10	3-9	
Buffy coat	300 µl	40-55	40-55	3-6
<i>E. coli</i>	3.5 x 10 <sup>6</sup> cells	2.5-2.8	1.3-1.6	1.6-1.8
Yeast*	2.2 x 10 <sup>6</sup> cells			11-18
( <i>S. cerevisiae</i> )	1.1 x 10 <sup>7</sup> cells			70-78

Many different, diverse sample types have been purified by MasterPure. Several are shown in Table 1. MasterPure is available for virtually any type of sample.

MasterPure may be used to purify total nucleic acid, DNA or RNA from any sample. Total nucleic acid purification permits you to compare DNA and RNA from the same sample to gain a deeper understanding of your sample.

Workflow



Figure 2. MasterPure easily purifies many different sample types.



MasterPure has been published showing excellent purification of nearly every type of sample. Here we see just a few of the types of samples that have been published using MasterPure.

#### Purify any sample

MasterPure has been shown to work for many types of human tissue and blood samples, plants, and bacteria.

MasterPure is safe and nontoxic. No dangerous chemicals, phenol or hazards are used in the method. MasterPure is a wise choice for safety and high yields of RNA, DNA, or total nucleic acids.

Stop stocking three different kits for small, moderate and abundant samples! MasterPure is designed to be used with small, moderate and abundant samples without the need for many kits. One MasterPure kit permits you to purify RNA, DNA or total nucleic acid from any amount of sample.

One kit to purify your choice of nucleic acids.

#### Suitable for Illumina® sequencing

Total nucleic acid, DNA or RNA purified by MasterPure is suitable for use with Illumina sequencing. All sequencing applications begin with MasterPure, including:

- ▶ Ribo-Zero
- ▶ RNA-Seq
- ▶ Bisulfite sequencing for epigenetics
- ▶ DNA-Seq
- ▶ Exome capture
- ▶ More...

Cat. # Quantity

#### MasterPure™ Complete DNA and RNA Purification kit

MC85200 200 DNA Purifications  
100 RNA Purifications

MC89010 10 DNA Purifications  
5 RNA Purifications

#### MasterPure™ DNA Purification Kit

MCD85201 200 Purifications



DNA-Sequencing	
Description	Catalog Number
MasterPure™ Complete DNA and RNA Purification Kit	MC85200
MasterPure™ DNA Purification Kit	MCD85201
TruSeq DNA PCR-Free LT Sample Preparation Kit - Set A	FC-121-3001
TruSeq DNA PCR-Free LT Sample Preparation Kit - Set B	FC-121-3002
TruSeq DNA PCR-Free HT Sample Preparation Kit	FC-121-3003
TruSeq Nano DNA LT Sample Preparation Kit - Set A	FC-121-4001
TruSeq Nano DNA LT Sample Preparation Kit - Set B	FC-121-4002
TruSeq Nano DNA HT Sample Preparation Kit	FC-121-4003
Nextera Rapid Capture Exome (8 rxn x 1 Plex)	FC-140-1000
Nextera Rapid Capture Exome (8 rxn x 3 Plex)	FC-140-1083
Nextera Rapid Capture Exome (8 rxn x 6 Plex)	FC-140-1086
Nextera Rapid Capture Exome (8 rxn x 9 Plex)	FC-140-1089
Nextera Rapid Capture Exome (2 rxn x 12 Plex)	FC-140-1001
Nextera Rapid Capture Exome (4 rxn x 12 Plex)	FC-140-1002
Nextera Rapid Capture Exome (8 rxn x 12 Plex)	FC-140-1003
Nextera Rapid Capture Expanded Exome (2 rxn x 12 Plex)	FC-140-1004
Nextera Rapid Capture Expanded Exome (4 rxn x 12 Plex)	FC-140-1005
Nextera Rapid Capture Expanded Exome (8 rxn x 12 Plex)	FC-140-1006
EpiGenome™ Methyl-Seq Kit	EGMK81312


ChIP	
Description	Catalog Number
TruSeq ChIP Sample Preparation Kit - Set A	IP-202-1012
TruSeq ChIP Sample Preparation Kit - Set B	IP-202-1024

Methylation Arrays	
Description	Catalog Number
HumanMethylation450 DNA Analysis BeadChip Kit (24 samples)	WG-314-1003
HumanMethylation450 DNA Analysis BeadChip Kit (48 samples)	WG-314-1001
HumanMethylation450 DNA Analysis BeadChip Kit (96 samples)	WG-314-1002





**Figure 4: Adapter Ligation Results in Sequence-Ready Constructs without PCR**



to the A-tailed fragmented DNA. These newly redesigned adapters contain the full complement of sequencing primer hybridization sites for single, paired-end, and multiplexed reads. This eliminates the need for additional PCR steps to add the index tag and multiplex primer sites (Figure 4D). Following the denaturation and amplification steps (Figure 4E), libraries can be pooled with up to 12 samples per lane (96 sample per flow cell) for cluster generation on either cBot or the Cluster Station.

Master-mixed reagents and an optimized protocol improve the library construction workflow, significantly decreasing hands-on time and reducing the number of clean-up steps when processing samples for large-scale studies (Table 1). The simple and scalable workflow allows for high-throughput and automation-friendly solutions, as well as simultaneous manual processing for up to 96 samples. In addition, enhanced troubleshooting features are incorporated into each step of the workflow, with quality control sequences supported by Illumina RTA software.

**Enhanced Quality Controls**

Specific Quality Control (QC) sequences, consisting of double-stranded DNA fragments, are present in each enzymatic reaction of the TruSeq sample preparation protocol: end repair, A-tailing, and ligation. During analysis, the QC sequences are recognized by the RTA software (versions 1.8 and later) and isolated from the sample data. The presence of these controls indicates that its corresponding step was successful. If a step was unsuccessful, the control sequences will be substantially reduced. QC controls assist in comparison between experiments and greatly facilitate troubleshooting.

**Designed For Automation**

The TruSeq Sample Preparation Kits are compatible with high-throughput, automated processing workflows. Sample preparation can be performed in standard 96-well microplates with master-mixed reagent pipetting volumes optimized for liquid-handling robots. Barcodes on reagents and plates allow end-to-end sample tracking and ensure that the correct reagents are used for the correct protocol, mitigating potential tracking errors.

**Part of an Integrated Sequencing Solution**

Samples processed with the TruSeq Sample Preparation Kits can be amplified on either the cBot Automated Cluster Generation System or the Cluster Station and used with any of Illumina's next-generation sequencing instruments, including HiSeq™ 2000, HiSeq 1000, HiScan™SQ, Genome Analyzer<sub>IIx</sub> (Figure 5).


**Summary**

Illumina's new TruSeq Sample Preparation Kits enable simplicity, convenience, and affordability for library preparation. Enhanced multiplexing with 24 unique indexes allows efficient high-throughput processing. The pre-configured reagents, streamlined workflow, and automation-friendly protocol save researchers time and effort in their next-generation sequencing pursuits, ultimately leading to faster discovery and publication.

Learn more about Illumina's next-generation sequencing solutions at [www.illumina.com/sequencing](http://www.illumina.com/sequencing).

AAAGAATGATAACAGTAAACACACTTCTGTTAACTTAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC  
 AATCAAGCTACCGTAAACGAACGATCAATTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC  
 AACGACGAAAGGATGATAACAGTAAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCT  
 TTAACCTAACATTAAAGAGCTACCGTAAACGAACGATCAACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC  
 AAAGATTAAAGAGCTACCGTAAACGAACGATCAACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC  
 AAGATTACTTGATCCACTGATTCACAGCTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC  
 AACCTATCAATTGAGACTAAATATTAAGTACCATTAAAGATTACTTGATCCACTGATTCACAGCTACCGTAAACGAACGATCAATTGAGACTAAATATTAAGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTC

Figure 5. Illumina's Complete Sequencing Solution



#### Ordering Information

Product	Catalog No.
<b>For RNA Preparation</b>	
TruSeq RNA Sample Preparation Kit v2, Set A (12 indexes, 48 samples)	RS-122-2001
TruSeq RNA Sample Preparation Kit v2, Set B (12 indexes, 48 samples)	RS-122-2002
<b>For DNA Preparation</b>	
TruSeq DNA Sample Preparation Kit v2, Set A (12 indexes, 48 samples)	FC-121-2001
TruSeq DNA Sample Preparation Kit v2, Set B (12 indexes, 48 samples)	FC-121-2002
<b>For Cluster Generation on cBot and Sequencing on the HiSeq 2000/1000 and HiScanSQ</b>	
TruSeq Paired-End Cluster Kit v3—cBot—HS (1 flow cell)	PE-401-3001
TruSeq Single-Read Cluster Kit v3—cBot—HS (1 flow cell)	GD-401-3001
<b>For Cluster Generation on cBot and Sequencing on the Genome Analyzer<sub>IIx</sub></b>	
TruSeq Paired-End Cluster Kit v2—cBot—GA (1 flow cell)	PE-300-2001
TruSeq Single-Read Cluster Kit v2—cBot—GA (1 flow cell)	GD-300-2001
<b>For Cluster Generation on the Cluster Station and Sequencing on the Genome Analyzer<sub>IIx</sub></b>	
TruSeq Paired-End Cluster Kit v5—CS—GA (1 flow cell)	PE-203-5001
TruSeq Single-Read Cluster Kit v5—CS—GA (1 flow cell)	GD-203-5001

Illumina, Inc. • 9885 Towne Centre Drive, San Diego, CA 92121 USA • 1.800.809.4566 toll-free • 1.858.202.4566 tel • techsupport@illumina.com • illumina.com

#### FOR RESEARCH USE ONLY

© 2011 Illumina, Inc. All rights reserved.  
Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CSPro, DASL, Eco, Genetic Energy, GAllx, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
Pub. No. 970-2009-039 Current as of 27 April 2011




AGAAATGATAACAGTAACACACTTTCTGTAAACCTTAAAGATTACCTTTCACACTGATTCACAGCTACCGTAAGAAACGTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCTTCTGTTAAACCTTAAAGATTACCTTGATCCACTGATTCAAC  
TGACGTACCGTACCGAAGCTATCATTAAGATTACTTTGATCCACTGATTGACCGTAAAGCGTAAACCGCTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCTTCTGTTAAACCTTAAAGATTACCTTGATCCACTGATTCAAC  
CGACGAAAAGAAATGATAACAGTAACACACTTTCTGTAAACCTTAAAGATTACCTTTCACACTGATTCACAGCTACCGTAAGAAACGTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCTTCTGTTAAACCTTAAAGATTACCT  
TACCTACCATTAAGAGACTAACCGTCCACAGCTAAACACTTTCTGTAAACCTTAAAGATTACCTTGATCCACTGATTCACAGCTACCGTAAAGCGTAAACCGCTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCCACAGCTAACAAAGATGATA  
AGATGATAAACAGTAACACACTTTCTGTAAACCTTAAAGATTACTTTGATCCACTGATTCACAGCTACCGTAAAGCGTAAACCGCTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCTTCTGTTAAACCTTAAAGATTACCTT  
GATTACTTTGATCCACTGATTCACAGCTTAAAGATTACTTTGATCCACTGATTCACAGCTACCGTAAAGCGTATGAATTGAGACTTTCTGTAAACCTTAAAGATTACTTTGATCCACTGATTCACAGCTACCGTAAAGCGTATGAATTGAGACTAAATAT  
CGTATGAATTGAGACTAAATATTAACGTAACCCATTAAAGGCTAACCGTCTTCTGTTAAACCTTAAAGATTACCTTGATCCACTGATTCAAC






Figure 5: Increased Coverage of Challenging Regions



When compared to libraries generated by PCR-based workflows, such as TruSeq DNA Sample Preparation, PCR-Free libraries show improved coverage for challenging regions of the genome. These regions include known human protein coding and non-protein coding exons and genes defined in the RefSeq Genes track in the UCSC Genome Browser.<sup>2</sup> G-Rich regions denote 30 bases with  $\geq 80\%$  G. High GC regions are defined as 100 bases with  $\geq 75\%$  GC content. Huge GC regions are defined as 100 bases with  $\geq 85\%$  GC content. "Difficult" promoters denote the set of 100 promoter regions that are insufficiently covered, which have been empirically defined by the Broad Institute of MIT and Harvard.<sup>3</sup> AT dinucleotides indicate 30 bases of repeated AT dinucleotide.

Figure 6: PCR-Free Protocol Eliminates Coverage Gaps in GC-Rich Content



Increased coverage of TruSeq DNA PCR-Free libraries results in fewer coverage gaps, demonstrated here in the GC-rich coding regions of the *RNPEPL1* promoter (A) and the *CREBBP* promoter (B). PCR-Free sequence information is shown in the top panels of A and B, while sequence data generated using TruSeq DNA protocol (with PCR) are shown in the lower panels.

The TruSeq LT kit includes up to 24 indices with two sets of 12each, and the TruSeq HT kit offers 96 indices for efficient experimental design.

Multi-sample studies can be conveniently managed using the Illumina Experiment Manager, a freely available software tool that provides easy reaction setup for plate-based processing. It allows researchers to quickly configure the index sample sheet (i.e., sample multiplexing matrix) for the instrument run, enabling automatic demultiplexing.

### Flexible and Inclusive Sample Preparation

The TruSeq family of sample preparation solutions offers several kits for sequencing applications, compatible with a range of research needs and study designs (Table 1). All TruSeq kits support high- and low-throughput studies. The TruSeq DNA PCR-Free kit provides superior coverage quality and drastically reduces library bias and coverage gaps, without requiring PCR amplification. These kits enhance the industry's most widely adopted DNA sample preparation method, empowering next-generation sequencing applications.

### Simplified Solution

The comprehensive solution includes sample preparation reagents, sample purification beads, and robust TruSeq barcodes for sample multiplexing, providing a complete preparation method optimized for the highest performance on all Illumina sequencing platforms. The TruSeq DNA PCR-Free kit leverages the flexibility of two kit options, 24-sample and 96-sample, for a scalable experimental approach. With a simplified workflow and multiplexing options, the TruSeq DNA PCR-Free protocol offers the fastest library preparation method for the highest data quality.

AAAGAATGATAACAGTAAACACACTTCTGTTAACTTAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
AATCAACGTACCGTAACGAAACGATATCAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
AACGACGAAAGGATGATAACAGTAAACACACTTCTGTTAACTTAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
TTAAGCTAACATTAAAGAGCTACCGTACAGCACTAACACACTTCTGTTAACTTAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
AAAGATTGATAACAGTAAACACACTTCTGTTAACTTAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
AAGATTACTTGATCCACTGATTCACAGGTAAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA  
AACCTATCAATTTGAGACTAAATTAAGGTACCATTAAAGATTACTTGATCCACTGATTCACAGGTACCGTAACGAAACGATATCAATTTGAGACTAAATTAAGGTACCATTAAAGAGCTACCGTCTTCTGTTAACTTAAGATTACTTGATCCACTGATTCOA



Data Sheet: Sequencing

**Table 1: TruSeq DNA Sample Preparation Kits**

Specification	TruSeq Nano DNA	TruSeq DNA PCR-Free	TruSeq DNA
Description	Based upon widely adopted TruSeq sample prep, with lower input and improved data quality	Superior genomic coverage with radically reduced library bias and gaps	Original TruSeq next-generation sequencing sample preparation method
Input quantity	100–200 ng	1–2 µg	1 µg
Includes PCR	Yes	No	Yes
Assay time	~6 hours	~5 hours	1–2 days
Hands-on time	~5 hours	~4 hours	~8 hours
Target insert size	350 bp or 550 bp	350 bp or 550 bp	300 bp
Gel-Free	Yes	Yes	No
Number of samples supported	24 (LT) or 96 (HT) samples	24 (LT) or 96 (HT) samples	48 (LT) or 96 (HT) samples
Supports enrichment	No*	No*	Yes
Size-selection beads	Included	Included	Not included
Applications	Whole-genome sequencing applications, including whole-genome resequencing, <i>de novo</i> assembly, and metagenomics studies		
Sample multiplexing	24 single indices or 96 dual-index combinations		
Compatible Illumina sequencers	HiSeq®, HiScanSQ™, Genome Analyzer™, and MiSeq® systems		

\*Nextera Rapid Capture products support a variety of enrichment applications. For more information, visit [www.illumina.com/NRC](http://www.illumina.com/NRC).

**Summary**

The TruSeq DNA PCR-Free Sample Preparation Kit optimizes the TruSeq workflow to deliver a faster sample preparation method for any species. The choice between protocol options provides greater flexibility to support a variety of applications and genomic studies. The PCR-Free kit also removes PCR-induced bias to facilitate detailed and accurate insight into the genome. By leveraging a faster workflow and superior data quality, the TruSeq DNA PCR-Free Sample Preparation Kit enables researchers to obtain high-quality genomic data, faster.

**References**

1. Saunders CJ, Miller NA, Soden SE, Dinwiddie DL, Noll A, et al. (2012) Rapid whole-genome sequencing for genetic disease diagnosis in neonatal intensive care units. *Science Translational Medicine* 4(154): 154ra135.
2. [genome.ucsc.edu](http://genome.ucsc.edu)
3. [www.broadinstitute.org](http://www.broadinstitute.org)

**Ordering Information**

Product	Catalog No.
TruSeq DNA PCR-Free LT Sample Preparation Kit Set A (24 samples)	FC-121-3001
TruSeq DNA PCR-Free LT Sample Preparation Kit Set B (24 samples)	FC-121-3002
TruSeq DNA PCR-Free HT Sample Preparation Kit (96 samples)	FC-121-3003

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • [techsupport@illumina.com](mailto:techsupport@illumina.com) • [www.illumina.com](http://www.illumina.com)

**FOR RESEARCH USE ONLY**

© 2013 Illumina, Inc. All rights reserved.  
 Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GalX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
 Pub. No. 770-2013-001 Current as of 16 May 2013




AGAATGATACAGTAAACACACTTCCTGTTAAACCTTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 TCAACGTACGTAACGAAACGTATCATTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 CGACGAAAGAAAGATGATAACAGCTAACACACTTCCTGTTAAACCTTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 TACGTACCATTAAAGGCTACCGTCTTCTGTTAAACCTTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 AAGATTGATAACAGCTAACACACTTCCTGTTAAACCTTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 GATTACTTGTATCCACTGATTC AACGTTAAAGTTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC  
 GGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTCTGTAAAGCTTAAGATTACTTGTATCCACTGATTC AACGTAACCGTAACGAAACGTATCAATTGAGACCTAAATTAACGTACCATTAAAGGCTACCGTCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTC AAC




## Data Sheet: Sequencing

**Figure 2: TruSeq Nano DNA Workflow**




The TruSeq Nano DNA LT indexing solution features a single-index adapter at Step E.

**Figure 3: Fewer Gaps in Coverage**



TruSeq Nano DNA libraries show significant reduction in the number and total size of gaps when compared to libraries prepared using the TruSeq DNA protocol. A gap is defined as a region  $\geq 10$  bp in length, where an accurate genotype cannot be determined due to low depth, low alignment scores, or low base quality.

**Figure 4: Greater Coverage Uniformity**



TruSeq Nano DNA libraries provide greater coverage uniformity across the genome when compared to those generated using the TruSeq DNA protocol.

```

AATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTACCGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
ACGTAACGTAACGAAAGCTATGATTAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
ACGAAAGAAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
GTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
ATGATAACGTAACGAAAGCTATGATTAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
TTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
TATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAAACGAAACGATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGGTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCACGGTAA
    
```



### Efficient Sample Multiplexing

Using a simple procedure, indices are added to sample genomic DNA fragments to provide an innovative solution for sample multiplexing. For the greatest operational efficiency, up to 96 pre-plated, uniquely indexed samples can be pooled and sequenced together in a single flow cell lane on any Illumina sequencing platform. After sequencing, the indices are used to demultiplex the data and accurately assign reads to the proper samples in the pool. The TruSeq LT kit uses a single index for demultiplexing, while the TruSeq HT kit employs a dual-indexing strategy, using a unique combination of two indices to demultiplex. The LT kit includes up to 24 indices with two sets of 12 each, and the HT kit offers 96 indices.

### Streamlined Solution

This inclusive kit contains sample preparation reagents, sample purification beads, and robust TruSeq indices for multiplexing, providing a complete preparation method optimized for the highest performance on all Illumina sequencing platforms. The TruSeq Nano DNA kit leverages the flexibility of two kit options, 24-sample and 96-sample, for scalable experimental design. With a simplified workflow and flexible multiplexing options, the TruSeq Nano DNA protocol offers a streamlined library preparation method that delivers high-quality sequencing data.


### Summary

The TruSeq Nano DNA Sample Preparation Kit optimizes the TruSeq workflow to deliver a low-input sample preparation method for any sequencing application. Low- and high-throughput options and varied insert sizes provide greater flexibility to support a variety of applications and genomic studies. Workflow innovations reduce PCR-induced bias to facilitate detailed and accurate insight into the genome. By leveraging a faster workflow and enhanced data quality, the TruSeq Nano DNA Sample Preparation Kit provides an all-inclusive sample preparation method for genome sequencing applications.

### References

1. genome.ucsc.edu
2. www.broadinstitute.org

**Figure 6: TruSeq Nano DNA Protocol Reduces Number of Coverage Gaps**



Increased coverage of TruSeq Nano DNA libraries results in fewer coverage gaps, demonstrated here in the GC-rich coding regions of the *RINPEPL1* promoter (A) and the *ZBTB34* promoter (B). Sequence information generated by TruSeq Nano DNA prep is shown in the top panels of A and B, while sequence data generated using TruSeq DNA protocol are shown in the lower panels.

### Ordering Information

Product	Catalog No.
TruSeq Nano DNA LT Sample Preparation Kit Set A (24 samples)	FC-121-4001
TruSeq Nano DNA LT Sample Preparation Kit Set B (24 samples)	FC-121-4002
TruSeq Nano DNA HT Sample Preparation Kit (96 samples)	FC-121-4003

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

#### FOR RESEARCH USE ONLY

© 2013 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAlix, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Pub. No. 770-2013-012 Current as of 16 May 2013

AGAATGATAACAGTAAACACACTTGTGTTAAACCTTAAGATTACTTGTATCCACTGATTCACACGTACCGTAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
TCAACGTACCGTAAACGACTATATCATTAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
CGACGAAAGAAATGATAACAGTAAACACACTTGTGTTAAACCTTAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
RACGTACCATTAAGAGCTAACCGTCTTCTGTTAAACCTTAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
AGATGATAACAGTAAACACACTTGTGTTAAACCTTAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
GATTACTTGTATCCACTGATTCACACGTAAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC  
GGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACACGTACCGTAAACGAAACGTATCAATTGAGACTAAATTAACGTAACACATTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTGTATCCACTGATTCACAC























# Nextera® Rapid Capture Exomes

A rapid workflow and comprehensive exome content, with unparalleled flexibility.

### Highlights

- Rapid exome preparation and enrichment**  
Prep and enrich 96 exomes in only 1.5 days with less than 5 hours hands-on time
- Comprehensive exome coverage**  
Two different exome designs are available to access core exonic content or expanded content
- Kit configurations designed to fit your needs**  
Choose the optimal fit for your system, samples, and study, with more flexible options than ever before
- Complete support for entire process from sample preparation to sequencing**  
All-in-one kit for prep and enrichment from the world's leading sequencing provider

### Overview

Nextera Rapid Capture Exomes are all-in-one kits for sample preparation and exome enrichment that allow researchers to identify coding variants up to 70% faster than other methods. Nextera Rapid Capture Exome delivers 37 Mb of expertly selected exonic content, including challenging regions excluded from other exome designs.

### Rapid Exome Prep and Enrichment

Nextera Rapid Capture Exomes provide sample prep and exome enrichment in only 1.5 days. Sequencing with the HiSeq® 2500 or NextSeq™ 500 system enables experiments to go from DNA sample to data in as little as 2.5 days. The speed of Nextera Rapid Capture Exomes enables you to complete projects faster, return results faster, and ultimately publish faster.

### Focused Exonic Content

Nextera Rapid Capture Exome has been optimized to provide uniform and specific coverage of 37 Mb of expert-selected exonic content. The probe set was designed to enrich 214,405 exons (Table 1). This focused design, paired with uniform and specific enrichment, enables the most comprehensive exome sequencing available and reliable identification of true, coding variants (Table 2).

Table 1: Coverage Details

	Nextera Rapid Capture Exome	Nextera Rapid Capture Expanded Exome
<b>Coverage Specifications</b>		
Number of target exons	214,405	201,121
Target content	Coding exons	Exons, UTRs, and miRNA
<b>Percent of Exome Covered (by Database)</b>		
Refseq	98.3%	95.3%
CCDS	98.6%	96.0%
ENSEMBL	97.8%	90.6%
GENCODE v12	98.1%	91.6%

Table 2: Comparison of Rapid Capture Exomes

Specification	Nextera Rapid Capture Exome	Nextera Rapid Capture Expanded Exome
Target size	37 Mb	62 Mb
Genomic DNA input		50 ng
Hands-on time		5 hours
Total time		1.5 days
Batch size		1–96 exomes

```

AGAATGATAACAGTAAACACCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTTAACGAAAGCTATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAAC
TCAACGTAACCGTAAACGATATCATTAAGATTACTTTGATCCACTGATTCAACGTAACCGTTAACGAAAGCTATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTAAACGAT
CGACGAAAAGGATGATAACAGTAAACACCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTTAACGAAAGCTATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAAC
TAGCTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACG
GATTACTTTGATCCACTGATTCAACGTTAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACGTAACCTT
CGTATCAATTGAGACTAAATTAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACGTAACCTTAAGAGCTAACCGTCTTCTGTTAACCTTAAGATTACTTTGATCCACTGATTCAACGTAACCGTAAACGATATCAATTGAGACTAAATTAACGTAACCTTAAACGAT

```





Table 3: Nextera Rapid Capture Throughput by Illumina Sequencing Systems

Pooling Plexity	Exome Samples per Run				
	MiSeq	NextSeq 500— Mid Output	NextSeq 500— High Output	HiSeq 2500— Rapid Run Mode	HiSeq 2500— High Output
1	Up to 1	–	–	–	–
3	–	Up to 3	–	–	–
6	–	–	Up to 6	Up to 24	Up to 96
9	–	–	Up to 9	Up to 24	Up to 115
12	–	–	Up to 12	Up to 24	Up to 115

Table 3 helps identify which options provide optimal alignment across three vital study design considerations: sequencing instrument, number of exome samples sequenced per run, and the number of exome samples pooled together before enrichment (pooling plexity).

Integrated sample barcodes then allow the pooling of up to 12 samples for a single exome Rapid Capture pull down. Next, libraries are denatured into single-stranded DNA (Figure 1B) and biotin-labeled probes specific to the targeted region are used for the Rapid Capture hybridization (Figure 1C).

The pool is enriched for the desired regions by adding streptavidin beads that bind to the biotinylated probes (Figure 1D). Biotinylated DNA fragments bound to the streptavidin beads are magnetically pulled down from the solution (Figure 1E). The enriched DNA fragments are then eluted from the beads and hybridized for a second Rapid Capture. This entire process is completed in only 1.5 days, enabling a single researcher to efficiently process up to 96 exomes at one time—all without automation.

## Summary

Nextera Rapid Capture Exomes provide a fully integrated, rapid solution for exome library prep and enrichment. Available in a wide range of kit configurations (Table 3), as well as two unique designs, Nextera Rapid Capture Exomes provide unparalleled flexibility to optimally align with your specific needs.

## References

- Manolio TA, Collins FS, Cox NJ, Goldstein DB, Hindorf LA, et al. (2009) Finding the missing heritability of complex diseases. *Nature* 4618: 747–753.

## Ordering Information

Kit Description	Catalog No.
Nextera Rapid Capture Exome (8 rxn x 1 plex)	FC-140-1000
Nextera Rapid Capture Exome (8 rxn x 3 plex)	FC-140-1083
Nextera Rapid Capture Exome (8 rxn x 6 plex)	FC-140-1086
Nextera Rapid Capture Exome (8 rxn x 9 plex)	FC-140-1089
Nextera Rapid Capture Exome (2 rxn x 12 plex)	FC-140-1001
Nextera Rapid Capture Exome (4 rxn x 12 plex)	FC-140-1002
Nextera Rapid Capture Exome (8 rxn x 12 plex)	FC-140-1003
Nextera Rapid Capture Expanded Exome (2 rxn x 12 plex)	FC-140-1004
Nextera Rapid Capture Expanded Exome (4 rxn x 12 plex)	FC-140-1005
Nextera Rapid Capture Expanded Exome (8 rxn x 12 plex)	FC-140-1006

AAAGAATGATAACAGTAAACAACACTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
AATCAACGTACCGTAAACGAACACTTATAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
AACGACAAAAGAATGATAACAGTAAACAACACTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
TTAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
AAAGATTGATAACAGTAAACAACACTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
AAGATTACTTGATCCACTGATTCAACGTAAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA  
AACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACGAACGTATCAATTGAGACTAAATATAACGTACCATTAAAGAGTACCCTCTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCA



# Nextera® Rapid Capture Custom Enrichment

Leverage a superior sample preparation and enrichment workflow for unparalleled access to your regions of interest.

### Highlights

- Integrated sample preparation and enrichment workflow**  
 Nextera tagmentation and optimized hybridization reduce workflow duration and generate data faster
- Target your regions of interest**  
 Choose 0.5-15 Mb of custom content, and pool up to 12 samples per enrichment reaction
- Evolve your design with add-on content**  
 Supplement existing panels and keep adding on as your research needs expand

### Introduction

Nextera Rapid Capture Custom Enrichment is an all-in-one assay for sample preparation and custom target enrichment. Nextera tagmentation coupled with optimized target capture ensures the fastest enrichment workflow time for your custom content. The flexible, fully customizable design accommodates up to 15 Mb of custom content so you can focus on the regions of the genome that you care about. The new add-on feature in DesignStudio allows you to iteratively expand your content as new discoveries are made.

### Custom Probe Design


The first step in developing any Nextera Rapid Capture Custom Enrichment assay is to design your custom probe set. DesignStudio is a free online user-friendly tool accessed through your Myllumina account. Designate your regions of interest, refine your custom probe set and place an order for your custom design. DesignStudio uses a complex algorithm to optimize probe set design and alert you to any potential coverage gaps or challenging regions. Desired targets can be added individually or in batches by chromosomal coordinate or gene name.

### Unmatched Ease of Workflow

Nextera Rapid Capture Enrichment allows researchers to maximize the productivity of their lab personnel and Illumina sequencing technology. The simplicity and speed of the Nextera Rapid Capture assay enables a single technician to prepare and enrich 12 samples in only 1.5 days.

Nextera-based sample preparation generates adapter-tagged libraries from 50 ng input genomic DNA (Figure 2A). Nextera tagmentation of DNA simultaneously fragments and tags DNA without the need for mechanical shearing. Integrated sample barcodes allow the pooling of up to 12 of these adapter ligated sample libraries into a single, hybridization-based, pull down reaction. The pooled libraries are then

Figure 1: Overview of Nextera Rapid Capture Custom Enrichment



The Nextera Rapid Capture Custom Enrichment Kit is an integral part of a complete and fully supported solution for targeted resequencing.

denatured into single-stranded DNA (Figure 2B) and biotin-labeled probes complementary to the targeted region are used for the Rapid Capture hybridization (Figure 2C). Streptavidin beads are added, which bind to the biotinylated probes that are hybridized to the targeted regions of interest (Figure 2D). Magnetic pull down of the streptavidin beads enriches the targeted regions that are hybridized to biotinylated probes. (Figure 2E). The enriched DNA fragments are then eluted from the beads and a second round of Rapid Capture is completed to increase enrichment specificity. The entire process is completed in only 1.5 days, enabling a single researcher to efficiently process up to 12 samples at one time—all without automation.

### Data Analysis


Sequence data generated from custom enrichment samples on HiSeq® and NextSeq™ systems are analyzed using the Enrichment Workflow in the HiSeq Analysis Software (HAS). HAS analysis can be accessed directly via a linux kernel or by using the optional Analysis Visual Controller (AVC) interface¹.

Custom pools sequenced on MiSeq® are analyzed using MiSeq Reporter (MSR). The Enrichment Workflow from both HAS and MSR generates aligned sequence reads in the .bam format using the BWA algorithm and performs indel realignment using the GATK indel realignment tool. Variant calling occurs in the target regions specified in the manifest file. The GATK variant caller generates .vcf

GAATGATAACAGTAACACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAAC  
 AACGTACCGTAAACGAACGATGATTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAAC  
 GACGAAAAGATGATAACAGTAACACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCT  
 ACGTACGATTAAGAGCTACCGTGAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCT  
 GAATGAACAGTAACAGACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAAC  
 ATTAATTGATCCACTGATTCAACGTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAAC  
 GATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCCTGAACGAAGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGTACCGTCTTCTGTTAAACCTTAAGATTACTTGATCCACTGATTCAAC



Figure 3: High Coverage Uniformity Across Custom 12-plex Pools



Nextera Rapid Capture Custom Enrichment provides uniform target enrichment across different custom probe sets and individual samples within a 12-plex pool. A. Coverage uniformity is shown as % of targeted bases that are represented by >0.2x mean coverage. Mean coverage for these custom probe sets can be found in Table 1. Error bars show SD of uniformity across the 12 pooled samples for each project. B. Coverage uniformity for each of 12 pooled samples within Project 3 is shown. Mean coverage for this run was 300x, and % of targeted bases that were covered at > 60x are shown.

Table 1: Sequencing Details for Example Projects


Project	Content	Mean Coverage	% On Target Bases*
1**	0.5 Mb	1500x	88.6
2**	0.5 Mb	146x	79.5
3†	3.5 Mb	300x	80.1
4**	7 Mb	152x	72.5

\*Calculated using Picard Hybrid Selection tool with 250 bp padding?  
 \*\*Sequenced on HiSeq  
 †Sequenced on MiSeq

**Summary**

Nextera Rapid Capture Custom Enrichment leverages a superior integrated sample prep and enrichment workflow to provide unparalleled access to your genomic regions of interest. Not only will you be able to perform targeted sequencing using only 50 ng of input DNA, you'll do so faster and more efficiently than ever before. Take advantage of robust add-on functionality to refine your content over time, or add regions of unique interest to established panels such as Nextera Rapid Capture Exome or other TruSight™ content sets.

Figure 4: Add-On Content Retains High Coverage



High coverage uniformity is maintained when 3.5 Mb of add-on content is added to the Nextera Rapid Capture Exome. All samples were run as 12-plex pools.

AAGAAATGATAACAGTAAACACACTTGTGTTAAACCTTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 ATCAACGTAACGTAACGAAAGTATCAITTAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 ACGACGAAAGATGATAACAGTAAACACACTTGTGTTAAACCTTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 TAAGGTACATTAAAGAGTACCSTGTTAAACAGTAAACACACTTGTGTTAAACCTTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 AAGATGTTAAACAGTAAACAGTATTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 AGATTACTTGATCCACTGATTC AACGTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC  
 ACGTATCAATTGAGACTAAATATTAAAGTACCATTAAAGATTACTTGATCCACTGATTC AACGTACCGTAAACGAACGATACATTGAGACTAAATATTAAAGTACCATTAAAGAGTACCSTGTTCTGTAAACCTTAAAGATTACTTGATCCACTGATTC





# Methyl-Seq

## EpiGnome™ Methyl-Seq Kit

Unlock limited samples (50-100ng DNA input) to discover methylation patterns of all CpG, CHH & CHG regions.

- ▶ Unlock small samples (50-100ng DNA input)
- ▶ Pre-library bisulfite conversion
- ▶ Comprehensive, whole genome results
- ▶ 5 hour method
- ▶ Informatics app note demystifies analysis
- ▶ **Capture full sample diversity**




### Sequence the entire sample—no loss of information!

The process of bisulfite treatment denatures genomic DNA into single stranded DNA. EpiGnome converts single stranded DNA into an Illumina® sequencing library. All ssDNA fragments are captured into an Illumina sequencing library during the EpiGnome procedure, therefore eliminating sample loss associated with other methods.

Calico cats are domestic cats with a spotted or parti-colored coat that is predominantly white, with patches of two other colors. Calico cats are almost always female because the X chromosome determines the coat color. During embryonic development, one X chromosome is hypermethylated and inactivated. The remaining X chromosome determines coat color.

Figure 1. EpiGnome is sensitive to CpG methylation patterns.




CpG methylation patterns across region of chr1 show variable CpG methylation (red) from 50 ng input of GM12878 lymphoblastoid gDNA treated with bisulfite. Comparison to coverage patterns from non-bisulfite treated (green) gDNA shows the methylated regions of chromosome 1.

### Workflow



Figure 2. Deep coverage of genes of interest.



**ACMG\_genes:** designated as medically relevant.

**Cancer genes:** protein coding genes known to be involved in cancer.

**fosmid\_promoters:** of high interest and difficult to sequence.

EpiGnome WGBS method yields high coverage of genes of interest for Cancer genes and those that have been defined as medically relevant by the American College of Medical Genetics.

### Deep coverage of critical genomic regions

Depth of coverage is enhanced in genomic areas with biological utility (Figure 2). EpiGnome captures full sample diversity of critical areas including:

- Coding region start and end for exons from the canonical transcript of protein coding genes for genes known to be involved in cancer, taken from SOMA and CRUK panels as well as literature derived Cancer genes.
- Genes defined by the American College of Medical Genetics as being medically relevant (ACMG\_genes)
- Exonic coding regions from Ensemble 70 (exons\_ensemble70)
- List of 100 promoters defined by the Broad Institute as being of high interest and difficult to sequence (fosmid\_promoters)

Coverage was obtained from 125.4 million reads in a single lane of a HiSeq. Increasing throughput of the HiSeq Systems enables complete methylation information to be captured from a growing number of samples.

### Success begins with purification

#### MasterPure™ DNA Purification Kit

Purification is an important step to prepare your sample. MasterPure safely removes unwanted material to give you pure DNA.

MasterPure offers unique benefits:

- ▶ Very high yields
- ▶ Recover <90% of theoretical yield
- ▶ Safe and nontoxic
- ▶ Available for all sample sizes

Cat. # Quantity

#### MasterPure™ Complete DNA and RNA Purification Kit

MC85200	200 Purifications
MC89010	10 Purifications

Cat. # Quantity

#### EpiGnome™ Methyl-Seq Kit

EGMK81312	12 reactions
EGMK91324	24 reactions
EGMK91396	96 reactions

#### EpiGnome™ Index PCR Primers

EGIDX81312	12 indexes, 10 reactions each
------------	----------------------------------

#### FailSafe™ PCR Enzyme Mix

FSE51100	100 units
----------	-----------


The FailSafe™ PCR Enzyme Mix is required for EpiGnome Methyl-Seq Kit.





Data Sheet: Sequencing

Figure 2: Bioanalyzer Trace of Mafk Library



Bioanalyzer trace data for a library generated for transcription factor target Mafk using the TruSeq ChIP Sample Preparation Kit with 5 ng of input DNA. The center peak indicates robust yield within the desired insert size range.

The TruSeq ChIP process begins with the enrichment of specific cross-linked DNA-protein complexes using an antibody against a protein of interest (Figure 1A-B). The stretches of DNA bound to the target protein are then isolated and used as input DNA for library generation. DNA fragments are end-repaired and an 'A'-base added to the blunt ends of each strand, preparing them for ligation to the sequencing adapters (Figure 1C-D). Each TruSeq adapter contains a 'T'-base overhang on the 3'-end providing a complementary overhang for ligating the adapter to the A-tailed fragmented DNA (Figure 1E). Final product is created (Figure 1F) and after size selection, all of the ChIP DNA fragments are simultaneously sequenced.

For maximum flexibility, TruSeq ChIP Sample Preparation Kits can be used to prepare samples for single-read or paired-end sequencing, and are compatible with any Illumina sequencing instrument, including MiSeq® and all instruments in the HiSeq® system family.

Table 1: Motif-Finder Analysis of Peaks Identified using TruSeq Sample Preparation Kits Compared to ENCODE Reference Peak Data

Name	% Top Peaks with Mafk Motif
TruSeq ChIP	95%
ENCODE HELA	92%
ENCODE HES	86%


TruSeq Data Quality

Proven TruSeq data quality delivers the most complete and accurate profile of target protein-DNA interactions, enabling an optimal percentage of passing filter reads, percent alignable reads, and coverage uniformity, as well as high sensitivity to detect low-abundance hits.

Robust Multiplex Performance

The TruSeq ChIP Sample Preparation Kits provide up to 24 total indexes to increase throughput and consistency without compromising results. The TruSeq universal adapters ligate to sample fragments during library construction, allowing samples to be pooled and individually identified during downstream analysis. This indexing capability improves workflow efficiency and enables robust multiplex sequencing. By enhancing study design flexibility, indexing aids researchers in deriving the most value from each run by efficiently distributing read output based on optimal per-sample read depth requirements.

Figure 3: Peak Finding Output for Mafk



TruSeq ChIP Sample Preparation Kits enable the generation of libraries across a broad range of study designs. Above is peak data for a negative Ig control, the transcription factor target Mafk, and a reference peak for Mafk from the ENCODE database.

AGAATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAGAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTAA  
 TCAACGTACCGTAAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 CGACGAAAGGATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 VACGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 AGATGATAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 GATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 TCAACGTACCGTAAACAGTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT  
 GATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATTCAACGTACCGTAAACCGTATCAATTGAGACTAAATATTAAGGTACCACTTAAGAGCTACCGTCTCTGTAAACCTTAAGATTACTTGATCCACTGATT



## Data Sheet: Sequencing

### Summary

TruSeq ChIP Sample Preparation Kits offer proven TruSeq accuracy, and a simple, streamlined workflow, enabling highly-multiplexed, cost-effective ChIP sequencing. Supporting analysis of a broad range of targets across the genome even from low sample input, the kits provide a complete, accurate profile of DNA-protein binding interactions and enhanced visibility to the mechanics of gene regulation.

### References

1. Johnson DS, Mortazavi A, Myers RM, Wold B (2007) Genome-wide mapping of *in vivo* protein-DNA interactions. *Science* 316: 1497–1502.
2. Barski A, Cuddapah S, Cui K, Roh TY, Schonnes DE et al. (2007) High-resolution profiling of histone methylations in the human genome. *Cell* 129: 823–837.
3. Marban C, Su T, Ferrari R, Li B, Vatakis D, et al. (2011) Genome-wide binding map of the HIV-1 Tat protein to the human genome. *PLoS One* 6: e26894.
4. Fujiki R, Hashiba W, Sekine H, Yokoyama A, Chikanishi T, et al. (2011) GlcNAcylation of histone H2B facilitates its monoubiquitination. *Nature* 480: 557–560.
5. Botti E, Spallone G, Moretti F, Marinari B, Pinetti V, et al. (2011) Developmental factor IRF6 exhibits tumor suppressor activity in squamous cell carcinomas. *Proc Natl Acad Sci U S A* 108: 13710–13715.
6. Bernt KM, Zhu N, Sinha AU, Vempati S, Faber J, et al. (2011) MLL-rearranged leukemia is dependent on aberrant H3K79 methylation by DOT1L. *Cancer Cell* 20: 66–78.
7. de Almeida SF, Grosso AR, Koch F, Fenouil R, Carvalho S, et al. (2011) Splicing enhances recruitment of methyltransferase HYPB/Setd2 and methylation of histone H3 Lys36. *Nat Struct Mol Biol* 18: 977–983.
8. Wu H, D'Alessio AC, Ito S, Xia K, Wang Z, et al. (2011) Dual functions of Tet1 in transcriptional regulation in mouse embryonic stem cells. *Nature* 473: 389–393.
9. ENCODE Project Consortium, Myers RM, Stamatoyannopoulos J, Snyder M, Dunham I, Hardison RC, et al. (2011) A user's guide to the encyclopedia of DNA elements (ENCODE). *PLoS Biol*. 9:e1001046. PMID: 21526222; PMCID: PMC3079585.

### Ordering Information

Product	Catalog No.
TruSeq ChIP Sample Preparation Kit, Set A (12 indexes, 48 samples)	IP-202-1012
TruSeq ChIP Sample Preparation Kit, Set B (12 indexes, 48 samples)	IP-202-1024

illumina, inc. • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

#### FOR RESEARCH USE ONLY

© 2012 Illumina, Inc. All rights reserved.  
Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GALx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
Pub. No. 770-2012-029 Current as of 22 August 2012

illumina®

AAGAATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AATGAGCTACCGTAAAGAGAGCTATGATTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AACGAGCGAAAGATGATACACAGTAAACCAACTTCTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
TTAAGCTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AAAGGATGATAACAGTAACACACTTCTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA  
AAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTAACGTAAGAACAGCTATCAATTGAGACTAAATATTAAAGTACCAATTAAGAGCTACCGTCTTGTGTTAACCTTAAGATTACTTGATCCACTGATCAACGTA

## DNA-Sequencing

Description	Catalog Number
MasterPure™ Complete DNA and RNA Purification Kit	MC85200
MasterPure™ DNA Purification Kit	MCD85201
TruSeq DNA PCR-Free LT Sample Preparation Kit - Set A	FC-121-3001
TruSeq DNA PCR-Free LT Sample Preparation Kit - Set B	FC-121-3002
TruSeq DNA PCR-Free HT Sample Preparation Kit	FC-121-3003
TruSeq Nano DNA LT Sample Preparation Kit - Set A	FC-121-4001
TruSeq Nano DNA LT Sample Preparation Kit - Set B	FC-121-4002
TruSeq Nano DNA HT Sample Preparation Kit	FC-121-4003
Nextera Rapid Capture Exome (8 rxn x 1 Plex)	FC-140-1000
Nextera Rapid Capture Exome (8 rxn x 3 Plex)	FC-140-1083
Nextera Rapid Capture Exome (8 rxn x 6 Plex)	FC-140-1086
Nextera Rapid Capture Exome (8 rxn x 9 Plex)	FC-140-1089
Nextera Rapid Capture Exome (2 rxn x 12 Plex)	FC-140-1001
Nextera Rapid Capture Exome (4 rxn x 12 Plex)	FC-140-1002
Nextera Rapid Capture Exome (8 rxn x 12 Plex)	FC-140-1003
Nextera Rapid Capture Expanded Exome (2 rxn x 12 Plex)	FC-140-1004
Nextera Rapid Capture Expanded Exome (4 rxn x 12 Plex)	FC-140-1005
Nextera Rapid Capture Expanded Exome (8 rxn x 12 Plex)	FC-140-1006
EpiGnome™ Methyl-Seq Kit	EGMK81312

## ChIP

Description	Catalog Number
TruSeq ChIP Sample Preparation Kit - Set A	IP-202-1012
TruSeq ChIP Sample Preparation Kit - Set B	IP-202-1024

## Methylation Arrays

Description	Catalog Number
HumanMethylation450 DNA Analysis BeadChip Kit (24 samples)	WG-314-1003
HumanMethylation450 DNA Analysis BeadChip Kit (48 samples)	WG-314-1001
HumanMethylation450 DNA Analysis BeadChip Kit (96 samples)	WG-314-1002

## RNA-sequencing

Description	Catalog Number
MasterPure™ Complete DNA and RNA Purification Kit	MC85200
TotalScript™ RNA-Seq Kit	TSRNA 12924
ScriptSeq™ Complete Gold Kit (Blood)	BGGB1306
ScriptSeq™ Complete Gold Kit (Blood) - Low Input	SCL24GBL
Ribo-Zero Magnetic Gold Kit (Yeast)	MRZY1324
ScriptSeq™ Complete Gold Kit (Yeast)	BGY1324
ScriptSeq™ Complete Gold Kit (Yeast) - Low Input	SCGL6Y
ARTseq™ Ribosome Profiling Kit - Mammalian	RPHMR12126
ARTseq™ Ribosome Profiling Kit - Yeast	RPYSC12116
<b>Any species</b>	
TruSeq® Stranded mRNA LT Set A	RS-122-2101
TruSeq® Stranded mRNA LT - Set B	RS-122-2102
TruSeq® Stranded mRNA HT	RS-122-2103
TruSeq™ RNA Sample Prep Kit v2 -Set A (48rxn)	RS-122-2001
TruSeq™ RNA Sample Prep Kit v2 -Set B (48rxn)	RS-122-2002
<b>Human/Mouse/Rat</b>	
TruSeq® Strnd Total RNA LT(w/Ribo-Zero™ Human/Mouse/Rat)Set A	RS-122-2201
TruSeq® Strnd Total RNA LT(w/Ribo-Zero™ Human/Mouse/Rat)Set B	RS-122-2202
TruSeq® StrndTotal RNA HT (w/ Ribo-Zero™ Human/Mouse/Rat)	RS-122-2203
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Gold) Set A	RS-122-2301
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Gold) Set B	RS-122-2302
TruSeq® Stranded Total RNA HT (w/ Ribo-Zero™ Gold)	RS-122-2303
<b>Human/Mouse/Rat (Blood-derived)</b>	
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Globin) Set A	RS-122-2501
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Globin) Set B	RS-122-2502
TruSeq® Stranded Total RNA HT (w/ Ribo-Zero™ Globin)	RS-122-2503
<b>Plant</b>	
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Plant) Set A	RS-122-2401
TruSeq® Stranded Total RNA LT (w/ Ribo-Zero™ Plant) Set B	RS-122-2402
TruSeq® Stranded Total RNA HT (w/ Ribo-Zero™ Plant)	RS-122-2403

### Small RNA-sequencing

Description	Catalog Number
TruSeq® Small RNA Sample Prep Kit -Set A	RS-200-0012
TruSeq® Small RNA Sample Prep Kit -Set B	RS-200-0024
TruSeq® Small RNA Sample Prep Kit -Set C	RS-200-0036
TruSeq® Small RNA Sample Prep Kit -Set D	RS-200-0048

### Targeted RNA-Sequencing


Description	Catalog Number
TruSeq Targeted RNA Expression Custom Components	
TruSeq Targeted RNA Custom Kit (48 Samples)	RT-101-1001
TruSeq Targeted RNA Custom Kit (96 Samples)	RT-102-1001
TruSeq Targeted RNA supplemental content (48 Samples)	RT-801-1001
TruSeq Targeted RNA supplemental content (96 Samples)	RT-802-1001
TruSeq Targeted RNA Index Kit	RT-401-1001
<b>TruSeq Targeted RNA Expression Fixed Panels</b>	
TruSeq Targeted RNA Apoptosis Panel Kit (48 Samples)	RT-201-1010
TruSeq Targeted RNA Apoptosis Panel Kit (96 Samples)	RT-202-1010
TruSeq Targeted RNA Cardiotoxicity Panel Kit (48 Samples)	RT-201-1009
TruSeq Targeted RNA Cardiotoxicity Panel Kit (96 Samples)	RT-202-1009
TruSeq Targeted RNA Cell Cycle Panel Kit (48 Samples)	RT-201-1003
TruSeq Targeted RNA Cell Cycle Panel Kit (96 Samples)	RT-202-1003
TruSeq Targeted RNA Cytochrome p450 Panel Kit (48 Samples)	RT-201-1006
TruSeq Targeted RNA Cytochrome p450 Panel Kit (96 Samples)	RT-202-1006
TruSeq Targeted RNA HedgeHog Panel Kit (48 Samples)	RT-201-1002
TruSeq Targeted RNA HedgeHog Panel Kit (96 Samples)	RT-202-1002
TruSeq Targeted RNA Neurodegeneration Panel Kit (48 Samples)	RT-201-1001
TruSeq Targeted RNA Neurodegeneration Panel Kit (96 Samples)	RT-202-1001
TruSeq Targeted RNA NFkB Panel Kit (48 Samples)	RT-201-1008
TruSeq Targeted RNA NFkB Panel Kit (96 Samples)	RT-202-1008
TruSeq Targeted RNA Stem Cell Panel Kit (48 Samples)	RT-201-1005
TruSeq Targeted RNA Stem Cell Panel Kit (96 Samples)	RT-202-1005
TruSeq Targeted RNA TP53 Pathway Panel Kit (48 Samples)	RT-201-1007
TruSeq Targeted RNA TP53 Pathway Panel Kit (96 Samples)	RT-202-1007
TruSeq Targeted RNA Wnt Pathway Panel Kit (48 Samples)	RT-201-1004
TruSeq Targeted RNA Wnt Pathway Panel Kit (96 Samples)	RT-202-1004



## RNA-Seq without rRNA Depletion

### TotalScript™ RNA-Seq Kits

- ▶ Powered by Nextera™
- ▶ 1-5 ng Input RNA
- ▶ Designed for precious samples
- ▶ Rapid method with 5 hr workflow
  
- ▶ rRNA removal not required, begin with total RNA
- ▶ Directional libraries
- ▶ 12 Indexes available




Powered by Nextera! TotalScript RNA-Seq Kit is designed for RNA-Seq of precious samples, and only 1-5 ng of intact total RNA is needed for each sample. No need for poly(A) enrichment or rRNA removal. Sequencing data is similar to data from libraries using much more RNA.

#### New kinds of samples can now be sequenced, including:

- ▶ Cancer samples
- ▶ Stem cells
- ▶ Other low input samples

Figure 1. Consistent results from TotalScript.



#### Retain more sample

#### Prevent transcript loss


TotalScript produces consistent results from small amounts of sample (Fig. 1). 1 ng or 5 ng of total RNA was prepared with TotalScript. Results show similar amounts of coding and non-coding coverage between samples. The sample was Universal Human Reference RNA (UHR) total RNA.

#### Workflow



Figure 2. You choose the coverage profile.

Method	Total RNA Input (ng)	% rRNA	Coverage
Random Priming	1-5	<40%	Even
Mixed Priming	1-5	<25%	Slight 3' Bias
dT Priming	1-5	<5%	3' Bias



**You choose the desired rRNA content and transcript coverage with TotalScript™**

Three options are included in every TotalScript kit (Fig. 2). All options produce directional libraries from very small amounts of total RNA.

1. Random Hexamer Primer option produces even transcript coverage with <40% of reads mapping to rRNA.
2. Mixed Primer option produces good transcript coverage with <25% rRNA mapped reads.
3. Oligo(dT) Primer option produces <5% rRNA reads with transcript coverage strongest at the 3' end.

Different sources of RNA may produce different levels of rRNA contamination.

TotalScript RNA-Seq libraries shown were made from 5 ng of total UHR RNA using the Optimized Buffer included with TotalScript (Fig 2).

**Success begins with purification**

**MasterPure™ RNA Purification Kit**

Purification is an important step to prepare your sample. MasterPure safely removes unwanted material to give you pure, intact total RNA.

MasterPure offers unique benefits:

- ▶ Keep RNA intact (does not degrade RNA)
- ▶ Retain RNA diversity (including small RNA)
- ▶ Maximize genes discovered
- ▶ Available for all sample sizes

Cat. #	Quantity
<b>MasterPure™ RNA Purification Kit (for isolating RNA only)</b>	
MCR85102	100 Purifications

Cat. #	Quantity
<b>TotalScript™ RNA-Seq Kit</b>	
TSRNA12924	24 Reactions
TSRNA1296	12 Reactions
<b>TotalScript™ Index Kit</b>	
TSIDX12910	11 indexes









Data Sheet: Illumina® Sequencing

Figure 5. Illumina's Complete Sequencing Solution



Ordering Information

Product	Catalog No.
<b>For RNA Preparation</b>	
TruSeq RNA Sample Preparation Kit v2, Set A (12 indexes, 48 samples)	RS-122-2001
TruSeq RNA Sample Preparation Kit v2, Set B (12 indexes, 48 samples)	RS-122-2002
<b>For DNA Preparation</b>	
TruSeq DNA Sample Preparation Kit v2, Set A (12 indexes, 48 samples)	FC-121-2001
TruSeq DNA Sample Preparation Kit v2, Set B (12 indexes, 48 samples)	FC-121-2002
<b>For Cluster Generation on cBot and Sequencing on the HiSeq 2000/1000 and HiScanSQ</b>	
TruSeq Paired-End Cluster Kit v3—cBot—HS (1 flow cell)	PE-401-3001
TruSeq Single-Read Cluster Kit v3—cBot—HS (1 flow cell)	GD-401-3001
<b>For Cluster Generation on cBot and Sequencing on the Genome Analyzer<sub>FLX</sub></b>	
TruSeq Paired-End Cluster Kit v2—cBot—GA (1 flow cell)	PE-300-2001
TruSeq Single-Read Cluster Kit v2—cBot—GA (1 flow cell)	GD-300-2001
<b>For Cluster Generation on the Cluster Station and Sequencing on the Genome Analyzer<sub>FLX</sub></b>	
TruSeq Paired-End Cluster Kit v5—CS—GA (1 flow cell)	PE-203-5001
TruSeq Single-Read Cluster Kit v5—CS—GA (1 flow cell)	GD-203-5001

Illumina, Inc. • 9885 Towne Centre Drive, San Diego, CA 92121 USA • 1.800.809.4566 toll-free • 1.858.202.4566 tel • techsupport@illumina.com • illumina.com

FOR RESEARCH USE ONLY

© 2011 Illumina, Inc. All rights reserved. Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CSPRO, DASL, Eco, Genetic Energy, GALiX, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. Pub. No. 970-2009-039 Current as of 27 April 2011



```

AGAATGATAACAGTAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
TCAGCGTACCGTAAACGAAAGGTTATCAATTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
CGACGAAAGGATTGATAGACTAAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
TAGCTACCATTAAGAGCTACCGTCCACAGATAACACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
AGATTGATAACAGTAAACACTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
GATTACTTGATCCACTGATTCACGTAACGTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG
GATTAACCTTAAGATTACTTGATCCACTGATTCACGTAACCGTAAACGAAACGATCAATTGAGACTAAATATAACGTACCATTAAAGACTACCGTCTTCTGTAAACCTTAAGATTACTTGATCCACTGATTCACG

```















Figure 5: Visualization of TruSeq Targeted RNA Expression Data using MiSeq Reporter



Data visualization with MiSeq Reporter allows easy comparison of data sets.

**Product Specifications**

Specification	Value
Database content	> 400,000 designs (mouse, human, rat)
Target types	Gene, transcript, exon, splice junction, cSNP, fusion
Dynamic range	5 orders of magnitude
Time to answer	1.5 days
Hands-on time	4 hours
RNA quality	> 200 bp unfixed or FFPE

**Simple Data Analysis**

After a sequencing run on the MiSeq system, data are automatically aligned and can be viewed using the MiSeq Reporter. As shown in Figure 5, pairwise comparisons for relative expression between samples or groups of samples is simple and intuitive. Customizable significance thresholds allow you to quickly identify differentially expressed targets. The TruSeq Targeted RNA Expression user experience is customized and streamlined, and keeps project data highly accessible.

**Summary**

Designed for the MiSeq system, TruSeq Targeted RNA Expression provides rapid and economical RNA profiling and validation for your gene expression studies. Go from sample to answer in less than two days with a simple, streamlined workflow and automated data visualization. Choose validated, pre-designed panels or add custom content to your existing assays for the ultimate flexibility to evolve your research.

**References**

1. <https://icm.illumina.com/>
2. Considerations for Designing a Successful TruSeq Targeted RNA Expression Experiment Technical Note, 2013.

ATGATACAGTAAACAACTTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
ACGTACCGTAAACAACTTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
GCGAAGGATGATACAGTAAACAACTTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
GTACCAATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
ATGATACAGTAAACAACTTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
TACTTGATCCACTGATTCAACGTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA  
ATGATACAGTAAACAACTTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTACCGTAAACGAAACGATCAATTGAGACTAAATATAAGCTACCATTAAAGAGCTACCGCTCTCTGTAACTTAAAGTACTTGATCCACTGATTCAACGTA

Data Sheet: Sequencing

Ordering Information

Product Name	Number of Samples	Catalog No.
<b>TruSeq Targeted RNA Expression Custom Components</b>		
TruSeq Targeted RNA Custom Kit	48	RT-101-1001
	96	RT-102-1001
TruSeq Targeted RNA Supplemental Content	48	RT-801-1001
	96	RT-802-1001
<b>TruSeq Targeted RNA Expression Fixed Panels</b>		
TruSeq Targeted RNA Apoptosis Panel Kit	48	RT-201-1010
	96	RT-202-1010
TruSeq Targeted RNA Cardiotoxicity Panel Kit	48	RT-201-1009
	96	RT-202-1009
TruSeq Targeted RNA Cell Cycle Panel Kit	48	RT-201-1003
	96	RT-202-1003
TruSeq Targeted RNA Cytochrome p450 Panel Kit	48	RT-201-1006
	96	RT-202-1006
TruSeq Targeted RNA Hedgehog Panel Kit	48	RT-201-1002
	96	RT-202-1002
TruSeq Targeted RNA Neurodegeneration Panel Kit	48	RT-201-1001
	96	RT-202-1001
TruSeq Targeted RNA NFκB Panel Kit	48	RT-201-1008
	96	RT-202-1008
TruSeq Targeted RNA Stem Cell Panel Kit	48	RT-201-1005
	96	RT-202-1005
TruSeq Targeted RNA TP53 Pathway Panel Kit	48	RT-201-1007
	96	RT-202-1007
TruSeq Targeted RNA Wnt Pathway Panel Kit	48	RT-201-1004
	96	RT-202-1004
<b>TruSeq Targeted RNA Expression Index Kits</b>		
TruSeq Targeted RNA Index Kit	48	RT-401-1001
TruSeq Targeted RNA Index Kit A	96	RT-402-1001
TruSeq Targeted RNA Index Kit B	96	RT-402-1002
TruSeq Targeted RNA Index Kit C	96	RT-402-1003
TruSeq Targeted RNA Index Kit D	96	RT-402-1004

illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

FOR RESEARCH USE ONLY

© 2012-2013 Illumina, Inc. All rights reserved.  
 Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAlx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
 Pub. No. 470-2012-004 Current as of 14 August 2013



AGAATGATAACAGTAAACACACTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCTTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC  
 TCAAGCTACCGTAACGAAAGGTATCATTAAGATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCTGCAACGCAAAAGAAATGATAACAGTAAACACTCTGTTAACCTT  
 CCAACGAAAGAAATGATAGAGCTAAACACACTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCTGCAACGCAAAAGAAATGATA  
 AACGTACCATTAAAGAGCTACCGCTGCAACGCAAAAGAAATGATAACAGTAAACACTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCT  
 AAGATTGATAACAGTAAACACTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCTGCAACGCAAAAGAAATGATA  
 GATTACTTGATCCACTGATTCAACGTAACCGTAACGAAACGTATCAATTGAGACTAAATATTAACGTACCATTAAAGAGCTACCGCTGCAACGCAAAAGAAATGATAACAGTAAACACTCTGTTAACCTTAAGATTACTTGATCCACTGATTCAAC



## RNA-Seq of Blood

### ScriptSeq™ Complete Gold (Blood)


- ▶ Removes globin mRNA and ribosomal RNA
- ▶ Creates an Illumina® sequencing library
- ▶ The data contains high amounts of coding and non-coding information
  
- ▶ Find more genes
- ▶ Find more coding and non-coding RNA's
- ▶ Good for small samples
- ▶ All phases of research



Blood is an important sample for research into 6,000 rare diseases and 12,000 disease groups. The data from samples treated with ScriptSeq Complete Gold (Blood) is focused on valuable RNA. Finding new genes, splice variants and isoforms is important to disease and health research.

ScriptSeq Complete (Blood) offers the most informative sequencing results by removing unwanted globin mRNA and ribosomal RNA prior to sequencing.

Figure 1. RNA-Seq libraries contain coding and non-coding RNA.



#### Find more coding and non-coding RNA


ScriptSeq Complete Gold (Blood) libraries were prepared from 5 µg of RNA isolated from human whole blood and sequenced on an Illumina® sequencer. Greater than 98% of all reads contain useful information.

RNA-Seq data is very useful to study disease (or health). Figure 1 shows an example in which 47% of the sequencing reads contain coding RNA and 52% contain non-coding RNA.

#### Workflow



Figure 2. Gene coverage from large (5 µg) or small (100 ng) of RNA.



#### Available for all sample sizes

ScriptSeq Complete Gold (Blood) is available for 100 ng + of total RNA. Results from small amounts of total RNA are very similar to results from high amounts of total RNA. Figure 2 shows coverage of the COX5B gene when either a small amount (100 ng) of total RNA or large amount (5 µg) of total RNA was treated with ScriptSeq Complete Gold (Blood).

#### Strong gene coverage

In figure 2, the height of the blue bars show how many reads align to that sequence. Taller bars show more reads and deeper (better) coverage. Coding (thick blue bars) regions in both the small and large input ranges is similar.

#### Success begins with purification

##### MasterPure RNA purification kit

Purification is an important step to prepare your sample. MasterPure safely removes unwanted material to give you pure, intact total RNA.

MasterPure offers unique benefits:

- ▶ Keep RNA intact (does not degrade RNA)
- ▶ Retain RNA diversity (including small RNA)
- ▶ Maximize genes discovered
- ▶ Available for all sample sizes

Cat. #	Quantity
<b>MasterPure™ RNA Purification Kit (for isolating RNA only)</b>	
MCR85102	100 Purifications

Cat. #	Quantity
<b>ScriptSeq™ Complete Gold Kit (Blood)—Low Input</b>	
SCL24GBL	24 Reactions
SCL6GBL	6 Reactions
For 100 ng – 1 µg total blood RNA.	
<b>ScriptSeq™ Complete Gold Kit (Blood)</b>	
BGGB1306	6 Reactions
BGGB1324	24 Reactions
For 1 µg – 5 µg total blood RNA.	
<b>FailSafe™ PCR Enzyme Mix</b>	
FSE51100	100 Units
Patents: <a href="http://www.illumina.com/patents">www.illumina.com/patents</a>	



## RNA-Seq of Yeast

### ScriptSeq™ Complete Gold Kit (Yeast)


- ▶ Removes ribosomal RNA with Ribo-Zero™
- ▶ Creates an Illumina® sequencing library with ScriptSeq v2
- ▶ Results contain coding and non-coding RNA
  
- ▶ One day method
- ▶ Find more genes
- ▶ Good for small samples



The yeast transcriptome is more complex than previously thought. RNA-Seq of yeast is a valuable approach for mapping the transcriptome and characterizing novel and

low abundance transcripts. The ScriptSeq Complete Gold Kit (Yeast) offers the most informative sequencing results by removing unwanted ribosomal RNA prior to sequencing.

Figure 1. RNA-Seq libraries contain coding and non-coding RNA.



Library composition of ScriptSeq™ Complete Gold Kit (Yeast) samples. ScriptSeq libraries were constructed from 1 µg of *S. cerevisiae* total RNA samples and sequenced on an Illumina® MiSeq™.

### Find more coding and non-coding RNA

#### ScriptSeq Complete Gold Kit (Yeast):


- ▶ Find more coding RNA
- ▶ Removes ribosomal RNA
- ▶ Creates Illumina® sequencing libraries
- ▶ Data contains high amounts of coding information

RNA-Seq data is very useful to study yeast gene expression. Figure 1 shows an example in which 95.6% of the sequencing reads contain coding RNA and 4.4% contain non-coding RNA.

### Workflow



Figure 2. Enhanced coverage with ScriptSeq™ Complete Gold Kit (Yeast).



**Enhanced coverage with ScriptSeq Complete Gold Kit (Yeast)**

ScriptSeq Complete Gold (Yeast) contains Ribo-Zero Gold (Yeast) for depletion of yeast rRNA. Gene coverage of the TEF2 gene (Figure 2) shows that rRNA depletion reveals more reads. In the figure, the height of the blue bars shows how many reads align to that sequence. Taller bars show more reads and deeper (better) coverage.

**ScriptSeq Yeast is a powerful tool to study:**

- ▶ Transcriptome mapping
- ▶ Gene structure
- ▶ Characterization of novel and low abundance transcripts

**Success begins with purification**

**MasterPure™ RNA Purification Kit**

Purification is the first critical step to prepare samples for sequencing. MasterPure produces sequencer-ready RNA safely and easily.

**MasterPure offers unique benefits:**

- ▶ Keep RNA intact (does not degrade RNA)
- ▶ Retain RNA diversity (including small RNA)
- ▶ Maximize genes discovered
- ▶ Available for all sample sizes

Cat. #	Quantity
<b>MasterPure™ RNA Purification Kit (for isolating RNA only)</b>	
MCR85102	100 Purifications

Cat. #	Quantity
--------	----------

**Ribo-Zero™ Magnetic Gold Kit (Yeast)**

Suitable for 1-5 µg of total RNA.

MRZY1306	6 Reactions
MRZY1324	24 Reactions

**ScriptSeq™ Complete Gold Kit (Yeast)**

Includes Ribo-Zero Gold (Yeast). Suitable for 1-5 µg of total RNA.

BGY1306	6 Reactions
BGY1324	24 Reactions

**ScriptSeq™ Complete Gold Kit (Yeast)- Low Input**

Includes Ribo-Zero Gold (Yeast). Suitable for 100 ng - 1 µg of total RNA.

SCGL6Y	6 Reactions
SCGL6Y	24 Reactions

**FailSafe™ PCR Enzyme Mix**

FSE51100	100 Units
----------	-----------

The FailSafe PCR Enzyme Mix is required for ScriptSeq Complete Gold (Yeast) RNA-Seq library preparation.





## Ribosome Profiling

### ARTseq™ Ribosome Profiling Kits

ARTseq (Active mRNA Translation) Ribosome profiling is a powerful technique to study translation.

- ▶ Sequence ribosome protected mRNA
- ▶ Rapid, scalable spin-column method
- ▶ No ultracentrifuge required!
- ▶ Compatible with yeast and mammalian samples
  
- ▶ Predict protein abundance
- ▶ Investigate translational control
- ▶ Measure gene expression



#### Sequencing actively translated transcripts

Sequence mRNA fragments undergoing translation by ribosomes. These mRNA fragments are called “footprinted” or ribosome protected mRNA fragments.

ARTseq sequences ribosome-protected mRNA fragments to provide a “snapshot” of the active ribosomes in a cell. You can identify proteins being actively translated from samples prepared with ARTseq. Samples collected at different times often show changes in translation. Samples treated with different drugs often show different translation patterns.

Figure 1. Identify actively translated RNA's with ARTseq.



#### Sequence only the protein coding regions

Samples prepared with ARTseq are enriched for ORF and devoid of UTR sequences. The start and stop codons are easily seen. Sequences are focused on protein coding regions.

#### Workflow









## Data Sheet: Epigenetics

### Ordering Information

Catalog No.	Product	Description
WG-314-1003	Infinium HumanMethylation450 BeadChip Kit (24 samples)	Each package contains two BeadChips and reagents for analyzing DNA methylation in 24 human DNA samples.
WG-314-1001	Infinium HumanMethylation450 BeadChip Kit (48 samples)	Each package contains four BeadChips and reagents for analyzing DNA methylation in 48 human DNA samples.
WG-314-1002	Infinium HumanMethylation450 BeadChip Kit (96 samples)	Each package contains eight BeadChips and reagents for analyzing DNA methylation in 96 human DNA samples.

Each HumanMethylation450 BeadChip can process 12 samples in parallel and assay >450,000 methylation sites per sample.

### Summary

The HumanMethylation450 BeadChip's unique combination of comprehensive, expert-selected coverage, high sample throughput capacity, and affordable price makes it an ideal solution for large sample-size, genome-wide DNA methylation studies.

### References

- Portela A, Esteller M (2010) Epigenetic modifications and human disease. *Nat Biotechnology* 28: 1057–1068.
- [http://www.illumina.com/products/humanht\\_12\\_expression\\_beadchip\\_kits\\_v4.ilmn](http://www.illumina.com/products/humanht_12_expression_beadchip_kits_v4.ilmn)
- Infinium HD FFPE DNA Restoration Protocol
- [http://www.illumina.com/products/infinium\\_ffpe\\_dna\\_restoration\\_solution.ilmn](http://www.illumina.com/products/infinium_ffpe_dna_restoration_solution.ilmn)
- Infinium HD FFPE Methylation Assay, Manual Protocol
- Infinium HD FFPE Methylation Assay, Automated Protocol
- Illumina FFPE QC Assay Protocol

Illumina • +1.800.809.4566 toll-free • 1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

#### FOR RESEARCH USE ONLY

© 2012 Illumina, Inc. All rights reserved.  
Illumina, IlluminaDx, BaseSpace, BeadArray, BeadExpress, cBot, CSPro, DASL, DesignStudio, Eco, GALix, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
Pub. No. 270-2010-001 Current as of 09 March 2012

AGAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAACGAACGATCAATTAAGACTAAATATAACGTACCAATTAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
TGAAGGTACCGTAAAGAGCTATGTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
CGAGGAAAGAAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
AGGTACCAATTAAGAGCTATGTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
AGAAATGATAACAGTAACACACTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
GATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC  
CGTATCAATTAAGAGCTATGTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAGGTACCGTAAAGAGCTACCGTCTTCTGTTAACTTAAAGATTACTTGATCCACTGATTCAAC

FailSafe™ PCR System



## PCR Optimization

### FailSafe™ PCR System


- ▶ PCR of difficult or high GC templates
- ▶ PCR Amplifications up to 20 kb
- ▶ Works the first time, every time
- ▶ 3-fold lower error rate than *Taq* DNA Polymerase



The FailSafe PCR System uses the patented Epicentre PCR enhancement technology to allow PCR reactions to work the first time and every time. Twelve buffer options run at the same time, allowing quick and easy optimization of your PCR.

FailSafe has enabled many difficult samples to be used successfully in PCR and published. FailSafe will ensure your PCR is successful.

Figure 1. Ensure successful PCR with FailSafe.




### Works the first time, every time

Optimizing PCR is easy with FailSafe. Create a master mix of FailSafe Enzyme blend, template DNA and primers. Add the FailSafe PCR PreMix Selection Kit buffers to test which buffer is optimal for your reaction.

### Workflow



Figure 2. The FailSafe™ PCR System will work with nearly all DNA – animal, bacterial, plant or viral.



FailSafe PCR will amplify DNA from a range of different sequences and sequence sizes. PCR products shown are up to 20 kb for lambda DNA, up to 21.5 kb for human DNA, and up to 18 kb for *E. coli* DNA.

#### Amplify any sample

FailSafe will amplify DNA sequences from almost any source, up to 20 kb in length in a single round.

FailSafe PCR products can be used in many applications, including TA cloning and blunt-end cloning.

FailSafe components (Polymerase and Premixes) are available in the FailSafe Premix Choice kits and as single products.

#### Success begins with MasterPure™ DNA Purification Kit

MasterPure produces DNA safely and easily to be used in many applications.

MasterPure offers unique benefits:

- ▶ Keep DNA intact (does not degrade DNA)
- ▶ Scalable reaction sizes
- ▶ Available for multiple sample types

Cat. #	Quantity
<b>MasterPure™ Complete DNA and RNA Purification Kit</b>	
MC85200	200 Purifications
MC89010	10 Purifications

#### FailSafe PCR PreMix Selection Kit

FS99060 (Contains all 12 Premixes and FailSafe PCR Polymerase) – sufficient reagent for 48 reactions (four full template and primer optimizations)

#### FailSafe PCR System

FS99100 (100 units of FailSafe Polymerase and 1 PreMix of choice)  
 FS99250 (250 units of FailSafe Polymerase and two Premixes of choice)  
 FS9901K (1000 U of FailSafe Polymerase and eight Premixes of choice)

#### FailSafe PCR Polymerase

FSE51100 (100 U)

#### FailSafe PCR Polymerase

FSE5101K (1000 U)

#### FailSafe PCR Premixes

FSP995A-L (A through L), 2.5 ml (100 reactions)

FailSafe™ PCR Premix Selection Kits: Purchase of this product includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research. No other patent rights (such as 5' Nuclease Process patent rights) are conveyed expressly, by implication, or by estoppel. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 944.



## What Will You Create Today?

### ***Innovative Enzyme Solutions***

- ▶ Stringent QC
- ▶ No affinity tags
- ▶ ISO 13485 compliant by end of 2013

#### **Example applications:**

- ▶ Reverse transcriptases
- ▶ RNA polymerases
- ▶ RNase-free DNase
- ▶ Much more...

Epicentre develops and manufactures the highest purity enzymes for life science research. Epicentre specializes in unique and bulk enzyme projects to meet your specific needs. Since opening in 1987, Epicentre has a proven track record in manufacturing high purity enzymes for life sciences in our state-of-the-art facility in Madison, Wisconsin, USA.

Standard enzymes for molecular biology research are available. Unique, hard-to-find enzymes are available to your specifications.

#### **OEM opportunities available**

- ▶ Custom manufacturing for alternate size requirements or specific concentrations
- ▶ In-house technical expertise with over 25 years of experience
- ▶ Competitive pricing to meet your budget constraints
- ▶ Flexible, quick turnaround time for OEM needs

#### **Bulk availability**

- ▶ Standard or custom offerings
- ▶ Flexible (custom) concentrations and package sizes
- ▶ Bulk capabilities





## Browse the possibilities...

### DNA Polymerases

Klenow DNA Polymerase  
Exo-Minus Klenow DNA Polymerase (D355A, E357A)  
RepliPhi™ Phi29 DNA Polymerase  
Terminal deoxynucleotidyl Transferase, Recombinant  
T4 DNA Polymerase

### RNA Exonucleases

RNase R  
Terminator™ 5'-Phosphate-  
Dependent Exonuclease

### DNA Ligases

T4 DNA Ligase, Cloned  
Ampligase® Thermostable DNA Ligase  
CircLigase™ ssDNA Ligase  
CircLigase™ II ssDNA Ligase  
*E. coli* DNA Ligase

### RNA Polymerases

T7 RNA Polymerase  
T7 R&DNA™ Polymerase

### RNA Ligases

T4 RNA Ligase  
T4 RNA Ligase 2, Deletion Mutant  
Thermostable RNA Ligase

### Phosphatases/Kinases

APex™ Heat-Labile Alkaline Phosphatase  
Tobacco Acid Pyrophosphatase (TAP)  
RNA 5' Polyphosphatase  
T4 Polynucleotide Kinase, Cloned

### DNA Endonucleases


Baseline-ZERO™ DNase  
Endonuclease IV, *E. coli*  
T4 Endonuclease V  
Lambda Terminase  
RNase-Free DNase I  
Pvu Rts11 Endonuclease

### RNA Endonucleases

RNase A  
RNase I, *E. coli*  
RNase III, *E. coli*  
RNase H, *E. coli*  
Hybridase™ Thermostable RNase H  
RNase T1, *Aspergillus oryzae*  
RiboShredder™ RNase Blend  
RNase A

### DNA Exonucleases

Exonuclease I, *E. coli*  
Exonuclease III, *E. coli*  
Exonuclease VII  
Plasmid-Safe™ ATP-Dependent  
DNase  
Lambda Exonuclease  
RecBCD Nuclease, *E. coli*  
Rec J Exonuclease  
T5 Exonuclease



From the power of the HiSeq X to the speed of MiSeq, Illumina has the sequencer that's just right for you.



**MiSeq**  
 Focused power. Speed and simplicity for targeted and small genome sequencing.



**NextSeq 500**  
 Flexible power. Speed and simplicity for everyday genomics.



**HiSeq 2500**  
 Production power. Power and efficiency for large-scale genomics.



**HiSeq X\***  
 Population power. \$1,000 human genome and extreme throughput for population-scale sequencing.

Key applications	Small genome, amplicon, and targeted gene panel sequencing.	Everyday genome, exome, transcriptome sequencing, and more.		Production-scale genome, exome, transcriptome sequencing, and more.		Population-scale human whole-genome sequencing.
Run mode	N/A	Mid-Output	High-Output	Rapid Run	High-Output	N/A
Flow cells processed per run	1	1	1	1 or 2	1 or 2	1 or 2
Output range	0.3-15 Gb	20-39 Gb	30-120 Gb	10-180 Gb	50-1000 Gb	1.6-1.8 Tb
Run time	5-65 hours	15-26 hours	12-30 hours	7-40 hours	< 1 day - 6 days	< 3 days
Reads per flow cell†	25 Million‡	130 Million	400 Million	300 Million	2 Billion	3 Billion
Maximum read length	2 x 300 bp	2 x 150 bp	2 x 150 bp	2 x 150 bp	2 x 125 bp	2 x 150 bp

\* Specifications shown for an individual HiSeq X System. HiSeq X is only available as part of the HiSeq X Ten.  
 † Clusters passing filter.  
 ‡ For MiSeq V3 Kits only.

Use our Sequencing Platform Comparison Tool to find the right sequencing system for your needs.

© 2014 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NextSeq, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. APex, ARTseq, EpiGnome, FailSafe, MasterPure, Ribo-Zero, ScriptSeq, and TotalScript are trademarks or registered trademarks of Epicentre (an Illumina company). All other brands and names contained herein are the property of their respective owners.

Illumina • 1.800.809.4566 toll-free (U.S.) • +1.858.202.4566 tel • techsupport@illumina.com • www.illumina.com

**FOR RESEARCH USE ONLY**

© 2014 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.  
Pub No. 073-2014-001 Current as of 29 May 2014

